



Amino Acids in Peas (*Pisum sativum*) by HPLC-HRAM-MS





Method Details

CONDITIONS

Column:	Avantor® ACE® 3 AQ						
Particle Size:	3 µm						
Dimensions:	150 x 3.0 mm						
Mobile Phase:	A: 0.1% formic acid in H ₂ O B: 0.1% formic acid in MeCN						
Gradient:	<table border="1"> <thead> <tr> <th>Time (mins)</th><th>% B</th></tr> </thead> <tbody> <tr> <td>0</td><td>0</td></tr> <tr> <td>10</td><td>100</td></tr> </tbody> </table>	Time (mins)	% B	0	0	10	100
Time (mins)	% B						
0	0						
10	100						
Flow Rate:	0.4 mL/min						
Injection:	5 µL						
Temperature:	30 °C						
Detection:	Exactive Orbitrap high resolution MS ESI positive ion mode Capillary temperature: 350 °C						

ORDERING TABLE

Product	Details	Size	Part Number
Avantor® ACE® 3 AQ	HPLC Column	150 x 3.0 mm	ACE-116-1503

Reproduced with permission of Fera Science Ltd, York, UK

Avantor® ACE®