

An Introduction to the Theory and Practice of LC Method Translations and LC Method Transfers

Gemma Lo



Method Translations: What Are You Trying to Achieve?

- ◆ **Faster separations with the same performance?**
 - ◆ eg HPLC → UHPLC methods for increased productivity / sample throughput with similar efficiency, selectivity and resolution?
- ◆ **Converting UHPLC → HPLC methods?**
 - ◆ eg for offshoring or third party labs or manufacturing?
- ◆ **Porous particle → solid core particle method change?**
 - ◆ Take advantage of solid core efficiency, speed or low backpressure?
- ◆ **Faster method development?**
 - ◆ eg reduce overall cycle time?
- ◆ **Higher resolution / peak capacity?**
 - ◆ eg for related substances or complex samples?
- ◆ **Same methods...transferred between instruments?**



Isocratic Translations

Isocratic Method Translations: General Principles

- Maintain a **constant length to particle size ratio, L/d_p** (for the same phase type and phase vendor)

$$N = \frac{L}{HETP}$$

- Will give ~ **similar** performance i.e. efficiency (selectivity, resolution)*
- Thus, **300 x 3.9 mm, 10 μm \approx 150 x 4.6 mm, 5 μm = 30,000**

| | Column Length (mm) | | | | | |
|-----|--------------------|---------------|--------|---------------|---------|---------------|
| | 50 | 75 | 100 | 150 | 250 | 300 |
| 1.7 | 29,412 | 44,118 | 58,824 | | | |
| 1.8 | 27,778 | 41,667 | 55,556 | | | |
| 1.9 | 26,316 | 39,474 | 52,632 | | | |
| 2 | 25,000 | 37,500 | 50,000 | 75,000 | | |
| 2.5 | 20,000 | 30,000 | 40,000 | 60,000 | 100,000 | |
| 2.6 | 19,231 | 28,846 | 38,462 | 57,692 | 96,154 | |
| 2.7 | 18,519 | 27,778 | 37,037 | 55,556 | 92,593 | |
| 3 | 16,667 | 25,000 | 33,333 | 50,000 | 83,333 | |
| 5 | 10,000 | 15,000 | 20,000 | 30,000 | 50,000 | |
| 10 | 5,000 | 7,500 | 10,000 | 15,000 | 25,000 | 30,000 |

 Columns meeting L/dp of 30,000

Isocratic Method Translations: General Principles

- Once the new column format has been selected, isocratic translations are **fairly straightforward**.

- **Step 1: Scale injection volume (V_i) to new column dead volume (V_M)**

- **Step 2: Scale flow rate (F)**

$$V_{i2} = \frac{V_{i1} \times V_{M2}}{V_{M1}}$$

- **If no change in d_p** $F_2 = \frac{F_1 \times d_{c2}^2}{d_{c1}^2}$

- **If d_p changes can use** $F_2 = \frac{F_1 \times d_{c2}^2 \times d_{p1}}{d_{c1}^2 \times d_{p2}}$

- **Step 3: determine new run time**


$$t_2 = \frac{t_1 F_1 V_{M2}}{F_2 V_{M1}}$$

- **Backpressure of new method can be estimated**

$$P_2 = \frac{P_1 \times F_2 \times L_2 \times d_{c1}^2 \times d_{p1}^2}{F_1 \times L_1 \times d_{c2}^2 \times d_{p2}^2}$$



Isocratic Method Translations: HPLC to UHPLC



| | HPLC Method | UHPLC Method |
|----------------|-------------------------|--------------------------|
| Column | 150 x 4.6 mm, 5 μ m | 50 x 3.0 mm, 1.7 μ m |
| V_M | 1.571 mL | 0.223 mL |
| Injection Vol. | 5 μ L | |
| Flow Rate | 1.0 mL/min | |
| Run Time | 35 mins. | |

➤ **ACE 5 C18-PFP 150 x 4.6 mm**

$V_M = 1.571$ mL

➤ **$L/d_p = 30,000$**

➤ **ACE 1.7 C18-PFP 50 x 3.0 mm**

$V_M = 0.223$ mL

➤ **$L/d_p = 29,412$**

➤ **V_M experimentally determined**


➤ **Or use**

$$V_M = \pi \left(\frac{d_c}{2} \right)^2 L \varepsilon$$

➤ **ε = column porosity = ~0.63 for fully porous, ~0.55 for solid core**



Isocratic Method Translations: HPLC to UHPLC



| | HPLC Method | UHPLC Method |
|----------------|-------------------------|--------------------------|
| Column | 150 x 4.6 mm, 5 μ m | 50 x 3.0 mm, 1.7 μ m |
| V_M | 1.571 mL | 0.223 mL |
| Injection Vol. | 5 μ L | 0.7 μ L |
| Flow Rate | 1.0 mL/min | 1.25 mL/min |
| Run Time | 35 mins. | |

➤ **Step 1: Scale injection volume**


$$V_{i2} = \frac{V_{i1} \times V_{M2}}{V_{M1}} = \frac{5 \times 0.223}{1.571} = 0.7 \mu\text{L}$$

➤ **Step 2: Scale flow rate**

$$F_2 = \frac{F_1 \times d_{c2}^2 \times d_{p1}}{d_{c1}^2 \times d_{p2}} = \frac{1.00 \times 3.0^2 \times 5}{4.6^2 \times 1.7} = 1.25 \text{ mL/min}$$



Isocratic Method Translations: HPLC to UHPLC



| | HPLC Method | UHPLC Method |
|----------------|-------------------------|--------------------------|
| Column | 150 x 4.6 mm, 5 μ m | 50 x 3.0 mm, 1.7 μ m |
| V_M | 1.571 mL | 0.223 mL |
| Injection Vol. | 5 μ L | 0.7 μ L |
| Flow Rate | 1.0 mL/min | 1.25 mL/min |
| Run Time | 35 mins. | 4 mins |


➤ **Step 3: Calculate UHPLC run time**

$$t_2 = \frac{t_1 F_1 V_{M2}}{F_2 V_{M1}} = \frac{35 \times 1.00 \times 0.223}{1.25 \times 1.571} = 4 \text{ minutes}$$


➤ **We can now estimate the backpressure of the translated method**

$$P_2 = \frac{P_1 \times F_2 \times L_2 \times d_{c1}^2 \times d_{p1}^2}{F_1 \times L_1 \times d_{c2}^2 \times d_{p2}^2} = \frac{66 \times 1.25 \times 50 \times 4.6^2 \times 5^2}{1.00 \times 150 \times 3^2 \times 1.7^2} = 559 \text{ bar}$$

Isocratic Method Translations: HPLC to UHPLC




| | HPLC Method | UHPLC Method |
|----------------|-------------------------|--------------------------|
| Column | 150 x 4.6 mm, 5 μ m | 50 x 3.0 mm, 1.7 μ m |
| V _M | 1.571 mL | 0.223 mL |
| Injection Vol. | 5 μ L | 0.7 μ L |
| Flow Rate | 1.0 mL/min | 1.25 mL/min |
| Run Time | 35 mins. | 4 mins |




| | Original method | Translated method | Difference |
|---------------------|-----------------|-------------------|---------------|
| Run time | 35 minutes | 4 minutes | -88.6% |
| Solvent consumption | 35.0 mL | 5.0 mL | -85.7% |



Isocratic Method Translations: HPLC to UHPLC



| | HPLC Method | UHPLC Method |
|----------------|-------------------------|--------------------------|
| Column | 150 x 4.6 mm, 5 μ m | 50 x 3.0 mm, 1.7 μ m |
| V_M | 1.571 mL | 0.223 mL |
| Injection Vol. | 5 μ L | 0.7 μ L |
| Flow Rate | 1.0 mL/min | 1.25 mL/min |
| Run Time | 35 mins. | 4 mins |




- **Note: the same translation principles are also applicable to UHPLC to HPLC method translations**

ACE LC Translator

- ◆ All of the **translation / transfer** discussed in this talk have been conveniently combined in the **ACE LC Translator**

V1.3

ACE LC Translator



This tool has been developed to aid the practicing chromatographer to efficiently translate and transfer LC methods between different format columns and different LC systems. Also included are a set of useful tools for the calculation of everyday chromatographic and method parameters.

The various tools can be accessed via the links below.

Tools

| | |
|--|---|
| Method Translation Method Transfer Dwell Volume Calculator | Column Porosity Calculator Extra Column Volume Calculator Column Equilibration Calculator Buffer Calculator Mobile Phase Consumption Calculator |
|--|---|

Free Literature

| | |
|---------------------------------------|--|
| ACE Application Notes | ACE Product Portfolio Contact Us / Feedback |
|---------------------------------------|--|

[*Register for Future Updates* - register to receive future updates to the ACE LC Translator](#)

© 2017 ACT Ltd. All rights reserved. All products are available worldwide.
 For more information contact: info@ace-hplc.com or www.ace-hplc.com

[Disclaimer](#)


ACE LC Translator

- ◆ All of the **translation / transfer** discussed in this talk have been conveniently combined in the **ACE LC Translator**

- ◆ Download **free** at:
mac-mod.com

V1.3

Method Translation


 Methods are frequently translated from one column format to another and among different brands and models of instruments (different system volumes, dwell volumes etc). In order to maintain chromatographic and method performance through the translation process, a number of method parameter changes are required, such as flow rate, injection volume and t_R . Separate tools for isocratic and gradient methods are included on this page. For details on how to determine system dwell volumes, please refer to the Dwell Volume tab.

Please complete all input fields

Isocratic

| Column Information | | Translated | |
|-------------------------------------|---|-------------------------------------|---|
| Current | | Translated | |
| Column Length (L) | <input type="text" value="150"/> mm | Column Length (L) | <input type="text" value="50"/> mm |
| Column i.d. (d _c) | <input type="text" value="4.6"/> mm | Column i.d. (d _c) | <input type="text" value="3.0"/> mm |
| Particle Diameter (d _p) | <input type="text" value="5.0"/> μm | Particle Diameter (d _p) | <input type="text" value="1.7"/> μm |
| L/d _p | <input type="text" value="30000"/> | L/d _p | <input type="text" value="29412"/> |
| Column Porosity | <input type="text" value="0.63"/> <small>mL/min</small> | Column Porosity | <input type="text" value="0.63"/> <small>mL/min</small> |
| Column Volume (V _H) | <input type="text" value="1570"/> mL | Column Volume (V _H) | <input type="text" value="0.223"/> mL |

| Method | | Translated | |
|-----------------------|--|---------------------------------------|--|
| Current | | Translated | |
| Injection Volume | <input type="text" value="10.0"/> μL | Injection Volume | <input type="text" value="1.4"/> μL |
| Flow Rate | <input type="text" value="1.00"/> mL/min | Flow Rate (scaled to linear velocity) | <input type="text" value="0.43"/> mL/min |
| | | Flow Rate (scaled to particle size) | <input type="text" value="1.25"/> mL/min |
| | | Input Flow Rate | <input type="text" value="1.25"/> mL/min |
| Run Time | <input type="text" value="25.0"/> mins | New Run Time | <input type="text" value="2.8"/> mins |
| Recorded Backpressure | <input type="text" value="60"/> bar | Estimated Backpressure | <input type="text" value="508"/> bar |
| Solvent Use | <input type="text" value="25"/> mL | Estimated Solvent Use Difference | <input type="text" value="-86"/> % |

Main Menu

Gradient

| Column Information | | Translated | |
|-------------------------------------|--|-------------------------------------|--|
| Current | | Translated | |
| Column Length (L) | <input type="text"/> mm | Column Length (L) | <input type="text"/> mm |
| Column i.d. (d _c) | <input type="text"/> mm | Column i.d. (d _c) | <input type="text"/> mm |
| Particle Diameter (d _p) | <input type="text"/> μm | Particle Diameter (d _p) | <input type="text"/> μm |
| L/d _p | <input type="text"/> | L/d _p | <input type="text"/> |
| Column Porosity | <input type="text"/> <small>mL/min</small> | Column Porosity | <input type="text"/> <small>mL/min</small> |
| Column Volume (V _H) | <input type="text"/> mL | Column Volume (V _H) | <input type="text"/> mL |

| Method | | Translated | |
|--------------------------------|-----------------------------|---------------------------------------|-----------------------------|
| Current | | Translated | |
| Injection Volume | <input type="text"/> μL | Injection Volume | <input type="text"/> μL |
| Flow Rate | <input type="text"/> mL/min | Flow Rate (scaled to linear velocity) | <input type="text"/> mL/min |
| | | Flow Rate (scaled to particle size) | <input type="text"/> mL/min |
| | | Input Flow Rate | <input type="text"/> mL/min |
| LC Name | <input type="text"/> | LC Name | <input type="text"/> |
| Dwell Volume (V _D) | <input type="text"/> mL | Dwell Volume (V _D) | <input type="text"/> mL |

Method Translation Tool Example - Isocratic

Fill in the grey input boxes to translate the method

Displays new run time and solvent saving

Optional input for flow rate – great for allowing greater analyst choice

V1.3

Method Translation

Methods are frequently translated from one column format to another and among different brands and models of instruments (different system volumes, dwell volumes etc). In order to maintain chromatographic and method performance through the translation process, a number of method parameter changes are required, such as flow rate, injection volume and t_R . Separate tools for isocratic and gradient methods are included on this page. For details on how to determine system dwell volumes, please refer to the Dwell Volume tab.

ACE

Isocratic

Column Information

| Current | Translated |
|-----------------------------|-----------------------------|
| Column Length (L) | Column Length (L) |
| Column i.d. (d_c) | Column i.d. (d_c) |
| Particle Diameter (d_p) | Particle Diameter (d_p) |
| L/ d_p | L/ d_p |
| Column Porosity | Column Porosity |
| Column Volume (V_M) | Column Volume (V_M) |

Method

| Current | Translated |
|-----------------------|---------------------------------------|
| Injection Volume | Injection Volume |
| Flow Rate | Flow Rate (scaled to linear velocity) |
| | Flow Rate (scaled to particle size) |
| | Input Flow Rate |
| Run Time | New Run Time |
| Recorded Backpressure | Estimated Backpressure |
| Solvent Use | Estimated Solvent Use Difference |

[Main Menu](#)

Gradient

Please complete all input fields

Please complete all input fields



Isocratic Method Translations: Pharmacopeial Methods


- Revised guidance on allowable parameter changes recently updated for **USP <621> chapter***
 - **Significant update** to the general chapter – reading recommended
 - Isocratic monograph LC methods can be **-25% to +50% of the original L/d_p** according to the USP guidance
 - Also, L/d_p ratios outside this range can be **-25% to +50% of the original peak efficiency, N** ← large flexibility
- **Very few** changes / flexibility for **gradient** monograph methods
- When translating to smaller columns dimensions (and therefore smaller column volume), **system dispersion effects are important**

* www.usp.org. New guidance in general chapter USP <621> issued August 2014.



Method Translations: Isocratic Estradiol Porous → Porous (I)

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm



1. Find suitable 'new' column: 300x3.9mm, 10 μ m = L/d_p of 30,000 = 150x4.6, 5 μ m

2. Geometrically scale flow rate: $F_2 = F_1 \times \frac{d_{c2}^2}{d_{c1}^2} = 1 \times 4.6^2 / 3.9^2 = \underline{1.39 \text{ mL/min}}$


3. Volumetrically scale injection: $Inj_2 = Inj_1 \times \frac{V_{m2}}{V_{m1}} = 25 \times 1.620 / 2.329 = \underline{17.4 \mu\text{L}}$

$$V_M \approx \pi \left(\frac{d}{2} \right)^2 L \varepsilon$$




Method Translations: Isocratic Estradiol Porous → Porous (II)

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm



First Translated Monograph Separation: ACE Excel 5 μ m, L1 (C18), 150 x 4.6 mm



Method Translations: Isocratic Methods

- **Isocratic monograph LC methods can be -25% to +50% of the original L/d_p according to the USP guidance***

| | |
|----------------------|--------|
| Column Length (mm) | 300 |
| Particle Size (um) | 10 |
| Lower L/d_p (-25%) | 22,500 |
| L/d_p | 30,000 |
| Upper L/d_p (+50%) | 45,000 |

| Particle Size (μm) | Column Length (mm) | | | | | |
|--------------------|--------------------|--------|--------|--------|---------|--------|
| | 50 | 75 | 100 | 150 | 250 | 300 |
| 1.7 | 29,412 | 44,118 | 58,824 | | | |
| 1.8 | 27,778 | 41,667 | 55,556 | | | |
| 1.9 | 26,316 | 39,474 | 52,632 | | | |
| 2 | 25,000 | 37,500 | 50,000 | 75,000 | | |
| 2.5 | 20,000 | 30,000 | 40,000 | 60,000 | 100,000 | |
| 2.6 | 19,231 | 28,846 | 38,462 | 57,692 | 96,154 | |
| 2.7 | 18,519 | 27,778 | 37,037 | 55,556 | 92,593 | |
| 3 | 16,667 | 25,000 | 33,333 | 50,000 | 83,333 | |
| 5 | 10,000 | 15,000 | 20,000 | 30,000 | 50,000 | |
| 10 | 5,000 | 7,500 | 10,000 | 15,000 | 25,000 | 30,000 |


 Columns meeting L/d_p of 22,400 to 45,000

* www.usp.org. New guidance in general chapter USP <621> issued August 2014.




Method Translations: Isocratic Estradiol Porous → Porous (II)

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm




Translated Monograph Separation: ACE Excel 2 μ m, L1 (C18), 50 x 3.0 mm






Method Translations: Isocratic Estradiol Porous → Solid Core

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm



Translated Monograph Separation: ACE UltraCore 2.5 μ m, L1 (SuperC18), 100 x 4.6 mm




| Requirement | Requirement met? |
|-------------|------------------|
| Rs >2.0 | ✓ |
| %RSD <2.0% | ✓ |

➤ Look out for negative effects of extra-column band broadening!

$$V_M \approx \pi \left(\frac{d}{2}\right)^2 L \varepsilon$$



Isocratic Method Translations: Aspirin Porous → Solid Core (I)



Translation to ACE UltraCore 5 μ m SuperC18, 150 x 4.6 mm:

1. Column dimensions and particle size similar so same flow rate used



2. Determine column volumes and scale injection volume $V_M \approx \pi \left(\frac{d}{2}\right)^2 L \varepsilon$

$$Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}}\right)$$

Conditions: 60:35:5:0.2 v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm, 25°C, 1 mL/min, 5 μ L injection




Isocratic Method Translations: Aspirin Porous → Solid Core (II)





- **Reduced hydrophobicity of solid core particles leads to ‘faster’ analysis**
- **Can we go faster...? (within allowable range)**


Conditions: (Top): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (2.5 Hz), 25°C, 1 mL/min, 5 μ L injection
(Bottom): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (20 Hz), 25°C, 1 mL/min, 3.9 μ L injection







Isocratic Method Translations: Aspirin Porous → Solid Core (III)



Conditions (Top): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm, 25°C, 1 mL/min, 5 μ L injection
 (Bottom): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (20 Hz), 25°C, 0.43 mL/min, 0.9 μ L injection



Isocratic Method Translations: Aspirin Porous → Solid Core (IV)



Conditions ACE UltraCore 2.5 Super C18 75 x 4.6 mm:
60:35:5:0.2 v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (20 Hz), 25°C, 1.0 mL/min, 2 µL injection



Gradient Translations



Translations of Gradient Methods – The Basics

1. Calculate column volumes

- Better to experimentally determine porosity for accuracy

$$V_M \approx \pi \left(\frac{d}{2} \right)^2 L \varepsilon$$

2. Translate injection volume

- To give similar response

$$Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}} \right)$$

Experimentally determined

3. Translate flow rate

- Constant linear velocity

$$F_2 = F_1 \times \frac{d_{c2}^2}{d_{c1}^2}$$

- Scaled to new d_p

$$F_2 = \frac{F_1 \times d_{c2}^2 \times d_{p1}}{d_{c1}^2 \times d_{p2}}$$

4. Translate gradient time


- To maintain constant k^*

$$t_{G2} = \frac{t_{G1} V_{M2} F_1}{V_{M1} F_2}$$


5. Consider V_D/V_M ratio

Gradient Method Translations: Porous → Porous


ACE Excel 5 C18-Amide 250 x 4.6 mm, 1 mL/min




| | Original Method |
|----------------|-------------------------|
| Column | 250 x 4.6 mm, 5 μ m |
| V_M | 2.668 mL |
| Injection Vol. | 10 μ L |
| Flow Rate | 1.0 mL/min |
| Gradient time | 45 min. |




1. Metronidazole




2. 4-acetamidophenol




3. amiloride




4. caffeine




5. Hydrochlorothiazide




6. 4-Hydroxybenzoic acid




7. Acetophenone




8. 2-nitrophenol




9. Nitrobenzene



10. 1,2-Dinitrobenzene



11. ethylbenzoate



12. 3,4-dichlorobenzoic acid


Gradient analysis, A: 20 mM KH₂PO₄ pH 2.7 (aq), B: 20 mM KH₂PO₄ pH 2.7 in MeCN:water 65:35 v/v, 5-100%B in 45.0 mins, hold 100%B for 5.0 mins, 40°C, 1.0 mL/min, 254 nm.

1. metronidazole, 2. 4-acetamidophenol, 3. amiloride, 4. caffeine, 5. hydrochlorothiazide, 6. 4-hydroxybenzoic acid, 7. acetophenone, 8. 2-nitrophenol, 9. nitrobenzene, 10. 1,2-dinitrobenzene, 11. ethylbenzoate, 12. 3,4-dichlorobenzoic acid.



Gradient Method Translations: Porous → Porous

ACE Excel 5 C18-Amide 250 x 4.6 mm, 1 mL/min



| | Original Method |
|----------------|-------------------------|
| Column | 250 x 4.6 mm, 5 μ m |
| V_M | 2.668 mL |
| Injection Vol. | 10 μ L |
| Flow Rate | 1.0 mL/min |
| Gradient time | 45 min. |

➤ ACE 5 C18-Amide 250 x 4.6 mm

$V_M = 2.668 \text{ mL}$

➤ $L/d_p = 50,000$


Gradient analysis, A: 20 mM KH₂PO₄ pH 2.7 (aq), B: 20 mM KH₂PO₄ pH 2.7 in MeCN:water 65:35 v/v, 5-100%B in 45.0 mins, hold 100%B for 5.0 mins, 40°C, 1.0 mL/min, 254 nm.

1. metronidazole, 2. 4-acetamidophenol, 3. amiloride, 4. caffeine, 5. hydrochlorothiazide, 6. 4-hydroxybenzoic acid, 7. acetophenone, 8. 2-nitrophenol, 9. nitrobenzene, 10. 1,2,-dinitrobenzene, 11. ethylbenzoate, 12. 3,4-dichlorobenzoic acid.



Gradient Method Translations: Porous → Porous

ACE Excel 5 C18-Amide 250 x 4.6 mm, 1 mL/min



➤ **ACE 5 C18-Amide 250 x 4.6 mm**

➤ **$L/d_p = 50,000$**

➤ **ACE 3 C18-Amide 150 x 4.6 mm**

➤ **$L/d_p = 50,000$**

➤ **1. V_M experimentally determined**

➤ **Or use**

| | Original Method | Translated Method |
|----------------|-------------------------------|---|
| Column | 250 x 4.6 mm, 5 μm | 150 x 4.6 mm, 3 μm |
| V_M | 2.668 mL | 1.571 mL |
| Injection Vol. | 10 μL | |
| Flow Rate | 1.0 mL/min | |
| Gradient time | 45 min. | |

$V_M = 2.668 \text{ mL}$

$V_M = 1.571 \text{ mL}$

$$V_M = \pi \left(\frac{d_c}{2} \right)^2 L \varepsilon$$


➤ **$\varepsilon = \text{column porosity} = \sim 0.63$ for fully porous, ~ 0.55 for solid core**

Gradient analysis, A: 20 mM KH₂PO₄ pH 2.7 (aq), B: 20 mM KH₂PO₄ pH 2.7 in MeCN:water 65:35 v/v, 5-100%B in 45.0 mins, hold 100%B for 5.0 mins, 40°C, 1.0 mL/min, 254 nm.
 1. metronidazole, 2. 4-acetamidophenol, 3. amiloride, 4. caffeine, 5. hydrochlorothiazide, 6. 4-hydroxybenzoic acid, 7. acetophenone, 8. 2-nitrophenol, 9. nitrobenzene, 10. 1,2,-dinitrobenzene, 11. ethylbenzoate, 12. 3,4-dichlorobenzoic acid.



Gradient Method Translations: Porous → Porous

ACE Excel 5 C18-Amide 250 x 4.6 mm, 1 mL/min



| | Original Method | Translated Method |
|----------------|-------------------------------|-------------------------------------|
| Column | 250 x 4.6 mm, 5 μm | 150 x 4.6 mm, 3 μm |
| V_M | 2.668 mL | 1.571 mL |
| Injection Vol. | 10 μL | 5.9 μL |
| Flow Rate | 1.0 mL/min | |
| Gradient time | 45 min. | |

2. Translate injection volume

$$V_{i2} = \frac{V_{i1} \times V_{M2}}{V_{M1}} = \frac{10 \times 1.571}{2.668} = 5.9 \mu\text{L}$$


Gradient analysis, A: 20 mM KH₂PO₄ pH 2.7 (aq), B: 20 mM KH₂PO₄ pH 2.7 in MeCN:water 65:35 v/v, 5-100%B in 45.0 mins, hold 100%B for 5.0 mins, 40°C, 1.0 mL/min, 254 nm.

1. metronidazole, 2. 4-acetamidophenol, 3. amiloride, 4. caffeine, 5. hydrochlorothiazide, 6. 4-hydroxybenzoic acid, 7. acetophenone, 8. 2-nitrophenol, 9. nitrobenzene, 10. 1,2-dinitrobenzene, 11. ethylbenzoate, 12. 3,4-dichlorobenzoic acid.



Gradient Method Translations: Porous → Porous

ACE Excel 5 C18-Amide 250 x 4.6 mm, 1 mL/min



| | Original Method | Translated Method |
|----------------|-------------------------|-------------------------|
| Column | 250 x 4.6 mm, 5 μ m | 150 x 4.6 mm, 3 μ m |
| V_M | 2.668 mL | 1.571 mL |
| Injection Vol. | 10 μ L | 5.9 μ L |
| Flow Rate | 1.0 mL/min | 1.67 mL/min |
| Gradient time | 45 min. | |

2. Translate injection volume

$$V_{i2} = \frac{V_{i1} \times V_{M2}}{V_{M1}} = \frac{10 \times 1.571}{2.668} = 5.9 \mu\text{L}$$

3. Translate flow rate

$$F_2 = \frac{F_1 \times d_{c2}^2 \times d_{p1}}{d_{c1}^2 \times d_{p2}} = \frac{1.00 \times 4.6^2 \times 5}{4.6^2 \times 3} = 1.67 \text{ mL/min}$$


Gradient analysis, A: 20 mM KH₂PO₄ pH 2.7 (aq), B: 20 mM KH₂PO₄ pH 2.7 in MeCN:water 65:35 v/v, 5-100%B in 45.0 mins, hold 100%B for 5.0 mins, 40°C, 1.0 mL/min, 254 nm.

1. metronidazole, 2. 4-acetamidophenol, 3. amiloride, 4. caffeine, 5. hydrochlorothiazide, 6. 4-hydroxybenzoic acid, 7. acetophenone, 8. 2-nitrophenol, 9. nitrobenzene, 10. 1,2,-dinitrobenzene, 11. ethylbenzoate, 12. 3,4-dichlorobenzoic acid.



Gradient Method Translations: Porous → Porous

ACE Excel 5 C18-Amide 250 x 4.6 mm, 1 mL/min



| | Original Method | Translated Method |
|----------------|--------------------|--------------------|
| Column | 250 x 4.6 mm, 5 μm | 150 x 4.6 mm, 3 μm |
| V_M | 2.668 mL | 1.571 mL |
| Injection Vol. | 10 μL | 5.9 μL |
| Flow Rate | 1.0 mL/min | 1.67 mL/min |
| Gradient time | 45 min. | 15.87 min. |

3. Translate gradient time



$$t_{G2} = \frac{t_{G1} V_{M2} F_1}{V_{M1} F_2} = \frac{45 \times 1.571 \times 1.0}{2.668 \times 1.67} = 15.87 \text{ min.}$$

Gradient analysis, A: 20 mM KH₂PO₄ pH 2.7 (aq), B: 20 mM KH₂PO₄ pH 2.7 in MeCN:water 65:35 v/v, 5-100%B in 45.0 mins, hold 100%B for 5.0 mins, 40°C, 1.0 mL/min, 254 nm.

1. metronidazole, 2. 4-acetamidophenol, 3. amiloride, 4. caffeine, 5. hydrochlorothiazide, 6. 4-hydroxybenzoic acid, 7. acetophenone, 8. 2-nitrophenol, 9. nitrobenzene, 10. 1,2-dinitrobenzene, 11. ethylbenzoate, 12. 3,4-dichlorobenzoic acid.



Gradient Method Translations: Porous → Porous

ACE Excel 5 C18-Amide 250 x 4.6 mm, 1 mL/min

ACE Excel 3 C18-Amide 150 x 4.6 mm, 1.67 mL/min


| | Original Method | Translated Method |
|----------------|-------------------------|-------------------------|
| Column | 250 x 4.6 mm, 5 μ m | 150 x 4.6 mm, 3 μ m |
| V_M | 2.668 mL | 1.571 mL |
| Injection Vol. | 10 μ L | 5.9 μ L |
| Flow Rate | 1.0 mL/min | 1.67 mL/min |
| Gradient time | 45 min. | 15.87 min. |

Total run time: -64.7%

Solvent consumption: -41%

Gradient analysis, A: 20 mM KH₂PO₄ pH 2.7 (aq), B: 20 mM KH₂PO₄ pH 2.7 in MeCN:water 65:35 v/v, 5-100%B in 45.0 mins, hold 100%B for 5.0 mins, 40°C, 1.0 mL/min, 254 nm.

1. metronidazole, 2. 4-acetamidophenol, 3. amiloride, 4. caffeine, 5. hydrochlorothiazide, 6. 4-hydroxybenzoic acid, 7. acetophenone, 8. 2-nitrophenol, 9. nitrobenzene, 10. 1,2,-dinitrobenzene, 11. ethylbenzoate, 12. 3,4-dichlorobenzoic acid.



Step 1. Calculate Column Volumes

- Due to column volume differences with porous and solid core columns, **errors can creep into gradient translation work**

| Porous (0.63) column volumes | | Column Length | | | | | | |
|----------------------------------|-------|---------------|-------|-------|-------|-------|-------|-------|
| | | 20mm | 30mm | 50mm | 75mm | 100mm | 150mm | 250mm |
| Col. i.d. | 1mm | 0.010 | 0.015 | 0.025 | 0.037 | 0.049 | 0.074 | 0.124 |
| | 2.1mm | 0.044 | 0.065 | 0.109 | 0.164 | 0.218 | 0.327 | 0.546 |
| | 3.0mm | 0.089 | 0.134 | 0.223 | 0.334 | 0.445 | 0.668 | 1.113 |
| | 4.6mm | 0.209 | 0.314 | 0.523 | 0.785 | 1.047 | 1.570 | 2.617 |
| Solid core (0.55) column volumes | | Column Length | | | | | | |
| | | 20mm | 30mm | 50mm | 75mm | 100mm | 150mm | 250mm |
| Col. i.d. | 1mm | 0.009 | 0.013 | 0.022 | 0.032 | 0.043 | 0.065 | 0.108 |
| | 2.1mm | 0.038 | 0.057 | 0.095 | 0.143 | 0.190 | 0.286 | 0.476 |
| | 3.0mm | 0.078 | 0.117 | 0.194 | 0.292 | 0.389 | 0.583 | 0.972 |
| | 4.6mm | 0.183 | 0.274 | 0.457 | 0.686 | 0.914 | 1.371 | 2.285 |

$$V_M \approx \pi \left(\frac{d}{2} \right)^2 L \varepsilon$$


Porosity values:
 porous = 0.63
 solid core = 0.55

- Solid core particle columns have **~13% LESS** column volume than porous particle columns
- Better to **experimentally determine** porosity for vendor column



Gradient Method Translations: Acetaminophen Porous → Solid Core

ACE Excel 5 μ m SuperC18 150 x 4.6 mm



| | |
|-----------|-----------|
| t_G | 20 min |
| Flow: | 1 mL/min |
| Inj. Vol. | 5 μ L |
| P_{MAX} | 138 bar |

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$

Translation to ACE UltraCore 5 μ m SuperC18 150 x 4.6 mm:

1. Calculate column volumes...solid core porosity differences
2. Translate gradient time $\frac{t_{G1}F_1}{V_{M1}} = \frac{t_{G2}F_2}{V_{M2}}$
3. No flow rate change needed $F_2 = F_1 \times \frac{d_{c2}^2}{d_{c1}^2}$
4. Translate injection volume $Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}}\right)$
5. Calculate whether an injection hold or pre-injection is needed

$$V_M \approx \pi \left(\frac{d}{2}\right)^2 L \epsilon$$

$$\Delta = \left(\frac{V_D}{V_M}\right)_{original} - \left(\frac{V_D}{V_M}\right)_{new}$$


150x4.6mm, 5 μ m, gradient analysis, A= 20mM ammonium acetate pH 6.0 (aq), B= 20mM ammonium acetate pH 6.0 in MeCN:water 9:1 v/v, 5-95%B in 20.0 mins, hold 95%B for 5.0 mins, 30°C, 1.0 mL/min, 230 nm.

1. 4-aminophenol 2. Hydroquinone 3. acetaminophen 4. 4-acetamidophenol 5. phenol 6. 4-nitrophenol 7. 2-nitrophenol 8. 4-chloroacetanilide



Gradient Method Translations: Acetaminophen Porous → Solid Core


ACE Excel 5 μ m SuperC18 150 x 4.6 mm



| | |
|-----------|-----------|
| t_G | 20 min |
| Flow: | 1 mL/min |
| Inj. Vol. | 5 μ L |
| P_{MAX} | 138 bar |

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$

ACE UltraCore 5 μ m SuperC18 150 x 4.6 mm



| | |
|-----------|-------------|
| t_G | 16.4 min |
| Flow: | 1 mL/min |
| Inj. Vol. | 3.9 μ L |
| P_{MAX} | 159 bar |

Corrected

- V_m (porosity)
- t_G adjusted for constant k^*
- Injection volume adjusted


Run time reduced whilst peak capacity and resolution are increased

$$\frac{t_{G1}F_1}{V_{M1}} = \frac{t_{G2}F_2}{V_{M2}}$$



Gradient Method Translations: Acetaminophen Porous → Solid Core


ACE Excel 5 μ m SuperC18 150 x 4.6 mm



| | |
|-----------|-----------|
| t_G | 20 min |
| Flow: | 1 mL/min |
| Inj. Vol. | 5 μ L |
| P_{MAX} | 138 bar |

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$

ACE UltraCore 2.5 μ m SuperC18 50 x 3.0 mm




| | |
|-----------|--------------|
| t_G | 5.92 min |
| Flow: | 0.43 mL/min |
| Inj. Vol. | 0.63 μ L |
| P_{MAX} | 162 bar |

Corrected
 - V_m (porosity)
 - Constant k^* (t_G & F)
 - Injection volume (V_m)


Significant reduction in run time whilst maintaining resolution

$$\frac{t_{G1}F_1}{V_{M1}} = \frac{t_{G2}F_2}{V_{M2}}$$



Gradient Method Translations: Acetaminophen Porous → Solid Core


ACE Excel 5µm SuperC18 150 x 4.6 mm



| | |
|-----------|----------|
| t_G | 20 min |
| Flow: | 1 mL/min |
| Inj. Vol. | 5 µL |
| P_{MAX} | 138 bar |

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$

ACE UltraCore 2.5µm SuperC18 50 x 3.0 mm



| | |
|-----------|-------------|
| t_G | 2.99 min |
| Flow: | 0.85 mL/min |
| Inj. Vol. | 0.63 µL |
| P_{MAX} | 315 bar |

Corrected
- Constant k^* (t_G & F)

Significant reduction in run time whilst still maintaining resolution

$$\frac{t_{G1}F_1}{V_{M1}} = \frac{t_{G2}F_2}{V_{M2}}$$



Gradient Method Translations: Correcting for V_D/V_M

- It is **important** when translating to **smaller i.d.** (therefore **small V_M**) UHPLC columns to consider the **V_D/V_M ratio**.

$$\Delta = \left(\frac{V_D}{V_M} \right)_{\text{original}} - \left(\frac{V_D}{V_M} \right)_{\text{new}} \text{ must approach zero}$$

- For **maximum accuracy** of gradient translations, it may be necessary to correct for the change in ratio of V_D/V_M between the **original** and **translated** method.
- The correction can be made using


$$x = \left[\left(\frac{V_{D1}}{V_{M1}} \right) - \left(\frac{V_{D2}}{V_{M2}} \right) \right] \times \frac{V_{M2}}{F_2}$$

- **Negative** value: **injection must be delayed** x minutes after gradient starts.
- **Positive** value: a **pre-gradient isocratic hold** of x minutes should be added to gradient program.



Gradient Method Translations: Why correct for V_D/V_M ?


ACE 5 C18-PFP 150 x 4.6 mm



$t_G = 15$ minutes
 $F = 1.0$ mL/min

ACE Excel 1.7 C18-PFP 50 x 3.0 mm -

Injection delayed by 0.62 mins.




$t_G = 1.70$ minutes
 $F = 1.25$ mL/min

- Calculations tell us to **delay injection until 0.62 minutes** after the gradient starts to correct for V_D/V_M

$$x = \left[\left(\frac{V_{D1}}{V_{M1}} \right) - \left(\frac{V_{D2}}{V_{M2}} \right) \right] \times \frac{V_{M2}}{F_2} = \left[\left(\frac{1.098}{1.570} \right) - \left(\frac{0.926}{0.223} \right) \right] \times \frac{0.223}{1.25} = \mathbf{-0.62 \text{ minutes}}$$


Gradient analysis, A= 20mM ammonium acetate pH 6.0 (aq), B= 20 mM ammonium acetate pH 6.0 in MeCN:water 80:20 v/v, 5-95%B in 15.0 mins, hold 95%B for 2.0 mins, 40°C, 1.0 mL/min, 230 nm.

1. 4-aminophenol, 2. hydroquinone, 3. 4-acetamidophenol (paracetamol), 4. 2-aminophenol, 5. 2-acetamidophenol, 6. phenol, 7. 4-nitrophenol, 8. 4-chloroacetanilide, 9. 2-nitrophenol. Impurities were spiked at 0.5% w/w.



Gradient Method Translations: Why correct for V_D/V_M ?


ACE 5 C18-PFP 150 x 4.6 mm



$t_G = 15$ minutes
 $F = 1.0$ mL/min

ACE Excel 1.7 C18-PFP 50 x 3.0 mm -


Injection delayed by 0.62 mins.



$t_G = 1.70$ minutes
 $F = 1.25$ mL/min

ACE Excel 1.7 C18-PFP 50 x 3.0 mm -

No delayed injection



$t_G = 1.70$ minutes
 $F = 1.25$ mL/min

Method Translator Tool Example - Gradient

Input data into the grey boxes

Input dwell volume of LC systems

Gradient

| Column Information | | Translated | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|-------------------------------------|---|---------------------------------------|----------------------------------|------|---|-------|----|-------|----|-------|---|--|--|--|--|----------|--|------|----|------|---|------|----|------|----|------|---|--|--|--|--|
| Current | | Translated | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Column Length (L) | 150 mm | Column Length (L) | 50 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Column i.d. (d _c) | 4.6 mm | Column i.d. (d _c) | 3.0 mm | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Particle Diameter (d _p) | 5.0 μm | Particle Diameter (d _p) | 1.7 μm | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| L/d _p | 30000 | L/d _p | 29412 | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Column Porosity | 0.63 <small>What's This?</small> | Column Porosity | 0.63 <small>What's This?</small> | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Column Volume (V _M) | 1.570 mL | Column Volume (V _M) | 0.223 mL | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Method | | Translated | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Current | | Translated | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Injection Volume | 10.0 μL | Injection Volume | 1.4 μL | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Flow Rate | 1.00 mL/min | Flow Rate (scaled to linear velocity) | 0.43 mL/min | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | Flow Rate (scaled to particle size) | 1.25 mL/min | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | Input Flow Rate | 1.25 mL/min | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| LC Name | HPLC 1 | LC Name | UHPLC 1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Dwell Volume (V _D) | 1.098 mL | Dwell Volume (V _D) | 0.926 mL | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Recorded Backpressure | 67 bar | Estimated Backpressure | 568 bar | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | Estimated Run Time Difference | -87 % | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | Estimated Solvent Use Difference | -83 % | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Gradient | <table border="1"> <thead> <tr> <th>Time</th> <th>%B</th> </tr> </thead> <tbody> <tr><td>0.00</td><td>5</td></tr> <tr><td>15.00</td><td>95</td></tr> <tr><td>17.00</td><td>95</td></tr> <tr><td>17.50</td><td>5</td></tr> <tr><td></td><td></td></tr> <tr><td></td><td></td></tr> </tbody> </table> | Time | %B | 0.00 | 5 | 15.00 | 95 | 17.00 | 95 | 17.50 | 5 | | | | | Gradient | <table border="1"> <thead> <tr> <th>Time</th> <th>%B</th> </tr> </thead> <tbody> <tr><td>0.00</td><td>5</td></tr> <tr><td>1.70</td><td>95</td></tr> <tr><td>1.93</td><td>95</td></tr> <tr><td>1.98</td><td>5</td></tr> <tr><td></td><td></td></tr> <tr><td></td><td></td></tr> </tbody> </table> | Time | %B | 0.00 | 5 | 1.70 | 95 | 1.93 | 95 | 1.98 | 5 | | | | |
| Time | %B | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.00 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 15.00 | 95 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 17.00 | 95 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 17.50 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Time | %B | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.00 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1.70 | 95 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1.93 | 95 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1.98 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Suggested Re-equilibration Time | 16.8 mins | Suggested Re-equilibration Time | 2.6 mins | | | | | | | | | | | | | | | | | | | | | | | | | | | | |



New Gradient

This section tells you if a delayed injection of X mins is required

Time to delay injection after the gradient begins 0.62 mins (770 μL)

Correct translation of this method requires that the injection is delayed until after the gradient begins by this time. A delayed injection can be added to a method in many LC instrument software packages.



Transferring Between Instruments



Instrument to Instrument Same Method Transfers

◆ **Isocratic**

- If column format **remains identical, no changes** necessary.
- (If column format changes, translation required.)


◆ **Gradient**

- If column format **remains identical**, need to **correct for influence of differing system dwell volumes** only. Flow rate, injection volume and gradient times **remain unchanged**.
- (If column format changes, translation plus differing system dwell effects need to be calculated.)




Instrument to instrument method transfer: Isocratic

Transfer of an isocratic related substances method on **ACE Excel 1.7 Super C18 50 x 3.0 mm** between UHPLC systems.




| | UHPLC system 1 | | UHPLC system 2 | |
|---------------------|----------------|----------------|----------------|----------------|
| Extra column volume | 13 μ L | | 15.5 μ L | |
| | N | R _s | N | R _s |
| Peak 1 | 4285 | | 3748 | |
| Peak 2 | 4975 | 2.6 | 4388 | 2.5 |
| Peak 3 | 2455 | 5.0 | 1831 | 4.8 |
| Peak 4 | 9330 | 5.6 | 8303 | 5.1 |
| Peak 5 | 11346 | 19.7 | 10561 | 19.4 |
| Peak 6 | 11461 | 10.3 | 10644 | 10.1 |




Isocratic analysis 0.2% H₃PO₄ in MeCN:MeOH:water 35:5:60 v/v/v, 40°C, 1.25 mL/min, 0.7 μ L, Injection vol., 254 nm.
 1. 4-hydroxybenzoic acid, 2. 4-hydroxyisophthalic acid, 3. acetylsalicylic acid (aspirin), 4. salicylic acid, 5. acetylsalicylsalicylic acid, 6. salsalate. Impurities were spiked at 0.5% w/w.

Binary Pump Configuration: High Pressure Mixing




Quarternary Pump Configuration: Low Pressure Mixing





Gradient Cycles: What Really Happens With Dwell Volume





Instrument Data Can Vary Depending Upon Configuration

| Instrument | Typical Dwell Volume (µL) |
|----------------------------|---------------------------|
| Agilent HP1100 Binary | 180-900 |
| Agilent HP1100 Quaternary | 800-1100 |
| Agilent 1200 RRLC | ~300 |
| Dionex P680A Quaternary | <400 |
| Thermoquest P4000 Quat | <600 |
| Waters Alliance 2695 Quat | 600 |
| Waters Varian 9012 Ternary | 1000 |
| Waters Acquity UPLC | ~100 |

| | ACQUITY H-Class | 1290B (Standard) | 1290B (Optimised) | 1290Q | Chromaster UltraRs | 1200RR | 1100 LC125 | 1200MD | 1100 LC94 | 1100 LC98 |
|-------------------|-----------------|------------------|-------------------|-------|--------------------|--------|------------|--------|-----------|-----------|
| ECV (µL) | 21.4 | 13.6 | 8.3 | 21.3 | 18.3 | 18.7 | 23.9 | 29.8 | 62.0 | 102.5 |
| Bandwidth | 9.96 | 9.82 | 7.40 | 12.78 | 10.94 | 11.45 | 15.48 | 14.28 | 25.62 | 42.87 |
| Dwell volume (µL) | 372 | 198 | 198 | 950 | 262 | 740 | 1082 | 1110 | 1104 | 1100 |

See ACE Knowledge Note #0001 for how to measure system dwell volume



Instrument to Instrument Same Method Transfers


- ◆ For **gradient method transfer**, t_G , F and V_M **all remain identical**.
- ◆ Only need to correct for the **system dwell volume**
 - Given that **V_M does not change** for method transfer, the V_D/V_M equation simplifies to:

$$x = \left(\frac{V_{D1} - V_{D2}}{F} \right)$$

- **Negative** value: **injection must be delayed** x minutes after gradient starts.
- **Positive** value: a **pre-gradient isocratic hold** of x minutes should be added to gradient program.




Method Translations: Dwell Volumes / Injection Times



Gradient starts $\Delta (V_D/V_M) = \text{negative}$

e.g. small column on a high dwell volume system. The dwell time is artificially long.

Injection is delayed until after the start of gradient




Gradient starts $\Delta (V_D/V_M) = \text{positive}$

e.g. large column on a low dwell volume system. The dwell time is artificially short.

Pre-gradient hold is added. Effectively extends dwell time

UHPLC Method Transfer - Vanillins



Original method (Agilent 1290 Quaternary)

Dwell = 0.574 mL

Column: **ACE Excel 1.7 C18-Amide 50 x 3.0 mm**


A: 0.1% formic acid

B: 0.1% formic acid in MeCN


Flow rate: 1.3 mL/min

| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |


| | |
|------------|-----------|
| t_G | 1.32 |
| Post time: | 2 min |
| Inj. Vol. | 1 μ L |
| P_{MAX} | 502 bar |




1. vanillic acid




2. 4-hydroxybenzoic acid




3. vanillin




4. 4-hydroxybenzaldehyde




5. guaiacol



6. o-vanillin



7. ethyl vanillin



8. eugenol





UHPLC Method Transfer - Vanillins

◆ **Aim**

- To translate the method from Agilent 1290 Quaternary to Agilent 1290 binary
- Need to correct for system dwell volume only
- $V_{D1} - V_{D2} = \Delta = 0.574 - 0.202 = 0.372$ mL
- **Positive** value, therefore a **pre-gradient hold (x)** is required
- A negative value would require a delayed injection
- $x = \frac{|\Delta|}{F_2} = 0.29$ min **pre-gradient hold**



Instrument to Instrument Method Transfer






| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |



**Method transfer –
no correction for
system dwell**

| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |

Instrument to Instrument Method Transfer



| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |



**Method transfer –
no correction for
system dwell**

| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |



**Method transfer –
Corrected for ΔV_D**

| t | %B |
|-------------|-----------|
| 0 | 25 |
| 0.29 | 25 |
| 1.61 | 75 |
| 1.78 | 75 |
| 1.89 | 25 |

Method Transfer Tool Example - Gradient

V1.3

Method Transfer



Method transfer involves moving a method from one LC instrument to another, whilst keeping the column format constant. Isocratic analysis is simple with no method changes required. Gradient methods however, should be adjusted to account for any change in system dwell volume in order to ensure accurate method transfer. This tool automatically determines any correction required to ease method transfer. See the Dwell Volume tab for details of how to determine system dwell volume.

Column Information

| | | |
|-----------------------------|-------|-----------------------------|
| Column Length (L) | 50 | mm |
| Column i.d. (d_c) | 3.0 | mm |
| Particle Diameter (d_p) | 1.7 | μm |
| L/ d_p | 29412 | |
| Column Porosity | 0.63 | <small>What's This?</small> |
| Column Volume (V_M) | 0.223 | mL |

Method

Current

Flow Rate mL/min

LC Name
 Dwell Volume (V_D) mL

Gradient

| Time | %B |
|------|------|
| 0.00 | 25.0 |
| 1.32 | 75.0 |
| 1.49 | 75.0 |
| 1.60 | 25.0 |
| 3.60 | 25.0 |
| | |
| | |

Transferred

LC Name
 Dwell Volume (V_D) mL

Gradient


| Time | %B |
|------|------|
| 0.00 | 25.0 |
| 0.29 | 25.0 |
| 1.61 | 75.0 |
| 1.78 | 75.0 |
| 1.89 | 25.0 |
| 3.89 | 25.0 |

New Gradient with
isocratic hold



Transferring Between Instruments Part II

UHPLC Method Transfer - Vanillins



Original method (Agilent 1290 Quaternary)

Dwell = 0.574 mL

Column: ACE Excel 1.7 C18-Amide 50 x 3.0 mm


A: 0.1% formic acid

B: 0.1% formic acid in MeCN


Flow rate: 1.3 mL/min

| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |


| | |
|-------------------|-----------|
| t_G | 1.32 |
| Post time: | 2 min |
| Inj. Vol. | 1 μ L |
| P_{MAX} | 502 bar |




1. vanillic acid




2. 4-hydroxybenzoic acid




3. vanillin




4. 4-hydroxybenzaldehyde




5. guaiacol



6. o-vanillin




7. ethyl vanillin



8. eugenol




Instrument to Instrument Method Transfer



| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |




no correction for system dwell



| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |


Instrument to Instrument Method Transfer



| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |




no correction for system dwell



| t | %B |
|------|----|
| 0 | 25 |
| 1.32 | 75 |
| 1.49 | 75 |
| 1.6 | 25 |









Method corrected for ΔV_D



| t | %B |
|------|----|
| 0 | 25 |
| 0.22 | 25 |
| 1.54 | 75 |
| 1.71 | 75 |
| 1.82 | 25 |

Transferring Methods - Success or Not?




Summary

- ◆ **Relatively accurate** method translations are achievable
- ◆ **Isocratic** translations
 - ◆ Use **L/dp** ratio translations
 - ◆ Translate to **columns that have scalable** bonded phases
 - ◆ Be aware of **system dispersion effects** on smaller columns
- ◆ **Gradient** translations
 - ◆ **More complex** but easy calculations described
 - ◆ Consider **V_D/V_m impact** upon translation accuracy ← **few** online calculators currently offer this
- ◆ **High peak capacities** are possible for **complex samples** using **column coupling**
- ◆ **Instruments transfers** are **simple for isocratic** methods but **gradient methods** must **consider dwell volumes**


Useful Resources

- ◆ ACE Translation Tool:
 - ◆ (+help file)
 - ◆ (+ AKN#0023)




- ◆ ACE Knowledge Notes (AKNs):

- ◆ **AKN0001** - How to Determine System Dwell Volume
- ◆ **AKN0006** – UHPLC Column Connections
- ◆ **AKN0011** – Practical UHPLC
- ◆ **AKN0012** - Understanding the Relationship between Particle Size, Performance and Pressure
- ◆ **AKN0017** - How to Determine Extra Column Dispersion and Extra Column Volume
- ◆ **AKN0023** - Gradient Method Translation Using the ACE LC Translator



Download Resources at mac-mod.com
 or contact info@mac-mod.com



Thank You For Your Attention

info@mac-mod.com

MAC-MOD.com