


Applications of Fused-Core[®], Superficially Porous Particles (SPP) in Environmental Analysis

Andrew F. Harron, Ph.D.
Applications Scientist

Advanced Materials Technology, Inc.
Wilmington, Delaware, USA
aharron@advanced-materials-tech.com




Outline

- How SPP works
- What are the advantages of using Fused-Core[®] particle technology
- Versatile columns available for application development



How SPP Works



SPP Benefits

- The benefits of superficially porous particles (SPPs) include highly efficient separations in less time and with reduced backpressure requirements.
- In the applied/small molecule market throughput is the one of the most important factors in method development.

Higher Overall Efficiency (Maximum N)

- Comparison of different columns and maximum plates (theoretical calculations for small molecule at 50% ACN/water)

Column type	Length (mm)	Flow rate (mL/min)	Plates (N)	Pressure (bar)
5 μm FPP	150	0.6	14,600	100
3 μm FPP	150	0.6	24,200	309
1.8 μm FPP	150	0.6	30,840	771
5 μm SPP	150	0.6	28,300	78
2.7 μm SPP	150	0.6	38,300	284



Applications

- High throughput-Environmental labs, drug screening labs-Screening methods for drugs and mycotoxins. PFP (Pentafluorophenylpropyl) and Biphenyl phase used. Multiple isobaric species
- Triglyceride profiling-Food and beverage industry-analysis of edible oils-C30 phase used. Multiple isobaric species

Mycotoxins

- Produced by mold actively growing on crop, or can also begin growing during crop storage
- Toxic-crop destruction
- Heat stable to canning
- Danger to humans and animals (feed crops)- secondary metabolites of fungus
- HALO[®] PFP

TEST CONDITIONS:

Columns: HALO[®] 90 Å PFP, 2 μm, 2.1 x 50mm

Part Number: 91812-409

Mobile Phase A: Water/2mM ammonium formate/0.1% Formic acid

Mobile Phase B: Methanol/2mM ammonium formate/0.1% Formic acid

Gradient:	Time	%B
	0.01	15
	1.0	25
	2.0	40
	2.50	41
	4.50	100
	5.50	100
	5.51	15
	6.50	Finished

Flow Rate: 0.4 mL/min

Initial Pressure: 485 bar

Temperature: 40 ° C


Injection Volume: 1 μL

Sample Solvent: 95/5 water/methanol


LC System: Shimadzu Nexera X2

Detection: +ESI MS/MS

MS: Orbitrap Exactive mass spectrometer



Mycotoxins




Peak Number	Compound	Retention Time	Precursor Ion	Product Ion
1	Nivalenol	0.71	313.1235	175.10
2	Deoxynivalenol	1.38	297.1335	249.09
3	Deoxynivalenol-3-glucoside	1.70	459.1850	193.10
4	Fusarenon X	2.37	355.1387	247.10
5	Neosolaniol	2.87	383.1702	365.16
6	15-Acetyldeoxynivalenol	3.33	339.1378	321.15
7	3-Acetyldeoxynivalenol	3.36	339.1378	231.15
8	Gliotoxin	3.97	327.0436	196.08
9	Aflatoxin G2	4.27	331.0759	312.97
10	Aflatoxin M1	4.39	329.0604	273.12
11	Aflatoxin G1	4.40	329.0601	242.90
12	Aflatoxin B2	4.44	315.0820	284.87
13	HT-2 + Na	4.47	447.1934	345.10
14	Diacetoxyscirpenol	4.49	367.2637	307.15
15	Aflatoxin B1	4.52	313.0662	286.99
16	Ochratoxin A	4.67	404.0855	238.99
17	T-2 + Na	4.72	489.2049	245.09
18	Ochratoxin B	4.88	370.1321	324.15
19	Citrinin	4.96	251.0860	233.09
20	Zearalenone	5.11	319.1491	283.08
21	Patulin +MEOH	5.11	187.0723	98.95
22	Fumonisin B1	5.24	722.3868	334.25
23	Fumonisin B3	5.41	706.3901	336.25
24	Fumonisin B2	5.44	704.3901	336.25

Mycotoxin Isomeric resolution

Notice the resolution of the isomers, 15-acetyl deoxynivalenol and 3-acetyl deoxynivalenol (Peaks 6 & 7). Also the isobaric compounds aflatoxin M1 & aflatoxin G1 (Peaks 10 & 11) are resolved.

Peak Number	Compound	Retention Time	Precursor Ion	Product Ion
1	Nivalenol	0.71	313.1235	175.10
2	Deoxynivalenol	1.38	297.1335	249.09
3	Deoxynivalenol-3-glucoside	1.70	459.1850	193.10
4	Fusarenon X	2.37	355.1387	247.10
5	Neosolaniol	2.87	383.1702	365.16
6	15-Acetyldeoxynivalenol	3.33	339.1378	321.15
7	3-Acetyldeoxynivalenol	3.36	339.1378	231.15
8	Gliotoxin	3.97	327.0436	196.08
9	Aflatoxin G2	4.27	331.0759	312.97
10	Aflatoxin M1	4.39	329.0604	273.12
11	Aflatoxin G1	4.40	329.0601	242.90
12	Aflatoxin B2	4.44	315.0820	284.87
13	HT-2 + Na	4.47	447.1934	345.10
14	Diacetoxyscirpenol	4.49	367.2637	307.15
15	Aflatoxin B1	4.52	313.0662	286.99
16	Ochratoxin A	4.67	404.0855	238.99
17	T-2 + Na	4.72	489.2049	245.09
18	Ochratoxin B	4.88	370.1321	324.15
19	Citrinin	4.96	251.0860	233.09
20	Zearalenone	5.11	319.1491	283.08
21	Patulin+MEOH	5.11	187.0723	98.95
22	Fumonisin B1	5.24	722.3868	334.25
23	Fumonisin B3	5.41	706.3901	336.25
24	Fumonisin B2	5.44	704.3901	336.25



SAMHSA-5 Panel

- SAMHSA-The Substance Abuse and Mental Health Services Administration
- Amphetamines (Amphetamines, Methamphetamines, Speed)
- Cocaine (Cocaine, Crack, Benzoyllecognine)
- Marijuana – THC Tetrahydrocannabinol (Cannabinoids, Hash)
- Opiates (Heroin, Opium, Codeine, Morphine)
- Phencyclidine (PCP)
- HALO® Biphenyl

Columns: HALO® 90 Å Biphenyl, 2 µm, 2.1 x 100

Part Number: 91812-611

Mobile Phase A: Water/0.1% Formic acid

Mobile Phase B: Methanol/0.1% Formic acid

Gradient:	Time	%B
	0.0	5
	4.00	98
	5.00	98
	5.01	5
	7.00	END

Flow Rate: 0.4 mL/min

Initial Pressure: 325 bar


Temperature: 40 ° C

Injection Volume: 2 µL


Sample Solvent: 95/5 MEOH/Water

LC System: Shimadzu Nexera X2

MS: Orbitrap Exactive mass spectrometer




SAMHSA-5 Panel



SAMHSA-5 Panel isomeric resolution

Peak number	Compound	Observed ion
1	Morphine	286.341 m/z
2	Amphetamine	136.206 m/z
3	Methamphetamine	150.237 m/z
4	MDA	180.221 m/z
5	Phentermine	150.233 m/z
6	Codeine	300.364 m/z
7	6-MAM	328.380 m/z
8	MDMA	194.246 m/z
9	MDEA	208.271 m/z
10	Benzocgonine	290.331 m/z
11	PCP	244.387 m/z
12	THC-COOH	345.415 m/z



Notice the resolution of the isomers, methamphetamine and phentermine (Peaks 3 and 5)

TAG profiles of Edible oils on a C30

- Long chain fatty acids and esters
- Triglycerides & diglycerides
- Hydrophobic



Columns: HALO® 160 Å C30, 2.7 µm, 2.1 x 150

Part Number: 92115-730

Mobile Phase A: MEOH 10mM ammonium formate/0.1% formic acid

Mobile Phase B: IPA/0.1% Formic acid

Gradient:	Time	%B
	0.0	5.0
	10.00	20
	15.00	80
	25.00	80
	25.50	5.0
	32.00	END

Flow Rate: 0.4 mL/min

Initial Pressure: 325 bar

Temperature: 40 ° C



Injection Volume: 2 µL

Sample Solvent: MEOH


LC System: Shimadzu Nexera X2

MS: Orbitrap Exactive mass spectrometer


C18 vs C30 comparison



Corn oil




Coconut oil




Major TAG components in Coconut oil	Observed Ions	Peak Number
TAG (C32)	516.506 m/z	1
TAG (C34)	544.554 m/z	2
TAG (C34)	572.668 m/z	3
TAG (C36)	600.502 m/z	4
TAG (C36)	628.616 m/z	5
TAG (C36)	656.713 m/z	6
TAG (C36)	684.774 m/z	7
TAG (C38)	712.598 m/z	8
TAG (C38)	740.518 m/z	9
TAG (C40)	768.566 m/z	10
TAG (C40)	796.614 m/z	11

Grape Seed oil



Major TAG components in Grape seed oil	Observed Ions	Peak Number
DAG	634.475 m/z	1
DAG	636.491 m/z	2
TAG (POL)	894.718 m/z	3
TAG (LLL)	896.358 m/z	4
TAG (OLL)	898.471 m/z	5
TAG (PLL)	900.385 m/z	6
TAG (LLL)	900.485 m/z	7

Rapeseed/Canola Oil



Major components in canola oil	Observed Molecular mass (m/z)	Peak Number
DAG	634.425 m/z	1
DAG	636.495 m/z	2
DAG	638.437 m/z	3
TAG (LnLO)	858.767 m/z	4
TAG (LLO)	872.755 m/z	5
TAG (POO)	874.355 m/z	6
TAG (LOO)	885.432 m/z	7
TAG (LnOO)	896.432 m/z	8
TAG (LOP)	900.734 m/z	9
TAG (OOO)	902.499 m/z	10

Conclusion

- HALO columns are an attractive option for applications requiring high throughput and high selectivity.
- In addition, the columns exhibit excellent cross market appeal (Including biologics) by enabling separation and resolution of hydrophobic, as well as isomeric compounds.

Acknowledgements

- **Scientists at Advanced Materials Technology**
Conner McHale, Stephanie Schuster, William Johnson