



Arsenolipids from Edible Seaweed by LC-ICP-MS and LC-ESI-MS

Separation of arsenic species from methanolic extract of the edible seaweed *Alaria esculenta*

Arsenic-containing hydrocarbon:


m/z 389 [M + H]⁺ for C₂₁H₄₆AsO

Arsenic-containing phospholipids:

m/z 959 [M + H]⁺ for C₄₅H₈₉AsO₁₄P (C16:0/C16:0)

m/z 987 [M + H]⁺ for C₄₇H₉₃AsO₁₄P (C18:0/C16:0)

m/z 1015 [M + H]⁺ for C₄₉H₉₇AsO₁₄P (C20:0/C16:0)



ACE C18, 3µm 150 x 4.6mm

Gradient analysis

A = 0.1% formic acid in H₂O

B = 0.1% formic acid in CH₃OH

Time (mins) %B

0 0

20 100

45 100

Flow rate: 1ml/min

Injection volume: 100µl

Split ratio: 75% ESI-MS: 25% ICP-MS

Thermo Scientific Element 2 ICP-MS

Mode: Organic mode

Medium resolution

Thermo Scientific Orbitrap Discovery

Positive ESI mode


Spray voltage: 4.5kV

Capillary temperature: 320°C

Capillary voltage: 42V

Microcystins From Blue/Green Algae In Drinking Water

ACE Excel 2 μm C18, 100 x 2.1 mm



0.05 ppb each

Variants	R	L
MC-LR	Leucine	Arginine
MC-RR	Arginine	Arginine
MC-YR	Tyrosine	Arginine

Bruker Advance UHPLC system

ACE Excel 2 μm C18, 100 x 2.1mm

Gradient elution

A = 0.1% formic acid in water

B = Acetonitrile

T (mins) %B T (mins) %B

0 30 7.1 30

1 30 10 30

7 95

Flow rate: 0.4mL/min

Column temperature: 40°C

Injection volume: 50 μL

Concentration each microcystin: 0.05ppb


Bruker EVOQ Elite triple quad MS

VIP heated-ESI temperature: 350°C


Cone gas temperature: 200°C

Spray voltage: 4500V (+)


Collision gas: argon 1.5mTorr



Polycyclic Tetracarboxylic Acids



C80-82 polycyclic tetracarboxylic acids isolated from oilfield deposits



Tetramethyl ester of C_{80:8} ring acid


ACE UltraCore SuperPhenylHexyl 2.5µm, 100 x 2.1mm
 Gradient analysis
 A: CH₃OH – H₂O (98:2) containing 10mM ammonium acetate
 B: IPA-H₂O (98:2) containing 10mM ammonium acetate

T (mins)	%B	T (mins)	%B
0	0	15	100
1	0	25	100

Flow rate: 0.15ml/min
 Column temperature: Ambient
 Injection volume: 5µl
 LCQ Ion trap MS
 LC-ESI-MS extracted ion chromatograms
 Compounds detected as ammoniated quasimolecular ions
 [M+NH₄]⁺
 Detection limit ~ 0.1ppm



Perfluoro acids by LC-MS/MS



Peak ID	Analyte
1	Heptafluorobutyric acid
2	Perfluorohexanoic acid
3	Perfluorobutylsulphonic acid
4	Perfluoroheptanoic acid
5	Perfluorooctanoic acid
6	Perfluorohexylsulphonic acid
7	Perfluorononanoic acid
8	Perfluorodecanoic acid
9	Perfluorooctanesulphonic acid

ACE Excel 2 C18 2 μ m, 50 x 2.1mm

Gradient analysis

A = 2mM NH₄OAc, 0.1% acetic acid/CH₃CN (95:5)

B = 2mM NH₄OAc, 0.1% acetic acid/CH₃CN (5:95)

T (mins)	% B	T (mins)	% B
0	25	7.5	95
0.5	25	8.0	25
5.5	95	10.0	25

Flow rate: 0.5ml/min

Column temperature: 40°C

Injection volume: 20 μ l





AB SCIEX triple quad 5500

Negative ESI MRM

Source temperature: 450°C

IonSpray voltage: -2400V

Pesticides by LC-MS/MS



Reproduced with permission of Barry Whatmore, Kent Scientific Services

ACE UltraCore SuperC18, 2.5µm, 50 x 2.1mm
Gradient analysis

A = 0.1% HCOOH + 5mM NH₄CO₂H in 9:1 v/v H₂O: MeOH
B = 0.1% HCOOH + 5mM NH₄CO₂H in 1:9 v/v H₂O: MeOH

Flow Rate: 0.4ml/min Gradient conditions
Temperature: 40°C Time (mins) 0 1 15 18 18.05 20
Injection volume: 20µl %B 0 0 100 100 0 0

Agilent 6420 Triple Quadrupole MS, +ve mode ESI
Dynamic MRM

Also analysed under same conditions:


- | | |
|----------------------------|--------------------|
| Acephate | Hexaconazole |
| Acetamiprid | Hexaflumuron |
| Aldicarb | Imidacloprid |
| Aldicarb sulphone | Indoxacarb |
| Aldicarb sulphoxide | Mandipropamid |
| Benomyl | Methamidophos |
| Carbendazim | Methomyl |
| Carbofuran | Monocrotophos |
| Clofentezine | Nicotine |
| Clothianidin | Omethoate |
| Cyfluthrin | Oxamyl |
| Demeton S-methylsulphone | Pencycuron |
| Demeton S-methylsulphoxide | Prochloraz |
| Dicrotophos | Propargite |
| Dimethoate | Thiabendazole |
| Dinotefuran | Thiacloprid |
| DMA | Thiamethoxam |
| DMPF | Thiodicarb |
| Flubendiamide | Thiophanate methyl |
| Folpet | Triforine |
| Formetanate | |



Organophosphorus Flame Retardants in Water by LC-MS/MS

Symbol	Compound Name	Q1 Mass	Q3 Mass
TMP	tri-methyl phosphate	141	109
TEP	tri-ethyl phosphate	183	127
TiPP	tri-iso-propyl phosphate	225	99
TPrP	tri-n-propyl phosphate	225	99
TBP	tri-n-butyl phosphate	267	211
TCEP	tris(2-chloroethyl) phosphate	285	223
TCPP	tris((2R)-1-chloro-2-propyl) phosphate	327	99
TDCPP	tris(1,3-dichloro-2-propyl) phosphate	431	99
BDCP	bis(1,3-dichloro-2-propyl) phosphate	321	99
TPP	triphenyl phosphate	327	215
TBEP	tris(2-butoxyethyl) phosphate	399	299
BEHP	bis(2-ethylhexyl) phosphate	323	99

ACE C18 3µm, 100 x 2.1mm
 Gradient analysis
 A = 0.05mM NH₄CO₂H + 0.005% HCO₂H in water
 B = CH₃OH/CH₃CN (95:5)
 Time (mins) %B Curve
 0.1 50 -3
 12.0 90
 13.0 100
 15.0 100
 15.1 50
 20.0 50
 Flow rate: 0.25ml/min
 Column temperature: 25°C
 Sample volume: 80µl
 Detection: MS/MS





ACE UltraCore SuperC18: Impurity Profile of a Herbicide

ACE UltraCore SuperC18, 2.5 μ m, 150 x 4.6mm

Gradient analysis

A = CH₃CN – H₂O – TFA (5:95:0.05 v/v/v)

B = CH₃CN – TFA (99.9:0.05 v/v/v)

Time (mins)	%B	Time (mins)	%B
0	10	55	100
3	10	56	10
35	100	60	10

Flow rate: 0.60ml/min

Column temperature: 25°C

Injection volume: 10 μ l

Detection: UV, 240nm

Sample: Technical Grade Herbicide

