INDUSTRIAL

HALO





PEAK IDENTITIES:

- 1. Benzene
- 2. Toluene
- 3. o-Xylene
- 4. Ethylbenzene
- 5. m-, p-Xylenes

TEST CONDITIONS:

Column: HALO 90 Å C18, 2.7 μm, 0.25 x 170mm **Mobile Phase:** A: 94.5/5.0/0.5 Water/ Acetonitrile/ formic acid B: Acetonitrile

Gradient: %В Time 0.0 3 9.0 5 10.0 50 22.0 62 30.0 90 Flow Rate: 1.5 µL/min Pressure: 276 bar **Temperature:** Ambient Detection: 255 nm, UV Delay Volume: 1 µL LC System: Axcend Focus LC Data Courtesy of: Axcend

Benzene, toluene, ethylbenzene, and xylene (BTEX) are compounds that occur naturally in crude oil and can be found in sea water close to natural gas and petroleum deposits where drilling occurs.

BTEX compounds are created and used during the processing of petroleum products. In this application, a separation of BTEX at 6 ppm is performed using a HALO® C18 capillary column on an Axcend Focus portable LC.

Liquid chromatography is an alternative method to gas chromatography, which eliminates the use of a flame and hydrogen gas which is ideal for off- shore oil rigs. Smaller column internal diameters, like capillaries, allow for higher sensitivity which can benefit those who are working with very small sample concentrations.



AMT_AN_Rev_0
advancedmaterialstechnology
Made in the USA

INNOVATION YOU CAN TRUST, PERFORMANCE YOU CAN RELY ON HALO® and Fused-Core® are registered trademarks of Advanced Materials Technology