# Bridging the Gap between UHPLC and HPLC: Easy Method Transfer Using Fused-Core® Columns

## Abstract

Many chromatographers are now exploiting the speed and efficiency advantages of UHPLC as part of their method development scheme. The use of low dispersion UHPLC icient columns allows them to screen different column ever before. Moreover, with this approach the resulting high speed and/or ethod will be more robust and efficient, and will be able to general analytical results that enable faster and better decisions with high productivity

transfer such methods to a quality control or productio od developers often must cope with the limited availability of UHPLC and operator expertise in such laboratories. Often, it is necessary to transfe the method to a longer and larger ID column with a much larger particle size to be able to transfer the method to conventional instrumentation (300–400 bar pressure limit).

Fused-Core® UHPLC columns with 2.7 µm particle size can deliver performance comparable to sub-2-µm columns at 40–50% of the back pressure. This benefit makes method transfer much easier between UHPLC and HPLC systems. It is only necessary to make modifications to the extracolumn volume and dispersion of the HPLC system, with some additional adjustments to the analysis conditions.

We will discuss the parameters that must be considered and adjusted when transferring adient methods from UHPLC instruments to conventional HPLC system Isocratic and gradient examples will be shown in which Fused-Core UHPLC separations are transferred to different HPLC instruments having different extracolumn dispersion and different delay volumes. Guidance will be offered to make such method transfers more successful

# Objectives

- Discuss important parameters that affect method transfer for isocratic and gradient methods.
- Demonstrate transfer of an isocratic and gradient method from UHPLC to HPLC, and summarize results
- Show importance of injection delay for gradient method transfer.

#### Mobile Phase A Phase B Impact on **Method Transfer Pumping Syste** High pressure mixing systems typicall have lower delay volumes, with less • High Pressu "rounding" of gradient profiles Mixing Some systems have pressure-depender • Low Pressure delay volume. Mixing • Affects delay volume Affects gradient shape and especially Mixer Volum baseline noise with some solvents and additives (e.g., TFA) Affects Extracolumn Dispersion **Tubing Volume** (aka Peak Variance) and delay volume Differences can change selectivity Delay Volume Affects retention times in gradients Can affect extracolumn dispersion, es Flow Rate for different instruments Can affect extracolumn dispersion, esp Injector Type for different instruments • Autosampler Fixed loop has lower dispersion that (Fixed loop, flo flow through needle type through needle • Manual injectors often least dispersion

Manual injector

# Guidance for Maximum Flow Rates for UHPLC Methods to be Transferred to HPLC with HALO



• Calculations are approximate.

• Pressure for maximum flow rate set for 80% of 400 bar.

• Estimates made with 0.005" ID tubing in typical length (~60-70 cm) to account for system pressure (not including flow cell and autosampler)

# General (U)HPLC Method Transfer Concerns

- Column length and diameter
- Instrument brand and model
- Instrument pressure limit
- Pumping system type: low pressure mixing vs. high pressure mixing • Gradient delay volume and flow
- rate/pressure dependence
- Mixer volume
- Flow Rate
- Column Oven type and Column Temperature
- Detector flow cell and volume
- Detector data rate and response time

•	Injector type	and injection volume
---	---------------	----------------------

ACN/Water						
	ACN viscosities	1.01	0.91	0.75	0.62	0.52
ID	Length					
(mm)	(mm)	25°C	30°C	40°C	50°C	60°C
3.0	20	3.12	3.46	4.20	5.00	5.00
3.0	30	2.26	2.51	3.05	3.69	4.40
3.0	50	1.46	1.62	1.97	2.38	2.84
3.0	75	1.01	1.12	1.36	1.65	2.00
3.0	100	0.78	0.86	1.04	1.26	1.50
3.0	150	0.53	0.59	0.71	0.86	1.00
ACN/Water						
	ACN viscosities	1.01	0.91	0.75	0.62	0.52
ID	Length					
(mm)	(mm)	25°C	30°C	40°C	50°C	60°C
2.1	20	1.53	1.70	2.06	2.45	2.45
2.1	30	1.11	1.23	1.49	1.81	2.16
2.1	50	0.72	0.79	0.97	1.17	1.39
2.1	75	0.49	0.55	0.67	0.81	0.98
0.4	100	0.38	0.42	0.51	0.62	0.74
2.1	100					

#### Original Shimadzu Inst. Method

Column:	4.6 x 50 mm, HALO C18
Mobile Phase:	43:57 A/B
Mobile Phase A:	0.020M sodium phospha
	(pH=2.5)
Mobile Phase B:	50:50 MeOH/ACN prem
	20. T/.
Flow Rate:	3.0  mL/min.
Pressure:	338 Bar
Temperature:	35°C
Detection:	UV 254 nm, VWD
<b>Injection Volume:</b>	2.0 μL
Sample Solvent:	50:50 MeOH/water
Data Rate:	50 Hz.
<b>Response Time:</b>	0.02 sec.
Flow Cell:	2.5 μL semi-micro
LC System:	Shimadzu Prominence
	UFLC XR
Pressure maximum:	600 bar

S. A. Schuster<sup>1</sup> and <u>T.J. Waeghe<sup>2</sup></u> <sup>1</sup>Advanced Materials Technology Inc., 3521 Silverside Road, Suite 1-K, Quillen Building, Wilmington, DE 19810; <sup>2</sup>MAC-MOD Analytical, Inc., 103 Commons Court, Chadds Ford, PA 19317

## **Considerations for Method Transfer**

osampler Heat E	umr atei Exchai	nger	Detector Flow Cell Waste		
Parameter	Ι	G	Impact on Method Transfer		
Column Heater • Forced air • Block heater (contact)	Y	Y	<ul> <li>Actual measured temperatures will vary among instruments.</li> <li>Measured temperature may not match method set point.</li> </ul>		
Detector Type • VWD • Diode Array (DAD)	Y	Y	<ul> <li>Each detector brand and model can vary in performance within model and within brand.</li> <li>Performance among brands can vary for noise, sensitivity, linearity</li> </ul>		
Flow Cell • Design • Volume • Temperature	Y	Y	<ul> <li>Flow cell design can affect dispersion significantly</li> <li>Flow cell pathlength affects noise and signal.</li> <li>Mismatch of column and detector temperatures increases noise.</li> </ul>		
<ul> <li>Detector Settings</li> <li>Data rate (pts/sec)</li> <li>Response Time</li> <li>Detection wavelength bandwidth</li> <li>Reference wavelength</li> </ul>		Y	<ul> <li>Inadequate data rate causes poor peak fidelity, decreases observed efficiency, and increases tailing.</li> <li>Bandwidth mismatch or use/absence of reference wavelength for DAD can cause peak area ratio changes and RI anomalies, respectively.</li> </ul>		

# **Isocratic Method Transfer Example**



scaled with ratio of column IDs squared and by volume ratio, respectively

$$F_{col2} = F_{col1} \times \left(\frac{ID_{col2}}{ID_{col1}}\right)^{2}$$
$$\left(V_{inj}\right)_{col2} = \left(V_{inj}\right)_{col1} \times \left(\frac{ID_{col2}}{ID_{col1}}\right)^{2} \times \frac{L_{col2}}{L_{col1}}$$

# Instruments, Configurations and Method Parameter Settings

Parameter	Agilent 1100	Agilent 1200	Shimadzu Prominence
Pressure Limit (bar)	400	600	600
Column Heater Type	Block	Block	Forced Air
Detector Type	VWD	DAD	VWD
Flow Cell Volume (µL)	5	2	2.5
Flow Cell Path length (mm)	10	6	5
Min. Response Time (sec.)	0.062	0.02	0.02
Maximum Data Rate	13.7 Hz	80 Hz (40 Hz used)	50 Hz

## Isocratic Method Transfer from Shimadzu to 1100 and 1200



## **Results for Isocratic Method Transfer from** Shimadzu to Agilent 1100 and 1200

#### <u>Results</u>

- Resolution results for respective peak pairs ranged from 91-117% for the Agilent 1100, and from 95-107% for the Agilent 1200.
- Theoretical plate counts increased with k for respective analytes as expected, up to  $\sim k = 8-9$
- Plate counts increased slightly with increasing ID (4.6 > 3.0 > 2.1).

#### **Observations**

- DAD reference signal had to be turned off to avoid RI disturbance near start of run.
- Data rate for Agilent 1200 was set at 40 Hz to more closely match 50 Hz rate of Shimadzu.
- Flow cell path length differences among instruments caused peak areas and heights to vary.
- Analytical wavelength bandwidth and slit width on Agilent 1200 DAD had to be adjusted lower to give comparable peak areas.

### **Gradient Method Transfer Example** Transfer method to Agilent 1100 and Shimadzu Prominence Systems

Original Agilent 1200 Inst. Method

Flow Rate:

Detection

Data Rate:

Flow Cell:

Pressure max.:

• Standard mixer was removed

• All tubing was minimal length and 0.005" II

• Automatic delay volume reduction was turned ON

Max. Pressure:

water with 0.1% HCOOH ACN with 0.1% HCOOH 3% ACN to 70% ACN in .7min 0.42 mL/min. 116 Bar DAD 275 nm, Bandwidth Injection Volume: 2.0 µL Sample Solvent: 50:50 MeOH/water 40 Hz. **Response Time:** 0.02 sec. 2 μL micro flow cell 600 bar

1 x 50 mm, HALO C18

Flow rates and injection volumes are scaled as with

$$F_{col2} = F_{col1} \times \left(\frac{ID_{col2}}{ID_{col1}}\right)$$

isocratic methods.

$$(V_{inj})_{col2} = (V_{inj})_{col1} \times \left(\frac{ID_{col2}}{ID_{col1}}\right)^2 \times \frac{L_{col2}}{L_{col1}}$$

Pelay volumes of each instrument should be measured, and, if available, injector delay should be used to correct the effective delay volume to that of the original method.

$$T_{injdelay} = rac{V_{D(HPLC)} - V_{D(UHPLC)}}{F}$$

An injector program can be set on some instruments so that the gradient starts, and the sample is injected after a defined time delay, T<sub>inidelay</sub>





# **Results for Gradient Method Transfer from** Agilent 1200 to Agilent 1100 and Shimadzu

Results

- Resolution results for respective peak pairs ranged from 97-120% for the Agilent 1100, and from 81-108% for the Shimadzu Prominence.
- for adjustment of effective delay volume using injection delay.
- with different injection delay times.
- changes in selectivity.

#### **Observations**

- Flow cell volume and path length differences among instruments caused peak areas and heights to vary.
- Analytical wavelength bandwidth and slit width on Agilent 1200 DAD
- had to be adjusted lower to give comparable peak areas.

# How Should One Measure Method Transfer Success?

- Comparable peak shape, efficiency, selectivity, and resolution
- Comparable injection-to-injection repeatability
- Comparable S/N Ratio (impact on LOD and LOQ, linearity, UV spectral quality, etc.)
- Comparable accuracy and precision and other figures of merit
- Not necessarily identical retention times!



# **Gradient Transfer to Agilent 1100** with Different HALO Column IDs



# **Gradient Method Comparison Chromatograms**

• Injector program feature of Agilent 1100 system and software allowed • Demonstrated changes in selectivity and resolution on Agilent 1100

• Gradients run on Shimadzu system w/o injection delay showed major

# **Recommendations for Method Transfer from UHPLC to HPLC with HALO Columns:**



# **Summary and Conclusions**

- Discussed important parameters for transfer of an isocratic and a gradient method using HALO columns with 3 different IDs in same length
- Presented results from method transfer experiments, and discussed issues that can arise.
- Provided guidance for successful transfer of methods from UHPLC to HPLC using HALO **Fused-Core columns.**