



# Method Details

## CONDITIONS

Column: Avantor® ACE® 3 C18-300  
 Particle Size: 3 µm  
 Dimensions: 50 x 2.1 mm  
 Mobile Phase: A: 0.1% formic acid in water  
 B: 0.1% formic acid in acetonitrile

Gradient:

Time (mins)	% B
0	10
1.5	33
1.6	10

Flow Rate: 0.4 mL/min  
 Injection Volume: 10 µL  
 Temperature: 40 °C  
 Detection: XEVO TQ-S triple quad MS  
 Desolvation temperature: 500 °C  
 Ion source temperature: 150 °C  
 Positive ion mode ESI, MRM

## ORDERING TABLE

Product	Details	Size	Part Number
Avantor® ACE® 3 C18-300	HPLC Column	50 x 2.1 mm	ACE-211-0502

Reproduced with permission of Biotage GB Ltd. For extraction conditions see Biotage Application Note AN882.