

# Development of a Method for the Separation of Major Cannabinoids Using Six Column Selectivity Screening

Geoffrey Faden<sup>1</sup>, Alan P McKeown<sup>2</sup>

<sup>1</sup>MACMOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317 USA, <sup>2</sup>Advanced Chromatography Technologies Ltd, 1 Berry Street, Aberdeen, Scotland, AB25 1HF UK

## 1. Introduction


- Cannabinoids are a major compound class and are therefore of primary importance.
- Cannabis samples are complex and contain several target compound classes that are of interest from analytical and regulatory perspectives.
- Regulatory requirements for testing are varied and may require determination of several key components, for example<sup>1</sup>:
  - $\Delta^9$ -Tetrahydrocannabinol ( $\Delta^9$ -THC)
  - Cannabidiol (CBD)
  - $\Delta^9$ -Tetrahydrocannabinolic acid A (THCA-A)
  - Cannabidiolic acid (CBDA)
  - Cannabigerol (CBG)
  - Cannabinol (CBN)
- This poster summarises work performed to develop an LC-MS compatible method for the analysis of an extended list of 10 cannabinoids of interest.

<sup>1</sup>State of California, AB 266 Medical marijuana, Article 9, Section 19344 (2015-2016)


## 2. Experimental

- The sample was screened on six reversed-phase columns, including five ACE novel chemistries to identify the most suitable stationary phase chemistry.
- The ACE novel chemistries have been engineered to provide unique selectivity and provide a comprehensive method development tool.


**ACE C18** (Hydrophobic interactions)




**ACE C18-AR** (Hydrophobic and phenyl interactions)




**ACE C18-Amide** (Embedded amide group to increase retention of polar components and alternative selectivity)




**ACE C18-PFP** (Hydrophobic and PFP interactions)



**ACE SuperC18** (Highly retentive with extended pH stability)



**ACE CN-ES** (Hydrophobic and polar interactions)



## 3. Six Column Screening

- Sample was screened on the six ACE phases using a 5-95% gradient.
- LC/MS compatible mobile phase.
- The six ACE phases provide alternative selectivity – ideal for method development.
- The **ACE SuperC18** was found to be the most promising phase.


**Conditions**  
Column Format: 100 x 3.0 mm, 2  $\mu$ m  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeCN

Gradient:

Time (mins)	%B
0.0	5
10.0	95
13.0	95
13.5	5


Flow Rate: 0.6 mL/min  
Injection: 5  $\mu$ L  
Temperature: 40 °C  
Detection: UV, 210 nm

Sample
1. (-)-11-Nor-9-carboxy- $\Delta^9$ -tetrahydrocannabinol (THC-COOH)
2. Cannabidiol (CBDV)
3. Cannabidiolic acid (CBDA)
4. Cannabigerol (CBG)
5. Cannabidiol (CBD)
6. Tetrahydrocannabinol (THCV)
7. Cannabinol (CBN)
8. (-)-trans- $\Delta^9$ -Tetrahydrocannabinol (THC)
9. Cannabichromene (CBC)
10. $\Delta^9$ -Tetrahydrocannabinolic acid A (THCA-A)




## 4. Optimised Separation

- The gradient conditions were optimised to reduce analyte retention and maintain resolution of the analytes.
- An initial isocratic hold was required to maintain separation of peaks 1&2.
- Full separation achieved in <8 minutes.



## 5. Rapid Analysis

- The method was translated to an **ACE UltraCore 2.5 SuperC18** solid core column using the ACE LC Translator Tool (Download free at [www.ace-hplc.com](http://www.ace-hplc.com)).



## 6. Conclusions

- Full separation of 10 cannabinoids was successfully achieved.
- Method development using six column screening rapidly identified the **ACE SuperC18** as the optimum stationary phase for the separation.
- Optimisation of the gradient resulted in full separation in less than 8 minutes using LC-MS compatible conditions.
- The final separation was translated to a 50 x 3.0 mm format **ACE UltraCore 2.5 SuperC18** column to further reduce the analysis time to < 3 minutes.