



# Exploiting Selectivity in HPLC and UHPLC With Rational Stationary Phase Design

Alan P McKeown  
amckeown@ace-hplc.com

Advanced Chromatography Technologies



[www.ace-hplc.com](http://www.ace-hplc.com)  
[www.mac-mod.com](http://www.mac-mod.com)




## Outline

- ◆ Chromatographic **selectivity**
- ◆ Stationary phase **design** concepts
- ◆ The unique **ACE<sup>®</sup> C18-AR** and **ACE<sup>®</sup> C18-PFP** phases
- ◆ Introducing the **NEW ACE<sup>®</sup> Excel<sup>™</sup>** UHPLC products
- ◆ Examples
- ◆ Conclusions

# Chromatographic Peak Resolution

Efficiency      Selectivity      Retention

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1+k}$$



# The Importance of $N$ , $k$ and $\alpha$ For Resolution

## Typical separation:

$N = 10,000$  plates

$k = 3.8 / 4.2$  (4.0 mean)


$\alpha = 1.1$

$$R_s = \frac{1}{4} \sqrt{10,000} \left( \frac{1.1 - 1}{1.1} \right) \left[ \frac{4}{1 + 4} \right]$$

$$R_s = 1.8$$

$$R_s = \frac{\sqrt{N}}{4} \frac{\alpha - 1}{\alpha} \frac{k}{1 + k}$$

Which looks like



# The Importance of $N$ , $k$ and $\alpha$ For Resolution


**Double Efficiency** (eg 5  $\mu\text{m}$   $\rightarrow$  2.5  $\mu\text{m}$ ):

$N = 10,000 \rightarrow 20,000$  plates

$$R_s = \frac{\sqrt{N}}{4} \frac{\alpha-1}{\alpha} \frac{k}{1+k}$$

$$R_s = \frac{1}{4} \sqrt{20,000} \left( \frac{1.1-1}{1.1} \right) \left[ \frac{4}{1+4} \right]$$

$$R_s = 2.6$$



Opportunity to optimise further eg reduce column length to speed up

# The Importance of $N$ , $k$ and $\alpha$ For Resolution


**Double Retention Factor (eg decrease solvent strength):**

$k = 4 \rightarrow 8$

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1+k}$$

$$R_s = \frac{1}{4} \sqrt{10,000} \left( \frac{1.1 - 1}{1.1} \right) \left[ \frac{8}{1+8} \right]$$

$$R_s = 2.0$$



$R_s = 1.8$



**~10% Increase**



$R_s = 2.0$

Slight improvement in resolution has led to increased analysis time

# The Importance of $N$ , $k$ and $\alpha$ For Resolution

**Increase Selectivity (eg change column):**

$\alpha = 1.1 \rightarrow 1.2$

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1+k}$$

$$R_s = \frac{1}{4} \sqrt{10,000} \left( \frac{1.2 - 1}{1.2} \right) \left[ \frac{4}{1 + 4} \right]$$

$$R_s = 3.3$$



$R_s = 1.8$



**~80% Increase**



$R_s = 3.3$

**Significant opportunity to speed up for modest change in selectivity**

# Selectivity: The Key to Chromatographic Peak Resolution


$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1+k}$$

Efficiency →  $\frac{\sqrt{N}}{4}$     
 Selectivity →  $\frac{\alpha - 1}{\alpha}$     
 Retention →  $\frac{k}{1+k}$

~40%     ~80%     ~10%



From the examples, selectivity has the greatest impact on increasing peak resolution





# Which Factors<sup>a</sup> Affect Selectivity?

**MOST  
Influence**



**LEAST  
Influence**

## Isocratic Separations


- ◆ **Column stationary phase**
- ◆ **Organic Modifier**
- ◆ **pH (ionised analytes only)**
- ◆ **% Organic modifier**
- ◆ **Buffer selection**
- ◆ **Column temperature**
- ◆ **Buffer concentration**

## Gradient Separations


- ◆ **All parameters for isocratic**
- ◆ **Gradient steepness**
- ◆  **$k^*$**
- ◆ **Dwell volume**
- ◆ **Column dimensions**

# Influencing Selectivity – Bonded Phase Effects / Basic Analytes


**ACE C18 – Increase Retention**




**ACE CN – Elution Order**




**ACE C8 (start point)**



**ACE C4 – Decrease Retention**



**ACE Phenyl – Elution Order**




Use **ultra high purity** silica for **good** chromatography and **reproducibility**

## HPLC End User Surveys<sup>a</sup> ...Listening To The Analyst

- ◆ **Column reproducibility** and **column lifetime** are major factors for analysts
  - Have been the **top 2** feedback points since **2007**
  - **Critical** in **pharmaceutical** and other **major industries** for **method transfers / consistency** and **long term** performance
  
- ◆ **Reversed-phase** is the **dominant** separation mode
  - **C18 & C8 = 60%**; **Phenyl = 16%**; **CN = 9.5%**; **Fluorinated = 5.9%**
  - **92% analysts** use **C18** at some time in their work...they typically meet the above criteria
  
  - **BUT limited** selectivity

# 16 Pharmaceutically Relevant Analytes – C18 Columns



**C18 phases show 'similar' selectivity...**

## The Challenge...

- ◆ To **engineer** new phases with **alternative** selectivity but with the robust properties of the **C18** ligand
  - **Reproducible** (column-to-column & batch-to-batch)
  - **Excellent column lifetime**
  - **Superb efficiency** provided by **ultra-inert, ultra-pure** silica particle
  - **Low MS bleed**
  - Usable in **100% aqueous** eluents
  
- ◆ **Available** for **HPLC & UHPLC** separations
  
- ◆ **Available** as a '**Phase III Ready**' product family
  - ◆ Globally available, supply chain, reproducible, multiple batches etc



## Aromatic Functionality – Engineering New Stationary Phases

- ◆ Phases with aromatic functionality include **phenyl** and pentafluorophenyl (**PFP**) based ligands
  
- ◆ Advantages
  - Aromatic functionality potentially offer **unique interactions** with analytes (c.f. C18) giving **alternative** selectivity
  - Provides **enhanced retention** of polar compounds
  - Many aromatic functionality-based phases can be used in **100% aqueous** eluents
  
- ◆ Disadvantages
  - Phenyl / PFP phases **may suffer** phase bleed
  - **Batch-to-batch reproducibility** & **robustness** may be weak


## Aromatic Functionality: $\pi$ – $\pi$ Interactions

- ◆ A type of electron **donor-acceptor** interaction
- ◆ Originates from  $\pi$  systems in **unsaturated functional groups** on analytes and the stationary phase
- ◆ Types of  $\pi$ - $\pi$  interaction can be **manipulated** for maximum effect (**orthogonality**) in phase design
  - eg phenyl: **electron rich** ring on the stationary phase also acts as  **$\pi$ -base** and interacts well with **electron poor analytes**
  - eg PFP: **electron poor** ring on the stationary phase also acts as  **$\pi$ -acid** and interacts well with **electron rich analytes**


# The Power of $\pi$ ...Scientific Led Stationary Phase Design



Electron Donating Groups  
eg  $NH_2$ ,  $NR_2$ , alkyl,  $OCH_3$   
 $OR$ ,  $CH_3$ , Ar etc




Electron **Rich** Ring




Activity:  $\pi$ -donor ( $\pi$ -base)

Electron Withdrawing Groups  
eg  $NO_2$ , halides,  $NR_3^+$ ,  $CO_2H$ ,  
 $CN$ ,  $CO_2R$ ,  $SO_3H$ ,  $COH$  etc




Electron **Deficient** Ring



Activity:  $\pi$ -acceptor ( $\pi$ -acid)

Classic  $\pi$ - $\pi$   
interaction



How do we exploit these properties for new stationary phases?



**C18+Phenyl = ACE® C18-AR**

**C18+PFP = ACE® C18-PFP**



## Uniquely Designed Stationary Phases

- ◆ **ACE<sup>®</sup> C18-AR (USP L1)**
  - Ligand has C18 **hydrophobic** element PLUS **phenyl** character
- ◆ **ACE<sup>®</sup> C18-PFP (USP L1)**
  - Ligand has C18 **hydrophobic** element PLUS **PFP** character
- ◆ **Ultra-inert, ultra-pure** silica particle technology as used in all ACE<sup>®</sup> products for **high** peak efficiency
- ◆ Available in **3, 5 & 10 $\mu$ m**, (ACE<sup>®</sup>) and **2 $\mu$ m** (ACE<sup>®</sup> *Exce/™*)

**Multi-mode interaction mechanisms** result in **enhanced** chromatographic **selectivity** giving the analyst **new options** for method development

## ACE<sup>®</sup> C18-AR: Multi-Mode Separation Mechanisms

- Combining the character of **C18+phenyl** into a single individual phase harnesses **the best** of both phases for **unique** selectivity


Separation mechanism	Typical C18	Typical Phenyl	ACE <sup>®</sup> C18-AR
Hydrophobicity	++++	+ / ++	++++
$\pi$ - $\pi$ Interaction	-	+++	+++
Dipole - Dipole	-	+	+
Hydrogen Bonding	-	++	++
Shape Selectivity	++	++	++ / +++

- The predominance of each retention mechanism will be dictated by the analyte's physicochemical properties, its structure and the chromatographic conditions applied


**Multi-Mode Interactions Offer the Chromatographer More**

# ACE® C18-AR Aromatic Selectivity

- Illustrating **hydrophobicity** and  **$\pi$ -base character / aromatic selectivity** with a **simple example** using substituted aromatics



$\pi$ -acidity (order)      1                      >                      2                      >                      3                      -



\* Predicted data from ACD Labs software, 30May12

## ACE<sup>®</sup> C18-PFP: Multi-Mode Separation Mechanism

- Combining the character of **C18+PFP** into a single individual phase harnesses **the best** of both phases for **unique** selectivity

Separation mechanism	Typical C18	Typical PFP	ACE <sup>®</sup> C18-PFP
Hydrophobicity	++++	+ / ++	++++
$\pi$ - $\pi$ Interaction	-	+++	+++
Dipole - Dipole	-	++++	++++
Hydrogen Bonding	-	+++	+++
Shape Selectivity	++	+++	++++


- The predominance of each retention mechanism will be dictated by the analyte's physicochemical properties, its structure and the chromatographic conditions applied

**Multi-Mode Interactions Offer the Chromatographer More**

# ACE<sup>®</sup> C18-PFP Selectivity\*

Peak Number:

1



1,2,3-  
TMB


Log P:

1.7

$\pi$ -basicity (order):

1

2




1,2,4-  
TMB

1.6

1

3




1,2-  
DMB

1.7

2

4




1,4-  
DMB

2.1

2

5




MB

2.2

3

6




1,3-  
DMB

2.2

2

7




1,3,5  
TMB

1.6

1

8




Tol

2.7

-

- ◆ Elution / retention **not simply** a function of  $\pi$ -basicity and Log P
- ◆ Retention mechanism for C18-PFP **multi-modal**

# ACE® C18-PFP Selectivity





- ◆ **C18 or PFP mechanisms alone not enough** to fully resolve the methoxybenzene **isomers**
- ◆ **ACE C18-PFP mechanism combines hydrophobicity, shape selectivity, dipole-dipole and  $\pi$ - $\pi$  interactions**
- ◆ **Elution order, retention and selectivity all seen to differ**
- ◆ **Powerful positional isomer and shape selectivity**


1) 1,2,3-trimethoxybenzene, 2) 1,2,4-trimethoxybenzene, 3) 1,2-dimethoxybenzene, 4) 1,4-dimethoxybenzene 5) methoxybenzene, 6) 1,3-dimethoxybenzene, 7) 1,3,5-trimethoxybenzene, 8) toluene (ref) Mobile phase 50:50 v/v MeOH / H<sub>2</sub>O; Column= 150 x 4.6 mm id; 1.00 ml/min; 40C; 254 nm

# ACE® Phase Comparisons With The Selectivity Descriptor\*

**Selectivity = 100 x  $\sqrt{1 - R^2}$**



**Selectivity = 100 x  $\sqrt{1 - R^2}$**   
**= 100 x  $\sqrt{1 - 0.9888}$**   
**= 10.6**



\* Neue, O’Gara, Méndez “Selectivity in Reversed-Phase Separations: Influence of the Stationary Phase”, J. Chromatogr. A 1127 (2006), 161-174

# Ranking ACE® Phase Orthogonality With MeOH and MeCN

- ◆ For the 102 acidic, basic and neutral analytes assessed

## MeOH

Column 1	Column 2	Selectivity 'S'
C18	C18-AR	12
C18	C18-PFP	11
C18-AR	C18-PFP	10

## MeCN

Column 1	Column 2	Selectivity 'S'
C18	C18-AR	8
C18-AR	C18-PFP	8
C18	C18-PFP	7

MeOH	MeCN	Selectivity Value
C18-PFP	C18	19
C18-AR	C18	18
C18-AR	C18-PFP	18
C18-PFP	C18-AR	18
C18-PFP	C18-PFP	18
C18	C18-AR	17
C18	C18-PFP	17
C18	C18	15
C18-AR	C18-AR	15

*Shows value of using the 3 phases in a 2 solvent screen for method development work*






## What Do I Use These Novel Phases For: ACE C18-PFP?

- ◆ Useful for analytes that contain **electron donating** moieties eg -NH<sub>2</sub>, -NR<sub>2</sub>, -OCH<sub>3</sub>, -OH, -alkyl, -Ar etc
- ◆ eg nucleotides, nucleosides, nucleobases, halogenated aryl / aromatics, catecholamines, tetracyclines, beta blockers, structural isomers, coumarins etc
- ◆ **Excellent** shape and **positional isomer** selectivity

# What Do I Use These Novel Phases For: ACE C18-PFP?

- ◆ **C18-PFP: chloroacetophenone halogenated isomers separation**




**Regioisomer  
Selectivity**

## What Do I Use These Novel Phases For: ACE C18-AR?

- ◆ Useful for analytes that contain **electron withdrawing** moieties eg -NO<sub>2</sub>, -halides, -NR<sub>3</sub><sup>+</sup>, -SO<sub>2</sub>, -CO<sub>2</sub>H, -SO<sub>3</sub>H, -CO<sub>2</sub>R, -CHO etc
- ◆ eg aromatic compounds, anthocyanins, steroids, analgesics, phenolics, water soluble vitamins, sulphur containing compounds, quinolones, positional isomers etc
- ◆ **Moderate** shape selectivity

# What Do I Use These Novel Phases For: ACE C18-AR?

## ◆ C18-AR: steroids separation




## Combining Selectivity With The NEW ACE<sup>®</sup> *Excel*<sup>™</sup> Format

- ◆ **NEW high efficiency, ultra-inert** 2 $\mu$ m silica particles suitable for UHPLC at 1000bar (15000psi)
- ◆ **Nine selectivities** – including the **unique C18-AR** and **C18-PFP**
- ◆ **High reproducibility**: column-to-column and batch-to-batch
- ◆ **Ultra-robust** phases: **NEW** low dispersion column **hardware** and **NEW** High Stability Column (**HSC<sup>™</sup>**) packing technology
- ◆ Engineered with **lower back pressures** compared to other <2 $\mu$ m phases due to **2 $\mu$ m** particle size and **frit technology**
- ◆ **Fully scalable** to ACE<sup>®</sup> 3 $\mu$ m, 5 $\mu$ m and 10 $\mu$ m phases
- ◆ **Fully compatible** with all commercial HPLC and **UHPLC** kit

# ACE® Excel™ UHPLC Column Robustness


- ◆ **1000 bar for ~2000 gradient runs**
  - ◆ Isocratic efficiency assessments every ~100 runs ← more demanding!



100x2.1mm; MPA 0.1% FA (aq); MPB: 0.1% FA in MeOH; 0.73mL/min; gradient: 20-90%B in 6 mins.

**NEW High Stability Column (HSC™) Packing Technology**  
**Significantly Improves UHPLC Column Robustness**

# ACE<sup>®</sup> Excel<sup>™</sup> UHPLC Columns – Scalability & Reproducibility



**UHPLC**

**HPLC**


**Prep LC**



# ACE<sup>®</sup> Excel<sup>™</sup> Has Typically Lower Back Pressure For UHPLC


◆ Specifically engineered for **lower UHPLC** backpressures

ACE Excel 2 C18-AR




**P<sub>max</sub>: 364 bar**

Waters Acquity 1.7um BEH C18




**P<sub>max</sub>: 581 bar**

Phenomenex Kinetex 1.7um C18



**P<sub>max</sub>: 540 bar**

Agilent Zorbax Eclipse 1.8 um XDB C18



**P<sub>max</sub>: 540 bar**

Conditions: A=5mM formic acid (aq); B=5mM formic acid in MeOH; tg= 3 to 100%B in 5 min; 0.6 ml/min; 40C; 254nm

1. Paracetamol 2. Hydrochlorothiazide 3. Methylphenylsulphoxide 4. Methylphenylsulphone 5. Aspirin 6. Phenacetin 7. 1,3-dinitrobenzene 8. 1,2,4-trimethoxybenzene 9. Ethylbenzoate 10. Nimesulide 11. Ibuprofen 12. Indomethacin 13. Mefenamic acid

All trademarks are recognised...comparative separations may not be representative of all applications




**Selectivity, Speed & Scaling**  
**Isocratic & Gradient HPLC → UHPLC**

# ACE<sup>®</sup> Excel<sup>™</sup> C18-PFP Selectivity & Throughput (Isocratic)


**Aim: obtain  $R_s \geq 1.7$  in shortest possible time for mixture**

**Waters XBridge 5 $\mu$ m C18**  
 150 x 4.6 mm  
 1.00 ml/min  
 163 bar



**HPLC: 5 $\mu$ m C18**  
13.6 mins  
 ( $L/d_p = 3$ )

**Waters Acquity 1.7 $\mu$ m BEH C18**  
 50 x 2.1 mm  
 0.21 ml/min  
 246 bar




**UPLC: <2 $\mu$ m C18**  
4.75 mins  
 ( $L/d_p = 2.9$ )

To maintain  $R_s$  and reduce run time, keep  $L / d_p$  ratio constant

**~ x23 Quicker**

**ACE Excel 2 $\mu$ m C18-PFP**  
 30 x 2.1 mm  
 1.30 ml/min  
 492 bar



**UHPLC: ACE Excel 2 $\mu$ m C18-PFP**  
< 1 min


Using UHPLC and selectivity, it is possible to dramatically improve resolution allowing shorter columns & increased flow rates

Sample: 1) 1,2-dimethoxybenzene, 2) 1,3-dimethoxybenzene, 3) 1,3,5-trimethoxybenzene, 4) toluene (reference).  
 Mobile phase 50:50 MeOH / H<sub>2</sub>O; Temperature 40°C; 254 nm


# ACE<sup>®</sup> Excel<sup>™</sup> C18-PFP Selectivity & Throughput (Gradient)

**Aim: obtain  $R_s \geq 1.7$  in shortest possible time for mixture**


**ACE 5 $\mu$ m C18**  
 100 x 4.6 mm  
 1 ml/min,  $t_G = 29$  min  
 max pressure: 92 bar  
 40 min cycle time



**ACE Excel 2 $\mu$ m C18**  
 50 x 2.1 mm  
 0.6 ml/min,  $t_G = 5$  min  
 max pressure: 367 bar  
 9 min cycle time




**ACE Excel 2 $\mu$ m C18-PFP**  
 30 x 2.1 mm  
 2.5 ml/min,  $t_G = 0.7$  min  
 max pressure: 914 bar  
 1 min cycle time



**~ x25 Quicker**

**$R_s \geq 1.7$**




1, aspirin; 2, phenacetin; 3, 1,3-dinitrobenzene ; 4, ethylbenzoate; 5, nimesulide; 6, ibuprofen; 7, indomethacin.

**Pressure Effects**  
**HPLC ↔ UHPLC**

## Background

- ◆ Pressure is a **complex** physical parameter that affects **many elements** of a chromatography system
- ◆ Chromatographic **selectivity** and **retention** changes at **elevated pressures** have been investigated and reported<sup>a</sup>
- ◆ Observations are **highly dependent** upon the analytes and may be seen with **any manufacturer** phases operated under UHPLC conditions
- ◆ Changes are typically **not helpful** for HPLC ↔ UHPLC activities

# Effect of Pressure on Selectivity and Retention Factor



- ◆ Initial 2 $\mu$ m and 3 $\mu$ m data are similar (A, B)
  - ◆ Scalability looks good
  
- ◆ Retention and selectivity seen to change with pressure (B $\rightarrow$ E)

Agilent 1290, 50 x 2.1 mm (constant flow and restrictor capillary used)  
 Mobile phase: A=0.1% FA in water; B=0.1% FA in MeOH (51:49 v/v)  
 Flow Rate: 0.21 ml/min, Temperature: 40 °C  
 K= Ketoprofen; S= Sulindac; N=Naproxen

## Summary: Unwanted Selectivity Changes

- ◆ Pressure **induced  $k$**  and  **$\alpha$**  changes may be seen for **any manufacturer** phases under UHPLC conditions
- ◆ Changes in selectivity and retention may be **significant** with **ionised analytes** and **large MW** analytes<sup>a</sup>, but the impact on neutral molecules is **typically smaller**
- ◆ **Current** discussions / theory focus on changes in analyte **molar volume** as the **principle** cause for **changes** in  **$k$**  and  **$\alpha$**  observed
- ◆ Successful HPLC ↔ UHPLC **possible**...the analyst just needs to be **vigilant**

**Connections : Losses in  $N$  and  $A_s$**   
**Peak Dispersion**





## Background

- ◆ **UHPLC / optimised HPLC** instruments are **very sensitive** to the introduction of **extra column volume**
- ◆ Any time you **install** a column (from **any manufacturer**) it is vital to ensure **good** connections
- ◆ **Aim** for a '**fresh connection**' every time to ensure a **snug fit** between tubing and column and **reduce** the likelihood of an **unwanted** gap and / or tubing slippage
- ◆ **Free movement** of the **ferrule** and **nut** when installing the column gives you a **fresh connection**

# Losses in Performance Due to Incorrect Column Fitting



## Correctly fitted column



$N_{5\%h} = 12,800$  plates  
 $PW_{5\%h} = 0.291$  min  
 $A_s = 0.94$

◆ **Correctly fitted columns make the most of your column and system**

## Incorrectly fitted column




$N_{5\%h} = 9,850$  plates  
 $PW_{5\%h} = 0.335$  min  
 $A_s = 0.84$

◆ **Incorrectly connected columns lead to reduced efficiency, reduced asymmetry and possibly leaks**

◆ **Loss of ~23% for  $N$**   
**Loss of ~11% for  $A_s$**

# Summary: Column Connections

- ◆ Extra column volume **reduces** peak efficiency and asymmetry
- ◆ Make a **fresh connection every time** you install **any** column
- ◆ ACE recommend **reusable fittings** for a **fresh connection every time**
- ◆ All ACE® *Exce/™* columns have a **FREE** ‘Making Great UHPLC Connections’ leaflet in **every** box



Also downloadable from the ACE website:

[www.ace-hplc.com](http://www.ace-hplc.com)

## Overall Summary & Conclusions


- ◆ Understanding the **properties** of building blocks in stationary phase design led to these **unique ACE<sup>®</sup>** products
- ◆ **ACE<sup>®</sup> C18-AR** and **ACE<sup>®</sup> C18-PFP** are powerful tools for method development due to **unique** but **complementary** selectivities
- ◆ These **unique** phases are available for HPLC as the **ACE<sup>®</sup>** range and also UHPLC as the **NEW ACE<sup>®</sup> Excel<sup>™</sup> 2 μm** format
- ◆ These phases **meet** analyst **demands** of **reproducibility**, **robustness** & **low** phase **bleed** with **excellent** peak efficiency
- ◆ Operating at **high pressures** can deliver **excellent results** but remain vigilant - **selectivity** and **retention** may be affected...and even **column connections** become critical!

# Full Information On All ACE Products Available

## Unique Selectivities



## Free Guides



**MACMOD Analytical = <http://www.mac-mod.com/>  
ACT = <http://www.ace-hplc.com>**



**Thank You For Your Attention**

[amckeown@ace-hplc.com](mailto:amckeown@ace-hplc.com)



[www.ace-hplc.com](http://www.ace-hplc.com)  
[www.mac-mod.com](http://www.mac-mod.com)



ACE® is a registered trademark of Advanced Chromatography Technologies Ltd. ACE Excel™ and HSC™ are trademarks of Advanced Chromatography Technologies Ltd.

UPLC, Xbridge, ACQUITY are trademarks of Waters Corporation; ZORBAX, Eclipse are trademarks of Agilent Technologies Inc.; Kinetex is a trademark of Phenomenex Inc.; GOLD is a trademark of Thermo Fisher Scientific.