



SMARTER CHROMATOGRAPHY

# FOOD & BEVERAGE GUIDE

INFORMATION & RESEARCH COLLECTED IN COLLABORATION WITH:



# CONTENTS

## ALCOHOLS

|  |   |
|--|---|
| Analysis of Whiskey via Direct Split Injection — GC.....                   | 4 |
| Beer Fermentation Analysis — HPLC/UHPLC .....                              | 5 |
| Amino Acids in Peas ( <i>Pisum sativum</i> ) by HPLC-HRAM-MS — LC-MS ..... | 6 |

## AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

|   |    |
|---|----|
| Free Acids and Lactic Acid — GC.....                      | 8  |
| Hop Acids — HPLC .....                                    | 9  |
| LC-MS Separation of Corn Oil — LC-MS .....                | 10 |
| Organic Acids (I) — HPLC .....                            | 12 |
| Phenolic Acids on a 2 µm RP-Amide Column — UHPLC .....    | 13 |
| Polar Organic Acids on an AQ-C18 Column — HPLC/UHPLC..... | 15 |
| Wine Acid Analysis — HPLC.....                            | 16 |
| Anthocyanins in Blueberries — HPLC/UHPLC .....            | 17 |

## ANTHOCYANINS

|                                       |    |
|---------------------------------------|----|
| Aminoglycosides in Eggs — UHPLC ..... | 18 |
|---------------------------------------|----|

## ANTIBIOTICS

|                                |    |
|--------------------------------|----|
| Green Tea Extract — UHPLC..... | 19 |
|--------------------------------|----|

## ANTIOXIDANTS

|   |    |
|---|----|
| HPLC Analysis of Chlorogenic Acid in Green Coffee Extract — HPLC/UHPLC..... | 20 |
| Flavor and Fragrance Compounds — GC .....                                   | 21 |

## AROMATIC COMPOUNDS

|  |    |
|--|----|
| LC-MS of PAHs in Grilled Meat — HPLC/UHPLC ..... | 22 |
| Artificial Food Colors — HPLC .....              | 24 |

## ARTIFICIAL COLORS, DYES AND PRESERVATIVES

|                                    |    |
|------------------------------------|----|
| Illegal Dyes in Spices — HPLC..... | 25 |
| Preservatives (I) — UHPLC.....     | 26 |
| Capsaicinoids — UHPLC.....         | 27 |

## CAPSAICINS

|   |    |
|---|----|
| Capsaicins in Chili Powder — HPLC/UHPLC .....                                       | 28 |
| Curcumin Analysis in Turmeric on a 1.5 mm Stainless Steel Column — HPLC/UHPLC ..... | 29 |

## CURCUMINS

|   |    |
|---|----|
| Common Catechins Found in Tea via LC-MS — LC-MS ..... | 30 |
|---|----|

## FLAVONOIDS AND OTHER POLYPHENOLS

|  |    |
|--|----|
| Cyclodextrin-Encapsulated Flavor Compounds in Beer — HPLC..... | 31 |
| Five Flavonoids — UHPLC .....                                  | 33 |
| Four Flavonoids — HPLC.....                                    | 34 |
| Vanillins — HPLC .....   | 35 |
| Isoflavones — HPLC.....  | 36 |
| Polyphenols in Wine — HPLC/UHPLC.....                          | 37 |
| Six Flavonoids — HPLC/UHPLC.....                               | 38 |
| Food Additives — HPLC/UHPLC .....                              | 39 |

## FOOD ADDITIVES

|   |    |
|---|----|
| Food Additives Assay — HPLC/UHPLC .....             | 40 |
| Determination of Sterols in Sunflower Oil — GC..... | 41 |

## LIPIDS

|   |    |
|---|----|
| FAMEs and Conjugated Linoleic Acid (OLA) isomers — GC ..... | 42 |
| Steroids Spiked in Ground Beef — HPLC/UHPLC .....           | 43 |
| 250 Pesticide Screen using LC-MS/MS — LC-MS ....            | 44 |

## MYCOTOXINS, PESTICIDES AND TOXINS

|   |    |
|---|----|
| LC-MS Screening Comparison of Mycotoxins in Craft and Home Brewed Beers — HPLC/UHPLC..... | 49 |
| Mycotoxins by LC-MS/MS — LC-MS .....  | 51 |
| Paralytic Shellfish Poisoning (PSP)<br>Toxins — HPLC.....                                 | 52 |
| Patulin and HMF — HPLC/UHPLC.....   | 53 |
| Pesticide Screening of Barley — UHPLC.....  | 54 |
| Whey Proteins from Whole Milk — HPLC .....  | 55 |

## PROTEINS

|                                 |    |
|---------------------------------|----|
| Garlic Analysis (I) — HPLC..... | 56 |
|---------------------------------|----|

## SULFUR COMPOUNDS: GARLIC

|   |    |
|---|----|
| Garlic Analysis (II) — HPLC.....              | 57 |
| Additives and Intense Sweeteners — HPLC ..... | 58 |

## SWEETENERS

|  |    |
|--|----|
| Analysis of Sucrose and Lactose in Chocolate Using a Penta-HILIC Column — HPLC/UHPLC ..... | 59 |
| Analysis of Sugars in Juice — HPLC/UHPLC .....   | 60 |
| Analysis of Sugars in Pure Honey — HPLC/UHPLC.....   | 61 |
| Caffeine, Sweeteners and Preservatives — HPLC ..   | 62 |
| Fast Separation of Oligosaccharides — HPLC/UHPLC .....                                     | 64 |
| High Resolution Separation of Oligosaccharides — HPLC/UHPLC .....                          | 65 |
| LC-MS Analysis of Stevia Extract — LC-MS .....   | 66 |
| Sugars – Orange Juice — HPLC .....   | 67 |
| Sugars Separation — HPLC.....  | 68 |
| Carotenoid Analysis in Pumpkin — HPLC/UHPLC .....  | 69 |

## VITAMINS

|  |    |
|--|----|
| Carotenoids Extracted From Carrot Juice — HPLC/UHPLC ..... | 70 |
| Carotenoids on a C30 Column — HPLC/UHPLC.....              | 71 |
| Tocopherols — HPLC .....                                   | 72 |
| Vitamin K — UHPLC .....                                    | 73 |
| Vitamins - Water Soluble (Gradient) C18 — HPLC..           | 74 |
| Vitamins - Water Soluble (Gradient) C8 — HPLC.....         | 75 |
| Vitamins – Fat Soluble — HPLC .....                        | 76 |
| Vitamins D2 and D3 — UHPLC .....                           | 78 |
| Water Soluble Vitamins (Polar Molecules) — HPLC.           | 79 |
| Caffeine and Metabolites — HPLC .....                      | 80 |

## XANTHINES

|   |    |
|---|----|
| Caffeine and Related Compounds Using HILIC Mode — HPLC..... | 81 |
| Caffeine in Soda — HPLC.....                                | 82 |


# ALCOHOLS

## Analysis of Whiskey via Direct Split Injection—GC

### TEST CONDITIONS

Column: Avantor® Hichrom HI-624, 1.4 µm, 0.25 mm x 60 m  
Part Number: [HI082514060](#)  
Oven Program: 35 °C (5 min), 1 °C /min, 50 °C, 5 °C/min  
Carrier Gas: Hydrogen, 120 kPa  
Injector: Split, 200 °C, 120 mL/min split flow, 0.5 µL injected volume  
Detector: FID, 250 °C

The HI-624 is a versatile general-purpose column for a wide range of applications, including residual solvents and volatile organic compounds. It is a mid-polarity phase with a 6% cyanopropylphenyl, 94% methyl polysiloxane composition and is classified as USP G43.



### ANALYTES

1. Acetaldehyde
2. Methanol
3. Ethanol
4. Acetone
5. n-Propanol
6. Ethyl Acetate
7. Isobutanol
8. Isopentanol
9. 2-Methyl-1-butanol

# ALCOHOLS

## Beer Fermentation Analysis — HPLC/UHPLC

### TEST CONDITIONS


Column: HALO® Penta-HILIC, 90 Å, 2.7 µm, 3.0 x 50 mm  
Part Number: 92813-405  
Mobile Phase A: Water  
Mobile Phase B: Acetonitrile

| Gradient: | Time (min) | %B |
|-----------|------------|----|
|           | 0.0        | 92 |
|           | 8.0        | 52 |

Flow Rate: 0.75 mL/min  
Temperature: 65 °C  
Detection: ELSD, 40 °C, 45 psi  
Injection Volume: 2 µL  
Data Rate: 10 Hz, 2 sec filter

Data Courtesy of Merlin K. L. Bicking, Ph. D. (ACCTA, Inc.)

A Belgian ale is analyzed with a HALO® Penta-HILIC column using an evaporative light scattering detector (ELSD). Sugars, oligosaccharides, and polysaccharide levels are monitored throughout the fermentation process in order to track yeast behavior. These levels will decrease over time as the yeast converts the sugars to ethanol. The Penta-HILIC/ ELSD combination is a great way to perform rapid sugar analysis providing high resolution and good peak shape at elevated temperatures.




# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Amino Acids in Peas (*Pisum sativum*) by HPLC-HRAM-MS—LC-MS


### TEST CONDITIONS

|               |  |              |                                      |
|---------------|--|--------------|--------------------------------------|
| Column:       | Avantor® ACE® AQ, 100 Å, 3 µm, 3.0 x 150 mm                            | Flow Rate:   | 0.4 mL/min                           |
| Part Number:  | ACE1161503   | Temperature: | 30 °C                                |
| Mobile Phase: | A: 0.1% Formic acid in H <sub>2</sub> O<br>B: 0.1% Formic acid in MeCN | Injection:   | 5 µL                                 |
| Gradient:     | Time (min)   | %B           | Detection:                           |
|               | 0  | 0            | Exactive Orbitrap high resolution MS |
|               | 10   | 100          | ESI positive ion mode                |
|               |  |              | Capillary temperature: 350 °C        |

Reproduced with permission of Fera Science Ltd, York, UK



# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS




# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Free Acids and Lactic Acid – GC

### TEST CONDITIONS

Column: Avantor® Hichrom HI-FFAP EXT, 0.5 µm, 0.32 mm x 30 m  
Part Number: [HI443205030](#)  
Oven Program: 70 °C (0.5 min), 5 °C/min, 242.5 °C  
Carrier Gas: Helium, 90 kPa  
Injector: Split, 250 °C  
Detector: FID, 250 °C

The HI-FFAP EXT is a nitroterephthalic acid modified polyethylene glycol (PEG) phase of high polarity, with an extended temperature range and aqueous sample compatibility. Useful for the analysis of free acids, flavor compounds, alcohols and polar compounds.



### ANALYTES

1. Acetic acid
2. Propionic acid
3. Iso-butyric acid
4. Pivalic acid
5. Butyric acid
6. Iso-valeric acid
7. Caproic acid
8. Caprylic acid
9. Nonanoic acid
10. Lactic acid
11. Decanoic acid
12. Dodecanoic acid

# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Hop Acids—HPLC

### TEST CONDITIONS


Columns: HALO® Biphenyl, 90 Å, 5 µm, 4.6 x 150 mm  
Part Number: 95814-711  
Mobile Phase A: Water, 0.1% Formic acid  
Mobile Phase B: Acetonitrile, 0.1% Formic acid  
Gradient:

| Time (min) | %B |
|------------|----|
| 0.0        | 60 |
| 3.0        | 60 |
| 6.0        | 80 |


Flow Rate: 2.0 mL/min  
Initial Pressure: 236 bar  
Temperature: 30 °C  
Detection: 270 nm, PDA  
Injection Volume: 5 µL  
Sample Solvent: Acetonitrile  
Data Rate: 100 Hz  
Response Time: 0.025 sec  
Flow Cell: 1 µL  
LC System: Shimadzu Nexera X2

Hops are primarily made up of essential oils and alpha and beta acids. They have many benefits in the beer brewing process, including their antiseptic nature and bitterness flavor they give to the beer. Alpha and beta acids from the International Calibration Standard Extract (ICE-3) are separated on a HALO® Biphenyl column.


### STRUCTURES




Cohumulone




Humulone




Adhumulone



Colupulone




Lupulone



Adlupulone

### ANALYTES

Alpha acids      Beta acids  
1. Cohumulone    4. Colupulone  
2. Humulone       5. Lupulone  
3. Adhumulone     6. Adlupulone



# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## LC-MS Separation of Corn Oil — LC-MS

### LC TEST CONDITIONS

Columns: HALO® C18, 90 Å, 2.7 µm, 2.1 x 150 mm  
Part Number: 92812-702  
Columns: HALO® C30, 160 Å, 2.7 µm, 2.1 x 150 mm  
Part Number: 92112-730  
Mobile Phase A: Methanol  
Mobile Phase B: IPA/0.1% Formic acid

| Gradient: | Time (min) | %B  |
|-----------|------------|-----|
|           | 0.00       | 10  |
|           | 10.00      | 10  |
|           | 14.00      | 40  |
|           | 22.00      | 40  |
|           | 22.01      | 10  |
|           | 24.00      | END |


Flow Rate: 0.3 mL/min  
Initial Pressure: 325 bar  
Temperature: Ambient  
Injection Volume: 2 µL  
Sample Solvent: MeOH  
LC System: Shimadzu Nexera X2

Corn oil, composed mainly of long chain fatty acids and esters, is an edible oil which comprises approximately 5-10% of edible oil consumption. In recent years, corn oil has been used in biodiesel, pharmaceutical, and cosmetic applications as well. The use of a C18 column for the analysis of edible oils is difficult due to the high concentration of hydrophobic triglycerides (TAGs); therefore, the C30 phase has seen increased application in this area. Here we show a comparison between the C18 and C30 phase, and demonstrate that the 2.7 µm HALO® C30 is an ideal choice for the separation and resolution of high mass triglycerides found in edible oils such as corn oil. C30 offers superior specificity compared to C18 columns by exhibiting higher shape selectivity, enabling better separation of hydrophobic, long-chain, structures.


### MS TEST CONDITIONS

MS system: Shimadzu LCMS-2020  
Ionization: +ESI  
Spray voltage: 4.50 kV  
Drying line temp: 300 °C  
Heat Block: 450 °C

### STRUCTURES



DAGs





TAGs

# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## LC-MS Separation of Corn Oil

(continued)



DAGs = diacylglycerols or diglycerides  
TAGs = triacylglycerols or triglycerides






# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Organic Acids (I) — HPLC

### TEST CONDITIONS


Column: Avantor® ACE® Excel C18-Amide, 100 Å, 3 µm, 2.1 x 250 mm  
Part Number: EXL11122502U  
Mobile Phase: 20 mM H<sub>3</sub>PO<sub>4</sub> in MeOH/H<sub>2</sub>O (40:60 v/v)  
Flow Rate: 0.21 mL/min  
Temperature: 20 °C  
Injection Volume: 5 µL  
Detection: UV, 210 nm

### STRUCTURES



### ANALYTES

1. Sinapic acid
2. Ferulic acid
3. 3,4-Dimethoxycinnamic acid
4. Cinnamic acid
5. 4-Methoxycinnamic acid



# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Phenolic Acids on a 2 µm RP-Amide Column — UHPLC

### TEST CONDITIONS

Column: HALO® RP-Amide, 90 Å, 2 µm, 2.1 x 100 mm

Part Number: 91812-607

Mobile Phase A: 20mM Phosphoric acid

Mobile Phase B: Methanol

Gradient:

| Time (min) | %B |
|------------|----|
| 0.00       | 30 |
| 3.75       | 60 |
| 4.25       | 60 |

Flow Rate: 0.5 mL/min

Initial Pressure: 716 bar

Temperature: 35 °C

Detection: UV 220 nm, PDA

Injection Volume: 0.5 µL

Sample Solvent: Methanol

Data Rate: 40 Hz

Response Time: 0.025 sec


Flow Cell: 1 µL

LC System: Shimadzu Nexera X2

Phenolic acids can be found in many plant-based foods and beverages. Fruits, vegetables, and even olive oils all contain different varieties of these acids. For example, sinapic acid can be found in wine and caffeic acid can be found in coffee. These compounds act as antioxidants and can also help anti-inflammatory conditions in the body. They also affect the flavors of the food or oil.

### ANALYTES


1. Homovanillic acid
2. Caffeic acid
3. Syringic acid
4. Vanillic acid
5. Chlorogenic acid
6. Sinapic acid
7. Ferulic acid
8. *p*- Coumaric acid
9. *trans*- Cinnamic acid




# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Phenolic Acids on a 2 µm RP-Amide Column (*continued*)


### STRUCTURES OF PHENOLIC ACIDS




Homovanillic acid




Caffeic acid




Syringic acid




Vanillic acid




Chlorogenic acid




Sinapic acid




Ferulic acid



p- Coumaric acid



trans- Cinnamic acid



Resveratrol


# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Polar Organic Acids on an AQ-C18 Column—HPLC/UHPLC


### TEST CONDITIONS

Column: HALO® AQ-C18, 90 Å, 2.7 µm, 4.6 x 250 mm  
Part Number: 92814-922  
Isocratic: 20 mM Potassium phosphate buffer pH: 2.7  
Flow Rate: 1.0 mL/min  
Pressure: 307 bar  
Temperature: 40 °C  
Detection: UV 214 nm, PDA  
Injection Volume: 20 µL  
Sample Solvent: Mobile phase  
Data Rate: 100 Hz  
Response Time: 0.025 sec  
Flow Cell: 1 µL  
LC System: Shimadzu Nexera X2


### STRUCTURES




Oxalic acid




Tartaric acid



Malic acid




Ascorbic acid




L-Lactic acid




Acetic acid



Citric acid




Succinic acid



Fumaric acid

Organic acids are common in the food and beverage industry and can be found in many sample types such as fruits, vegetables, and wines. This separation of nine polar organic acids is performed on a HALO® AQ-C18 column using 100% aqueous mobile phase at low pH. The 250 mm column length was chosen to provide excellent resolution with reasonable run time for this polar mixture.



### ANALYTES

1. Oxalic acid
2. Tartaric acid
3. Malic acid
4. Ascorbic acid
5. L-Lactic acid
6. Acetic acid
7. Citric acid
8. Succinic acid
9. Fumaric acid


# AMINO, BITTER, FATTY, ORGANIC, PHENOLIC AND POLAR ACIDS

## Wine Acid Analysis—HPLC


### TEST CONDITIONS

Column: Avantor® ACE® Excel C18-Amide, 100 Å, 3 µm, 2.1 x 250 mm  
Part Number: EXL11122502U  
Mobile Phase: 40 mM Ammonium phosphate pH 2.5 in H<sub>2</sub>O  
Flow Rate: 0.21 mL/min  
Temperature: 25 °C  
Injection: 5 µL  
Detection: UV, 214 nm


### STRUCTURES




Oxalic acid




Tartaric acid




Malic acid




Lactic acid



Ascorbic acid



Citric acid



### ANALYTES

1. Oxalic acid
2. Tartaric acid
3. Malic acid
4. Lactic acid
5. Ascorbic acid
6. Citric acid

# ANTHOCYANINS

## Anthocyanins in Blueberries — HPLC/UHPLC

### TEST CONDITIONS

Column: HALO® LPH-C18, 90 Å, 2.7 µm, 2.1 x 100 mm

Part Number: 92822-616

Mobile Phase A: Water/ 3% Phosphoric acid (pH: 1.4)

Mobile Phase B: Methanol

| Gradient: | Time (min) | %B |
|-----------|------------|----|
|           | 0.0        | 23 |
|           | 10.8       | 26 |
|           | 29.8       | 60 |

Flow Rate: 0.27 mL/min

Pressure: 144 bar

Temperature: 30 °C

Detection: UV 525 nm, PDA

Injection Volume: 4.5 µL

Sample Solvent: Water


Data Rate: 100 Hz

Response Time: 0.025 sec

Flow Cell: 1 µL

LC System: Shimadzu Nexera X2

Anthocyanins, a category of polyphenols, are a type of pigment found in plants that offer several health benefits. These flavonoids have antioxidant effects that can be found in a variety of different fruits and vegetables, including blueberries. A separation of anthocyanins is performed on a HALO® 90 Å LPH-C18 column, which is ideal for the low pH conditions of this method. Blueberries were mixed with methanol, water, and formic acid followed by being spun down and filtered. Due to the sterically protected ligand, the LPH-C18 column can withstand these low pH conditions and maintain stable retention times while standard C18 columns will show retention loss over time.




# ANTIBIOTICS

## Aminoglycosides in Eggs — UHPLC

### TEST CONDITIONS

| Columns:      | Avantor® ACE® Excel C18-PFP, 100 Å, 2 µm, 2.1 x 100 mm   | Flow Rate:   | 0.4 mL/min                              |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
|---------------|--|--------------|---|------------|----|--|-----|-----|--|-----|------|--|-----|------|--|-----|------|--|-----|------|--|-----|------|--|-----|------|--|-----|-----|--|------|-----|--|
| Part Number:  | EXL10101002U   | Temperature: | 40 °C                                   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| Mobile Phase: | A: 20 mM HFBA in H <sub>2</sub> O/MeCN (98:2 v/v)<br>B: 20 mM HFBA in MeCN/H <sub>2</sub> O (98:2 v/v)   | Detection:   | Positive ESI MRM (transitions as shown) |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| Gradient:     | <table><thead><tr><th>Time (min)</th><th>%B</th><th></th></tr></thead><tbody><tr><td>0.0</td><td>5.0</td><td></td></tr><tr><td>2.0</td><td>15.0</td><td></td></tr><tr><td>4.5</td><td>19.0</td><td></td></tr><tr><td>5.5</td><td>19.5</td><td></td></tr><tr><td>6.0</td><td>22.0</td><td></td></tr><tr><td>7.0</td><td>35.0</td><td></td></tr><tr><td>9.0</td><td>48.0</td><td></td></tr><tr><td>9.5</td><td>5.0</td><td></td></tr><tr><td>12.0</td><td>5.0</td><td></td></tr></tbody></table> |              |   | Time (min) | %B |  | 0.0 | 5.0 |  | 2.0 | 15.0 |  | 4.5 | 19.0 |  | 5.5 | 19.5 |  | 6.0 | 22.0 |  | 7.0 | 35.0 |  | 9.0 | 48.0 |  | 9.5 | 5.0 |  | 12.0 | 5.0 |  |
| Time (min)    | %B   |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 0.0           | 5.0  |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 2.0           | 15.0   |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 4.5           | 19.0   |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 5.5           | 19.5   |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 6.0           | 22.0   |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 7.0           | 35.0   |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 9.0           | 48.0   |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 9.5           | 5.0  |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| 12.0          | 5.0  |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |
| Sample:       | Extraction at low pH, clean up with WCX SPE cartridge<br>Egg sample spiked at 100 ug/kg (CCa)  |              |   |            |    |  |     |     |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |      |  |     |     |  |      |     |  |

Reproduced with permission of Veterinary Drug Residues Department, Centro Nacional de Alimentación (CNA, AECOSAN), Madrid, Spain




GNT Gentamicin  
TBM Tobramycin  
KNM Kanamycin  
DHS Dihydrostreptomycin  
STR Streptomycin  
HGB Hygromycin-B  
SPT Spectinomycin

## Green Tea Extract — UHPLC

## TEST CONDITIONS

| Column:       | Avantor® ACE® Excel SuperC18, 100 Å, 1.7 µm, 3.0 x 100 mm  |  | Flow Rate:   | 0.8 mL/min   |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
|---------------|--|--|--------------|--|---|---|----|----|----|----|----|----|----|---|----|---|------------|------------|
| Part Number:  | <u>EXL17111003U</u>  |  | Temperature: | 80 °C  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| Mobile Phase: | A: 0.1% Formic acid in H <sub>2</sub> O<br>B: 0.1% Formic acid in MeCN   |  | Injection:   | 2 µL   |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| Gradient:     | <table border="1"> <thead> <tr> <th>Time (min)</th> <th>%B</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>5</td> </tr> <tr> <td>15</td> <td>20</td> </tr> <tr> <td>17</td> <td>95</td> </tr> <tr> <td>18</td> <td>95</td> </tr> <tr> <td>20</td> <td>5</td> </tr> <tr> <td>27</td> <td>5</td> </tr> </tbody> </table> |  | Time (min)   | %B   | 0 | 5 | 15 | 20 | 17 | 95 | 18 | 95 | 20 | 5 | 27 | 5 | Detection: | UV, 260 nm |
| Time (min)    | %B   |  |              |  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| 0             | 5  |  |              |  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| 15            | 20   |  |              |  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| 17            | 95   |  |              |  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| 18            | 95   |  |              |  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| 20            | 5  |  |              |  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
| 27            | 5  |  |              |  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
|               |  |  | Sample:      | Tablet ground to fine powder and extracted with MeCN/H <sub>2</sub> O (1:1 v/v) with ultrasonication. Supernatant diluted with H <sub>2</sub> O and filtered using Whatman Mini-Uniprep syringless filter. |   |   |    |    |    |    |    |    |    |   |    |   |            |            |
|               |  |  | System:      | Chromaster Ultra Rs  |   |   |    |    |    |    |    |    |    |   |    |   |            |            |




## HPLC Analysis of Chlorogenic Acid in Green Coffee Extract—HPLC/UHPLC

### TEST CONDITIONS

| Column:           | HALO® C18, 90 Å, 2.7 µm, 3.0 x 100 mm  |            |    |     |    |     |    |     |    |      |     |      |     |
|-------------------|--|------------|----|-----|----|-----|----|-----|----|------|-----|------|-----|
| Part Number:      | 92813-602  |            |    |     |    |     |    |     |    |      |     |      |     |
| Mobile Phase:     | A/B  |            |    |     |    |     |    |     |    |      |     |      |     |
|                   | A= Water + 0.1% Formic acid  |            |    |     |    |     |    |     |    |      |     |      |     |
|                   | B= Acetonitrile + 0.1% Formic acid   |            |    |     |    |     |    |     |    |      |     |      |     |
| Gradient:         | <table border="1"> <thead> <tr> <th>Time (min)</th> <th>%B</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>10</td> </tr> <tr> <td>4.0</td> <td>10</td> </tr> <tr> <td>9.0</td> <td>50</td> </tr> <tr> <td>11.0</td> <td>100</td> </tr> <tr> <td>13.0</td> <td>100</td> </tr> </tbody> </table> | Time (min) | %B | 0.0 | 10 | 4.0 | 10 | 9.0 | 50 | 11.0 | 100 | 13.0 | 100 |
| Time (min)        | %B   |            |    |     |    |     |    |     |    |      |     |      |     |
| 0.0               | 10   |            |    |     |    |     |    |     |    |      |     |      |     |
| 4.0               | 10   |            |    |     |    |     |    |     |    |      |     |      |     |
| 9.0               | 50   |            |    |     |    |     |    |     |    |      |     |      |     |
| 11.0              | 100  |            |    |     |    |     |    |     |    |      |     |      |     |
| 13.0              | 100  |            |    |     |    |     |    |     |    |      |     |      |     |
| Flow Rate:        | 0.75 mL/min  |            |    |     |    |     |    |     |    |      |     |      |     |
| Pressure:         | 250 Bar (initial pressure)   |            |    |     |    |     |    |     |    |      |     |      |     |
| Temperature:      | 30 °C  |            |    |     |    |     |    |     |    |      |     |      |     |
| Detection:        | UV 254, 325 nm, VWD  |            |    |     |    |     |    |     |    |      |     |      |     |
| Injection Volume: | 1.0 µL   |            |    |     |    |     |    |     |    |      |     |      |     |
| Sample Solvent:   | 50/50 water/acetonitrile   |            |    |     |    |     |    |     |    |      |     |      |     |
| Response Time:    | 0.02 sec   |            |    |     |    |     |    |     |    |      |     |      |     |
| Data rate:        | 25 Hz  |            |    |     |    |     |    |     |    |      |     |      |     |
| Flow Cell:        | 2.5 µL semi-micro  |            |    |     |    |     |    |     |    |      |     |      |     |
| LC System:        | Shimadzu Prominence UFC XR   |            |    |     |    |     |    |     |    |      |     |      |     |
| ECV:              | ~14 µL   |            |    |     |    |     |    |     |    |      |     |      |     |

Green coffee extract is sold as a dietary supplement to aid in weight loss. Chlorogenic acid is the active ingredient. The chlorogenic acid and caffeine can be easily analyzed using HPLC. A bottle of green coffee extract capsules was purchased at a local pharmacy. The powder from one capsule was placed into a vial containing 20 mL of 50/50: acetonitrile/methanol and vortexed to mix. The vial was sonicated for 10 minutes and then allowed to settle. Five mL of the supernate was filtered through a 13 mm, 0.45 µm porosity Nylon syringe filter. A portion of this solution was diluted 1:10 using a 50/50 mixture of water and acetonitrile. The diluted sample was injected into the HPLC. Caffeine was only detectable using the 254 nm wavelength, while chlorogenic acid was detected with increased sensitivity at the 325 nm wavelength. The solvent gradient extends past 10 minutes to remove some minor late eluting peaks.

### STRUCTURES




Chlorogenic Acid

Caffeine

### ANALYTES

1. Chlorogenic acid
2. Caffeine




# AROMATIC COMPOUNDS

## Flavor and Fragrance Compounds — GC

### TEST CONDITIONS

Column: Avantor® Hichrom HI-ACID (FFAP), 0.25 µm, 0.25 mm x 30 m  
Part Number: [HI102502530](#)  
Oven Program: 75 °C, 2 °C/min, 140 °C  
Carrier Gas: Helium, 1.5 mL/min  
Injector: Split, 250 °C  
Detector: FID, 250 °C

The HI-ACID (FFAP) column is a nitroterephthalic acid modified polyethylene glycol (PEG) phase of high polarity. Useful for the analysis of free fatty acids, flavors and fragrances, amides, alcohols and glycols.



### ANALYTES

1. Isoamyl Acetate
2. alpha-Phellandrene
3. Cumene
4. 1-Pentanol
5. 6-Methyl-5-hepten-2-one
6. 2-Nonalone
7. Ethyl Octanoate
8. Furfural
9. Benzaldehyde
10. Linalool
11. Isobutyric acid
12. n-Butyric acid
13. Ethyl Decanoate
14. Furfuryl alcohol
15. alpha and beta-Terpineol
16. Carvone
17. Methyl Salicylate

# AROMATIC COMPOUNDS

## LC-MS of PAHs in Grilled Meat — HPLC/UHPLC

### TEST CONDITIONS


Column: HALO® EnviroClass PAH, 90 Å, 2.7 µm, 2.1 x 100 mm  
Part Number: 92842-612  
Flow Rate: 0.4 mL/min  
Pressure: 289 bar  
Column Temperature: 30 °C  
Injection Volume: 1 µL  
Sample Solvent: Methanol  
LC System: Shimadzu Nexera  
Mobile Phase  
A: Water/0.1% Formic acid  
B: Acetonitrile/0.1% Formic acid

Gradient:

|  | Time (min) | %B  |
|--|------------|-----|
|  | 0.0        | 40  |
|  | 5.0        | 100 |
|  | 8.0        | 100 |
|  | 8.01       | 40  |

### MASS SPECTROMETRY CONDITIONS

MS System: Thermo Scientific™ Q Exactive™ HF  
ESI voltage: 5.5 kV  
Heater Temp: 400 °C  
Sheath gas: 35 (arbitrary units)  
Aux gas: 8 (arbitrary units)  
Tube lens voltage: 40 Vt




# AROMATIC COMPOUNDS

## LC-MS of PAHs in Grilled Meat (continued)

The HALO® PAH column continues in the tradition of HALO® products by offering high resolution separations, in high throughput time frames. 18 PAH compounds with 6 sets of isomeric compounds were able to be quickly and efficiently resolved in under 8 minutes. In addition, the high resolution separation of the HALO® PAH column, enabled chrysene and benzo[a]pyrene to be resolved from a complex meat matrix, enabling quantitation of PAH contamination present in barbequed steak. The concentration of PAHs in the sample, were below those established by the EU, and demonstrates that not only can the HALO® PAH column be used in the stringent regulatory testing of current established methods, but also be relied upon as future regulations dictate the establishment of new methods, requiring lower limits of detection. The HALO® PAH column offers a rugged and reproducible particle design meeting the needs of complex matrix testing. Fused-Core® technology is ideal for PAH analysis in particular, enabling customers to achieve analytical goals of speed, accuracy, and precision LC separations.

## ANALYTES AND ELUTION ORDER

| Peak # | Compound               | Precursor Ion | Fragment 1 | Fragment 2 |
|--------|------------------------|---------------|------------|------------|
| 1      | Naphthalene            | 128           | 78         | 102        |
| 2      | Acenaphthylene         | 152           | 126        | 151        |
| 3      | 1-Methylnaphthalene    | 142           | 89         | 115        |
| 4      | 2-Methylnaphthalene    | 142           | 115        | 141        |
| 5      | Acenaphthene           | 154           | 126        | 153        |
| 6      | Fluorene               | 166           | 115        | 165        |
| 7      | Phenanthrene           | 178           | 151        | 176        |
| 8      | Anthracene             | 178           | 152        | 176        |
| 9      | Fluoranthene           | 202           | 150        | 200        |
| 10     | Pyrene                 | 202           | 150        | 200        |
| 11     | Benzo[a]anthracene     | 228           | 150        | 226        |
| 12     | Chrysene               | 228           | 200        | 226        |
| 13     | Benzo[b]fluoranthene   | 252           | 224        | 250        |
| 14     | Benzo[k]fluoranthene   | 252           | 224        | 250        |
| 15     | Benzo[a]pyrene         | 252           | 224        | 250        |
| 16     | Dibenz[a,h]anthracene  | 278           | 248        | 276        |
| 17     | Benzo[ghi]perylene     | 276           | 248        | 274        |
| 18     | Indeno[1,2,3-cd]pyrene | 276           | 246        | 274        |



# ARTIFICIAL COLORS, DYES AND PRESERVATIVES

## Artificial Food Colors — HPLC

### TEST CONDITIONS


Column: Avantor® ACE® C18, 100 Å, 3 µm, 4.6 x 100 mm  
Part Number: [ACE1111046](#)  
Mobile Phase: A: 1.0120 g Tetrabutylammonium bromide and .6804 g  $\text{KH}_2\text{PO}_4$  in 1 L  $\text{H}_2\text{O}$   
B: 1.0120 g Tetrabutylammonium bromide in 1 L MeOH

Gradient:

| Time (min) | %B |
|------------|----|
| 0          | 45 |
| 12         | 60 |
| 25         | 45 |

Flow Rate: 0.8 mL/min  
Temperature: Ambient  
Injection: 10 µL  
Detection: UV-Vis, 480 nm

Reproduced with permission of Lincoln Sutton & Wood, Norwich, UK



# ARTIFICIAL COLORS, DYES AND PRESERVATIVES

## Illegal Dyes in Spices — HPLC


### TEST CONDITIONS

Column: Avantor® ACE® C18, 100 Å, 3 µm, 4.6 x 100 mm  
Part Number: [ACE1111046](#)  
Mobile Phase:  
A: H<sub>2</sub>O  
B: MeOH  
C: 0.06 M Tetrabutylammonium bromide and 0.5 M KH<sub>2</sub>PO<sub>4</sub> in H<sub>2</sub>O pH 2.55

| Gradient: | Time (min) | %A | %B | %C |
|-----------|------------|----|----|----|
|           | 0          | 45 | 50 | 5  |
|           | 45         | 3  | 92 | 5  |
|           | 65         | 3  | 92 | 5  |
|           | 66         | 45 | 50 | 5  |
|           | 75         | 45 | 50 | 5  |

Flow Rate: 1 mL/min  
Temperature: Ambient  
Injection: 10 µL  
Detection: UV-Vis, 500 nm

Reproduced with permission of Lincoln Sutton & Wood, Norwich, UK



### ANALYTES

1. Rhodamine B
2. Orange II
3. Metanil Yellow
4. Butter Yellow
5. Para Red
6. Sudan Orange G
7. Toluidine Red
8. Sudan I
9. Sudan Red G
10. Sudan II
11. Sudan Black B
12. Sudan III
13. Sudan Red 7B
14. Sudan Red B
15. Sudan IV


# ARTIFICIAL COLORS, DYES AND PRESERVATIVES

## Preservatives (I) – UHPLC


## TEST CONDITIONS

Column: Avantor® ACE® Excel C18, 100 Å, 1.7 µm, 3.0 x 50 mm  
Part Number: [EXL1710503U](#)  
Mobile Phase: 20 mM Potassium phosphate pH 2.5 in MeCN/H<sub>2</sub>O (30:70 v/v)  
Flow Rate: 0.43 mL/min  
Temperature: 20 °C  
Injection Volume: 0.7 µL  
Detection: UV, 230 nm


## STRUCTURES




### Phthalic acid




#### 4-Hydroxybenzoic acid



### Sorbic acid



#### Methyl paraben



## ANALYTES

1. Phthalic acid
  2. 4-Hydroxybenzoic acid
  3. Sorbic acid
  4. Methyl acid
  5. Ethyl acid
  6. Propyl paraben


# CAPSAICINS

## Capsaicins — UHPLC


### TEST CONDITIONS

Column: Avantor® ACE® Excel C18, 100 Å, 2 µm, 3.0 x 100 mm  
Part Number: EXL1011003U  
Mobile Phase: 0.1% H<sub>3</sub>PO<sub>4</sub> in MeCN/H<sub>2</sub>O (60:40 v/v)  
Flow Rate: 0.6 mL/min  
Temperature: 40 °C  
Injection Volume: 1 µL  
Detection: UV, 220 nm  
Sample Preparation: Extraction with 10 mL ethanol  
Sonication for 20 minutes  
Centrifugation for 10 minutes at 6,000 RPM  
Filtration with 0.45 µm Whatman Mini-UniPrep  
syringeless filter vial


### STRUCTURES



Capsaicin



Dihydrocapsaicin



### ANALYTES

1. Capsaicin
2. Dihydrocapsaicin

# CAPSAICINS



## Capsaicins in Chili Powder—HPLC/UHPLC

### TEST CONDITIONS

| Column:       | HALO® C18, 90 Å, 2.7 µm, 3.0 x 100 mm  | Flow Rate:      | 0.8 mL/min                 |     |    |     |    |     |     |      |     |                   |        |
|---------------|--|-----------------|----------------------------|-----|----|-----|----|-----|-----|------|-----|-------------------|--------|
| Part Number:  | <u>92813-602</u>   | Pressure:       | 223 bar starting pressure  |     |    |     |    |     |     |      |     |                   |        |
| Mobile Phase: | A/B<br>A= water<br>B= acetonitrile   | Temperature:    | 40 °C                      |     |    |     |    |     |     |      |     |                   |        |
| Gradient:     | <table border="1"><thead><tr><th>Time (min)</th><th>%B</th></tr></thead><tbody><tr><td>0.0</td><td>40</td></tr><tr><td>5.0</td><td>60</td></tr><tr><td>7.0</td><td>100</td></tr><tr><td>20.0</td><td>100</td></tr></tbody></table> | Time (min)      | %B                         | 0.0 | 40 | 5.0 | 60 | 7.0 | 100 | 20.0 | 100 | Injection Volume: | 1.0 µL |
| Time (min)    | %B   |                 |                            |     |    |     |    |     |     |      |     |                   |        |
| 0.0           | 40   |                 |                            |     |    |     |    |     |     |      |     |                   |        |
| 5.0           | 60   |                 |                            |     |    |     |    |     |     |      |     |                   |        |
| 7.0           | 100  |                 |                            |     |    |     |    |     |     |      |     |                   |        |
| 20.0          | 100  |                 |                            |     |    |     |    |     |     |      |     |                   |        |
|               |  | Sample Solvent: | Acetonitrile               |     |    |     |    |     |     |      |     |                   |        |
|               |  | Detection:      | UV 230 nm, VWD             |     |    |     |    |     |     |      |     |                   |        |
|               |  | Response Time:  | 0.02 sec                   |     |    |     |    |     |     |      |     |                   |        |
|               |  | Data rate:      | 25 Hz                      |     |    |     |    |     |     |      |     |                   |        |
|               |  | Flow Cell:      | 2.5 µL semi-micro          |     |    |     |    |     |     |      |     |                   |        |
|               |  | LC System:      | Shimadzu Prominence UFC XR |     |    |     |    |     |     |      |     |                   |        |
|               |  | ECV:            | ~14 µL                     |     |    |     |    |     |     |      |     |                   |        |

### ANALYTES

1. Capsaicin 1
2. Capsaicin 2
3. Dihydrocapsaicin 1
4. Dihydrocapsaicin 2



Capsaicin and dihydrocapsaicin are two of the main components of chili powder that give it the “heat” when making a batch of chili. The amount of heat is often measured by a subjective test and then rated in terms of Scoville units that are a dilution factor beyond which the capsaicins and other hot compounds cannot be detected. One can also use HPLC to measure these compounds more objectively. Here these two ingredients are separated from an acetonitrile extract using a HALO® C18 column.

## Curcumin Analysis in Turmeric on a 1.5 mm Stainless Steel Column—HPLC/UHPLC


### TEST CONDITIONS

|                   |                                       |
|-------------------|---------------------------------------|
| Column:           | HALO® C18, 90 Å, 2.7 µm, 1.5 x 100 mm |
| Part Number:      | <a href="#">9281X-602</a>             |
| Mobile Phase A:   | Water/ 0.1% Formic acid               |
| Mobile Phase B:   | ACN/ 0.1% Formic acid                 |
| Isocratic:        | 50 %B                                 |
| Flow Rate:        | 0.2 mL/min                            |
| Back Pressure:    | 337 bar                               |
| Temperature:      | 30 °C                                 |
| Detection:        | UV 425 nm, PDA                        |
| Injection Volume: | 1 µL                                  |
| Sample Solvent:   | Ethanol                               |
| Data Rate:        | 100Hz                                 |
| Response Time:    | 0.025 sec                             |
| Flow Cell:        | 1 µL                                  |
| Instrument:       | Shimadzu Nexera X2                    |

Turmeric - along with its most active compound, curcumin have many proven health benefits such as improved heart health, antioxidant properties, and a potent anti-inflammatory. A separation of curcumins is performed on a HALO® 1.5 mm C18 column showing excellent resolution and sensitivity. Turmeric capsules were dissolved in ethanol at 0.5 mg/ml and filtered through a 0.45 um syringe filter before injection.

### ANALYTES

1. Bisdemethoxycurcumin
2. Demethoxycurcumin
3. Curcumin



# FLAVONOIDS AND OTHER POLYPHENOLS

## Common Catechins Found in Tea via LC-MS — LC-MS

### TEST CONDITIONS

Column: HALO® LPH-C18, 90 Å, 2 µm, 2.1 x 100 mm  
Part Number: 91822-616  
Mobile Phase A: Water, 0.2% Formic acid (pH 2.45)  
Mobile Phase B: Acetonitrile, 0.2% Formic acid

| Gradient: | Time (min) | %B |
|-----------|------------|----|
|           | 0.0        | 10 |
|           | 1.0        | 10 |
|           | 6.0        | 21 |
|           | 7.0        | 21 |


Flow Rate: 0.3 mL/min  
Pressure: 438 bar  
Temperature: 40 °C  
Detection: +/- ESI MS/MS  
Injection Volume: 2 µL  
Sample Solvent: Water  
MS System: Shimadzu 8040  
LC System: Shimadzu Nexera X2

### MS CONDITIONS:

Nebulizer Gas Flow: 2 L/min DL  
Temperature: 250 °C  
Heat Block Temperature: 400 °C  
Drying Gas Flow: 10 L/min

Catechins belong to the subgroup of polyphenols called flavonoids. These compounds contain antioxidant properties and exist in food and medicinal plants, including tea. An LC-MS separation of catechins and caffeine is demonstrated on a HALO® LPH-C18 column showing excellent resolution using purified standards. This column is ideal for low pH separations due to its sterically protected ligand, preventing acid hydrolysis and reducing retention drift over time.

| Peak # | Compound                 | M/Z |
|--------|--------------------------|-----|
| 1      | Gallic acid              | 169 |
| 2      | Gallocatechin            | 305 |
| 3      | Epigallocatechin         | 305 |
| 4      | Catechin                 | 289 |
| 5      | Caffeine                 | 195 |
| 6      | Epicatechin              | 289 |
| 7      | Epigallocatechin gallate | 457 |
| 8      | Gallocatechin gallate    | 457 |
| 9      | Epicatechin gallate      | 441 |
| 10     | Catechin gallate         | 441 |



# FLAVONOIDS AND OTHER POLYPHENOLS

## Cyclodextrin-Encapsulated Flavor Compounds in Beer — HPLC


### TEST CONDITIONS

Column: Avantor® ACE® C18, 100 Å, 3 µm, 4.0 x 150 mm  
Part Number: ACE1111504  
Mobile Phase: 0.1% Phosphoric acid in MeOH/H<sub>2</sub>O (53:47 v/v)  
Flow Rate: 0.5 mL/min  
Temperature: 35 °C  
Injection Volume: 20 µL  
Detection: UV, 225 nm (261 nm for 2-Methoxy-4-vinylphenol)

Reproduced with permission of FlavorActiV Ltd, UK

### ANALYTES


1. 2-Methylphenol
2. 2-Bromophenol
3. 2-Methoxy-4-vinylphenol
4. 4-Ethylphenol
5. 2,4-Dichlorophenol
6. 2,4,6-Trimethylphenol
7. 2-Isopropylphenol



# FLAVONOIDS AND OTHER POLYPHENOLS

## Cyclodextrin-Encapsulated Flavor Compounds in Beer (continued)

- 1 ■ 2-Methylphenol
- 2 ▲ 2-Methoxy-4-vinylphenol
- 3 □ 2,6-Dichlorophenol
- 4 ◆ 2,4,6-Trimethylphenol
- 5 ▨ 2-Bromophenol
- 6 ▲ 4-Ethylphenol
- 7 ▲ 2-Isopropylphenol




# FLAVONOIDS AND OTHER POLYPHENOLS

## Five Flavonoids—UHPLC


### TEST CONDITIONS

Column: HALO® C8, 90 Å, 2 µm, 2.1 x 100 mm  
Part Number: 91812-608  
Mobile Phase: 75/25: A/B  
A= 0.025 M Ammonium formate, pH=3  
B= Acetonitrile  
Flow Rate: 0.5 mL/min  
Pressure: 473 bar  
Temperature: 40 °C  
Detection: UV 276 nm, PDA  
Injection Volume: 0.1 µL  
Sample Solvent: Methanol  
Response Time: 0.025 sec  
Flow Cell: 1 µL  
LC System: Shimadzu Nexera  
Data rate: 100 Hz  
ECV: ~ 7 µL


### STRUCTURES




Naringin




Quercetin



Myricetin




Naringenin



Hesperetin

### ANALYTES

1. Naringin
2. Myricetin
3. Quercetin
4. Naringenin
5. Hesperetin



This separation of five flavonoids shows symmetrical peak shape on a HALO® C8 column using a mass-spec compatible mobile phase.


# FLAVONOIDS AND OTHER POLYPHENOLS

## Four Flavonoids—HPLC


### TEST CONDITIONS

Column: Avantor® ACE® C18, 100 Å,  
5 µm, 4.6 x 150 mm  
Part Number: [ACE1211546](#)  
Mobile Phase: MeCN/0.1% Formic acid  
in H<sub>2</sub>O (40:60 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Injection Volume: 1 µL  
Detection: UV, 254 nm


### STRUCTURES




Rutin




Quercetin



Quercitrin



Kaempferol



# FLAVONOIDS AND OTHER POLYPHENOLS

## Vanillins — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® Excel C18-Amide, 100 Å, 3 µm, 4.6 x 150 mm

Part Number: EXL11121546U

Mobile Phase: A: 0.1% Formic acid in H<sub>2</sub>O  
B: 0.1% Formic acid in MeCN

| Gradient:           | Time (min) | %B |
|---------------------|------------|----|
|                     | 0          | 30 |
|                     | 10         | 55 |
|                     | 10         | 55 |
|                     | 15         | 30 |
| Post time 5 minutes |            |    |


Flow Rate: 1 mL/min

Temperature: 40 °C


Injection: 5 µL

Detection: UV, 260 nm


### STRUCTURES




Vanillic acid




4-Hydroxybenzoic acid




Vanillin




4-Hydroxybenzaldehyde



Guaiacol




Ethyl Vanillin



Eugenol

### ANALYTES

1. Vanillic acid
2. 4-Hydroxybenzoic acid
3. Vanillin
4. 4-Hydroxybenzaldehyde
5. Guaiacol
6. Ethyl Vanillin
7. Eugenol




# FLAVONOIDS AND OTHER POLYPHENOLS

## Isoflavones — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® C18, 100 Å, 5 µm, 4.6 x 150 mm  
Part Number: ACE1211546  
Mobile Phase: MeCN/0.1% Formic acid in H<sub>2</sub>O (35:65 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Injection Volume: 1 µL  
Detection: UV, 254 nm



# FLAVONOIDS AND OTHER POLYPHENOLS

## Polyphenols in Wine—HPLC/UHPLC

### TEST CONDITIONS


Column: HALO® LPH-C18, 90 Å, 2.7 µm, 2.1 x 100 mm  
Part Number: 92822-616  
Mobile Phase A: Water/ 0.1% Formic acid  
Mobile Phase B: Acetonitrile/ 0.1% Formic acid

| Gradient: | Time (min) | %B  |
|-----------|------------|-----|
|           | 0.0        | 0   |
|           | 3.5        | 8   |
|           | 7.1        | 10  |
|           | 25.0       | 30  |
|           | 26.0       | 40  |
|           | 27.0       | 100 |
|           | 29.0       | 100 |
|           | 30.0       | 0   |
|           | 35.0       | 0   |

Flow Rate: 0.3 mL/min  
Pressure: 159 bar  
Temperature: 30 °C  
Detection: UV 280 nm, PDA

Injection Volume: 0.7 µL  
Sample Solvent: Water  
Data Rate: 100 Hz  
Response Time: 0.025 sec  
Flow Cell: 1 µL  
LC System: Shimadzu Nexera X2

Polyphenols can be found in a wide variety of plant-based foods and are packed with antioxidants and potential health benefits. There are more than 8,000 of these types of compounds which contain multiples of phenol units. Common polyphenols found in wine are separated using a HALO® 90 Å LPH-C18 column using analytical standards. This stationary phase contains a sterically protected ligand which is ideal for high stability under low pH conditions.



### ANALYTES

- |                             |                     |
|-----------------------------|---------------------|
| 1. Gallic acid              | 10. o-Coumaric acid |
| 2. Epigallocatechin         | 11. Quercitrin      |
| 3. Chlorogenic acid         | 12. Myricetin       |
| 4. Catechin                 | 13. Resveratrol     |
| 5. Caffeic acid             | 14. Morin           |
| 6. Epicatechin              | 15. Quercetin       |
| 7. Epigallocatechin gallate | 16. Naringenin      |
| 8. p-Coumaric acid          | 17. Apigenin        |
| 9. Ferulic acid             | 18. Kaempferol      |

# FLAVONOIDS AND OTHER POLYPHENOLS


## Six Flavonoids – HPLC/UHPLC

## TEST CONDITIONS


|                   |   |
|-------------------|---|
| Column:           | HALO® C18, 90 Å, 2.7 µm, 4.6 x 50 mm  |
| Part Number:      | <u>92814-402</u>  |
| Mobile Phase:     | 70/30: A/B<br>A= 0.02 M Phosphate buffer, pH=2.9, (adj.)<br>B= Acetonitrile |
| Flow Rate:        | 2.0 mL/min  |
| Pressure:         | 224 bar   |
| Temperature:      | 30 °C   |
| Detection:        | UV 240 nm, VWD  |
| Injection Volume: | 1.0 µL  |
| Sample Solvent:   | Methanol  |
| Response Time:    | 0.02 sec  |
| Flow Cell:        | 2.5 µL semi-micro   |
| LC System:        | Shimadzu Prominence UFC-XR  |
| ECV:              | ~14µL   |
| Data rate:        | 25 Hz   |

## ANALYTES


1. Catechin
  2. Naringin
  3. Myricetin
  4. Quercetin
  5. Naringenin
  6. Hesperetin




STRUCTURES




## Catechin




## Quercetin




Naringin



### Naringenin



## Myricetin



## Hesperetin

Flavonoids are naturally occurring polyphenols that are found in plant leaves, flowers and seeds. They have beneficial health effects and are often taken as dietary supplements. Analysis of this flavonoids mixture can be carried out in less than 2 minutes using a short HALO® Fused Core C18 column.

# FOOD ADDITIVES


## Food Additives—HPLC/UHPLC

### TEST CONDITIONS


Column: HALO® Phenyl-Hexyl, 90 Å, 2.7 µm, 4.6 x 50 mm  
HALO® RP-Amide, 90 Å, 2.7 µm, 4.6 x 50 mm  
Part Numbers: 92814-406, 92814-407, respectively  
Mobile Phase: 70/30—A/B  
A= 0.025 M phosphate buffer, pH=2.5  
B= Methanol  
Flow Rate: 1.5 mL/min  
Pressure: approximately 220 Bar  
Temperature: 40 °C  
Detection: UV 220 nm, VWD  
Injection Volume: 2.0 µL  
Sample Solvent: 50/50-Water/methanol  
Response Time: 0.02 sec  
Flow Cell: 2.5 µL semi-micro  
LC System: Shimadzu Prominence UFC-XR  
Extra column volume: ~14µL

These compounds are often added to foods to sweeten or preserve them. They can be rapidly analyzed using HALO® Phenyl-Hexyl or RP-Amide phases. Note the difference in retention and selectivity of the two phases when run under the same conditions. This allows for flexibility in method development and optimization of the separation.


### STRUCTURES




Ascorbic acid




Saccharin




Benzoic acid




Aspartame



Methyl paraben





Sorbic acid



Dehydroacetic acid

### ANALYTES

1. Ascorbic acid
2. Saccharin
3. Aspartame
4. Sorbic acid
5. Benzoic acid
6. Methyl paraben
7. Dehydroacetic acid



# FOOD ADDITIVES

## Food Additives Assay — HPLC/UHPLC


### TEST CONDITIONS

|                       |  |
|-----------------------|--|
| Column:               | HALO® AQ-C18, 90 Å, 5 µm, 4.6 x 150 mm |
| Part Number:          | <u>95814-722</u>                       |
| Mobile Phase A:       | 20 mM Ammonium acetate                 |
| Mobile Phase B:       | Methanol                               |
| Isocratic:            | 90/10 A/B                              |
| Flow Rate:            | 2 mL/min                               |
| Pressure:             | 336 bar                                |
| Temperature:          | 30 °C                                  |
| Detection wavelength: | 230 nm                                 |
| Injection Volume:     | 10 µL                                  |
| Sample Solvent:       | Mobile phase                           |
| Data Rate:            | 100 hz                                 |
| Response Time:        | 0.025 sec                              |
| Flow Cell:            | 1 µL                                   |
| LC System:            | Shimadzu Nexera X2                     |

A rapid and highly efficient assay <400 bar for food security and safety measurements is demonstrated with a HALO® 90 Å AQ-C18 5 um, 4.6 x 150 mm column. Determination of acesulfame, benzoic acid, sorbic acid and saccharin sodium food additives are specified in China's national standard regulation methods GB 5009.28-2016 and GB 5009.140-2016. These compounds are used as anti-septic/anti-microbial agents to prevent spoilage of food products by microorganisms. A baseline resolution separation is completed <1.7 min; modernization of this method is as easy as exploiting the 5 micron HALO® column - compatible with HPLC and UHPLC instruments.

### ANALYTES

1. Acesulfame
2. Benzoic acid
3. Sorbic acid
4. Saccharin sodium




## Determination of Sterols (TMS Derivatives) in Sunflower Oil – GC

### TEST CONDITIONS

Column: Avantor® Hichrom HI-LAP, 0.1 µm, 0.32 mm x 25 m  
Part Number: HI463201025  
Oven Program: 260 °C Isothermal  
Carrier Gas: Hydrogen, 2 mL/min  
Injector: Split, 300 °C, 1 µL, 1:80 Split Ratio  
Detector: FID, 300 °C

The HI-LAP (Lipid Analysis Phase) is a dedicated and unique stationary phase especially developed for the analysis of lipids, sterols and saturated and unsaturated triglyceride separations.

Acknowledgement: Dr C. Mariani, Stazione Sperimentale Oli e Grassi, Via Giuseppe Colombo 79, Milan, Italy




# FAMEs and Conjugated Linoleic Acid (CLA) Isomers — GC

## TEST CONDITIONS

Column: Avantor® Hichrom HI-10, 0.2 µm, 0.25 mm x 100 m  
 Part Number: [HI2625020100](#)  
 Oven Program: 40 °C (4 min), 10 °C/min, 120 °C (1min), 5 °C/min, 180 °C  
 (18 min), 2 °C /min, 200 °C (1 min), 2 °C/min, 240 °C  
 Carrier Gas: Hydrogen, 1.2 mL/min, constant flow  
 Injector: Split, 250 °C, 1:30 split ratio  
 Detector: FID, 250 °C

The HI-10 is a high polarity phase with a 100% cyanopropyl polysiloxane composition. It is well suited for the analysis of FAMEs, especially for the separation of cis/trans isomers of FAMEs.




# Steroids Spiked in Ground Beef—HPLC/UHPLC

## TEST CONDITIONS

Column: HALO® C18, 90 Å, 2.7 µm, 2.1 x 100 mm  
 Part Number: 92812-602  
 Mobile Phase A: Water, 5 mM Ammonium formate, 0.1 % Formic acid pH 4.0  
 Mobile Phase B: Methanol  
 Flow Rate: 0.3 mL/min  
 Pressure: 190 bar  
 Temperature: 50 °C  
 Injection Volume: 2.0 µL  
 Sample Solvent: 45/55/ MEOH/H<sub>2</sub>O  
 Detection: +ESI/ -ESI MS/MS  
 LC System: Shimadzu Nexera X2  
 ESI LCMS system: Shimadzu LCMS-8040

| Time (min) | %B   |
|------------|------|
| 2.0        | 14   |
| 3.0        | 60   |
| 3.5        | 60   |
| 8.0        | 100  |
| 10.0       | 100  |
| 10.5       | 0    |
| 12.5       | stop |



For over fifty years, the Food and Drug Administration (FDA) has approved the use of a number of steroids in beef cattle, including natural estrogen, progesterone, testosterone, and their synthetic versions such as trenbolone acetate (TBA). The function of these drugs is to increase growth rate and the efficiency by which the animals convert the feed they eat into muscle/meat. The drugs are usually administered as implants (dosing of 100-200 days), which are placed under the skin on the back side of the animal's ear. The implants dissolve slowly under the skin and are not re-moved. Although cooking the meat does have some effect on the stability of the steroids in beef, it does not eliminate the exposure, as many steroids are stable at elevated temperatures. A standard panel of steroids spiked into ground beef, and then run on the HALO® 90 Å C18, shows a highly resolved separation of all compounds. The panel consisted of common growth promoters and those used for therapeutic purposes, and was chosen to represent the most common steroids that can be expected to be found in beef, through therapeutic or growth promotion utilization.

## MS SOURCE CONDITIONS:

Spray Voltage: 3.0 kV  
 Nebulizing gas: 2 L/min  
 Drying gas: 15 L/min  
 DL temp: 250 °C  
 Heat Block: 400 °C

| Peak ID | Compound                | Transition         | RT (Min) |
|---------|-------------------------|--------------------|----------|
| 1       | ALDOSTERONE             | 361.0000>343.1000  | 1.154    |
| 2       | CORTICOSTERONE          | 347.6000>109.0000  | 1.965    |
| 3       | ZERANOL                 | 321.0000>277.0000  | 2.355    |
| 4       | MGA                     | 395.0000> 325.1000 | 3.100    |
| 5       | TESTOSTERONE            | 289.0000>109.0000  | 3.366    |
| 6       | 17A-METHYL-TESTOSTERONE | 303.1000>97.0000   | 3.839    |
| 7       | PROGESTERONE            | 315.0000>109.1000  | 5.085    |
| 8       | ESTRADIOL 17B           | 272.4000>159.1000  | 7.501    |

# MYCOTOXINS, PESTICIDES AND TOXINS

## 250 Pesticide Screen using LC-MS/MS — LC-MS


### TEST CONDITIONS

Column: Avantor® ACE® Excel C18, 100 Å, 2 µm, 2.1 x 100 mm  
Part Number: EXL1011002U  
Mobile Phase:  
A: 10 mM Ammonium formate + 0.05% Formic acid in H<sub>2</sub>O  
B: 10 mM Ammonium formate + 0.05% Formic acid in MeOH

| Gradient: | Time (min) | %B  |
|-----------|------------|-----|
|           | 0.00       | 2   |
|           | 0.25       | 30  |
|           | 10.00      | 100 |
|           | 12.00      | 100 |
|           | 12.50      | 2   |
|           | 14.50      | 2   |

Flow Rate: 0.5 mL/min  
Temperature: 50 °C  
Detection: TSQ Quantiva triple quad MS  
Positive mode HESI  
Spray voltage: 3500 V  
Ion transfer tube temperature: 350 °C  
Vaporizer temperature: 300 °C

Reproduced with permission of Centre for Analytical Chemistry, California Department of Food and Agriculture, USA



# MYCOTOXINS, PESTICIDES AND TOXINS

## 250 Pesticide Screen using LC-MS/MS (continued)

| Analyte             | $t_R$<br>mins | Adduct                            | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|---------------------|---------------|-----------------------------------|----------------------|------------------|-----------------|
| 3-OH Carbofuran     | 2.25          | [M+H] <sup>+</sup>                | 238.1                | 181.2            | 163.1           |
| 5-OH Thiabendazole  | 1.66          | [M+H] <sup>+</sup>                | 218.0                | 147.2            | 191.1           |
| Abamectin           | 9.45          | [M+NH <sub>4</sub> ] <sup>+</sup> | 890.5                | 305.3            | 567.5           |
| Acephate            | 1.26          | [M+H] <sup>+</sup>                | 184.0                | 143.1            | 125.1           |
| Acetamiprid         | 2.24          | [M+H] <sup>+</sup>                | 223.1                | 126.1            | 90.1            |
| Aldicarb            | 2.95          | [M+NH <sub>4</sub> ] <sup>+</sup> | 208.1                | 116.1            | 89.0            |
| Aldicarb sulfone    | 1.44          | [M+NH <sub>4</sub> ] <sup>+</sup> | 240.1                | 148.0            | 86.0            |
| Aldicarb sulfoxide  | 1.37          | [M+NH <sub>4</sub> ] <sup>+</sup> | 224.1                | 132.0            | 89.1            |
| Allethrin           | 8.33          | [M+H] <sup>+</sup>                | 303.2                | 135.1            | 123.1           |
| Ametoctradin        | 7.64          | [M+H] <sup>+</sup>                | 276.2                | 149.1            | 176.2           |
| Atrazine            | 4.64          | [M+H] <sup>+</sup>                | 216.1                | 174.0            | 104.0           |
| Azinphos ethyl      | 6.30          | [M+H] <sup>+</sup>                | 346.0                | 132.1            | 233.0           |
| Azinphos methyl     | 5.14          | [M+H] <sup>+</sup>                | 318.0                | 132.0            | 124.9           |
| Azinphos methyl OA  | 2.98          | [M+H] <sup>+</sup>                | 302.0                | 132.2            | 160.0           |
| Azoxystrobin        | 5.59          | [M+H] <sup>+</sup>                | 404.1                | 372.1            | 344.1           |
| Bendiocarb          | 3.72          | [M+H] <sup>+</sup>                | 224.1                | 167.1            | 109.1           |
| Benoxacor           | 5.23          | [M+H] <sup>+</sup>                | 260.1                | 134.1            | 120.1           |
| Bifenazate          | 6.27          | [M+H] <sup>+</sup>                | 301.1                | 198.0            | 170.1           |
| Bitertanol          | 7.41          | [M+H] <sup>+</sup>                | 338.2                | 269.3            | 99.1            |
| Boscalid            | 5.85          | [M+H] <sup>+</sup>                | 343.0                | 307.0            | 140.0           |
| Bupirimate          | 6.68          | [M+H] <sup>+</sup>                | 317.2                | 210.2            | 237.3           |
| Buprofezin          | 8.24          | [M+H] <sup>+</sup>                | 306.1                | 201.1            | 106.1           |
| Cadusafos           | 7.58          | [M+H] <sup>+</sup>                | 271.1                | 159.0            | 131.0           |
| Carbaryl            | 4.07          | [M+NH <sub>4</sub> ] <sup>+</sup> | 219.1                | 145.1            | 127.0           |
| Carbendazim         | 2.10          | [M+H] <sup>+</sup>                | 192.1                | 160.1            | 132.1           |
| Carbofuran          | 3.77          | [M+H] <sup>+</sup>                | 222.1                | 165.2            | 123.2           |
| Carboxin            | 3.97          | [M+H] <sup>+</sup>                | 236.1                | 143.0            | 93.0            |
| Carfentrazone ethyl | 6.88          | [M+H] <sup>+</sup>                | 412.0                | 346.1            | 366.0           |
| Chlorantraniliprole | 5.24          | [M+H] <sup>+</sup>                | 484.0                | 286.0            | 194.0           |
| Chlорfenvinphos     | 7.21          | [M+H] <sup>+</sup>                | 359.0                | 170.0            | 99.1            |
| Chlorimuron ethyl   | 5.73          | [M+H] <sup>+</sup>                | 415.1                | 186.0            | 83.0            |
| Chlorpyrifos        | 8.47          | [M+H] <sup>+</sup>                | 349.9                | 198.0            | 97.0            |

| Analyte            | $t_R$<br>mins | Adduct   | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|--------------------|---------------|--|----------------------|------------------|-----------------|
| Chlorpyrifos OA    | 6.65          | [M+H] <sup>+</sup>   | 334.0                | 278.0            | 197.9           |
| Clethodim          | 7.71          | [M+H] <sup>+</sup>   | 360.3                | 164.1            | 136.1           |
| Clofentezine       | 7.38          | [M+H] <sup>+</sup>   | 303.0                | 138.1            | 102.0           |
| Cloransulam methyl | 4.13          | [M+H] <sup>+</sup>   | 430.0                | 398.1            | 370.0           |
| Clothianidin       | 1.99          | [M+H] <sup>+</sup>   | 250.0                | 169.1            | 132.0           |
| Coumaphos          | 7.07          | [M+H] <sup>+</sup>   | 363.0                | 227.1            | 307.1           |
| crotoxyphos        | 5.86          | [M+NH <sub>4</sub> ] <sup>+</sup>  | 332.1                | 127.1            | 193.1           |
| Crufomate          | 6.77          | [M+H] <sup>+</sup>   | 292.1                | 236.1            | 108.1           |
| Cyantraniliprole   | 4.33          | [M+2H] <sup>+</sup>  | 475.0                | 286.0            | 444.1           |
| Cyazofamid         | 6.52          | [M+H] <sup>+</sup>   | 325.1                | 108.1            | 261.2           |
| Cyflufenamid       | 7.42          | [M+H] <sup>+</sup>   | 413.1                | 295.1            | 203.0           |
| Cymoxanil          | 2.48          | [M+H] <sup>+</sup>   | 199.1                | 128.1            | 111.1           |
| Cyphenothenrin     | 9.27          | [M+NH <sub>4</sub> ] <sup>+</sup>  | 393.2                | 151.2            | 123.2           |
| Cyprosulfamide     | 3.30          | [M+H] <sup>+</sup>   | 375.1                | 135.1            | 254.1           |
| Cyromazine         | 1.15          | [M+H] <sup>+</sup>   | 167.1                | 125.2            | 68.2            |
| DEF                | 9.20          | [M+H] <sup>+</sup>   | 315.1                | 169.0            | 113.0           |
| Demeton-S sulfone  | 2.55          | [M+H] <sup>+</sup>   | 291.1                | 235.1            | 263.1           |
| Dialifos           | 7.46          | [M+H] <sup>+</sup>   | 394.0                | 208.1            | 181.0           |
| Diazinon           | 7.12          | [M+H] <sup>+</sup>   | 305.1                | 169.1            | 153.2           |
| Diazinon OA        | 5.32          | [M+H] <sup>+</sup>   | 289.1                | 153.2            | 233.1           |
| Dichlormid         | 3.85          | [M+H] <sup>+</sup>   | 208.0                | 140.0            | 81.2            |
| Dichlorvos         | 3.63          | [M+H] <sup>+</sup>   | 221.0                | 109.1            | 127.0           |
| Dicrotophos        | 1.87          | [M+H] <sup>+</sup>   | 238.1                | 112.2            | 193.1           |
| Diethofencarb      | 5.53          | [M+H] <sup>+</sup>   | 268.2                | 124.1            | 180.2           |
| Diffubenzuron      | 6.66          | [M+H] <sup>+</sup>   | 311.0                | 158.0            | 141.0           |
| Dimethenamid       | 5.70          | [M+H] <sup>+</sup>   | 276.1                | 244.1            | 168.2           |
| Dimethoate         | 2.23          | [M+H] <sup>+</sup>   | 230.1                | 199.0            | 125.0           |
| Dimethomorph       | 5.76,<br>6.07 | [M+H] <sup>+</sup>   | 388.1                | 301.0            | 165.1           |
| Dimotefuran        | 1.36          | [M+H] <sup>+</sup>   | 203.1                | 129.1            | 114.2           |
| Dioxacarb          | 2.26          | [M+H] <sup>+</sup>   | 224.1                | 123.1            | 167.1           |
| Dioxathion         | 8.10          | [M-C <sub>2</sub> H <sub>5</sub> O <sub>2</sub> PS <sub>2</sub> ] <sup>-</sup> | 271.1                | 97.0             | 125.0           |
| Disulfoton sulfone | 4.59          | [M+H] <sup>+</sup>   | 307.0                | 261.1            | 125.0           |

# MYCOTOXINS, PESTICIDES AND TOXINS

## 250 Pesticide Screen using LC-MS/MS (continued)

| Analyte                | t <sub>R</sub><br>mins | Adduct                            | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|------------------------|------------------------|-----------------------------------|----------------------|------------------|-----------------|
| Disulfoton sulfoxide   | 4.49                   | [M+H] <sup>+</sup>                | 291.0                | 185.1            | 213.1           |
| Diuron                 | 4.82                   | [M+H] <sup>+</sup>                | 233.0                | 72.1             | 160.0           |
| DMST                   | 3.90                   | [M+H] <sup>+</sup>                | 215.1                | 106.1            | 151.0           |
| Dodine                 | 7.56                   | [M+H] <sup>+</sup>                | 228.3                | 186.3            | 60.1            |
| Emamectin              | 8.57                   | [M+H] <sup>+</sup>                | 886.5                | 158.1            | 126.1           |
| Ethiofencarb           | 4.27                   | [M+H] <sup>+</sup>                | 226.1                | 107.1            | 169.1           |
| Ethiofencarb sulfone   | 1.90                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 275.1                | 107.1            | 201.1           |
| Ethiofencarb sulfoxide | 1.98                   | [M+H] <sup>+</sup>                | 242.1                | 107.1            | 185.0           |
| Éthion                 | 8.31                   | [M+H] <sup>+</sup>                | 385.0                | 199.1            | 143.0           |
| Ethion monoxon         | 6.73                   | [M+H] <sup>+</sup>                | 369.0                | 199.0            | 143.0           |
| Ethiprole              | 5.77                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 413.9                | 351.0            | 255.0           |
| Ethofumesate           | 5.55                   | [M+H] <sup>+</sup>                | 287.1                | 121.1            | 241.1           |
| Ethoprop               | 6.46                   | [M+H] <sup>+</sup>                | 243.1                | 173.0            | 131.0           |
| Etofenprox             | 9.75                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 394.2                | 177.2            | 107.1           |
| Etoxazole              | 8.73                   | [M+H] <sup>+</sup>                | 360.2                | 141.0            | 304.2           |
| Famoxadone             | 7.24                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 392.2                | 331.1            | 238.0           |
| Fenamidone             | 5.76                   | [M+H] <sup>+</sup>                | 312.1                | 236.1            | 92.2            |
| Fenamiphos             | 6.71                   | [M+H] <sup>+</sup>                | 304.1                | 217.1            | 202.0           |
| Fenamiphos sulfone     | 4.10                   | [M+H] <sup>+</sup>                | 336.1                | 266.1            | 188.1           |
| Fenamiphos sulfoxide   | 3.96                   | [M+H] <sup>+</sup>                | 320.1                | 233.1            | 171.1           |
| Fenazaquin             | 9.21                   | [M+H] <sup>+</sup>                | 307.2                | 161.2            | 57.2            |
| Fenhexamid             | 6.39                   | [M+H] <sup>+</sup>                | 302.1                | 178.0            | 97.2            |
| Fenobucarb             | 5.49                   | [M+H] <sup>+</sup>                | 208.1                | 95.0             | 152.0           |
| Fenoxyprop ethyl       | 8.04                   | [M+H] <sup>+</sup>                | 362.1                | 288.1            | 91.1            |
| Fenoxy carb            | 6.80                   | [M+H] <sup>+</sup>                | 302.1                | 88.1             | 116.1           |
| Fenpropimorph          | 6.42                   | [M+H] <sup>+</sup>                | 304.3                | 147.2            | 119.1           |
| Fenpyroximate          | 8.90                   | [M+H] <sup>+</sup>                | 422.2                | 366.1            | 214.2           |
| Fensulfothion          | 4.89                   | [M+H] <sup>+</sup>                | 309.0                | 235.0            | 281.1           |
| Fenuron                | 2.17                   | [M+H] <sup>+</sup>                | 165.1                | 72.1             | 77.1            |
| Flonicamid             | 1.66                   | [M+H] <sup>+</sup>                | 230.1                | 203.0            | 98.0            |
| Fluazifop P butyl      | 8.12                   | [M+H] <sup>+</sup>                | 384.1                | 282.2            | 328.2           |
| Fludioxonil            | 5.76                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 266.1                | 158.1            | 131.0           |

| Analyte            | t <sub>R</sub><br>mins | Adduct                            | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|--------------------|------------------------|-----------------------------------|----------------------|------------------|-----------------|
| Flufenoxuron       | 8.79                   | [M+H] <sup>+</sup>                | 489.0                | 158.1            | 141.1           |
| Flufenpyr ethyl    | 6.72                   | [M+H] <sup>+</sup>                | 409.1                | 335.0            | 307.0           |
| Flumetsulam        | 2.03                   | [M+H] <sup>+</sup>                | 326.1                | 129.1            | 109.0           |
| Flumiclorac pentyl | 8.13                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 441.1                | 308.1            | 354.1           |
| Fluometuron        | 4.31                   | [M+H] <sup>+</sup>                | 233.1                | 72.2             | 46.3            |
| Fluopicolide       | 6.00                   | [M+H] <sup>+</sup>                | 383.0                | 173.0            | 145.0           |
| Fluopyram          | 6.33                   | [M+H] <sup>+</sup>                | 397.1                | 173.0            | 208.0           |
| Fluoxastrobin      | 6.40                   | [M+H] <sup>+</sup>                | 459.1                | 427.2            | 188.1           |
| Fluridone          | 5.32                   | [M+H] <sup>+</sup>                | 330.1                | 309.1            | 290.0           |
| Flusilazole        | 6.77                   | [M+H] <sup>+</sup>                | 316.1                | 247.2            | 165.1           |
| Fluthiacet methyl  | 6.88                   | [M+H] <sup>+</sup>                | 404.0                | 344.0            | 273.9           |
| Flutolanil         | 5.95                   | [M+H] <sup>+</sup>                | 324.1                | 262.0            | 282.0           |
| Flutriafol         | 4.74                   | [M+H] <sup>+</sup>                | 302.1                | 70.1             | 123.1           |
| Fluxapyroxad       | 6.02                   | [M+H] <sup>+</sup>                | 382.1                | 342.1            | 314.1           |
| Forchlorfenuron    | 4.78                   | [M+H] <sup>+</sup>                | 248.1                | 129.1            | 93.1            |
| Formetanate HCl    | 1.26                   | [M+H] <sup>+</sup>                | 222.0                | 165.1            | 120.0           |
| Fosthiazate        | 4.40                   | [M+H] <sup>+</sup>                | 284.1                | 104.1            | 228.1           |
| Hexaconazole       | 7.29                   | [M+H] <sup>+</sup>                | 314.1                | 158.9            | 70.0            |
| Hexythiazox        | 8.51                   | [M+H] <sup>+</sup>                | 353.1                | 228.0            | 168.0           |
| Imazalil           | 5.14                   | [M+H] <sup>+</sup>                | 297.1                | 159.1            | 255.1           |
| Imazosulfuron      | 5.28                   | [M+H] <sup>+</sup>                | 413.0                | 153.0            | 156.1           |
| Imidacloprid       | 1.96                   | [M+H] <sup>+</sup>                | 256.1                | 209.1            | 175.0           |
| Imiprothrin        | 6.34                   | [M+H] <sup>+</sup>                | 319.2                | 151.1            | 123.1           |
| Indoziflam         | 6.58                   | [M+H] <sup>+</sup>                | 302.2                | 158.1            | 145.1           |
| Indoxacarb         | 7.75                   | [M+H] <sup>+</sup>                | 528.1                | 249.0            | 150.1           |
| Ipcconazole        | 7.81                   | [M+H] <sup>+</sup>                | 334.2                | 70.1             | 125.0           |
| Iprovalicarb       | 6.31                   | [M+H] <sup>+</sup>                | 321.2                | 119.1            | 186.2           |
| Isofenphos         | 7.39                   | [M+H] <sup>+</sup>                | 346.1                | 217.0            | 245.1           |
| Isoprocarb         | 4.67                   | [M+H] <sup>+</sup>                | 194.1                | 95.1             | 152.2           |
| Isoproturon        | 4.79                   | [M+H] <sup>+</sup>                | 207.2                | 72.2             | 165.2           |
| Kresoxim methyl    | 6.90                   | [M+H] <sup>+</sup>                | 314.1                | 267.2            | 222.1           |
| Lactofen           | 8.22                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 479.1                | 344.1            | 223.0           |

# MYCOTOXINS, PESTICIDES AND TOXINS

## 250 Pesticide Screen using LC-MS/MS (continued)

| Analyte               | t <sub>R</sub><br>mins | Adduct                            | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|-----------------------|------------------------|-----------------------------------|----------------------|------------------|-----------------|
| Lenacil               | 4.67                   | [M+H] <sup>+</sup>                | 235.1                | 153.1            | 136.1           |
| Leptophos OA          | 7.75                   | [M+2+H] <sup>+</sup>              | 396.9                | 155.1            | 364.9           |
| Linuron               | 5.46                   | [M+H] <sup>+</sup>                | 249.0                | 182.1            | 160.1           |
| Malathion             | 5.92                   | [M+H] <sup>+</sup>                | 331.0                | 127.1            | 285.1           |
| Malathion OA          | 3.89                   | [M+H] <sup>+</sup>                | 315.1                | 127.1            | 99.0            |
| Mandipropamid         | 5.94                   | [M+H] <sup>+</sup>                | 412.1                | 328.2            | 356.2           |
| Mefenpyrdiethyl       | 7.26                   | [M+H] <sup>+</sup>                | 373.1                | 327.1            | 160.0           |
| Mepronipyrim          | 6.21                   | [M+H] <sup>+</sup>                | 224.1                | 106.2            | 77.1            |
| Mesotripone           | 2.01                   | [M+H] <sup>+</sup>                | 340.1                | 228.1            | 104.1           |
| Metaflumizone         | 8.30                   | [M+H] <sup>+</sup>                | 507.1                | 178.0            | 287.1           |
| Metalaxyl             | 4.91                   | [M+H] <sup>+</sup>                | 280.1                | 220.1            | 192.1           |
| Metaldehyde           | 2.02                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 194.1                | 62.2             | 45.3            |
| Metconazole           | 7.32                   | [M+H] <sup>+</sup>                | 320.2                | 70.1             | 125.0           |
| Methamidophos         | 1.16                   | [M+H] <sup>+</sup>                | 142.0                | 94.2             | 125.1           |
| Methidathion          | 4.97                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 320.0                | 145.1            | 85.1            |
| Methiocarb            | 5.64                   | [M+H] <sup>+</sup>                | 226.1                | 169.2            | 121.1           |
| Methiocarb sulfone    | 2.35                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 275.0                | 122.1            | 201.1           |
| Methiocarb sulfoxide  | 2.10                   | [M+H] <sup>+</sup>                | 242.1                | 185.1            | 122.1           |
| Methomyl              | 1.61                   | [M+H] <sup>+</sup>                | 163.1                | 106.1            | 88.1            |
| Methoxyfenozide       | 6.04                   | [M+H] <sup>+</sup>                | 369.2                | 149.1            | 313.1           |
| Metolcarb             | 3.28                   | [M+H] <sup>+</sup>                | 166.1                | 109.1            | 94.1            |
| Metribuzin            | 3.59                   | [M+H] <sup>+</sup>                | 215.1                | 187.1            | 131.1           |
| Mevinphos             | 2.70                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 242.1                | 193.1            | 127.1           |
| Monocrotophos         | 1.71                   | [M+H] <sup>+</sup>                | 224.1                | 193.0            | 127.0           |
| Monolinuron           | 4.16                   | [M+H] <sup>+</sup>                | 215.1                | 126.1            | 148.1           |
| Myclobutanil          | 6.15                   | [M+H] <sup>+</sup>                | 289.1                | 125.0            | 70.1            |
| Nicosulfuron          | 3.45                   | [M+H] <sup>+</sup>                | 411.1                | 182.0            | 213.0           |
| Norflurazon           | 4.98                   | [M+H] <sup>+</sup>                | 304.0                | 160.0            | 140.0           |
| Norflurazon desmethyl | 4.43                   | [M+H] <sup>+</sup>                | 290.0                | 179.0            | 140.0           |
| Omethoate             | 1.33                   | [M+H] <sup>+</sup>                | 214.0                | 183.0            | 125.0           |
| Oxamyl                | 1.48                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 237.1                | 72.0             | 90.0            |
| Oxamyl oxime          | 1.34                   | [M+H] <sup>+</sup>                | 163.1                | 72.1             | 90.1            |

| Analyte                   | t <sub>R</sub><br>mins | Adduct                            | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|---------------------------|------------------------|-----------------------------------|----------------------|------------------|-----------------|
| Oxydemeton methyl         | 1.57                   | [M+H] <sup>+</sup>                | 247.0                | 169.1            | 109.1           |
| Oxydemeton methyl sulfone | 1.62                   | [M+H] <sup>+</sup>                | 263.0                | 169.0            | 109.0           |
| Parathion methyl OA       | 3.10                   | [M+H] <sup>+</sup>                | 248.0                | 202.0            | 109.1           |
| Parathion OA              | 4.61                   | [M+H] <sup>+</sup>                | 276.1                | 220.1            | 248.1           |
| Pencycuron                | 7.50                   | [M+H] <sup>+</sup>                | 329.1                | 125.1            | 89.1            |
| Penflufen                 | 6.95                   | [M+H] <sup>+</sup>                | 318.2                | 234.1            | 141.0           |
| Penthiopyrad              | 7.05                   | [M+H] <sup>+</sup>                | 360.1                | 177.1            | 276.1           |
| Phenothrin                | 9.56                   | [M+H] <sup>+</sup>                | 351.2                | 183.1            | 168.0           |
| Phenthroate               | 6.81                   | [M+H] <sup>+</sup>                | 321.0                | 247.1            | 79.1            |
| Phorate OA                | 5.10                   | [M+H] <sup>+</sup>                | 245.0                | 75.2             | 47.2            |
| Phorate OA Sulfone        | 2.51                   | [M+H] <sup>+</sup>                | 277.0                | 155.0            | 127.0           |
| Phorate OA Sulfoxide      | 2.31                   | [M+H] <sup>+</sup>                | 261.0                | 153.0            | 81.0            |
| Phorate Sulfone           | 4.61                   | [M+H] <sup>+</sup>                | 293.0                | 114.9            | 171.0           |
| Phorate Sulfoxide         | 4.49                   | [M+H] <sup>+</sup>                | 277.0                | 170.9            | 199.0           |
| Phosalone                 | 7.35                   | [M+H] <sup>+</sup>                | 368.0                | 182.0            | 111.1           |
| Phosmet                   | 5.21                   | [M+H] <sup>+</sup>                | 318.0                | 160.1            | 133.1           |
| Phosmet OA                | 3.12                   | [M+H] <sup>+</sup>                | 302.0                | 160.0            | 133.0           |
| Phosphamidon              | 3.43                   | [M+H] <sup>+</sup>                | 300.1                | 127.1            | 174.1           |
| Phoxim                    | 7.25                   | [M+H] <sup>+</sup>                | 299.1                | 77.2             | 129.1           |
| Picoxystrobin             | 6.79                   | [M+H] <sup>+</sup>                | 368.1                | 145.0            | 115.0           |
| Pirimicarb                | 4.24                   | [M+H] <sup>+</sup>                | 239.2                | 182.1            | 72.0            |
| Pirimicarb Desmethyl      | 2.71                   | [M+H] <sup>+</sup>                | 225.1                | 168.2            | 72.1            |
| Pirimiphos Methyl         | 7.34                   | [M+H] <sup>+</sup>                | 306.1                | 164.2            | 108.1           |
| Prallethrin               | 7.69                   | [M+H] <sup>+</sup>                | 301.2                | 133.0            | 151.2           |
| Prochloraz                | 7.39                   | [M+H] <sup>+</sup>                | 376.0                | 308.1            | 70.1            |
| Profoxydim                | 7.71,<br>9.00          | [M+H] <sup>+</sup>                | 466.2                | 280.0            | 180.0           |
| Promecarb                 | 5.88                   | [M+H] <sup>+</sup>                | 208.1                | 109.0            | 151.1           |
| Propamocarb               | 1.41                   | [M+H] <sup>+</sup>                | 189.1                | 102.0            | 144.0           |
| Propaquizafop             | 8.21                   | [M+H] <sup>+</sup>                | 444.1                | 299.2            | 371.2           |
| Propargite                | 8.74                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 368.2                | 231.2            | 175.1           |
| Propetamphos              | 6.13                   | [M+H] <sup>+</sup>                | 282.1                | 138.1            | 156.1           |
| Propoxur (S)              | 3.69                   | [M+H] <sup>+</sup>                | 210.1                | 168.2            | 111.1           |

# MYCOTOXINS, PESTICIDES AND TOXINS

## 250 Pesticide Screen using LC-MS/MS (continued)

| Analyte              | t <sub>R</sub><br>mins | Adduct                            | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|----------------------|------------------------|-----------------------------------|----------------------|------------------|-----------------|
| Prosulfuron          | 5.29                   | [M+H] <sup>+</sup>                | 420.1                | 167.1            | 141.1           |
| Pymetrozine          | 1.44                   | [M+H] <sup>+</sup>                | 218.1                | 105.1            | 78.1            |
| Pyraclostrobin       | 7.30                   | [M+H] <sup>+</sup>                | 388.1                | 163.1            | 194.1           |
| Pyraoflufen Ethyl    | 7.13                   | [M+H] <sup>+</sup>                | 413.0                | 339.0            | 253.1           |
| Pyrozophos           | 7.31                   | [M+H] <sup>+</sup>                | 374.1                | 222.2            | 194.1           |
| Pyridaben            | 9.22                   | [M+H] <sup>+</sup>                | 365.1                | 309.0            | 147.1           |
| Pyridalyl            | 10.21                  | [M+2+H] <sup>+</sup>              | 492.0                | 110.9            | 164.0           |
| Pyrimethanil         | 5.45                   | [M+H] <sup>+</sup>                | 200.1                | 107.1            | 168.1           |
| Pyriproxyfen         | 8.39                   | [M+H] <sup>+</sup>                | 322.1                | 96.0             | 227.1           |
| Quinalphos           | 6.78                   | [M+H] <sup>+</sup>                | 299.1                | 163.1            | 147.1           |
| Quinoxifen           | 8.50                   | [M+H] <sup>+</sup>                | 308.0                | 197.1            | 214.1           |
| Quizalofop Ethyl     | 8.01                   | [M+H] <sup>+</sup>                | 373.1                | 299.2            | 255.1           |
| Resmethrin           | 9.40                   | [M+H] <sup>+</sup>                | 339.2                | 128.1            | 171.1           |
| Rimsulfuron          | 3.94                   | [M+H] <sup>+</sup>                | 432.1                | 182.1            | 139.0           |
| Rotenone             | 6.71                   | [M+H] <sup>+</sup>                | 395.2                | 213.2            | 192.1           |
| Saflufenacil         | 5.32                   | [M+H] <sup>+</sup>                | 501.1                | 349.1            | 198.0           |
| Sedaxane             | 6.20,<br>6.54          | [M+H] <sup>+</sup>                | 332.2                | 159.0            | 139.0           |
| Sethoxydim           | 8.03                   | [M+H] <sup>+</sup>                | 328.2                | 178.0            | 220.1           |
| Simazine             | 3.66                   | [M+H] <sup>+</sup>                | 202.1                | 104.1            | 132.1           |
| Spinetoram           | 8.14                   | [M+H] <sup>+</sup>                | 748.5                | 142.1            | 203.1           |
| Spinosad A           | 7.69                   | [M+H] <sup>+</sup>                | 732.5                | 142.1            | 98.0            |
| Spinosad D           | 8.10                   | [M+H] <sup>+</sup>                | 746.5                | 142.1            | 98.0            |
| Spirodiclofen        | 8.91                   | [M+H] <sup>+</sup>                | 411.1                | 313.1            | 71.1            |
| Spiromesifen         | 8.66                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 388.1                | 273.1            | 187.0           |
| Spiromesifen Alcohol | 5.01                   | [M+H] <sup>+</sup>                | 273.2                | 187.1            | 179.1           |
| Spirotetramat        | 6.38                   | [M+H] <sup>+</sup>                | 374.2                | 302.3            | 216.2           |
| Spiroxamine          | 5.95                   | [M+H] <sup>+</sup>                | 298.3                | 144.2            | 100.2           |
| Sulfoxaflor          | 2.39                   | [M+NH <sub>4</sub> ] <sup>+</sup> | 295.2                | 174.1            | 154.1           |
| Sulprofos            | 8.56                   | [M+H] <sup>+</sup>                | 323.0                | 219.1            | 139.1           |

| Analyte               | t <sub>R</sub><br>mins | Adduct               | Precursor<br>Ion m/z | Quant<br>Ion m/z | Conf<br>Ion m/z |
|-----------------------|------------------------|----------------------|----------------------|------------------|-----------------|
| TCMTB                 | 5.48                   | [M+H] <sup>+</sup>   | 239.0                | 180.0            | 136.0           |
| Tebufenozide          | 6.78                   | [M+H] <sup>+</sup>   | 353.2                | 133.0            | 104.8           |
| Tebufenpyrad          | 8.19                   | [M+H] <sup>+</sup>   | 334.2                | 117.1            | 145.1           |
| Tebuthiuron           | 3.89                   | [M+H] <sup>+</sup>   | 229.1                | 172.0            | 116.0           |
| Tepraloxydim          | 4.10,<br>6.19          | [M+H] <sup>+</sup>   | 342.2                | 250.1            | 166.1           |
| Terbufos Sulfone      | 5.46                   | [M+H] <sup>+</sup>   | 321.0                | 115.0            | 143.0           |
| Terbufos Sulfoxide    | 5.49                   | [M+H] <sup>+</sup>   | 305.1                | 97.0             | 187.0           |
| Terbutylazine         | 5.71                   | [M+H] <sup>+</sup>   | 230.1                | 174.1            | 104.1           |
| Tetrachlorvinphos     | 6.86                   | [M+2+H] <sup>+</sup> | 366.9                | 127.1            | 206.0           |
| Tetramethrin          | 7.91,<br>8.10          | [M+H] <sup>+</sup>   | 332.2                | 164.1            | 135.1           |
| Thiabendazole         | 2.48                   | [M+H] <sup>+</sup>   | 202.0                | 175.0            | 131.1           |
| Thiacloprid           | 2.55                   | [M+H] <sup>+</sup>   | 253.0                | 126.1            | 99.1            |
| Thiamethoxam          | 1.65                   | [M+H] <sup>+</sup>   | 292.0                | 211.1            | 181.1           |
| Thifensulfuron Methyl | 3.28                   | [M+H] <sup>+</sup>   | 388.0                | 167.1            | 205.0           |
| Thiobencarb           | 7.46                   | [M+H] <sup>+</sup>   | 258.1                | 125.0            | 89.0            |
| Thiodicarb            | 4.34                   | [M+H] <sup>+</sup>   | 355.1                | 163.2            | 88.1            |
| Thionazin             | 4.74                   | [M+H] <sup>+</sup>   | 249.1                | 193.1            | 97.0            |
| Topramezone           | 1.63                   | [M+H] <sup>+</sup>   | 364.1                | 334.1            | 125.1           |
| Triadimefon           | 6.07                   | [M+H] <sup>+</sup>   | 294.1                | 197.0            | 225.0           |
| Triadimenol           | 6.25                   | [M+H] <sup>+</sup>   | 296.1                | 70.2             | 99.0            |
| Triazophos            | 6.19                   | [M+H] <sup>+</sup>   | 314.1                | 162.1            | 119.1           |
| Tribenuron Methyl     | 4.59                   | [M+H] <sup>+</sup>   | 396.1                | 155.1            | 181.1           |
| Trichlorfon           | 2.26                   | [M+H] <sup>+</sup>   | 256.9                | 109.0            | 221.0           |
| Tricyclazole          | 2.80                   | [M+H] <sup>+</sup>   | 190.0                | 163.1            | 136.1           |
| Trifloxystrobin       | 7.78                   | [M+H] <sup>+</sup>   | 409.1                | 186.2            | 206.2           |
| Triflumizole          | 7.87                   | [M+H] <sup>+</sup>   | 346.1                | 278.0            | 73.0            |
| Triforine             | 5.23                   | [M+2+H] <sup>+</sup> | 434.9                | 213.0            | 98.2            |
| Zoxamide              | 7.09                   | [M+H] <sup>+</sup>   | 336.0                | 187.0            | 159.0           |

# MYCOTOXINS, PESTICIDES AND TOXINS

## LC-MS Screening Comparison of Mycotoxins in Craft and Home Brewed Beers — HPLC/UHPLC

### TEST CONDITIONS

Column: HALO® PFP, 90 Å, 2.7 µm, 2.1 x 100 mm

Part Number: 92812-609

Mobile Phase A: Water, 5 mM Ammonium formate, 0.1 % Formic acid

Mobile Phase B: Methanol, 0.1% Formic acid

Gradient:

| Time (min) | %B  |
|------------|-----|
| 0.0        | 0   |
| 0.5        | 14  |
| 2.0        | 14  |
| 3.0        | 60  |
| 3.5        | 60  |
| 8.0        | 100 |
| 10.0       | 100 |
| 10.5       | 0   |
| 12.5       | End |

Flow Rate: 0.4 mL/min

Pressure: 290 bar

Temperature: 40 °C

Injection Volume: 7.0 µL

Sample Solvent: 49/50/1 ACN/H<sub>2</sub>O/Acetic acid

+ESI MS/MS


LC System: Shimadzu Nexera X2

ESI LCMS System: Shimadzu LCMS-8040

Mycotoxin contamination can have serious health implications. Although there are no set regulatory limits for mycotoxins in beer, most governments have clear levels for mycotoxins in various types of grain and animal feed. For example, in the United States, most levels are in the mid to high ppb range. Despite relatively low levels of mycotoxin activity in the beer, given the propensity for people to indulge in excessive drinking, and the cumulative effects of the toxicity of these compounds, excessive consumption would lead to a cumulative toxic effect, which warrants further analysis and regulation. Beer analysis can be challenging due to matrix effects and interference, often resulting in low sensitivity and ambiguous results; therefore, it is critical to have a column that has superior performance. The HALO® 90 Å PFP can not only meet these challenges, but exceed them by demonstrating superior performance and sensitivity, making it an ideal column to be used in environmental, and, specifically, mycotoxin analysis.

# MYCOTOXINS, PESTICIDES AND TOXINS

## LC-MS Screening Comparison of Mycotoxins in Craft and Home (*continued*)




**CRAFT-BREWED BEER**

| Peak ID | Mycotoxin               | Retention Time (min) | Precursor Ion | Product Ion |
|---------|-------------------------|----------------------|---------------|-------------|
| 1       | T-2 Toxin               | 3.95                 | 489.2         | 245.1       |
| 2       | Aflatoxin G2            | 4.65                 | 331.1         | 189.2       |
| 3       | 15-acetyldeoxynivalenol | 4.88                 | 339.1         | 321.1       |
| 4       | Aflatoxin B2            | 5.52                 | 315.1         | 287.1       |
| 5       | Aflatoxin M1            | 5.75                 | 329.1         | 273.3       |
| 6       | Zearalenone             | 9.55                 | 319.1         | 283.2       |

**HOME BREWED BEER**

| Peak ID | Mycotoxin               | Retention Time (min) | Precursor Ion | Product Ion |
|---------|-------------------------|----------------------|---------------|-------------|
| 1       | T-2 Toxin               | 3.95                 | 489.2         | 245.1       |
| 2       | 15-acetyldeoxynivalenol | 4.88                 | 339.1         | 321.1       |
| 3       | Aflatoxin M1            | 5.75                 | 329.1         | 273.3       |
| 4       | Zearalenone             | 9.55                 | 319.1         | 283.2       |




# MYCOTOXINS, PESTICIDES AND TOXINS

## Mycotoxins by LC-MS/MS — LC-MS

### TEST CONDITIONS

| Column:       | Avantor® ACE® Excel C18-AR, 100 Å, 2 µm, 2.1 x 50 mm   |  | Flow Rate:   | 0.6 mL/min |     |    |     |    |     |    |     |    |            |  |
|---------------|--|--|--------------|------------|-----|----|-----|----|-----|----|-----|----|------------|--|
| Part Number:  | <u>EXL1090502U</u>   |  | Temperature: | 40 °C      |     |    |     |    |     |    |     |    |            |  |
| Mobile Phase: | A: 1 mM Ammonium acetate, 0.5% Acetic acid in H <sub>2</sub> O<br>B: 1 mM Ammonium acetate, 0.5% Acetic acid in 95% MeOH   |  | Injection:   | 2 µL       |     |    |     |    |     |    |     |    |            |  |
| Gradient:     | <table><thead><tr><th>Time (min)</th><th>%B</th></tr></thead><tbody><tr><td>0.0</td><td>40</td></tr><tr><td>1.0</td><td>40</td></tr><tr><td>2.4</td><td>60</td></tr><tr><td>6.8</td><td>87</td></tr></tbody></table> |  | Time (min)   | %B         | 0.0 | 40 | 1.0 | 40 | 2.4 | 60 | 6.8 | 87 | Detection: | AB SCIEX triple quad 5500<br>Positive ESI mode<br>Source temperature: 500 °C<br>IonSpray voltage: 5500 V |
| Time (min)    | %B   |  |              |            |     |    |     |    |     |    |     |    |            |  |
| 0.0           | 40   |  |              |            |     |    |     |    |     |    |     |    |            |  |
| 1.0           | 40   |  |              |            |     |    |     |    |     |    |     |    |            |  |
| 2.4           | 60   |  |              |            |     |    |     |    |     |    |     |    |            |  |
| 6.8           | 87   |  |              |            |     |    |     |    |     |    |     |    |            |  |

Reproduced with permission of Biotage GB Ltd



### ANALYTES

1. Aflatoxin G2 (m/z 331.1 - 313.1)
2. Aflatoxin G1 (m/z 329.0 - 243.1)
3. Aflatoxin B2 (m/z 315.1 - 287.0)
4. Aflatoxin B1 (m/z 313.1 - 285.0)
5. HT-2 toxin (m/z 442.2 - 263.1)
6. T-2 toxin (m/z 484.2 - 305.1)
7. Ochratoxin A (m/z 404.1 - 239.0)

# MYCOTOXINS, PESTICIDES AND TOXINS

## Paralytic Shellfish Poisoning (PSP) Toxins — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® UltraCore SuperC18, 95 Å, 5 µm, 4.6 x 150 mm  
 Part Number: CORE5A1546U  
 Mobile Phase: A: 0.1 M Ammonium formate in H<sub>2</sub>O  
 B: 0.1 M Ammonium formate in H<sub>2</sub>O/MeOH (95:5 v/v)



| Gradient: | Time (min) | %B |
|-----------|------------|----|
|           | 0.00       | 0  |
|           | 2.00       | 0  |
|           | 4.00       | 80 |
|           | 5.50       | 80 |
|           | 5.51       | 0  |
|           | 7.00       | 0  |

Flow Rate: 2 mL/min  
 Temperature: 20 °C  
 Injection Volume: 30 µL  
 Detection: Fluorescence  $\lambda_{\text{Ex}}$  340 nm,  $\lambda_{\text{Em}}$  395nm  
 Sample: Prechromatographic oxidation with hydrogen peroxide and periodate

Reproduced with permission of CEFAS, Weymouth, UK

### STRUCTURES

| PST Variant | R1 | R2                            | R3                            | R4                      |
|-------------|----|-------------------------------|-------------------------------|-------------------------|
| STX         | H  | H                             | H                             | H <sub>2</sub> N-COO    |
| NEO         | OH | H                             | H                             | H <sub>2</sub> N-COO    |
| GTX1        | OH | H                             | OSO <sub>3</sub> <sup>-</sup> | H <sub>2</sub> N-COO    |
| GTX2        | H  | H                             | OSO <sub>3</sub> <sup>-</sup> | H <sub>2</sub> N-COO    |
| GTX3        | H  | OSO <sub>3</sub> <sup>-</sup> | H                             | H <sub>2</sub> N-COO    |
| GTX4        | OH | OSO <sub>3</sub> <sup>-</sup> | H                             | H <sub>2</sub> N-COO    |
| GTX5 (B1)   | H  | H                             | H                             | O <sub>3</sub> S-NH-COO |
| C1          | H  | H                             | OSO <sub>3</sub> <sup>-</sup> | O <sub>3</sub> S-NH-COO |
| C2          | H  | OSO <sub>3</sub> <sup>-</sup> | H                             | O <sub>3</sub> S-NH-COO |
| dcSTX       | H  | H                             | H                             | OH                      |
| dcNEO       | OH | H                             | H                             | OH                      |
| dcGTX1      | OH | H                             | OSO <sub>3</sub> <sup>-</sup> | OH                      |
| dcGTX2      | H  | H                             | OSO <sub>3</sub> <sup>-</sup> | OH                      |
| dcGTX3      | H  | OSO <sub>3</sub> <sup>-</sup> | H                             | OH                      |



### ANALYTES

1. dcGTX2,3
2. GTX1/4 + dcGTX2,3
3. C1,2
4. dcSTX + dcNEO
5. dcSTX + NEO
6. GTX2/3 + GTX1/4
7. GTX5
8. STX + NEO

# MYCOTOXINS, PESTICIDES AND TOXINS

## Patulin and HMF — HPLC/UHPLC

### TEST CONDITIONS


Column: HALO® Biphenyl, 90 Å, 2.7 µm, 2.1 x 100 mm  
Part Number: 92812-611  
Mobile Phase A: Water with 0.1% Acetic acid  
Mobile Phase B: Acetonitrile with 0.1% Acetic acid

| Gradient: | Time (min) | %B |
|-----------|------------|----|
|           | 0.0        | 5  |
|           | 2.6        | 90 |


Flow Rate: 0.6 mL/min  
Initial Pressure: 285 bar  
Temperature: 40 °C  
Detection: UV 276 nm, PDA  
Injection Volume: 1.0 µL  
Sample: Apple Juice spiked with HMF and 50 ng/mL Patulin  
Data Rate: 100 Hz  
Response Time: 0.025 sec  
Flow Cell: 1 µL  
LC System: Shimadzu Nexera X2

In the United States the FDA maintains different limits for mycotoxins in many foods and beverages. Patulin, a mycotoxin that is produced from mold on a variety of fruits has a limit of 50 ug/kg. For analysis, patulin was spiked into apple juice and the sample was cleaned up using solid phase extraction. Interfering analytes such as 5-(Hydroxymethyl) furfural (HMF) can make analysis more challenging. This separation shows the two compounds separated on a HALO® Biphenyl column with enough resolution to easily check for sample recovery.


### STRUCTURES



5-(Hydroxymethyl) furfural



Patulin



# MYCOTOXINS, PESTICIDES AND TOXINS

## Pesticide Screening of Barley — UHPLC

### TEST CONDITIONS

Column: HALO® LPH-C18, 90 Å, 2 µm, 2.1 x 100 mm  
Part Number: 91822-616  
Mobile Phase A: Water, 0.1% Formic acid  
Mobile Phase B: Acetonitrile, 0.1% Formic acid

| Gradient: | Time (min) | %B  |
|-----------|------------|-----|
|           | 0.0        | 30  |
|           | 1.0        | 30  |
|           | 12.0       | 100 |
|           | 16.0       | 100 |


Flow Rate: 0.2 mL/min  
Pressure: 235 bar  
Temperature: 30 °C  
Detection: +ESI MS/MS  
Injection Volume: 2 µL  
Sample Solvent: Methanol  
MS System: Shimadzu 8040  
LC System: Shimadzu Nexera X2

### MS CONDITIONS

Nebulizer Gas Flow: 3 L/min  
DL Temperature: 250 °C  
Heat Block Temperature: 400 °C  
Drying Gas Flow: 18 L/min

Pesticide screening methods can help show whether there is a concern with your soil, crops, and even water supply. A pesticide screening is performed on a sample of barley using a HALO® 90 Å LPH-C18 column. This column is ideal for low pH testing conditions based on its sterically protected ligand which helps reduce acid hydrolysis of the stationary phase leading to an increase in column lifetime.

| Peak # | Compound        | Transition   | CE  |
|--------|-----------------|--------------|-----|
| 1      | Carbendazim     | 192>160.1    | -21 |
| 2      | Dicrotophos     | 238>112      | -22 |
| 3      | Azamethiphos    | 324.9>183    | -17 |
| 4      | Pyrimethanil    | 200.10>107.2 | -25 |
| 5      | Carbofuran      | 222>123      | -22 |
| 6      | Dodemorph       | 282.2>116.1  | -25 |
| 7      | Atrazine        | 216.03>174.1 | -17 |
| 8      | Diuron          | 232.94>72    | -17 |
| 9      | Iprovalicarb    | 321.1>119    | -30 |
| 10     | Azoxystrobin    | 404.04>372.1 | -14 |
| 11     | Fluopram        | 396.98>208   | -25 |
| 12     | Methoxyfenozide | 369.1>149.1  | -25 |
| 13     | Flutolanil      | 324>242.1    | -28 |
| 14     | Picoxystrobin   | 368>145.1    | -25 |



# PROTEINS

## Whey Proteins from Whole Milk — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® C4-300, 300 Å, 3 µm, 2.1 x 150 mm  
Part Number: ACE2131502  
Mobile Phase: A: 0.5% Formic acid in H<sub>2</sub>O  
B: 0.5% Formic acid in MeCN


| Gradient: | Time (min) | %B |
|-----------|------------|----|
|           | 0          | 35 |
|           | 16         | 43 |
|           | 17         | 80 |
|           | 20         | 80 |
|           | 21         | 35 |
|           | 31         | 35 |

Flow Rate: 0.4 mL/min  
Temperature: 40 °C  
Injection: 10 µL  
Detection: ESI-MS (+ve)

*Reproduced with permission of University of Vienna, Austria*

### ANALYTES

1. α-Lactalbumin
2. β-Lactoglobulin B
3. β-Lactoglobulin A



# SULFUR COMPOUNDS: GARLIC


## Garlic Analysis (I) — HPLC

### TEST CONDITIONS


Column: Avantor® ACE® C18, 100 Å, 5 µm, 4.6 x 250 mm  
Part Number: ACE1212546  
Mobile Phase: H<sub>2</sub>O/MeOH (50:50 v/v)  
Flow Rate: 1 mL/min  
Temperature: 30 °C  
Injection Volume: 20 µL  
Detection: UV, 210 nm

Reproduced with permission of Neem Biotech, Cardiff, UK


### STRUCTURES





Allicin



Alliin



Deoxyallin



### ANALYTES

1. Allicin
2. Alliin
3. Deoxyallin

# SULFUR COMPOUNDS: GARLIC

## Garlic Analysis (II) — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® C18, 100 Å, 5 µm, 4.6 x 250 mm

Part Number: ACE1212546

Mobile Phase: A: H<sub>2</sub>O

B: MeCN

Gradient:

|  | Time (min) | %B  |
|--|------------|-----|
|  | 0          | 40  |
|  | 20         | 100 |
|  | 25         | 100 |


Flow Rate: 1 mL/min

Temperature: 30 °C


Injection Volume: 20 µL

Detection: UV, 254 nm

### STRUCTURES




Allicin



Ajoene

Reproduced with permission of Neem Biotech, Cardiff, UK



# SWEETENERS

## Additives and Intense Sweeteners — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® C18, 100 Å, 5 µm, 4.0 x 250 mm

Part Number: ACE1212504

Mobile Phase: A: H<sub>2</sub>O

B: MeCN

C: 1% TFA in H<sub>2</sub>O

Gradient:

| Time (min) | %A | %B | %C |
|------------|----|----|----|
| 0          | 88 | 2  | 10 |
| 25         | 50 | 40 | 10 |
| 30         | 50 | 60 | 10 |
| 35         | 88 | 2  | 10 |

Flow Rate: 1 mL/min


Temperature: 30 °C

Detection: ELSD

Reproduced with permission of Durham County Council, UK

### ANALYTES

1. Acesulfame K
2. Theobromine
3. Theophylline
4. Cyclamate
5. Saccharin
6. Caffeine
7. Sucratose
8. Quinine sulphate
9. Aspartame
10. Neohesperidin dihydrochalcone



## Analysis of Sucrose and Lactose in Chocolate Using a Penta-HILIC Column — HPLC/UHPLC

### TEST CONDITIONS



|                   |  |
|-------------------|--|
| Column:           | HALO® Penta-HILIC 90 Å, 2.7 µm, 4.6 x 150 mm |
| Part Number:      | <u>92814-705</u>                             |
| Mobile Phase A:   | Water  |
| Mobile Phase B:   | Acetonitrile                                 |
| Flow Rate:        | 1.4 mL/min                                   |
| Pressure:         | 213 bar                                      |
| Temperature:      | 65 °C  |
| Detection:        | ELSD, 40 °C, 3.3 bar                         |
| Injection Volume: | 15 µL  |
| Sample Solvent:   | 80/20 ACN/ Water                             |
| Response Time:    | 0.10 sec                                     |
| DataRate:         | 10 Hz  |
| LC System:        | Shimadzu NexeraX2                            |

Chocolate is a very well-known, popular, food type worldwide. It is used for all occasions and can even have some health benefits as well, which include improved blood flow and brain function. There are four main types of chocolate to choose from- milk, white, dark, and raw.

Analysis of three different types of chocolate (milk, white, and dark) was carried out (or performed) in HILIC mode using an ELSD detector. The compounds of interest were sucrose and lactose. The HALO® Penta-HILIC column was used, which has a polar ligand with 5 hydroxyl groups tethered via novel proprietary linkage chemistry to Fused-Core® silica particles.

### ANALYTES

1. Sucrose
2. D-(+) Lactose monohydrate



# SWEETENERS

## Analysis of Sugars in Juice — HPLC/UHPLC


### TEST CONDITIONS

|                   |  |
|-------------------|--|
| Column:           | HALO® Penta-HILIC 90 Å, 2.7 µm, 4.6 x 150 mm |
| Part Number:      | <u>92814-705</u>                             |
| Mobile Phase A:   | Water  |
| Mobile Phase B:   | Acetonitrile                                 |
| Isocratic:        | 80 %B  |
| Flow Rate:        | 1.4 mL/min                                   |
| Pressure:         | 213 bar                                      |
| Temperature:      | 65 °C  |
| Detection:        | ELSD, 40 °C, 3.3 bar                         |
| Injection Volume: | 0.2 µL                                       |
| Sample Solvent:   | Water  |
| Data Rate:        | 10 Hz  |
| Response Time:    | 0.10 sec                                     |
| Flow Cell:        | 1 µL   |
| LC System:        | Shimadzu Nexera                              |

The main sugars in natural fruit juice are fructose, glucose, and sucrose. Each type of juice will contain different ratios of these sugars. Juices obtained from concentrate can also be found to have various amounts of artificial sweeteners. Analysis of sugars is performed on a HALO® Penta-HILIC column with excellent speed and resolution. A comparison of the different sugars in apple juice and orange juice is observed using an ELSD detector.

### ANALYTES

1. Sucrose
2. Glucose
3. Fructose



# SWEETENERS

## Analysis of Sugars in Pure Honey — HPLC/UHPLC

### TEST CONDITIONS


|                   |  |
|-------------------|--|
| Column:           | HALO® Penta-HILIC 90 Å, 2.7 µm, 4.6 x 150 mm |
| Part Number:      | <u>92814-705</u>                             |
| Mobile Phase A:   | Water  |
| Mobile Phase B:   | Acetonitrile                                 |
| Isocratic:        | 80% B  |
| Flow Rate:        | 1.4 mL/min                                   |
| Initial Pressure: | 213 bar                                      |
| Temperature:      | 65 °C  |
| Detection:        | ELSD, 40 °C, 3.3 bar                         |
| Injection Volume: | 15 µL  |
| Sample Solvent:   | 80/20 ACN/ Water                             |
| Data Rate:        | 10 Hz  |
| Response Time:    | 0.10 sec                                     |
| LC System:        | Shimadzu Nexera X2                           |

Honey can significantly range in quality depending on its purity and levels of sucrose and maltose. Natural honey primarily consists of fructose and glucose, while adulterated honey can contain high levels of sucrose and maltose.

A HALO® Penta-HILIC column separates the primary monosaccharides in pure honey clover showing no signs of adulteration.

### ANALYTES

1. D-(-) Fructose
2. D-(+) Glucose
3. Sucrose
4. D-(+) Maltose



# SWEETENERS

## Caffeine, Sweeteners and Preservatives — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® Excel C18-Amide, 100 Å, 3 µm, 4.6 x 150 mm

Part Number: EXL11121546U

Mobile Phase: A: 0.1% Phosphoric acid in H<sub>2</sub>O

B: 0.1% Phosphoric acid in MeCN

Gradient:

|  | Time (min) | %B |
|--|------------|----|
|  | 0          | 5  |
|  | 12         | 65 |
|  | 15         | 95 |
|  | 16         | 95 |
|  | 16.5       | 5  |
|  | 30         | 5  |

Flow Rate: 1 mL/min

Temperature: 30 °C


Injection: 5 µL

Detection: UV, 214 nm

Instrument: VWR Hitachi Chromaster


### ANALYTES

1. Pyridoxine (B6)
2. Citric acid
3. Pantothenic acid (B5)
4. Acesulfame K
5. Caffeine
6. Aspartame
7. Sorbate
8. Benzoate





# SWEETENERS

## Caffeine, Sweeteners and Preservatives (continued)



The method was used to test two brands of energy drink:



# SWEETENERS

## Fast Separation of Oligosaccharides — HPLC/UHPLC


### TEST CONDITIONS

Column: HALO® Penta-HILIC 90 Å, 2.7 µm, 4.6 x 50 mm  
Part Number: 92814-405  
Mobile Phase A: Water  
Mobile Phase B: ACN  
Gradient: 75-55% B in 10 min  
Flow Rate: 2.0 mL/min  
Pressure: 105 bar  
Temperature: 65 °C  
Detection: ELSD, 40 °C, 3.3 bar  
Injection Volume: 20 µL  
Sample Solvent: 70/30 ACN/Water  
Data Rate: 10 Hz  
Response Time: 0.10 sec  
LC System: Shimadzu Nexera X2

The combination of evaporative light scattering detection (ELSD) and a short 50 mm HALO® Penta-HILIC column enables a fast analysis of oligosaccharides in under 10 minutes whereas traditional columns could have analysis times as long as 30 minutes to more than an hour. Using ELSD eliminates the need to label the sugar with either a UV or fluorescent tag, which simplifies the analysis. Analytes are labeled by degree of polymerization (DP).

### ANALYTES

1. DP 3
2. DP 4
3. DP 5
4. DP 6
5. DP 7
6. DP 8
7. DP 9
8. DP 10
9. DP 11
10. DP 12
11. DP 13
12. DP 14
13. DP 15
14. DP 16
15. DP 17
16. DP 18
17. DP 19
18. DP 20
19. DP 21




# SWEETENERS

## High Resolution Separation of Oligosaccharides – HPLC/UHPLC

### TEST CONDITIONS

|                   |  |
|-------------------|--|
| Column:           | HALO® Penta-HILIC 90 Å, 2.7 µm, 2.1 x 150 mm |
| Part Number:      | 92812-705                                    |
| Mobile Phase A:   | Water  |
| Mobile Phase B:   | ACN  |
| Gradient:         | 75-55% B in 25 min                           |
| Flow Rate:        | 0.5 mL/min                                   |
| Pressure:         | 168 bar                                      |
| Temperature:      | 65 °C  |
| Detection:        | ELSD, 40 °C, 3.3 bar                         |
| Injection Volume: | 20 µL  |
| Sample Solvent:   | 70/30 ACN/Water                              |
| Data Rate:        | 10 Hz  |
| Response Time:    | 0.10 sec                                     |
| LC System:        | Shimadzu Nexera X2                           |

High resolution of oligosaccharides is demonstrated using a dextran ladder on a HALO® Penta-HILIC column with the simple mobile phases of acetonitrile and water. The use of the evaporative light scattering detector (ELSD) eliminates the need to label the sugars with either a UV or fluorescent tag, reducing the time required for sample preparation. Analytes are labeled by degree of polymerization (DP).



# SWEETENERS

## LC-MS Analysis of Stevia Extract — LC-MS


### TEST CONDITIONS

|                      |   |
|----------------------|---|
| Column:              | HALO® Penta-HILIC 90 Å, 5 µm, 3.0 x 250 mm  |
| Part Number:         | 95813-905   |
| Mobile Phase:        | A= 50/50 Water/Acetonitrile with 5 mM Ammonium formate, pH 3<br>B= 5/95 Water/Acetonitrile with 5 mM Ammonium formate, pH 3 |
| Gradient :           | 90% B to 67% B over 30 min  |
| Flow Rate:           | 0.5 mL/min  |
| Pressure:            | 60 bar  |
| Temperature:         | Ambient   |
| Injection Volume:    | 5 µL  |
| Sample Solvent:      | 80/20 Acetonitrile/Water  |
| LC System:           | Shimadzu Nexera   |
| MS:                  | Shimadzu LCMS 2020(single quadrupole)   |
| ESI:                 | +4.5 kV,  |
| Scan range:          | 200-1200 m/z  |
| Scan rate:           | 2 pps   |
| Capillary:           | 250 °C  |
| Heat block:          | 350 °C  |
| Nebulizing gas flow: | 1.5 L/min   |
| Drying gas flow:     | 15 L/min  |

### EXTRACTION PROCEDURE

1. Weigh 400 mg of Stevia rebaudiana leaves (Sigma S5381)
2. Crush leaves with mortar and pestle and transfer to vial
3. Add 8 mL of 50/50 (v/v) acetonitrile/water
4. Sonicate vial contents for 15 minutes
5. Filter sample using 25 mm syringe filter having 0.2 µm PTFE membrane (VWR 28145-495)
6. Centrifuge @ 10K rpm (5 min) and collect supernate
7. Dilute 400 µL of extract in 600 µL of acetonitrile for overall concentration of 80/20 acetonitrile/water
8. Centrifuge diluted sample @ 10K (5 min) rpm and inject the supernate

Stevia is a natural sweetener and is used as a substitute for sugar. LC/MS analysis of stevia glycosides from a Stevia extract is easily accomplished using HALO® Penta-HILIC column due to its unique bonded phase containing five OH groups and the high efficiency of the 5-micron Fused-Core particles.




# SWEETENERS

## Sugars – Orange Juice — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® Excel NH<sub>2</sub>, 100 Å, 3 µm, 4.6 x 150 mm  
Part Number: EXL11141546U  
Mobile Phase: MeCN/H<sub>2</sub>O (75:25 v/v)  
Flow Rate: 1 mL/min  
Temperature: 35 °C  
Injection Volume: 10 µL  
Detection: RI, 35 °C  
Sample: 100 µL orange juice in 900 µL mobile phase  
System: Chromaster 600




# SWEETENERS

## Sugars Separation — HPLC

### TEST CONDITIONS


Column: Avantor® ACE® Excel NH<sub>2</sub>, 100 Å, 5 µm, 4.6 x 150 mm  
Part Number: EXL12141546U  
Mobile Phase: MeCN/H<sub>2</sub>O (70:30 v/v)  
Flow Rate: 1 mL/min  
Temperature: 35 °C  
Injection Volume: 10 µL  
Detection: RI, 35 °C  
Sample: 5 mg/mL in mobile phase  
System: Chromaster 600



## Carotenoid Analysis in Pumpkin—HPLC/UHPLC


### TEST CONDITIONS

|                         |  |                   |                 |
|-------------------------|--|-------------------|-----------------|
| Column:                 | HALO® C30, 160 Å, 2.7 µm, 4.6 x 150 mm | Detection:        | 450 nm          |
| Part Number:            | 92114-730                              | Injection Volume: | 20.0 µL         |
| Competitor:             | FPP C30, 3.0 µm, 4.6 x 150 mm          | Sample Solvent:   | Methanol        |
| Isocratic:              | 100% Methanol                          | Data Rate:        | 14 Hz           |
| Flow Rate:              | 1.5 mL/min                             | Response Time:    | 0.12 sec        |
| Initial HALO® Pressure: | 277 bar                                | Flow Cell:        | 5 µL semi-micro |
| Temperature:            | 15 °C                                  | LC System:        | Agilent 1100    |



### ANALYTES

- 1. Astaxanthin
- 2. Lutein
- 3. Zeaxanthin
- 4. α-Carotene
- 5. β-Carotene
- i. unidentified isomers



Pumpkins contain high amounts of carotenoids, especially beta carotene. Carotenoids are fat-soluble compounds that can be split into two main groups called xanthophylls and carotenes. These compounds both contain anti-oxidant properties and some can be converted into vitamin A when released into the body. A liquid-liquid extraction is performed with 0.2g of pumpkin pulp. Carotenoids are extracted from the pumpkin and analyzed on a HALO® C30 column. The HPLC oven set at sub-ambient temperature enables optimum resolution of early eluting peaks.




## Carotenoids Extracted From Carrot Juice — HPLC/UHPLC

### TEST CONDITIONS

Column: HALO® C30, 160 Å, 2.7 µm, 2.1 x 50 mm  
Part Number: 92112-430  
Isocratic: 100 % Methanol  
Flow Rate: 0.4 mL/min  
Pressure: 100 bar  
Temperature: 30 °C  
Detection: UV 450 nm, PDA  
Injection Volume: 2.5 µL  
Sample Solvent: Methanol/ Isopropyl alcohol  
Data Rate: 40 Hz  
Response Time: 0.025 sec  
Flow Cell: 1 µL  
LC System: Shimadzu Nexera X2


The carotenoids lutein, α- carotene, and β- carotene were isolated from a commercially available carrot juice using liquid liquid extraction. Carotenoids are responsible for the orange color in vegetables such as carrots and are considered antioxidants. The separation was performed on a HALO® C30 column with high resolution between the α- and β-carotene peaks.

### STRUCTURES



### ANALYTES

1. Lutein
2. α-carotene
3. β-carotene
- i. unidentified isomers



## Carotenoids on a C30 Column — HPLC/UHPLC

### TEST CONDITIONS

Column: HALO® C30, 160 Å, 2.7 µm, 3.0 x 150 mm  
Part Number: 92113-730  
Mobile Phase A: Methanol  
Mobile Phase B: Ethanol


| Gradient: | Time (min) | %B |
|-----------|------------|----|
|           | 0.0        | 0  |
|           | 20.0       | 40 |

Flow Rate: 0.65 mL/min  
Temperature: 38 °C  
Detection: UV 471 nm, PDA  
Injection Volume: 0.6 µL  
Data Rate: 2.5 Hz  
Response Time: 2 sec  
Flow Cell: 13 µL  
LC System: Agilent 1100

Data Courtesy of Nature's Sunshine Products

Carotenoids can be split into two main classes called xanthophylls and carotenes. They are responsible for absorbing light for photosynthesis and protecting chlorophyll from photodamage. A separation done by Nature's Sunshine Products shows excellent resolution of carotenoids on a HALO® C30 column.


### STRUCTURES



General structure of a carotenoid

### ANALYTES

1. Lutein
2. *cis*- carotenoid 1
3. *cis*- carotenoid 2
4. α- Carotene
5. β- Carotene
6. *cis*- Lycopene
7. Lycopene




# VITAMINS

## Tocopherols — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® SIL, 100 Å, 5 µm, 4.6 x 250 mm  
Part Number: ACE1272546  
Mobile Phase: Hexane/IPA (98:2 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Injection Volume: 1 µL  
Detection: UV-Vis, 450 nm




# VITAMINS

## Vitamin K—UHPLC


### TEST CONDITIONS

Column: Avantor® ACE® Excel CN-ES, 100 Å, 2 µm, 3.0 x 100 mm  
Part Number: EXL10131003U  
Mobile Phase: 95:5 MeOH/H<sub>2</sub>O  
Flow Rate: 1 mL/min  
Temperature: 30 °C  
Injection: 5 µL  
Detection: UV, 270 nm


### STRUCTURES




Vitamin K3 (Menadione)




cis-Vitamin K1 (2,3-cis-phyloquinone)



trans-Vitamin K1 (2,3-trans-phyloquinone)



Vitamin K2 (MK-7) - Menaquinone 7



### ANALYTES

1. Vitamin K3 (Menadione)
2. cis-Vitamin K1 (2,3-cis-phyloquinone)
3. trans-Vitamin K1 (2,3-trans-phyloquinone)
4. Vitamin K2 (MK-7) - Menaquinone 7

# VITAMINS

## Vitamins - Water Soluble (Gradient) C18 — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® Excel C18, 100 Å, 5 µm, 4.6 x 250 mm

Part Number: EXL1212546U

Mobile Phase: A: 50 mM KH<sub>2</sub>PO<sub>4</sub> pH 3.0 in H<sub>2</sub>O

B: MeOH

Gradient:


| Time (min) | %B |
|------------|----|
| 0          | 3  |
| 5          | 3  |
| 45         | 45 |
| 50         | 80 |

Flow Rate: 1 mL/min

Temperature: Ambient

Injection Volume: 1 µL

Detection: UV, 205 nm



### ANALYTES

1. Pyridoxamine
2. Thiamine (Vitamin B1)
3. L-Ascorbic acid (Vitamin C)
4. Nicotinic acid
5. Pyridoxal
6. Impurity
7. Pyridoxine
8. p-Aminobenzoic acid
9. Pantothenic acid (Vitamin B5)
10. Folic acid (Vitamin B9)
11. Cyanocobalamin (Vitamin B12)
12. d-Biotin (Vitamin B7)
13. Riboflavin (Vitamin B2)

# VITAMINS

## Vitamins - Water Soluble (Gradient) C8—HPLC

### TEST CONDITIONS

Column: Avantor® ACE® C8, 100 Å, 5 µm, 4.6 x 250 mm

Part Number: [ACE1222546](#)

Mobile Phase: A: 50 mM KH<sub>2</sub>PO<sub>4</sub> pH 2.5 in H<sub>2</sub>O

B: MeOH


Gradient:

| Time (min) | %B |
|------------|----|
| 0          | 0  |
| 3          | 0  |
| 16.5       | 45 |
| 19.5       | 80 |

Flow Rate: 1 mL/min

Temperature: Ambient

Detection: UV, 205 nm



# VITAMINS


## Vitamins – Fat Soluble — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® Excel C18-Amide, 100 Å, 3 µm, 4.6 x 150 mm  
Part Number: EXL11121546U  
Mobile Phase: MeOH/MeCN (90:10 v/v)  
Flow Rate: 1 mL/min  
Temperature: 20 °C  
Detection: UV, 280 nm

### ANALYTES


1. Menadione (Vitamin K3)
2. Retinol (Vitamin A)
3. Vitamin A acetate
4. Ergocalciferol (Vitamin D2)
5. Cholecalciferol (Vitamin D3)
6. Vitamin E acetate
7.  $\alpha$ -Tocopherol (Vitamin E)
8. Vitamin K1




# VITAMINS

## Vitamins – Fat Soluble *(continued)*


## STRUCTURES




#### Mengdione (Vitamin K3)




#### **Retinol (Vitamin A)**




#### Vitamin A acetate




#### Ergocalciferol (Vitamin D<sub>2</sub>)




#### Cholecalciferol (Vitamin D3)



#### Vitamin E acetate



### **α-Tocopherol (Vitamin E)**



## Vitamin K1

## Vitamins D2 and D3 — UHPLC




### TEST CONDITIONS

Column: Avantor® ACE® Excel C18-Amide, 100 Å, 2 µm, 3.0 x 50 mm  
Part Number: EXL10120503U  
Mobile Phase: 100% MeCN  
Flow Rate: 0.43 mL/min  
Temperature: 20 °C  
Injection: 2 µL  
Detection: UV, 265 nm

### ANALYTES

1. Ergocalciferol (D2)
2. Cholecalciferol (D3)

### STRUCTURES




## Water Soluble Vitamins (Polar Molecules) — HPLC

### TEST CONDITIONS

Column: Avantor® ACE® C18-AR, 100 Å, 3 µm, 4.6 x 150 mm  
Part Number: ACE1191546  
Mobile Phase: 0.1% Phosphoric acid in H<sub>2</sub>O/MeOH (96.5:3.5 v/v)  
Flow Rate: 1 mL/min  
Temperature: 22 °C  
Injection: 2 µL  
Detection: UV, 260 nm

### ANALYTES

1. Pyridoxamine (Vitamin B6)
2. Thiamine (Vitamin B1)
3. Isonicotinamide
4. Nicotinamide
5. L-Ascorbic acid (Vitamin C)
6. Orotic acid
7. Hypoxanthine
8. Pyridoxal (Vitamin B6)
9. Pyridoxine (Vitamin B6)
10. p-Aminobenzoic acid



# XANTHINES

## Caffeine and Metabolites — HPLC







### TEST CONDITIONS

Column: Avantor® ACE® Excel SuperC18, 100 Å, 5 µm, 4.6 x 150 mm  
Part Number: EXL12111546U  
Mobile Phase: A: 20 mM Ammonium acetate pH 7.0 in H<sub>2</sub>O  
B: 20 mM Ammonium acetate pH 7.0 in MeCN/H<sub>2</sub>O (90:10 v/v)

| Gradient:            | Time (min) | %B |
|----------------------|------------|----|
|                      | 0          | 2  |
|                      | 45         | 15 |
|                      | 48         | 15 |
|                      | 49         | 2  |
| Post time 10 minutes |            |    |

Flow Rate: 1 mL/min  
Temperature: 60 °C  
Injection: 1 µL  
Detection: UV, 273 nm

### STRUCTURES



### ANALYTES

1. Xanthine
2. Theobromine
3. Paraxanthine
4. Theophylline
5. Caffeine

# XANTHINES

## Caffeine and Related Compounds Using HILIC Mode—HPLC


### TEST CONDITIONS

Column: Avantor® ACE® HILIC-N, 100 Å, 5 µm, 4.6 x 150 mm  
Part Number: HILN51546U  
Mobile Phase: A: 10 mM Ammonium formate pH 3.0 in MeCN/H<sub>2</sub>O (96:4 v/v)  
B: 10 mM Ammonium formate pH 3.0 in MeCN/H<sub>2</sub>O (50:50 v/v)


| Gradient: | Time (min) | %B  |
|-----------|------------|-----|
|           | 0          | 0   |
|           | 15         | 100 |

Flow Rate: 1.5 mL/min  
Temperature: 15 °C  
Injection: 2 µL  
Detection: UV, 275 nm  
Sample: Caffeine (25 mg/mL) with related substances theophylline, theobromine, xanthine and hypoxanthine at 0.5% w/w dissolved in MeCN/H<sub>2</sub>O (90:10 v/v)


### STRUCTURES




Caffeine




Theophylline




Theobromine



Xanthine



Hypoxanthine





### ANALYTES

1. Caffeine
2. Theophylline
3. Theobromine
4. Xanthine
5. Hypoxanthine

## Caffeine in Soda — HPLC

### TEST CONDITIONS


|                 |   |                   |                              |
|-----------------|---|-------------------|------------------------------|
| Column:         | HALO® C18, 90 Å, 5 µm, 3.0 x 50 mm with guard column      | Temperature:      | 30 °C                        |
| Part Number:    | <a href="#">95813-402</a>                                 | Injection Volume: | 1.0 µL                       |
| GC Part Number: | <a href="#">95813-102</a>                                 | Sample Solvent:   | (caffeine std.) mobile phase |
| Mobile Phase:   | 75/25: A/B<br>A= 0.1% Formic acid in water<br>B= Methanol | Detection:        | UV 254 nm, VWD               |
| Flow Rate:      | 0.8 mL/min  | Response Time:    | 0.02 sec                     |
| Pressure:       | 120 Bar   | Date rate:        | 25 Hz                        |
|                 |   | Flow Cell:        | 2.5 µL semi-micro            |
|                 |   | LC System:        | Shimadzu Prominence UFC XR   |
|                 |   | ECV:              | ~14 µL                       |



# XANTHINES

## Caffeine in Soda (*continued*)

### STRUCTURE:



Caffeine

| Sample                      | Caffeine tested mg/(355 mL) | Can value mg/(355 mL) |
|-----------------------------|-----------------------------|-----------------------|
| Store brand cola 1          | 12                          | N/A                   |
| Cola 2                      | 53                          | 54                    |
| Cola 3                      | 43                          | 43                    |
| Cola 4                      | 36                          | 38                    |
| Cola 5                      | 38                          | 38                    |
| Store brand diet cola 1     | 12                          | N/A                   |
| Diet cola 2                 | 45                          | 46                    |
| Diet cola 3                 | 34                          | 34                    |
| Diet cola 4                 | 36                          | 35                    |
| Energy drink 1*             | 160                         | 160                   |
| Energy drink 2**            | 79                          | 80                    |
| Diet Energy drink**         | 79                          | 80                    |
| Non-cola drink 1            | 53.3                        | 54                    |
| Non-cola drink 2            | 22                          | 22                    |
| Diet non-cola drink         | 43                          | 41                    |
| Diet cola 1 non caffeinated | 0                           | N/A                   |
| Diet cola 2 non-caffeinated | 0                           | N/A                   |
| Diet cola 3 non-caffeinated | 0                           | N/A                   |

355 mL = 12 oz

\*amount in 16 oz (473 mL) cans

\*\*amount in 8.4 oz (248 mL) cans

### DISCUSSION:

A selection of sodas was purchased along with several energy drinks at a local grocery store. An approximate 10 mL sample from each soda was placed in a separate 20 mL scintillation vial and capped. The vials were shaken several times and then the cap was loosened and the vial sonicated for 10 minutes to remove CO<sub>2</sub>. Then a 1 mL aliquot was placed into a 1.5 mL HPLC sample vial. A one microliter quantity was injected into the HPLC under conditions tabulated elsewhere on this page. A guard column was used to prevent the buildup of a brown material on the analytical column packing. The material was likely caramel coloring.

The chromatogram shown is from a regular cola drink.

A standard curve of peak area vs. caffeine concentration was made over the range of 0.11–0.44 mg/mL. These values were used to calculate the line fit of Y=mX + b. From the peak response the concentration of (mg of caffeine)/mL was calculated and then multiplied by the number of mL in the beverage can (usually 355 mL). Results are shown on the table to the left.



SMARTER CHROMATOGRAPHY

Phone: (610) 358-9696 | Toll-free: 1 (800) 441-7508

103 Commons Court | Chadds Ford, PA 19317

<https://www.mac-mod.com>

[technical@mac-mod.com](mailto:technical@mac-mod.com)

Avantor® and ACE® are registered trademarks of Avantor Inc.

HALO® and Fused-Core® are registered trademarks of Advanced Materials Technology, Inc.

MAC-MOD is an approved distributor of Avantor® ACE®, and HALO® products in the US.