

Application note # C-13032

Fast Separation of Fragrance Allergens using the Avantor® Hichrom HI-1701 Phase




Figure 1: Fast separation of fragrance allergens using the Avantor® Hichrom HI-1701 phase.

Method Details

CONDITIONS

Oven Program:	50 °C (0.1 min), 50 °C/min, 250 °C (5 min)
Carrier Gas:	Hydrogen, 0.5 mL/min
Injector:	Split 230 °C, 0.5 µL 1:300 split ratio
Detector:	FID, 250 °C
Sample Dilution:	1% in cyclohexane

The HI-1701 is a mid-polarity phase with a crossbonded 14% cyanopropylphenyl, 86% methyl polysiloxane composition. It is well suited for a variety of applications including pesticides, residual solvents, allergens and essential oils.

PEAK IDENTIFICATION

1. Limonene	11. Eugenol	21. Farnesol 2
2. Linalool	12. Anisole and Cinnamol	22. Farnesol 3
3. Benzyl alcohol	13. Alpha-Isomethyl ionone	23. Amyl cinnamol
4. Veratrol	14. Isoeugenol	24. Hexyl cinnamal
5. Methyl 2-octanoate	15. Lilial	25. Cinnamal
6. Citronellol	16. Coumarin	26. Benzyl salicylate
7. Unknown	17. Amyl cinnamal	27. Benzyl cinnamate
8. Neral	18. Lyrall 1	
9. Geranial	19. Lyrall 2	
10. Hydroxycitronellal	20. Farnesol 1	

ORDERING TABLE

Product	Details	Dimensions	Part Number
Avantor® Hichrom HI-1701	FAST GC Column	0.10 mm, 0.10 µm, 5 m	HI06-10-010-5