

## Conditions

Column: ACE Excel 1.7 C18-Amide  
 Dimensions: 100 x 2.1 mm  
 Part Number: EXL-1712-1002U  
 Mobile Phase: A: 0.01% formic acid in H<sub>2</sub>O  
 B: 0.01% formic acid in MeCN

Time (mins)	%B
0.0	3
2.5	10
8.0	100
8.5	3
10.0	3

Flow Rate: 0.5 mL/min  
 Detection: Exactive accurate mass MS system  
 ESI in negative ion mode

Sample: Analytes between *m/z* 70-800 monitored  
 Metabolites from green tea extracted into cold water by vortexing for 20 mins. Samples filtered prior to injection onto column and modular Accela LC system

