# High-Speed HPLC Separations with Polar Fused-Core® Packings

#### Abstract

Separations by reversed-phase HPLC are most often performed using C18 bonded phase columns. The reasons for this popularity are 1) the C18 phase was the first reversedphase packing in the commercial market, 2) the C18 phase has a reputation for being a very robust phase, and 3) there is a vast library of applications on the C18 phase on which to base starting conditions for new applications. Other reversed phase HPLC columns such as cyanopropyl, phenyl, amino, etc. are used significantly less often in published applications, but they do provide chromatographers with alternate selectivities for compounds relative to C18 phases. These other phases are used when C18 columns cannot separate compounds of interest or when orthogonal separations are required. This presentation introduces two new columns based on fused-core® particles that will provide alternate selectivities to the existing C8 and C18 fused-core columns. The advantages of the combination of the high-speed characteristics of the fused-core particles and the alternate selectivities of HILIC and embedded polar stationary phases will be demonstrated.

# Fused-Core Particle Technology

#### Cartoon of Particle



#### SEM Photo of Sliced Particles



#### Specifications for Fused-core Particles

- 2.7 µm overall diameter
- 1.7 µm solid core
- 0.5-µm-thick porous outer shell
- 9 nm pores in porous shell
- Highly purified Type B silica
- Surface area ~150 m2/g
- Pressures to at least 600 bar

#### Unusually Narrow Particle Size Distribution



# J. J. DeStefano, W. L. Johnson, and J. J. Kirkland Advanced Materials Technology, Inc. Wilmington, DE 19810

## Theory

Van Deemter equation

#### $\mathbf{H} = \mathbf{A} + \mathbf{B}/\mathbf{u} + \mathbf{C}\mathbf{u}$

- A-term dominated by packed bed homogeneity
- C-term related to resistance to mass transfer in the porous structure
- Small A and C-terms preferred
- Short diffusion path due to thin shell results in smaller C-term
- High efficiencies (h<2!) due to uniform particle size leading to smaller A-term
- Moderate backpressure due to particle size being larger than sub-2-µm particles

#### van Deemter Plots of Totally Porous vs. Fused-Core Particles



Columns 50 x 4.6mm; Mobile phase: 60% ACN/40% water Bonded phase: CI8;Temperature: 24 °C; Solute: naphthalene;

## Fused-Core vs. Sub-2-Micron



## Back Pressure Plots for High-speed Columns



#### Small-particle Column Comparison



#### Separation of Explosives



Column: 4.6 x 50 mm Halo C18; Mobile phase: 27% Methanol/73% Water; Temperature: 40°C; Flowrate: 3.3 ml/min; Detector: UV @ 254 nm; Pressure: 343 bar; Agilent 1100

## New Polar Phases Provide Alternate Selectivity to CI8

HALO RP-Amide (Polar-embedded Phase) HALO HILIC (Bare Fused-Core Silica)











# Polar Solutes





New HILIC Phase Provides Good Retention for



## Conclusions

- Efficiency and separation speed rivals that of sub-2-micron columns, but with about one-half the pressure drop
- Produce unusual column efficiencies; reduced plate heights *h* of ~1.5 for small molecules
- Thin porous outer shell results in superior mass transfer kinetics and better efficiency at high mobile phase velocities, especially for larger molecules
- New polar phases offer alternatives for separation selectivity and easier method development