


High Temperature/ Low pH Stability of HALO 1000 Å ES-C18, 2.7 μm



PEAK IDENTITY:

1. Trastuzumab

STRUCTURES:

TEST CONDITIONS:

Column: HALO 1000 Å ES-C18, 2.7 μm, 2.1 x 50mm

Part Number: 92712-402

Mobile Phase A: Water/ 0.1% TFA

Mobile Phase B: Acetonitrile/ 0.1% TFA

Gradient:	<u>Time</u>	<u>%B</u>
	0.0	32
	4.0	38

Flow Rate: 0.4 mL/min

Initial Pressure: 81 bar

Temperature: 80 °C

Detection: UV 280 nm, PDA

Injection Volume: 1.2 μL

Sample Solvent: Water

Data Rate: 40 Hz


Response Time: 0.025 sec

Flow Cell: 1 μL

LC System: Shimadzu Nexera X2



1000 Å 2.7μm particle



ES-C18 bonded phase

Trastuzumab (MW ~148 kDa) is a monoclonal antibody used to treat breast cancer. A stability experiment using a HALO 1000 Å ES-C18 column shows excellent reproducibility for 500 injections of trastuzumab. The sterically protected C18 bonded phase enables rugged stability at the elevated temperature and low pH conditions that are typically used for protein analysis.