



## 1. Introduction

- Exploring **chromatographic selectivity** is a **powerful** approach for LC method development.
- Efficient **method development** requires a **logical exploration** of the key chromatographic parameters affecting selectivity.
- Rationally designed method development strategies **assess key parameters** and allow **well informed decision making** leading to **robust stationary phase / mobile phase selection**.
- Method development strategies based on **column / mobile phase screening** and optimisation are commonly utilised for reversed phase.
- For **HILIC**, such strategies are **less common / less well defined**.
- This poster demonstrates a **simple, step-by-step approach** to HILIC method development, based on the **concept of exploring column selectivity**.



## 2. Selectivity in HILIC

- In HILIC, the **column stationary phase** has a significant effect on chromatographic **selectivity**.
  - The **ACE HILIC range** consists of three complementary phases **specifically designed to offer maximum selectivity differences** – ideal for method development:
- 

**ACE HILIC-B**

$S = 100 \times \sqrt{(1 - r^2)}$


$S = 0$ , Identical selectivity

$S = 100$ , Completely orthogonal

**ACE HILIC-A**      **ACE HILIC-N**


\*Generated from comparison of retention data for 54 compounds at pH 3.0
- ACE HILIC-A**
  - ACE HILIC-B**
  - ACE HILIC-N**
- Mobile phase pH** is also a powerful parameter and can affect ionisation of analytes and the **stationary phase itself**.
  - Method development strategies based on **screening different stationary phases and mobile phase pH** are therefore the optimum choice.
  - Buffer concentration** and **temperature** are **less influential**, however can be used to **fine-tune methods**.

## 3. ACE MD Protocol




## 5. Caffeine and Related Substances: Optimisation

- Isocratic conditions** were assessed. **Retention of later eluting compounds increased** but **resolution of the critical pairs was not improved**.
- Reducing the temperature improved the resolution** of both impurity peak pairs.
- Lower temperature** was used to achieve separation on the **ACE 5 HILIC-N** at pH 3.0.
- A **small increase in acetonitrile** in the **gradient starting conditions** was also found to be beneficial



## 4. Example: Caffeine and Related Substances

- A caffeine and related substances **HILIC method** was developed.
- All analytes are **polar neutral** with **negative logP values**.
- The **ACE HILIC-A, HILIC-B and HILIC-N** phases were screened using gradient conditions at pH 3.0 and 4.7 (Table 1).



First published in Chromatography Today, Volume 8, Issue 4, Nov/Dec 2015

## 6. Summary

- A **systematic and rationally designed method development strategy** can aid in **streamlining the method development process**.
- In **HILIC**, **column stationary phase** and **mobile phase pH** are the most **critical parameters** affecting selectivity.
- The **step-by-step method development strategy** proposed in this poster therefore provides a powerful means by which to **probe selectivity** of a new application.
- The **ACE HILIC-A, HILIC-B and HILIC-N** phases provide **complimentary selectivity**, ideal for **method development**.
- Screening** an analyte mixture on these three phases has been demonstrated as an effective method development strategy for **selecting an appropriate stationary phase/mobile phase combination**.
- Optimisation can be achieved** by altering parameters such as **ionic strength, % organic and temperature**.