



# INTACT PROTEIN AND GLYCOPROTEIN SEPARATIONS BY HILIC

Barry Boyes<sup>1,2</sup>, Ron Orlando<sup>2</sup>, and Joseph DeStefano<sup>1</sup>

<sup>1</sup>Advanced Materials Technology, Inc.

Wilmington, DE USA

<sup>2</sup>Complex Carbohydrate Research Center

University of Georgia, Athens, GA USA

bboyes@advanced-materials-tech.com

# Agenda

- HILIC separations have utility in glycoscience applications
- Polar interactions for glycans, glycopeptides and glycoproteins lead to useful separations
- Materials and methods of use for such separations are improving, but use remains complex
- Will show examples of protein separations based on Glycoform structure


# HILIC for Glycoscience

## Utility for Oligosaccharide Profiling


### Penta-HILIC Separations of Abundant bov. *a1*-AGP N-glycans

2.1 mm ID x 150 mm; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600 mL/min.


MALDI of Permethylated Glycans



HILIC of End-labeled Glycans



*trisialated-biantennary *a1*-AGP N-glycans*





# Peptide Separations by HILIC

## Generalities:

- acidic conditions (TFA, FA, HOAc, Phosphoric, etc.)
- can exhibit orthogonality in retention, relative to RP
  - M. Gilar, et al., Anal. Chem. 82(2010) 265, and others
  - Boyes, et al., ISPPP 2013
- does not have to be at a disadvantage for efficiency
- can be very MS friendly
- can show some solubility problems
- additional “polar modes” of operation are possible, eg., ERLIC, HILIC/SALT (Mant, et al., J. Chromatogr. 1277 (2013))
- *can allow separations that are very difficult using RP*

# Penta-HILIC Strongly Retains N-linked Glycopeptides

2.1 mm ID x 100 mm, 0.35 mL/min, 40 C, MS: SQ TIC (+ 300-2000 m/z) @ 0.35/s  
20 µg Bovine Ribonuclease B tryptic digest (CAM)




# LC/MS of O-GlcNAcylated Peptides

2.1 x 100 mm Halo Penta-HILIC, 0.4 mL/min at 60°C.

Gradient conditions: A – 0.1% formic acid/10 mM ammonium formate; B – 90% Acetonitrile in A.

500 pmol, MS operated +4.5 kV, with 0.45 s scan from 500 to 2000 m/z (TIC)

RP Gradient - 4% to 34% AcN/30 min (1%/min); HILIC Gradient - 90% to 60% AcN/30 min (-1%/min).



APP695-14GPeP

APP695-14Pep

VPTT(O GlcNAc)AASTPDAVK

VPTTAASSTPDAVK

See Poster P-W-2304

# Protein Separations by HILIC

- Limited use relative to RP
  - Variety of examples by Alpert and collaborators
  - Tetaz, et al., J. Chromatogr. 1218 (2011) 5892.
  - Carroll, et al., P.N.A.S. U.S.A. 103 (2006) 16170.
  - Focused on membrane proteins, histone and lipoprotein PTM modifications
- Expected to show very different selectivity than Reversed Phase separations of proteins
- Challenge of using high organic mobile phases for proteins
- Not clear what features are needed to improve intact protein separations – mobile phase, additives, stationary phases, or combinations of all

This could be a longer term project!

# Operational Considerations


- Mobile phases components:
  - TFA, Formic Acid, FA/AF (0.05%-1%)
  - Acetonitrile, *i*-propanol, *n*-propanol
  - HFIP, TFE, others

Solubility is always a consideration for both salts and proteins
- Stationary phase characteristics:
  - 2.7 µm Halo Penta-HILIC
    - hydroxylated neutral HILIC bonded phase
    - Fused-Core particle with 0.5 um SPP layer
  - 90 Å and 160 Å pore size
- Aqueous soluble proteins and glycoproteins
  - Samples adjusted to 1-3% Formic Acid in High Organic MP
- Combination of absorbance and MS detection, where possible


# Effect of Acid Modifier

2.1 x 100 mm Halo Penta-HILIC 90 Å pore size, 0.5 mL/min at 50°C.


Gradient conditions: A – 0.1% acid; B – 90% Acetonitrile in A. Gradient - 90% to 50% AcN/20 min (-2%/min).



0.1% TFA



0.1% FA




0.1% FA/10 mM AF


# Current Consensus MP & Pore Size

- FA/AF buffer (0.2%/20 mM) with 0.5% HFIP (does not HURT!)
- Mixture of Acetonitrile/*i*-Propanol (40%/35%) to 10% iPA
- Elevated column temperature: 40-60°C
- Pore Size is not a clear factor for proteins
- Adjusting sample solvent is crucial – addition of FA to 1-2% resolves many solubility problems

2.1 x 100 mm Halo Penta-HILIC, 0.5 mL/min at 50°C; Gradient conditions: B – 0.2% FA/20 mM AF/40% AcN/35% IPA/ 0.5% HFIP; A – 0.2% FA/20 mM AF/10% IPA Gradient - 90% to 55% B in 20 min.




160 Å pore  
- Lower Rt  
- Slightly Lower PW



90 Å pore

# Bovine Ribonuclease as a Glycoprotein



- Composition:
  - RNase A - 124 aa residues, 4 disulfide bridges, 13.7 kDa
  - RNase B – Asn<sup>34</sup> is the site of glycosylation, with 5 high mannose structures M5-M9 , leading to enzyme heterogeneity, 14.9-15.5 kDa
  - In addition to compositional heterogeneity, branch variants of the high mannose oligosaccharides M5, M7, M8 are indicated



Prien, et al., J. A.S.M.S. 20 (2009) 539

# RNase B N-Linked Glycans

Release of protein N-linked glycans using PNGase F yields oligosaccharides with a free reducing terminus (alditol) that is readily labeled by amines via formation of a Schiff's base, which can be reduced readily.




Procainamide  
Abs: 300 nm  
Fl: Ex 330/Em 380 nm  
Mass: glycan + 219.32

Many amines have been applied to labeling glycans, (Harvey, 2011, J. Chrom. 879).


Procainamide is favored due to improvements in ESI-MS detection (Klapoetke, et al., 2010, J. Pharm. Biomed. Anal. 53), and favorable fragmentation conditions in CID.

*2.1 mm ID x 150 mm Halo Penta-HILIC; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600 mL/min. Detection: 300 nm Abs; ESI-MS (MS-2020, (+) 4.2 kV, 400-2000 with SIM)*



# Separation of RNases by RP




2.1 x 100 mm Halo Peptide ES-C18, 0.5 mL/min at 50°C; Gradient conditions: A – 0.1% TFA/water, B – 0.1% TFA in 80% AcN; Gradient - 30% to 40% B in 20 min (0.4% AcN/min), 0.5 µg each.



RP allows Resolution of A/B, little resolution of B glycoforms

# HILIC Separation of RNase B Glycoforms


2.1 x 100 mm Halo Penta-HILIC, 0.35 mL/min at 50°C; Gradient conditions: B – 0.2% FA/20 mM AF/40% AcN/35% IPA/ 0.5% HFIP; A – 0.2% FA/20 mM AF/10% IPA Gradient - 85% to 75% B in 20 min.



ID	Mass (Da) Theory	Mass (Da) Observed
RNA A	13,681	13,686
RNA B M5	14,898	14,903
RNA B M6	15,060	15,068
RNA B M7	15,222	15,227
RNA B M8	15,384	15,390
RNA B M9	15,564	15,550


# Isolated RNase B Glycoforms

Separation Conditions: 4.6 x 250 mm Penta HILIC 2.7  $\mu$ m 90 $\text{\AA}$ , 60 C, 0.4 mL/min, 215 nm, 78-75 % B over 60 min, B – 40% AcN/35% IPA/0.2%FA 20 mM AF/0.5% HFIP




# Characterization of Resolved RNase B Glycoforms

*High Resolution ESI/MS (Orbitrap) Infused in 0.5% FA/50% AcN*




*Analysis of PNGase F released glycans from M5, M6 and M9 Proteins*



# Isolated RNase B Glycoforms

Assay Conditions: 50 mM NaOAc, pH 5.0, 25°C, 0.5 mg/mL RNA, ΔA300



# More on Glycoproteins

- Protein structure contributes to retention
- Mass and nature of glycans will contribute to retention
- Not all glycoproteins are as simple as RNase B
- LC/MS resolution requirements increase a lot with size
- MS Sensitivity is a problem as size and glycosylation rise

116

VEROSTEK, LUBOWSKI, AND TRIMBLE

TABLE 3  
Characteristics of Glycoproteins Chosen to Test the Organic Solvent Precipitation/Extraction Methodology


Glycoprotein	Source	Monomer (kDa)	Native form	CH <sub>2</sub> O <sup>a</sup> Content (g/100 g)	Oligosaccharides/monomer			Amount treated (mg)
					N-linked	High mannose	Complex	
RNase B	Bovine	13	Monomer	11		1	0	1, 20
Ovalbumin	Hen	43	Monomer	3		1	0	5 to 60
Invertase	Wt yeast	60	Octamer	50	9-10	0	0	20
	Sec18 yeast			28	9-10	0	0	10
Thyroglobulin	Bovine	330	Dimer	9	5-6	13-14	0	20, 40
Fetuin	Bovine	48	Monomer	23	0	3	3	40
BSSL	Human/ <i>P. pastoris</i>	76.3	Dimer	16	0.7	0	27	65, 80, 100

<sup>a</sup> Carbohydrate.


Verostek, et al., Anal. Biochem. 278 (2000) 111.

# Real Glycoprotein HILIC HUMPograms


2.1 x 100 mm Halo Penta-HILIC, 0.4 mL/min at 50°C; Gradient conditions: B – 0.2% FA/20 mM AF/40% AcN/35% IPA/ 0.5% HFIP; A – 0.2% FA/20 mM AF/10% IPA. Gradient - 90% to 50% B in 15 min.



Intact Fetuin



Asialo-Fetuin



PNGase F  
Fetuin  
(O-glycosylated)

# Conclusions

- HILIC separations methods have impact in glycobiology at the glycan, glycopeptide and glycoprotein structure levels.
- Results with intact protein and glycoprotein LC/MS separations of smaller proteins performed well, but expanding to a broader range of proteins requires further effort on mobile phase and use conditions (sample solvent dependencies).
- Selectivity for glycan separations, when attached to peptide and protein structures, appears excellent.
- RNase B glycoforms can be well resolved, under conditions that retain activity.




# Thank you for your attention!

## Also thanks to:


- AMT – Jack Kirkland, Stephanie Schuster, William Miles and Bob Moran for advice and technical assistance.
- NIH for Financial Support - GM093747

# Protein and Glycoprotein Surfaces


Glycan positions are shown for the ten best (red) and poorest (blue) fitting models.




RNaseB



a1AGP



Fetuin



# Separation of RNase B PNGase F Released N-Linked Glycans

## Penta-HILIC Separations of Labeled Oligosaccharides and Glycans

2.1 mm ID x 150 mm; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600 mL/min.

Detection: 300 nm Abs; ESI-MS (MS-2020, (+) 4.2 kV, 400-2000 with SIM)

