

# Investigating the Effects of Chromatographic Parameters on Column Equilibration in Isocratic and Gradient HILIC Separations


Alan P McKeown<sup>1</sup>, Ed Faden<sup>2</sup> and Geoff Faden<sup>2</sup>

<sup>1</sup>Advanced Chromatography Technologies Ltd, 1 Berry Street, Aberdeen, Scotland, AB25 1HF UK

<sup>2</sup>MACMOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317 USA

## 1. Introduction

- Hydrophilic Interaction Liquid Chromatography (HILIC) is a powerful technique for the separation of **hydrophilic** and **polar compounds**.
- HILIC utilises a **polar stationary phase** and **high-organic** containing mobile phase.
- The HILIC retention mechanism is **complex** and includes **electrostatic**, **polar** and **partitioning interactions**.
- The key to HILIC retention is formation of a **water-enriched layer** at the particle surface.
- Column equilibration** to form a stable water-enriched layer is **vitaly important** for **reproducible chromatography**.
- Long and variable** equilibration times are often cited as a negative aspect of the HILIC technique.
- This work aims to assess **key isocratic** and **gradient** method parameters and how they **affect column equilibration**.



## 2. Experimental

- The **three ACE HILIC phases** (100 x 3.0 mm, 3 µm) were assessed in this study:
- HILIC-A** (An acidic stationary phase showing high cation exchange capability)
- HILIC-B** (A basic character phase with reasonable anion exchange capacity)
- HILIC-N** (A neutral bonded phase with low anion and cation exchange capabilities)
- Isocratic equilibration:** The following parameters were assessed:

Parameter	Values assessed
Fresh vs used column	Freshly packed column vs previously used column
Mobile phase ionic strength	5, 10 and 20 mM ammonium formate pH 3.0
Percentage water in the mobile phase	6, 10, 15 and 20 % water

- Gradient equilibration:** The **effect of post gradient re-equilibration** time was investigated (re-equilibration times equivalent to 10, 20, 30 and 200 column volumes).

For all experiments, a test mix containing 4-hydroxybenzoic acid (acid), salbutamol (base) and 2'-deoxyguanosine (polar neutral) was injected.


## 3. Fresh vs Used Columns

- Freshly packed columns were equilibrated with mobile phase for 400 column volumes. Sample injected every 10 mins.
- The buffer was removed from the column (MeCN:H<sub>2</sub>O 1:1) and stored for 48 h in IPA.
- Columns were then re-equilibrated with mobile phase for 400 column volumes and sample injected every 10 mins.

$$\% \text{Equilibrated}_{t=x} = (k_x/k_{final}) \times 100$$


- Used columns** were found to **equilibrate faster** than **freshly packed columns**, potentially due to the presence of **residual non-aqueous packing solvents** in pores of freshly packed columns.
- Equilibration times** for used **ACE HILIC** columns were **similar to the ACE C18**.
- Under **these conditions**, **fresh bonded phases** appear to equilibrate faster than **bare silica**.

Mobile phase: 10 mM Ammonium formate pH 3.0 in MeCN/H<sub>2</sub>O 90:10 (v/v)  
Flow rate: 0.6 mL/min Temp: 35 °C Detector: 214 nm Injection volume: 5 µL



## 4. Ionic Strength


- All three phases were equilibrated with 5, 10 and 20 mM ammonium formate pH 3.0 in MeCN:H<sub>2</sub>O 9:1 (v/v) for 400 column volumes. Sample injected every 10 mins.
- All three phases were >95% equilibrated after 20 column volumes with 5 and 10 mM buffer concentrations.
- The **20 mM** ionic strength mobile phase required **slightly longer** to equilibrate (up to 50 column volumes on the HILIC-A).



Mobile phase: 5, 10 or 20 mM Ammonium formate pH 3.0 in MeCN/H<sub>2</sub>O 90:10 (v/v)  
Flow rate: 0.6 mL/min Temp: 35 °C Detector: 214 nm Injection volume: 5 µL


## 5. Percentage Organic

- All three ACE phases were equilibrated with 10 mM ammonium formate pH 3.0 in MeCN:H<sub>2</sub>O at different volume fractions for 400 column volumes. Sample injected every 10 mins.
- 20 column volumes** was sufficient to equilibrate (>95%) all ACE phases under most conditions.
- 94% MeCN** took **>100 column volumes** to equilibrate on the ACE HILIC-B and HILIC-N.
- This implies that formation of a **stable water layer** takes longer when using **low-aqueous containing mobile phases** in HILIC mode for **bonded phases**, possibly due to the presence of a more substantial water-enriched region at the particle surface.




## 6. Isocratic: Observations

- Similar trends** were observed for acidic (4-hydroxybenzoic acid), neutral (2'-deoxyguanosine) and basic (salbutamol) analytes.
- Stable retention** was achieved on **all ACE phases** once equilibrated.
- Fresh columns** took **longer** to equilibrate than used columns
  - Minimum **20 column volumes** for **used** columns
  - Minimum **60-80 column volumes** for **fresh** columns
- Higher ionic strength** mobile phases may require **longer to equilibrate**.
- Formation of a stable water-enriched layer may be **substantially slower** when using mobile phases containing **~5-6% water**.



## 7. Gradient Equilibration

- A **HILIC gradient** was run on each phase and the **post-gradient equilibration time** varied:
  - 10, 20, 30 and 200 column volumes
  - 1 column volume = 0.75 minutes (100 x 3.0 mm column, 0.6 mL/min)
  - 10 injections performed, 1<sup>st</sup> injection discarded
- The **length of the re-equilibration** stage was found to affect **analyte retention**.



## 8. Gradient Equilibration

- Replicate injections** (n=9) for each re-equilibration experiment showed **excellent reproducibility** (%RSD typically <0.3).
- This implies that **re-equilibration times of 10-20 column volumes** should be **generally applicable** to obtain **reproducible chromatography**.
- It is **essential** to **accurately control gradient re-equilibration times** in HILIC mode. Care should also be taken to account for differing **system dwell volumes**.
- Important to carry out a **blank run before injecting** or **discard the first injection**.

		Retention time (%RSD*)		
	Re-equilibration	4-Hydroxybenzoic acid	2'-Deoxyguanosine	Salbutamol
		10 column volumes	10 column volumes	10 column volumes
HILIC-A	10 column volumes	1.178 (0.09)	4.134 (0.22)	4.547 (0.25)
	20 column volumes	1.300 (0.06)	4.574 (0.18)	4.429 (0.12)
	30 column volumes	1.330 (0.11)	4.683 (0.21)	4.446 (0.13)
	200 column volumes	1.343 (0.61)	4.788 (0.61)	4.514 (0.46)
HILIC-B	10 column volumes	1.351 (0.18)	5.559 (0.09)	4.244 (0.14)
	20 column volumes	1.398 (0.30)	5.712 (0.07)	4.441 (0.13)
	30 column volumes	1.409 (0.31)	5.810 (0.09)	4.508 (0.15)
	200 column volumes	2.222 (0.14)	6.003 (0.41)	4.639 (0.61)
HILIC-N	10 column volumes	1.367 (0.08)	5.636 (0.05)	4.349 (0.12)
	20 column volumes	1.488 (0.08)	5.857 (0.04)	4.666 (0.08)
	30 column volumes	1.602 (0.08)	6.028 (0.05)	4.767 (0.13)
	200 column volumes	1.816 (0.11)	5.765 (0.76)	4.917 (0.43)

## 9. Conclusions

- It is **essential** to **fully and consistently equilibrate HILIC columns** to obtain **reproducible chromatography**.
- With adequate equilibration, **ACE HILIC** columns **show similar reproducibility** to **reversed-phase columns**.
- Fresh columns** should be equilibrated for at least **60-80 column volumes** before injecting.
- Equilibration times for **used columns** are shorter, **20 column volumes** is usually sufficient.
- Under **certain conditions**, column **equilibration may take longer** than 20 column volumes.
  - High buffer concentrations
  - Low percentage of water in the mobile phase
- For **HILIC gradients**, varying the post-gradient **re-equilibration time** affects **analyte retention** and **potentially selectivity**. It is crucial to **accurately control** and **record re-equilibration times**.
- When **transferring HILIC gradients**, **instrument dwell times** may need compensating for to maintain consistent re-equilibration times.