

Low Level Determination of Catecholamines: Epinephrine, Norepinephrine and Dopamine from Plasma By UHPLC-MS/MS Using a Novel C18-PFP Column

Edward Faden¹, Alan P McKeown²


¹MACMOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317 USA

²Advanced Chromatography Technologies Ltd, 1 Berry Street, Aberdeen, Scotland, AB25 1HF UK


1. BACKGROUND

- Catecholamines are small polar molecules which can provide helpful information for monitoring various human disease states and endogenous processes. Commonly determined catecholamines within the clinical setting include dopamine, epinephrine and norepinephrine.
- Traditional analytical methods routinely include the use of electrochemical detection which can bring its own challenges. This work describes the extraction from plasma of 3 key catecholamines with subsequent quantification using UHPLC-MS/MS armed with the novel ACE Excel 1.7 μ m C18-PFP.
- This phase is designed to maximise selectivity and separation of closely related polar species by mechanisms that include hydrophobicity, and pi-pi interactions.

2. RESOLUTION, SELECTIVITY, EFFICIENCY & RETENTION



3. RATIONAL PHASE DESIGN TO MAXIMISE SELECTIVITY



Robust C18-Based Phases Designed With Ultra-Low Bleed For MS

4. PLASMA SAMPLE EXTRACTION & INSTRUMENT CONDITIONS

Format: EVOLUTE EXPRESS® WCX 10mg fixed well plate, part number 602-0010-PX01.

MS Conditions: AB SCIEX 5500 triple quadrupole MS. Curtain gas = 35; collision gas = 7; ionspray voltage = 5500; temperature = 700C; ion source gas 1 = 50; ion source gas 2 = 50; setting time = 40 ms. Positive ions acquired in MRM mode.

| Compound | MRM Transition | DP (V) | EP (V) | CE (V) | CE (V) |
|-------------------|----------------|--------|--------|--------|--------|
| Epinephrine | 166.1 > 107.1 | 148 | 8 | 24 | 16 |
| Norepinephrine | 152.1 > 107.1 | 25 | 2 | 22 | 25 |
| Dopamine | 154.1 > 91.1 | 50 | 9 | 29 | 13 |
| d6-Epinephrine | 172.1 > 112.1 | 148 | 8 | 24 | 16 |
| d6-Norepinephrine | 158.1 > 111.1 | 25 | 2 | 22 | 25 |
| d4-Dopamine | 158.1 > 95.1 | 50 | 9 | 29 | 13 |

Column: ACE Excel 1.7 μ m C18-PFP, 2.1 x 100 mm, part number EXL-1710-1002U


Mobile Phase: 0.25mM ammonium acetate / 0.1% HCOOH (aq)
Mobile Phase B: 0.25mM ammonium acetate / 0.1% HCOOH in MeOH

Flow Rate: 0.4 mL/min
Gradient: Start at 5 %B and hold until 1.2mins then step gradient to 95 %B and hold to 3.2mins. Switch back to 5 %B and hold for 4 mins equilibration before next injection.

Internal Standards: d6-Epinephrine, d6-Norepinephrine, d4-Dopamine purchased from LGC (Teddington, UK)
Sample: Mix plasma sample (250 μ L) containing internal standards with 250 μ L 0.05% w/w HCOOH

Injection Volume: 20 μ L
Temperature: 40C

5. RESULTS: CATECHOLAMINES IN PLASMA BY UHPLC-MS/MS



6. SUMMARY AND CONCLUSIONS

- A rapid UHPLC-MS method for low level determination of 3 key catecholamines using the novel ACE Excel 1.7 μ m C18-PFP has been established.
- An extraction protocol using ISOLUTE® PLD+ Protein and Phospholipid Removal plates was developed to enable satisfactory recoveries and low level detection of the analytes from the serum matrix.
- Recoveries were found to be >70% with RSDs <10% for both analytes when spiked into serum samples across a wide range of 2 – 100 ng/mL. The method was also linear for both from 1 – 100 ng/mL in PBS/BSA with r^2 >0.99.

ACE® is a registered trademark of Advanced Chromatography Technologies Limited. Excel™ is a trademark of Advanced Chromatography Technologies Limited. Advanced Chromatography Technologies Limited acknowledges the registered and unregistered trademarks of Biotage® and has no affiliation with the company.

All trademarks are recognised. Extraction and MS data reproduced with permission of Dr Lee Williams, Biotage® GB Ltd. See Biotage Application note AN844 for further details.