

Metabolomic Analysis of Extracted Jurkat T Cells by LC-HRMS


Application #AN3980

Conditions

Column: ACE Excel 2 C18-PFP
Dimensions: 100 x 2.1 mm
Part Number: EXL-1010-1002U
Mobile Phase: A: 0.1% formic acid in H₂O
B: 0.1% formic acid in MeCN


Time (mins)	%B
0	0
1	0
11	65
13	65
18	95
20	95

Flow Rate: 0.35 mL/min
Injection: 5 µL
Temperature: 35 °C
Detection: Thermo Scientific Q Exactive Orbitrap MS
Heated electrospray ionisation in positive mode
Spray Voltage: 3.3 kV
Capillary Temperature: 300 °C
Heater Temperature: 350 °C
Mass Scan Range: *m/z* 70-1000
Resolution: 70,000




TIC overlay for Jurkat T-lymphocyte cells rinsed with either 0.3% ammonium formate (darker line) or 0.3% ammonium acetate

1. Caffeine-d3 (IS)
2. Tryptophan-d3 (IS)



3. 1,3,7-Octanetriol



4. Pimelylcarnitine

Ulmer CZ, Yost RA, Chen J, Mathews CE, Garrett TJ. Liquid-Chromatography-Mass Spectrometry Metabolic and Lipidomic Sample Preparation Workflow for Suspension-Cultured Mammalian Cells using Jurkat T Lymphocyte Cells, *J. Proteomics Bioinform.* (2015), 8(6), 126-132

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.