


Application note #7940

LC-MS/MS separation of 10 nitrosamines and DMF in valsartan API on an Avantor[®] ACE[®] Excel C18 column













Method Details

CONDITIONS

Column: Avantor® ACE® Excel C18
 Particle Size: 2 µm
 Dimensions: 100 x 2.1 mm
 Mobile Phases: A: 0.1% formic acid in H₂O:MeOH (98:2, v/v)
 B: 0.1% formic acid in H₂O:MeOH (2:98, v/v)

Time (mins)	% B
0	0
1	0
2.5	50
5	86.5
5.5	86.5
5.6	0

Flow Rate: 0.5 ml/min
 Temperature: 40 °C
 Injection volume: 40 µl
 Detection: Sciex QTRAP® 6500+ LC-MS/MS system.
 Ionisation mode: APCI, positive mode; Source temperature: 300 °C; Curtain gas: 33 psig;
 Ion source gas: 30 psig; Needle current: 2 µA
 Sample: Valsartan (66.67 mg/ml) spiked with nitrosamines at 1 ng/ml. NDMA-d6, NMBA-d3, NDEA-d10 and NDBA-d18 were used as internal standards.



MRM TRANSITIONS

Analyte	MRM	Optimised MS Parameters		
		Declustering potential (V)	Collision energy (V)	Cell exit potential (V)
NDMA	+75.0 amu → +43.0 amu	11	19	10
	+75.0 amu → +58.0 amu	11	17	28
NDMA-d6	+81.2 amu → +46.0 amu	40	22	11
	+81.2 amu → +64.1 amu	40	17	12
DMF	+74.0 amu → +42.0 amu	1	46	12
	+74.0 amu → +30.0 amu	1	33	14
NMOR	+117.0 amu → +87.0 amu	11	15	4
	+117.0 amu → +44.9 amu	11	24	4
NMBA	+147.1 amu → +117.1 amu	11	11	12
	+147.1 amu → +87.1 amu	11	17	10
NMBA-d3	+150.1 amu → +120.2 amu	16	11	8
	+150.1 amu → +47.1 amu	21	17	8
NDEA	+103.1 amu → +75.1 amu	16	21	10
	+103.1 amu → +47.1 amu	16	23	22
NDEA-d10	+113.2 amu → +34.2 amu	21	33	6
	+113.2 amu → +49.1 amu	6	23	6
NPIP	+115.0 amu → +69.1 amu	1	21	8
	+115.0 amu → +41.0 amu	1	31	10
NEIPA	+117.1 amu → +75.1 amu	26	17	10
	+117.1 amu → +47.1 amu	21	23	10
NDIPA	+131.1 amu → +89.1 amu	76	15	10
	+131.1 amu → +47.1 amu	71	23	10
NMPA	+137.1 amu → +66.0 amu	21	23	8
	+137.1 amu → +107.1 amu	16	21	12
NDPA	+131.1 amu → +89.1 amu	16	17	10
	+131.1 amu → +43.1 amu	16	21	10
NDBA	+159.2 amu → +57.1 amu	46	17	10
	+159.2 amu → +103.2 amu	51	15	10
NDBA-d18	+177.3 amu → +66.2 amu	46	23	8
	+177.3 amu → +46.2 amu	41	37	22

ORDERING TABLE

Product	Details	Size	Part Number
Avantor® ACE® Excel 2 C18	HPLC Column	100 x 2.1 mm	EXL-101-1002U
Methanol	VWR HiPerSolv CHROMANORM® for LC-MS	2.5 L	83638.320
Water	VWR HiPerSolv CHROMANORM® for LC-MS	2.5 L	83645.320
Formic acid	VWR HiPerSolv CHROMANORM® for LC-MS	10 x 1 mL	85048.001

ORDERING TABLE (US)

Product	Details	Size	Catalog Number
Avantor® ACE® Excel 2 C18	HPLC Column	100 x 2.1 mm	76392-538
Methanol ≥99.9%	ULTRA for LC-MS, tested for UHPLC, J.T.Baker®	1 L	JT9863-2
Water	ULTRA for LC-MS, tested for UHPLC, J.T.Baker®	1 L	JT9823-1