

Outperforming Sub-2- μm Totally Porous Particles using Fused-Core Technology

Robert E. Moran, Stephanie A. Schuster, Barry E. Boyes, William L. Miles

Advanced Materials Technology Inc., Wilmington, DE 19810

Presented at EAS 2014


Abstract

Fused-core particles have shown distinct advantages over comparable totally porous particles (TPP) for separating molecules of all sizes and characteristics. The introduction of a new 2.0 μm particle size to the Fused-core product line continues to advance high speed, high efficiency separations into the future. These 2.0 μm Fused-core particles permit separation speeds and resolution improvements greater than sub-2- μm totally porous particles, but with about 20% lower column back pressure. Based on chromatographic theory, the trend of using smaller particle sizes has been accepted as a way to increase efficiency. This report describes the effect of these new particles on multiple factors of separation performance, including reduced plate height, column efficiency, back pressure, and resolution.

Engineered for high performance

Comparing van Deemter Plots (h)


Reduced Plate Height vs. Mobile Phase Velocity Plots
Columns: 50 x 2.1 mm, Instrument: Shimadzu Nexera, Solute: naphthalene
Mobile phase: Halo: 50/50 ACN/water, $k = 8.3$
1.6 μm SPP: 48.5/51.5 ACN/water, $k = 8.3$, 1.7 μm TPP: 47/53 ACN/water, $k = 8.2$
1.7 μm TPP: 48.5/51.5 ACN/water, $k = 8.3$, Temperature: 35 °C, Injection volume: 0.2 μL



Industry leading efficiency over broad range of flow rates


Ballistic Separation of Anticoagulants

Column: 2.1 x 30 mm HALO 2 C18, 2 μm , 90 Å
Mobile phase A: 30 mM formic acid
Mobile phase B: 50/50 ACN/MeOH
Gradient: Time 0-1.06 min 20-75%
Flow rate: 1.1 mL/min
Temperature: 45 °C
Detection: 254 nm
Injection: 0.2 μL
Max Pressure: 430 bar



Steroids Separation : HALO vs. TPP


Column: 2.1 x 50 mm TPP
Mobile Phase: 57% H₂O/43% MeOH
Flow rate: 0.3 mL/min
Temperature: 30 °C
Detection: 280 nm
Injection volume: 1.0 μL



Faster separation with high efficiency


Explosive Separation : HALO vs. TPP

Column: 2.1 x 100 mm
Mobile Phase: 72/28 H₂O/MeOH
Instrument: Shimadzu Nexera
Flow rate: 0.4 mL/min
Temperature: 42 °C
Detection: 254 nm
Injection volume: 1.0 μL



Narrower and taller peaks at 15% lower pressure


HALO® Fused-Core Particle Designs



*Available in 90 Å and 160 Å pore sizes
†160 Å available by June 2015

Pressure vs. Flow Rate


Mobile Phase/Column Pressure Plots
Columns: 50 x 2.1 mm, C18, Instrument: Shimadzu Nexera, Solute: naphthalene
Mobile phase: Halo: 50/50 ACN/water, $k = 8.3$, SPP: 1.6 μm , 48.5/51.5 ACN/water, $k = 8.3$
SPP: 1.7 μm , 47/53 ACN/water, $k = 8.2$, Temperature: 35 °C, Injection volume: 0.2 μL



Incredible high pressure stability


Halo 2 μm High Pressure Stability

Column: 3.0 x 150 mm Penta-HLIC
Instrument: Shimadzu Nexera
Injection Volume: 1.0 μL
Detection: 254 nm
Temperature: 25 °C
Flow Rate: 0.6 mL/min



Plates per Bar

Column: 2.1 x 50 mm C18, Instrument: Shimadzu Nexera, Solute: naphthalene
Mobile phase: 50/50 - 47/53 ACN/water, Flow rate: 0.5 mL/min, Temperature: 35 °C




Peptide Separation : Peak Capacity

Column: 2.1 x 150 mm
Instrument: Agilent 1200
Injection Volume: 2.5 μL
Detection: 215 nm
Temperature: 60 °C

Flow rate: 0.5 mL/min
Mobile Phase A: H₂O + 0.1% TFA
Mobile Phase B: 80/20 ACN/H₂O + 0.1% TFA
Gradient: 15-70% B in 90 min

HALO 2 Peptide ES-C18, 2 μm , 160 Å
Peak Capacity: 480
Max Pressure: 901 bar

HALO 2 Peptide ES-C18 available June 2015



Outstanding peak capacity


Peptide Peak Capacity : HALO vs. TPP

Column: 4.6 x 150 mm
Instrument: Agilent 1100
Injection Volume: 5 μL
Detection: 215 nm
Temperature: 40 °C

Flow rate: 4.0 mL/min
Mobile Phase A: H₂O/0.1% TFA
Mobile Phase B: ACN/0.1% TFA
Gradient: 5-60% B in 5 min

HALO 5 Peptide ES-C18
Pressure: 333 bar
Peak capacity: 94

TPP 5 μm C18
Pressure: 331 bar
Peak capacity: 37



Peak capacity 2.5 times TPP-column

Conclusions

- Fused-Core columns consistently provide sharper peaks, lower reduced plate heights, and lower operating pressures when compared to totally porous particles of equivalent particle size
- Sharper peaks result in increased peak capacity and sensitivity for all separations
- 2 μm Fused-Core columns deliver similar performance to sub-2- μm columns at lower operating pressures
- Multiple particle sizes to fill the various objectives required in HPLC/UHPLC separations
 - Small particles for high speed and high resolution R&D separations
 - Larger particles for routine QA and QC separations
- Various pore sizes to accommodate the broad range of analyte molecular weights
 - 90 Å - small molecules
 - 160 Å - peptides
 - 400 Å - proteins

Acknowledgments

Special thanks to Joseph DeStefano for expert knowledge and insight.

