




HALO® PAH Stability at 600 bar

278-P



PEAK IDENTITIES:

- | | |
|-----------|------------------------|
| 1. Uracil | 3. 1-Cl-4-Nitrobenzene |
| 2. Phenol | 4. Naphthalene |

TEST CONDITIONS:

Column: HALO 90 Å PAH, 2.7 μ m, 2.1 x 150 mm

Part Number: 92842-712

Mobile Phase A: Water

B: Acetonitrile

Isocratic: 50% B

Flow Rate: 0.6 mL/min

Back Pressure: 597 bar

Temperature: 30 °C

Detection: 254 nm, PDA

Injection Volume: 0.5 μ L

Sample Solvent: 60/40 ACN/ Water

Data Rate: 100 Hz

Response Time: 0.025 sec.

Flow Cell: 1 μ L

LC System: Shimadzu Nexera X2

Polycyclic Aromatic Hydrocarbons (PAHs) are a group of more than 100 chemicals released from the combustion of coal, oil, gasoline, tobacco, and wood. They can also be found in cooked food. PAHs are persistent chemicals and must be closely regulated for early detection/monitoring to minimize hazardous exposure in the environment and/or use of contaminated raw materials in different industries. The HALO® PAH column shows excellent stability at elevated back pressure making it an excellent choice for polycyclic aromatic hydrocarbon analysis.

