




A Reversed Phase Separation of Polar Metabolites Using the HALO® 1.5 mm ID AQ-C18 Column

324



TEST CONDITIONS:

Column: HALO 90 AQ-C18, 2.7 μ m 1.5 x 150 mm

Part Number: 9281X-722

Mobile Phase A: 8 mM ammonium formate, pH 4.0 (aq.), in 50:50 acetonitrile:water

Mobile Phase B: 8 mM ammonium formate, pH 4.0 (aq.), in 95:5 acetonitrile:water

Gradient: Time	%B
0.0	0
1.5	0
12.0	95
14.0	95

Flow Rate: 0.2 mL/min

Temperature: 35 °C

Injection Volume: 1 μ L

Sample Solvent: 98/2 5mM ammonium acetate/methanol

LC System: Shimadzu Nexera X2

MS CONDITIONS:

System: ThermoFisher Q Exactive HF Hybrid Orbitrap

Spray Voltage (kV): 3.5

Capillary Temperature: 350 °C

Sheath gas: 40

Aux gas: 15

RF lens: 40

Metabolites from a yeast extract were separated using a HALO® 1.5 mm ID 90 Å AQ-C18, 2.7 μ m column. The isomers leucine and isoleucine are baseline resolved by the use of an isocratic hold at 0% B enabled by the 100% aqueous compatibility of the HALO® AQ-C18 phase. By using a 1.5 mm ID column, 50% less solvent is used compared to running on a 2.1 mm ID column.

