# PRACTICAL APPLICATION OF FUSED-CORE® COLUMNS WITH CONVENTIONAL HPLC INSTRUMENTATION: Understanding the Importance of Extracolumn Volume

### **Objectives of the Study**

- <sup>,</sup> Determine the effects of extracolumn volume and dispersion on performance of HALO<sup>®</sup> Fused-Core<sup>®</sup> columns when using conventional instrumentation.
- Determine how much extracolumn volume and dispersion is allowed to achieve acceptable performance from various HALO column geometries
- Provide instrument configuration and parameter recommendations to chromatographers so that Fused-Core columns can deliver acceptable and optimum performance levels for both high throughput and high resolution separations

### HALO Fused-Core Particle



### HALO Fused-Core Columns Exhibit UHPLC-like **Performance at Conventional HPLC Pressure**



### **Extracolumn dispersion:** Experimental parameters



- signal filtering
- Impact of extra-column dispersion varies inversely with column volume, peak volume, and retention



## Thomas J. Waeghe, Robert T. Moody, Carl L. Zimmerman, MAC-MOD Analytical, Inc., Chadds Ford, PA

### Estimating maximum acceptable extracolumn volume, w



### Maximum Acceptable ECV

- Lower for higher efficiency columns (N) (smaller d\_ lower H/h)
- Lower for smaller volume columns  $(V_m)$
- Higher for longer retention (k)

Example: Calculating maximum ECV to achieve 90% of resolution (81% of N) when k = 3.

<u>4.6 x 250 mm, 5 µm column</u>  $N_{obs}/N_{theor} = 0.81, N_{theor} = 20,000; V_m 2.5 mL$  $W_{m} = 125 \,\mu L$ 

<u>3 x 50 mm, 2.7 µm HALO</u>  $_{\rm m} \sim 177 \,\mu L$ , h ~1.45, N<sub>theor</sub> = 12770  $W_{ec} = 12 \,\mu L$ 

Most conventional HPLCs have  $30 < w_{cc} < 100 \,\mu L$ 

### Affect of sources of extracolumn dispersion on isocratic and gradient separations

Extracolumn dispersion before the column has much more effect on isocratic separations than gradient separations

- Injection volume
- Sample solvent composition
- Flow path in autosampler and valve
- Tubing from the autosampler to mobile phase preheater (if present)
- Precolumn heat exchanger
- Tubing from precolumn heat exchanger to column
- Volume due to any in-line filters, unions, guard columns, etc.

### Extracolumn dispersion after the column effects both isocratic and gradient separations

- Tubing from the column to flow cell
- Flow cell volume and design
- Signal filtering due to detector response time

### **Experimental: Instrumentation and Configurations**

### Agilent 1100 Quaternary System

- Popular HPLC system in many laboratories
- Quaternary pump with low pressure mixing,
- Standard autosampler with variable volume injector
- Standard tubing ID (0.007") and lengths in sample flow path
- Variable wavelength UV-VIS detector Flow cells
  - -1 μL, 5 μL and 14 μL
- \* Response times
- -0.0625, 0.125, 0.25, 0.5, 1.0, 2.0, 4.0, 8.0 sec

### Evaluation of 3 Instrument Configurations with Different **ECVs**

- ECVs chosen for acceptable performance with respective column sizes
- Standard configuration
- Low ECV configuration
- Ultra-Low ECV configuration

### **Experimental: Column dimensions, test mixture** and criterion selected

### Columns

- 4.6 mm ID: 50, 100, and 150 mm
- 3.0 mm ID: 50 and 100 mm
- 2.1 mm ID: 50 mm Column geometries selected based on (1) column volume, (2) greater

### Test Mixture

- Probe analytes chosen from among mono-substituted and di-substituted benzenes
- uracil, benzyl alcohol, benzonitrile, nitrobenzene, anisole,
- 1-chloro-4-nitrobenzene, toluene
- k values from ~0.35 to ~3.5 vs. uracil void marker

### Criterion

- Aim for  $N_{obs} \ge 81\%$  of  $N_{theor}$  ( $R_{obs} \ge 90\%$  of maximum theoretical R)
- Reduce ECV and decrease response times to improve N for smaller columns

likely use by more chromatographers and (3) historical sales data

### HALO Fused-Core columns deliver significantly better performance when ECV is reduced







### Standard Configured Agilent 1100

- ECV ~36 µL
- Standard flow cell; 14 μL 0.5 sec. response time Standard length and ID tubing (0.007" ID x 750 mm total)

### HALO C18 4.6 x 100 mm 2.7 µm, 2 mL/min



### Ultra-Low ECV Configured Agilent 1100

- ECV ~10 µL
- Semi-micro flow cell; 5 μL 0.5 sec. response time
- Reduced length and ID tubing (0.005" ID x 460 mm total)

### Theoretical plates vs. ECV for conventional and HALO Fused-Core columns



### How to improve resolving power with high efficiency, smaller volume columns

### **Reduce Extracolumn Dispersion**

- Use smaller volume flow cell ( $\leq 5 \,\mu$ L)
- Use fastest practical data rate and shortest detector response time - Set detector response time to fastest setting that provides acceptable S/N. (< 0.2 seconds recommended for HALO columns. • Set data rate to collect at least 20 data points across the narrowest peak of interest. (> 5 Hz recommended for HALO columns.) • Reduce tubing ID and length between injector and flow cell

- Connection from column to flow cell is more important than from autosampler to column - use a smaller volume pre-column heat exchanger  $(1.6 \,\mu\text{L})$
- (if necessary)
- Use smallest practical injection volume (repeatability) • Choose sample solvent composition to be weaker than mobile phase,
- solubility permitting

### Configuring your instrument for ultrafast and high resolution HPLC

	Standard Configuration	Low Volume Configuration	Ultralow Volume Configuration
Total ECV	~36 µL	~17 µL	~12 µL
Flow Cell	Standard, 14 or 13 µL	Semi-micro, 5 µL	Micro, 1 or 2 µL
Tubing	0.007" ID, Standard Lengths	0.005" ID, Minimal lengths	0.005″ ID, Minimal lengths
Heat Exchanger	3 µL, Standard	3 μL, Standard	1.6 µL, Micro
Acceptable Fused- Core Columns	4.6 mm ID columns > 50 mm lengths	All 4.6 mm ID columns and 3.0 mm ID columns ≥75 mm lengths	All 4.6, 3.0, and 2.1 mm ID columns

### Ultrafast and high resolution separations with low **ECV configured HPLC system**



### **Ultra-fast gradient separation of 9 component** sample in ~1 minute with HALO Fused-Core column on a low ECV configured HPLC system



6. valerophenone, 7. hexanonphenone, 8. heptanophenone, 9. octanophenone

### Conclusions

Extracolumn dispersion in conventional HPLC systems can limit the benefits of columns that generate fast eluting, small volume peaks.

Achieving "UPLC-like" performance is possible with conventional HPLC equipment by reducing extracolumn band broadening and using an appropriate HALO Fused-Core column.