

Recombinant hGMCSF Purified from *Escherichia Coli*

ACE[®]
Ultra-inert
UHPLC & HPLC Columns

Application #AN3840


Conditions

Column: ACE 5 C18
Dimensions: 150 x 4.6 mm
Part Number: ACE-121-1546
Mobile Phase: A: 0.1% TFA in H₂O/MeCN (90:10 v/v)
B: 0.1% TFA in H₂O/MeCN (10:90 v/v)

Time (mins)	%B
0	10
20	65
23	100

Flow Rate: 1 mL/min
Temperature: 30 °C
Detection: UV, 215 nm


hGMCSF = human Granulocyte Macrophage Colony Stimulating Factor, a 127 amino acid residue cytokine with a molecular weight of 14,477 Da



Das KMP, Banerjee S, Shekhar N, Damodaran K, Nair R, Somani S, Raiker VP, Jain S, Padmanabhan S. Cloning, Soluble Expression and Purification of High Yield Recombinant hGMCSF in *Escherichia coli*. Int. J. Mol. Sci. 2011, 12, 2064-2076; doi:10.3390/ijms12032064

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit
www.ace-hplc.com or email: info@ace-hplc.com


ACE[®]
UHPLC & HPLC
Columns