


# Separation of Remdesivir and Impurity



# Method Details

## CONDITIONS

Column: ACE Excel 2 C18-PFP  
 Dimensions: 100 x 3.0 mm  
 Mobile Phases: A: 0.1% Formic acid in H<sub>2</sub>O  
 B: 0.1% Formic acid in MeCN

Gradient:

Time (mins)	%B
0	5
0.5	5
2.5	95
3.5	95
4	5
11	5

Flow Rate: 0.6 mL/min  
 Injection: 1 µL  
 Temperature: 40 °C  
 Detection: UV, 254 nm  
 Instrument: VWR ChromasterUltra Rs  
 System Dwell Volume: 284 µL

## ORDERING TABLE

Product	Details	Size	Part number
ACE Excel 2 C18-PFP	UHPLC column	100 x 3.0 mm	EXL-1010-1003U
Acetonitrile	VWR CHROMANORM® gradient grade	2.5 L	83639.320
Water	VWR CHROMANORM® gradient grade	2.5 L	83650.320
Formic acid	VWR NORMAPUR® analytical reagent	100 mL	20318.180
ChromasterUltra Rs UHPLC system	Pre-configured complete system with DAD		903-0400
UHPLC guard cartridge	For 3.0 mm i.d. C18-PFP columns	3 pk	ACE-1P10-GD2U
UHPLC guard holder	Required for use of UHPLC guard cartridges	1	H0011