

Selective Detection HPLC Assays via In-Column Derivatisation

Colin Pipe¹, Thirada Suktham², Arianne Soliven¹, Tony Edge¹, Matt James¹, Andrew Jones³, Andrew Shalliker^{2,4} and Geoff Faden⁵


¹Avantor, VWR International Ltd., 1-3 The Markham Centre, Station Road, Theale, Reading, Berkshire RG7 4PE, UK. ²Western Sydney University, South Parramatta, NSW 2151, AUSTRALIA. ³Chemika 119 Magowar Rd, Girraween, NSW 2145 AUSTRALIA. ⁴Chromaspeed Pty Ltd. Tonsley SA 5042 AUSTRALIA. ⁵MAC-MOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317 USA

Background

- Selective detection assays via HPLC post-column derivatization (PCD) have many applications: amino acid, antioxidant and phenolic analysis.
- Advantages over pre-column derivatization include:
 - Reduced sample manipulation.
 - Ability to work with less stable derivatization products.
- For amino acid analysis using OPA, PCD requires an in-line reactor which adds significant dead volume:
 - This significantly reduces the chromatographic performance.
- This poster investigates an alternative approach, In-column derivatisation (ICD) to overcome this for amino acid analyses.
 - In ICD, reagent/eluent mixing is initiated at the outlet frit, providing enhanced mixing.
 - Use of reactor eliminated, improving chromatographic performance.
- Application of ICD is demonstrated for authentication of coffee via the antioxidant fingerprint, which is complex and difficult to imitate.


ICD Approach

- PCD Setup:**
 - More complex
 - Reactor, tee and tubing introduce excessive dead column
 - Degrades peak efficiency
- ICD Setup:**
 - Simplified setup
 - Reaction occurs at column frit
 - Dead volume significantly reduced







ICD Amino Acid Assays with O-phthalaldehyde (OPA)

- For the OPA reagent, a small amount of post column dead volume was found to be beneficial for ICD:



ICD (+80 µL) provided significant improvements vs. PCD (500 µL):

	Glycine	Aspartic acid	Histidine	Lysine
Plates	4620	1731	4534	1719
Peak Width (10%) (min)	0.16	0.24	0.17	0.25
Area (mAU)	3986	3886	1011	594
Height (mAU)	795	465	194	69



References


- Using HPLC with In-Column Derivatization to Authenticate Coffee Samples. *Molecules* 2023 **28** 1651.
- Determination of antioxidants by a novel on-line HPLC-cupric reducing antioxidant capacity (CUPRAC) assay with post-column detection. *Anal. Chim. Acta* 2010 **674** 79-88.

In collaboration and with permission on behalf of Dr. Andrew Jones, Chemika.


Authentication via ICD: Antioxidant Profiling¹

An ICD selective detection assay for antioxidants using a CUPRAC reagent² was used to authenticate coffee:

- Antioxidant profiles are complex and unique for different coffees.
- 15 antioxidant indicators varied in detection response.



- Authentication workflow:



- 5 samples matched to a library of 32 coffees.
- Compilation of a larger library for improved classification.

Highest match	Coffee				
	Unknown 1	Unknown 2	Unknown 3	Unknown 4	Unknown 5
	Starbucks Columbia	Guatemala	Vittoria Espresso	Profondo	Or Absolu
Slope	0.9765	0.9999	0.9782	0.9903	0.9711
R ²	0.9748	0.9887	0.9956	0.9914	0.9879
P.F.	14.3692	1.2538	1.5741	1.8529	4.0843
2 nd Highest Match	No Match	Inca Peru	Arabica Catuai	Long Black	No Match
Slope		1.0043	1.0022	1.0007	
R ²		0.9958	0.9831	0.9899	
P.F.		1.1441	1.5687	1.8319	

Conclusions

- ICD provides a simple solution to complex/challenging separations requiring PCD.
- Since the application of ICD requires minimal addition of post-column dead volume, post-column dispersion is reduced, conserving the columns theoretical separation efficiency (55 to 167% efficiency improvement and up to 1/3 reduction in peak width compared to PCD).
- Due to the high efficiency, complex samples can be chemically fingerprinted using compounds that are only visible through chemical reactions.
- As the chromatographic efficiency is higher than in PCD, sensitivity is often improved.
- High efficiency also means that shorter columns can be employed, and this increases the analytical throughput.