


1. Introduction

- Cyanobacteria / blue-green algae are common. Microcystins are cyclic heptapeptides produced by cyanobacteria and associated with human hepatotoxicity.
- Under nutrient-rich conditions, cyanobacteria can rapidly accumulate and form blooms.
- Due to their hepatotoxic nature, the WHO has established a provisional limit of 1mg/L in water supplies for total Microcystin-LR as a marker for cyanobacteria toxin levels.
- This poster describes a method for low ppb quantification of Microcystin-LR, Microcystin-RR and Microcystin-YR by UHPLC-MS/MS using an **ACE Excel 2mm C18** UHPLC column.

2. Algae Blooms and Microcystins



Variants	R	L
Microcystin-LR	Leucine	Arginine
Microcystin-RR	Arginine	Arginine
Microcystin-YR	Tyrosine	Arginine

Over 80 known microcystin variants have been reported

Int. J. Environ. Res. Public Health 2012, 9

3. UHPLC and MS Conditions

Bruker Advance™ UHPLC

Column: **ACE Excel 2mm C18**, 100 x 2.1 mm
 Part Number: EXL-101-1002U

Injection volume: 50 µL
 Flow rate: 0.4 mL/min
 Column temperature: 40 °C
 Mobile phase A: 0.1% formic acid (aq)
 Mobile phase B: MeCN

Gradient conditions:

Time (mins)	%B
0	30
1.0	30
7.0	95
7.1	30
10.0	30


Bruker EVOQ™

VIP Heated-ESI Temp: 350 °C
 Heated gas: 50 units
 Nebulizer gas: 50 units
 Cone gas temp: 200 °C
 Cone gas: 15 units
 Spray voltage: 4500 V (Positive)
 Active exhaust: On
 Collision gas: Argon 1.5 mTorr

Optimized MRM transitions:
 MC RR (MW: 1038) transition: 520 → 135 (CE at 24V)
 MC LR (MW: 995) transition: 498 → 135 (CE at 11V)
 MC YR (MW: 1045) transition: 523 → 135 (CE at 9V)




4. MS Scans and Sensitivity Investigations at 100 ppb




m/z 520.0 [M+2H]²⁺ m/z 1,045.5[M+H]⁺
 m/z 523 [M+2H]²⁺ m/z 523 [M+2H]²⁺

m/z 996 [M+H]⁺
 m/z 498 [M+2H]²⁺


Microcystin-YR:
 x 194 signal sensitivity increase with
 m/z 523 [M+2H]²⁺



Microcystin-LR:
 x 128 signal sensitivity increase with
 m/z 498 [M+2H]²⁺




5. Reproducibility and Calibration Data




Peak area reproducibility at 0.05 ppb (n=7)

	MC RR	MC YR	MC LR
1	1246	300	478
2	1239	364	584
3	1396	340	620
4	1157	345	585
5	1166	339	527
6	1189	318	585
7	1147	384	487
mean	1220.00	341.43	552.29
std dev.	86.78	27.71	55.00
(%RSD)	7.11	8.12	9.96



1. RR $r^2 = 0.9975$
 2. YR $r^2 = 0.9993$
 3. LR $r^2 = 0.9992$



5. Conclusions

- A new method was established using an **ACE Excel 2µm C18** UHPLC column and Bruker Advance / EVOQ UHPLC-MS/MS instrumentation.
- The method was found to be reproducible and sensitive at very low levels (0.05 ppb) for microcystins in water.
- Rapid separation and elution of Microcystin-LR, Microcystin-RR and Microcystin-RR was achieved with an excellent linear range of 0.05 ppb – 50 ppb.
- MS optimization work indicated significantly enhanced sensitivity was possible using doubly charged ions.