

Separation and Low-Level Determination of Catecholamines and Metanephrines from Urine by UHPLC-MS/MS using a Novel C18-Based Stationary Phase

Alan P McKeown

Advanced Chromatography Technologies Ltd, 1 Berry Street, Aberdeen, Scotland, AB25 1HF UK

1. Introduction

- Determination of urinary catecholamines and their O-methylated metabolites (metanephrines) as biomarkers for various disorders and disease states is important in the clinical setting.
- The polar nature of these compounds makes them challenging to retain by reversed phase chromatography.
- In this work, analysis was successfully achieved using a novel stationary phase: **ACE Excel 2 C18-PFP**.
- The **ACE C18-PFP** has been designed to maximise selectivity and separation with multiple interaction modes including hydrophobic and π - π mechanisms.
- The **ACE C18-PFP** can provide enhanced retention for electron-rich aromatic analyte moieties, such as those found in catecholamines and metanephrines.
- This poster demonstrates a UHPLC-MS/MS method for the extraction, separation and quantification of both compound classes from urine.

2. UHPLC-MS/MS Conditions

Column: ACE Excel 2 C18-PFP, 100 x 2.1 mm (EXL-1010-1002U)
Instrument: Shimadzu Nexera UHPLC with AB Sciex 5500 MS
Mobile phase: A: 0.25 mM ammonium formate + formic acid
 B: 0.25 mM ammonium formate + formic acid in Methanol
Flow rate: 0.4 mL/min
Temperature: 40 °C
Injection volume: 7.5 μ L
Gradient:







t	%B
0	5
1.3	5
1.4	95
4.4	95

MS parameters (positive mode):
Curtain gas: 35 psi
Ion spray voltage: 5500 V
Temperature: 700 °C
lon source gas 1: 50 psi
lon source gas 2: 50 psi

Sample: Analyte standards and deuterated internal standards spiked into urine at 2-20 ng/mL, (3-methoxytyramine 50 ng/mL).


Optimised SPE protocol:

Step	Volume	Standard SPE	Load-Wash-Elute SPE
Condition	500 μ L	MeOH	-
Equilibration	500 μ L	10 mM ammonium acetate	-
Sample Load	150 μ L	Urine:250 mM ammonium acetate (75 μ L:150 μ L)	-
Wash 1	500 μ L	10 mM ammonium acetate	-
Wash 2	500 μ L	IPA	-
Elution	125 μ L	IPA:H ₂ O (85:15 v/v) + 0.1% formic acid	-

Catecholamines	Metanephrines
 Epinephrine LogP: -0.43 LogD _{pH3} : -2.95	 Metanephrine LogP: 0.00 LogD _{pH3} : -2.80
 Dopamine LogP: 0.03 LogD _{pH3} : -2.25	 3-Methoxytyramine LogP: 0.53 LogD _{pH3} : -2.11
 Norepinephrine LogP: -0.68 LogD _{pH3} : -3.17	 Normetanephrine LogP: -0.39 LogD _{pH3} : -3.03


LogP and LogD data obtained from www.chemicalize.com

3. Rational Phase Design to Maximise Selectivity



4. UHPLC-MS/MS Results

- The separation was developed to maximise retention of norepinephrine and epinephrine. This compromised the peak shape of 3-methoxytyramine.
- Direct injection of SPE eluate – no evaporation step.
- Analysis of all six analytes was achieved in <3.5 minutes.




MRM Parameters:

Analyte	Transition
Norepinephrine	152.1 > 107.1
D ₆ -Norepinephrine	158.1 > 111.1
Epinephrine	166.1 > 107.1
D ₆ -Epinephrine	172.1 > 112.1
Normetanephrine	166.1 > 134.0
D ₃ -Normetanephrine	169.1 > 137.0
Dopamine	154.1 > 91.1
D ₄ -Dopamine	158.1 > 95.1
Metanephrine	180.1 > 148.0
D ₃ -Metanephrine	183.1 > 151.0
3-Methoxytyramine	151.2 > 90.9

5. Recovery and Linearity

- Excellent linearity determined across varying clinical ranges:
 - 0.1 to 25 ng/mL for epinephrine
 - 0.5 to 125 ng/mL for metanephrine and normetanephrine
 - 1 to 250 ng/mL for norepinephrine and 3-methoxytyramine
 - 2.5 to 625 ng/mL for dopamine
- Full SPE procedure compared to Load-Wash-Elute approach
 - Similar recovery and reproducibility for both approaches



6. Conclusions

- Retention and separation of analytes (full separation of isobaric species) was achieved by UHPLC using the novel **ACE Excel 2 C18-PFP**.
- The use of novel LC stationary phases can help optimise selectivity and retention.
- Enhanced aromatic π - π interactions between the electron-deficient PFP ring and electron-rich analyte rings, means the **ACE Excel C18-PFP** phase overcomes the challenge of poor retention in reversed phase due to catecholamine and metanephrine polarity.
- Chromatographic separation achieved in <3.5 minutes.
- Simple optimised SPE protocol with direct injection of SPE eluate.
- Excellent linearity and recoveries demonstrated across varied clinical ranges.