

Separation and Low Level Determination of Thyroid Hormones From Human Serum by UHPLC-MS/MS Using a Novel C18-Based Stationary Phase


Alan P McKeown¹ Edward Faden² Geoffrey Faden²

¹Advanced Chromatography Technologies Ltd, 1 Berry Street, Aberdeen, Scotland, AB25 1HF UK, ²MACMOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317 USA

1. BACKGROUND

- Determination of **thyroid hormones** helps the clinical understanding of a variety of **conditions and disease** states that include **metabolic disorders** and **depression**.
- Analytes of interest include **thyroxine (T4)** which is the dominant **prohormone** in humans. The bioactive **tri-iodothyronine (T3)** results from **deiodination of T4** whilst the **minor isobaric** metabolite **reversed tri-iodothyronine (rT3)** may also be present.
- A **semi-automated solid phase** extraction workflow for **serum** samples was devised. The **novel ACE Excel 2µm C18-AR** column provided separation with quantification by UHPLC-MS/MS.

2. RATIONAL PHASE DESIGN TO MAXIMISE SELECTIVITY



3. SERUM SAMPLE EXTRACTION & INSTRUMENT CONDITIONS

Format: EVOLUTE® EXPRESS AX 30 mg 96 Well Plate, part number 603-0030-PX01.

MS Conditions: XEVO TQS triple quadrupole (Waters Inc., USA). Desolvation temperature = 500C, Ion source temperature = 150C, positive ions acquired in MRM mode.



Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
T3	651.8 > 605.8	50	22
	(651.8 > 507.8)	(50)	(22)
rT3	651.8 > 478.9	(50)	(35)
	(651.8 > 605.8)	(50)	(22)
T3/T3-d6 ISTD	651.8 > 507.8	(50)	(22)
	(651.8 > 478.9)	(50)	(35)
T4	657.8 > 611.8	50	22
	777.7 > 731.7	50	25
	(777.7 > 351.0)	(50)	(45)
	(777.7 > 633.8)	(50)	(23)

Column: ACE Excel 2µm C18-AR, 2.1 x 100 mm
Mobile Phase A: 2mM ammonium acetate / 0.1% HCOOH (aq)
Mobile Phase B: 2mM ammonium acetate / 0.1% HCOOH in MeOH
Flow Rate: 0.4 mL/min
Gradient:



Time	% A	% B
0	40	60
3	23	77
3.1	40	60

Injection Volume: 10 µL
Temperature: 40C

4. THYROXINE & METABOLITES: LOW LEVEL LC-MS/MS ANALYSIS



5. EXAMPLE LINEARITY AND RECOVERY DATA



6. SUMMARY AND CONCLUSIONS

- Quantification** of thyroid analytes is helpful for **clinical diagnostics** and achievable using **UHPLC-MS/MS**.
- An extraction protocol using **EVOLUTE® EXPRESS AX** and a separation method using the novel **ACE® Excel C18-AR** column were developed to enable **low level detection** of the thyroid hormones.
- Recoveries** were excellent and the method was found to be **linear** for each thyroid analyte across a **wide concentration range** of **25-1000pg/mL**.

ACE® is a registered trademark of Advanced Chromatography Technologies Limited. Excel™ is a trademark of Advanced Chromatography Technologies Limited. Advanced Chromatography Technologies Limited acknowledges the registered and unregistered trademarks of Biotage® and has no affiliation with the company.

All trademarks are recognised. Extraction, UHPLC and MS data reproduced with kind permission of Dr Lee Williams, Biotage®, Cardiff, UK.