

Separation of Solvents using the Avantor® Hichrom HI-VOC 1 Phase




Figure 1: Separation of solvents using the Avantor® Hichrom HI-VOC 1 phase.

Method Details

CONDITIONS

Oven Program:	40 °C (12 min), 12 °C/min, 260 °C
Carrier Gas:	Helium, 0.7 mL/min
Injector:	Split, 250 °C, 0.5 µL, 150 mL/min Split flow
Detector:	FID, 290 °C

The HI-VOC 1 is a unique application specific column for the analysis of volatile organic compounds (VOCs), solvents and purgeable compounds.

ORDERING TABLE

Product	Details	Dimensions	Part Number
Avantor® Hichrom HI-VOC 1	GC Column	0.25 mm, 1.50 µm, 30 m	HI64-25-150-30