

# **Superficially Porous Particles for Peptide and Protein Analysis**

Barry E. Boyes, Joseph J. Kirkland, Stephanie A. Schuster, Brian M. Wagner, and  
Joseph J. DeStefano

Advanced Materials Technology, Inc., 3521 Silverside Rd., Wilmington, DE 19810

# Objective

The objective of this work is to demonstrate fast, efficient separations of proteins using superficially porous particles that have been optimized for such applications.

## Abstract

The original 2.7  $\mu\text{m}$  superficially porous particles (SPP) introduced in 2006 were created with an average pore size of 90  $\text{\AA}$ , which was suitable for small molecule analytical separations. This SPP particle technology now has been expanded to include wider pore sizes and larger particle sizes that are specifically designed for larger biomolecules. Novel particle designs with specially selected bonded phases for peptide and protein separations are described. This presentation includes fast separations of peptides and intact protein mixtures, as well as examples of very high resolution separations of larger proteins and associated variants and contaminants. Columns with specially engineered bonded phases for these particles demonstrate high temperature stability, which is ideally suited for the conditions that are often used for analytical and small scale preparative biomolecular separations. Protein recovery and sample loading investigations are included. The optimized shell thickness of the new 400  $\text{\AA}$  SPP particles represents a compromise between a short diffusion path versus adequate retention and mass load tolerance. Examples of high molecular weight protein separations highlight the advantages of using columns of superficially porous particles with wider pores. Some comparisons with conventional totally porous particles are also shown.


# Sample Loading Study for HALO Protein C4

Columns: 4.6 x 100 mm; Gradient: 39-49% B in 10 min.

Mobile phase: A – 0.1% trifluoroacetic acid



B – acetonitrile with 0.1% trifluoroacetic acid; Temperature: 60 °C

Flow rate: 0.5 mL/min.; Injection volume: 5  $\mu$ L




The particle with 0.20  $\mu$ m shell thickness offers a compromise between sample loadability and retention, as well as an optimized diffusion path for large MW biomolecules.

# HALO<sup>®</sup> Protein Particles




**SEM image of a HALO Protein particle**



**FIB (Focused Ion Beam) sliced particle**

# Column Stability Study

Column: 2.1 x 100 mm; Mobile phase gradient: 25-40% acetonitrile/0.1% aqueous trifluoroacetic acid in 10 min; Temperature: 90 °C; Flow rate: 0.5 mL/min; Detection: 215 nm




- The HALO Protein C4 bonded phase is stable up to 90 °C, showing very little loss of retention.

# Protein Separations: Effect of Temperature

Column: 2.1 x 100 mm HALO Protein C4  
Instrument: Agilent 1200 SL  
Injection Volume: 2  $\mu$ L  
Detection: 215 nm  
Temperature: as indicated


Mobile Phase A: water/0.1% TFA  
Mobile Phase B: acetonitrile/0.1% TFA  
Gradient: 28-58% B in 10 min.  
Flow rate: 0.45 mL/min

Peak Identities (in order):  
1. Lysozyme 14.3 kDa  
2. BSA 66.4 kDa  
3.  $\alpha$ -Chymotrypsinogen A 25.0 kDa  
4. Enolase 46.7 kDa  
5. Ovalbumin 44.0 kDa



- Protein peak shape and recovery improve with increased temperature of analysis.

# Protein Separations: Fused-Core compared to Totally Porous



- Separation is 3 times faster at the same back pressure on the HALO Protein column compared to the same sample run on a sub-2-μm totally porous particle column

# Protein Recovery Studies

Protein	Recovery
Cytochrome c	100 (5.8 SD)
Catalase	92 (18 SD)

- Proteins were fraction collected from a 4.6 x 100 mm HALO Protein C4 column run at 60 °C under gradient conditions with water/ACN/0.1% TFA mobile phase. Blanks were obtained by replacing the column with a union
- Lyophilized proteins were reconstituted using 3 M Urea/1% Triton X-100/0.25% acetic acid
- Protein recoveries were measured using QuantiPro™ BCA Assay Kit for 0.5-30 µg/mL protein (Sigma-Aldrich, St. Louis, MO)
- Samples were incubated at 37 °C for 100 min.
- Each sample was run in duplicate
- Absorbance values were measured at 562 nm
- HALO Protein C4 shows good recovery of proteins




# Protein Separations: Effect of Pore Size

Columns: 2.1 x 100 mm  
Instrument: Shimadzu Nexera  
Injection Volume: 1  $\mu$ L  
Detection: 280 nm  
Temperature: 60  $^{\circ}$ C

Mobile Phase A: water/0.1% TFA  
Mobile Phase B: 80/20 ACN/water/0.1% TFA  
Gradient: 40-47% ACN in 10 min.  
Flow rate: 0.3 mL/min

Peak Identities (in order):

- |                                 |                           |
|---------------------------------|---------------------------|
| 1. Catalase                     | 250 kDa [~60 kDa subunit] |
| 2. $\alpha$ -Chymotrypsinogen A | 25.0 kDa                  |
| 3. $\beta$ -Galactosidase       | 465 kDa [116 kDa subunit] |
| 4. $\beta$ -Amylase             | 200 kDa [~50 kDa subunit] |




- Peak widths in minutes provided above each peak.
- The 400 Å pores of HALO Protein enable sharp peaks for high MW biomolecules.

# Tryptic Digest using HALO-5 Peptide ES-C18


Column: 4.6 x 100 mm, HALO-5 Peptide ES-C18  
Instrument: Agilent 1100  
Injection Volume: 10  $\mu$ L  
Detection: 215 nm  
Temperature: 45  $^{\circ}$ C  
Pressure: 54 bar initial

Mobile Phase A: water/0.1% TFA  
Mobile Phase B: ACN/0.1% TFA  
Gradient: 5-40% B in 45 min.  
Flow rate: 1.0 mL/min  
Sample: Apomyoglobin Tryptic Digest [2 mg/mL]



The extremely low back pressure of the HALO-5 Peptide ES-C18 column enables fast, efficient proteomic separations with a low potential for plugging.


# Separation of Reduced IgG1 using TFA Mobile Phase



Column: 2.1 x 100 mm HALO Protein C4  
 Instrument: Shimadzu Nexera  
 Injection Volume: 1  $\mu$ L  
 Detection: 280 nm  
 Temperature: 80  $^{\circ}$ C

Mobile Phase A: water/0.1% TFA  
 Mobile Phase B: 80/20 ACN/water/0.1% TFA  
 Gradient: 33-40% B in 10 min.  
 Flow rate: 0.25 mL/min



Sample: 0.5 mg/mL IgG1 treated with 100 mM DTT  
 in 8 M Guanidine HCl at 50  $^{\circ}$ C for 35 min.  
 Sample reduced only.



Antibody structure: Afaneh, C, Aull, MJ, Kapur, S, *Modern Immunosuppression Regimens in Kidney Transplantation*, 2012.

# High Resolution Analysis of mAb IgG1 Light and Heavy Chains with LC/MS

2.1 mm ID x 100 mm HALO Protein C4; 0.4 mL/min.; A: 0.5 % formic acid with 20 mM Ammonium Formate  
B: 45% AcN/45% IPA/ 0.5 % formic acid with 20 mM Ammonium Formate; Gradient: 29-32% B in 20 min.; 80 °C  
Detection: 280 nm Abs; Shimadzu LCMS-2020, ESI +4.5 kV, 2 pps, 500-2000 m/z



ID	Mass (Da)
LC1	23,204
LC2	23,192
LC3	23,203
HC1	50,539
HC2	50,424
HC3	50,668
HC4	50,680
HC5	28,862

Masses  
deconvoluted  
using MagTran

# Conclusions

- Fused-core particles with 400 Å pores are effective for efficiently separating proteins without restricted diffusion
- Protein separations can be run approximately 3 times faster on columns of Fused-core particles compared to columns of sub-2-µm particles at the same back pressure
- Fused-core particles have performance advantages over totally porous particles for separating peptides and proteins
- Columns of 400 Å particles are both efficient and stable up to 90 °C
- With the low back pressure afforded by 5-µm 160 Å Fused-core particles, columns of these particles are less prone to overpressurizing due to plugging and longer columns can be run for high resolution separations of proteomic samples
- With the correct choice of mobile phase, high resolution LC-MS data can be obtained for mAb separations using 400 Å Fused-core particles

# Acknowledgment

Special thanks to Robert Moran for assistance with chromatographic measurements.