

# Superficially Porous Silica Particle Technology Developments for Pharma and Biopharma Applications

Barry Boyes<sup>1,2</sup>, Joseph DeStefano<sup>1</sup>, Tim Langlois<sup>1</sup>, and  
Joseph J. Kirkland<sup>1</sup>

<sup>1</sup>Advanced Materials Technology, Inc.  
Wilmington, DE USA

<sup>2</sup>Complex Carbohydrate Research Center  
University of Georgia, Athens, GA USA  
bboyes@advanced-materials-tech.com



# Agenda

- The Brief Logic Fused-Core (SPP) Particles
- Recent Developments for Small Molecules
  - Particle Size Developments: Smaller and Larger
  - Expansion of Bonded Phase Choices
- Recent Developments for Larger Molecules
  - 160 Å pore RP for peptides
  - 400 Å pore RP for proteins
  - Penta-HILIC for glycans and glycopeptides

# Faster HPLC Separations

- Smaller Particle Packed Beds
  - Totally Porous (including flow through)
  - Not Porous (Pellicular)
  - Partly Porous (Superficially Porous)
- Monolithic Materials
- Open Tubular Columns (channels)

# Superficially Porous (Fused-Core<sup>®</sup>) Particles



*Why  $2.7\ \mu\text{m}$ ?*


$$P \cong 2500L\eta F / (d_p^2 d_c^2)$$

$$d_p \sim (2500L\eta F / P d_c^2)^{1/2}$$

- Not as necessary to filter samples and mobile phase since frits are not as small as needed for sub- $2\text{-}\mu\text{m}$ : practical compromise.
- With  $0.5\ \mu\text{m}$  shell high resolution is maintained at high flow rates (flat C-term in van Deemter plot).

# Separation of Explosives (c. 2007 EAS)

Column: 4.6 x 50 mm Halo C18; Mobile phase: 27% Methanol/73% Water; Temperature: 40°C;  
Flowrate: 3.3 ml/min; Detector: UV @ 254 nm; Pressure: 343 bar; Agilent 1100



# Recent Developments for Small Molecules

- Particle Size Developments: Smaller and Larger

2.7  $\mu\text{m}$  Halo: high efficiency, preferred analytical choice

5  $\mu\text{m}$  Halo: lower back pressures,  $\sim 3 \mu\text{m}$  porous particle efficiency

2  $\mu\text{m}$  Halo: higher efficiency (higher pressures)

- Expansion of Bonded Phase Options

C18, C8, Phenyl-hexyl, RP-Amide, Pentafluoropropyl (PFP), ES-CN

Penta-HILIC

# 5 $\mu\text{m}$ Halo Particle Utility

## Effect of Particle Size on Reduced Plate Height


Columns: 4.6 x 150 mm; Temperature: 30 °C

Mobile phase: 50% acetonitrile/50% water

Solute: 1-Cl-4-nitrobenzene; Injection: 1  $\mu\text{L}$

Instruments: <400 bar, Agilent 1100; >600 bar, Agilent 1200

DeStefano, et al., (2012) J. Chromatogr. 1258, 76-83




h values lower for 5  $\mu\text{m}$  HALO (h = 1.2) - more homogeneously packed bed structure

# 5 μm Halo Particle Utility

## Effect of Flow Rate on Column Pressure

Columns: 4.6 x 150 mm; Temperature: 30 °C

Mobile phase: 50% acetonitrile/50% water




Pressure increases with  $\left(\frac{1}{d_p}\right)^2$ ; small particles require high pressure



# 5 $\mu\text{m}$ Halo Particle Utility

## Plates/Pressure for Various Particle Sizes



The 5  $\mu\text{m}$  HALO fused-core particle has more than double the number of plates/pressure of the 5  $\mu\text{m}$  totally porous particles and four times the number of plates/pressure of the 3  $\mu\text{m}$  totally porous particles. Data on 4.6 x 150 mm columns at the plate height minimum, except for 1.8  $\mu\text{m}$  particle (estimated).

# 5 $\mu\text{m}$ Halo Particle Utility

## 5 $\mu\text{m}$ HALO Fused-core vs. 5 $\mu\text{m}$ Totally Porous: NSAIDs


Columns: 4.6 x 150 mm Instrument: Shimadzu Prominence UFLC XR

Flow rate: 2.0 mL/min, Injection Volume: 2  $\mu\text{L}$ ,

Detection: 254 nm; Temperature = 35  $^{\circ}\text{C}$

Mobile Phase: A: 20 mM pH 2.5 Potassium Phosphate


B: 50/50 ACN/MeOH; A:B = 48% A:52% B



# Separation of TCAs on PFP Phases

Column: 4.6 x 100 mm; Temperature: 35°C; Flow: 0.75 ml/min; Detector: UV @ 254 nm; 1 µL inj.

Mobile Phase: A- 25 mM Ammonium Acetate, pH 5.5, B- AcN; 36-65% B in 7 min.




# 2 $\mu\text{m}$ Halo Particle Utility


Column: 4.6 x 50 mm Halo C18; Mobile phase: 50%AcN/50% Water; Temperature: 25°C; Detector: UV @ 260 nm; Nexera, PDA, Measurements for 1-Cl-4-nitrobenzene

1.4  $\mu\text{m}$  Core, 0.33  $\mu\text{m}$  Radial Shell, 90 Å pore size

### Column Efficiency Comparisons Varying Particles



### Column Reduced Efficiencies Fitted to Knox




# Halo Penta-HILIC

## Opportunities

- HILIC may be an attractive alternative for difficult polar analyte separations
- HILIC separation speed and robustness may be improved
- New bonded phases may broaden applications served by HILIC

## Development and Characterization of Fused-Core HILIC BP


- HILIC retention is complex
- Highly hydroxylic bonded phase is appealing
- Polymeric bonded phases are likely to limit mass transfer kinetics – prefer use of a monolayer silane



# Effect of Linear Velocity on Penta-HILIC Column Efficiency



4.6 mm ID x 50 mm; 90% AcN/10 mM NH<sub>4</sub>Form 3.0, 25 °C; 1 μL, 50 ng Adenosine

Data fitted to Knox Equation




# Effect of Buffer Concentration on HILIC Separations

90% Acetonitrile/ $\text{NH}_4\text{Form}$ , pH 3.0: 2.1 mm ID x 100 mm, 25 °C, 0.5 mL/min



# High Speed HILIC Separation of Catecholamines and Amino Acids

4.6 mm ID x 50 mm; 2 mL/min., 85% AcN/10 mM NH<sub>4</sub>Form 3.0, 25 °C; 3 μL inj



**Note: at 85% AcN all  $k'$  lower on Silica than on Penta-HILIC.**






# Comparative Retention of HILIC Columns

2.1 mm ID x 150 mm


Test: 90% AcN/10 mM NH<sub>4</sub>Form pH 3.0, 0.5 mL/min, 23 °C

Toluene/Amitriptyline/Salycic Acid/Thiamine




# Halo Penta-HILIC Separation of Nucleosides and Nucleobases

## LC-MS Compatible Gradient with Ammonium Formate pH=6



Compounds	
1	Thymine
2	Uracil
3	Thymidine
4	2'-Deoxyadenosine
5	Adenine
6	Uridine
7	Adenosine
8	Hypoxanthine
9	Xanthine
10	Cytosine
11	2'-Deoxycytidine
12	Guanine
13	2'-Deoxyguanosine
14	Cytidine
15	Guanosine




4.6 mm ID x 100 mm Penta-HILIC

# Halo Penta-HILIC Separation of Nucleosides and Nucleobases


## HILIC Separation with Coupled UV and MS Detection of Nucleosides and Nucleobases

Shimadzu Nexera LC: 2.1 mm ID x 100 mm Halo Penta HILIC column; 35°C; Abs. 260 nm  
Gradient: 6% B for 3.5 min, 65-100% B 3.5 to 9.5 min., hold 100% B for 0.5 min;  
A – 10mM Ammonium Formate in 93% AcN/7% H<sub>2</sub>O, pH 6.0; B – A at 85% AcN  
Flow 0.31 mL/min 0-10 min, 1.0 mL/min (6%B) 10-16 min, return to initial  
Sample: 2  $\mu$ L of 15 purine and pyrimidine bases and –sides (as above)



## HILIC Separation with LC-ESI/MS Detection of Nucleoside and Nucleobase Molecular Ions


Shimadzu LCMS-2020 Single Quadrupole: ESI (+) 4.5 kV, (-) 3.5 kV; SIM 1 amu; 0.25 sec each SIM Event (2 points/sec sample rate); 1.2 LPM N<sub>2</sub> neb; 12.5 LPM N<sub>2</sub> dry gas; Cap DL 300°C; Heat Block, 400°C



# Halo Penta-HILIC Separation of Nucleotides

Temperature: 35°C Flow Rate: 1.5 mL/min Injection Volume: 10.0µL UV Detection: 260nm  
 Mobile Phase: Pump A= 100mM Ammonium Formate pH=6.0  
 Pump B= Acetonitrile


ATP



% B	Tailing Factor 10%	Width Half-Height (min)
60	3.062	0.576
50	2.782	0.523
40	2.417	0.480

## Ammonium Phosphate Mobile Phase

4.6 mm ID x 100 mm; Flow Rate: 1.5 mL/min; Temp. 35°C; Abs (260nm); 1.0 uL Inject  
 A- 12.5 mM Amm Phos (pH 6.0) in 50% AcN, B- 6.25 mM Amm Phos (pH 6.0) in 25% AcN  
 Gradient 90%-40% B over 8 min.



Elution Order	Tailing Factor 10%
1 Adenosine Monophosphate	1.265
2 Guanosine Monophosphate	1.225
3 Adenosine Diphosphate	1.244
4 Guanosine Diphosphate	1.198
5 Adenosine Triphosphate	1.294
6 Guanosine Triphosphate	1.397

Retention of nDP's and nTP's overlap. Above is an example for limiting resolution for GDP and ATP ( $R_s=1.5$ ).  
 Overall the ammonium phosphate buffer shows excellent peak shape and allows resolution between the purine mono-, di-, and trinucleotides.

# Recent Developments for Larger Molecules

- 160 Å pore size Reversed-Phase for peptides and small proteins
  - Typical TFA Conditions
  - LC/MS for Proteomics
    - Ammonium Formate Conditions
    - Utility for Protein ID Workflow
- 400 Å pore size RP for Intact Proteins
  - Definition
  - Advantage for Larger Proteins
- Penta-HILIC for N-Linked Glycan and Glycopeptide Analysis


# Halo Peptide ES-C18

## Protein and Peptide Separations

Column: 4.6 x 100 mm; Flow rate: 1.5 mL/min; Temperature: 30° C


A: 0.1% TFA/10% ACN, B: 0.1% TFA/70% ACN

Gradient: 0% to 50% B in 15 min.; Injection volume: 5  $\mu$ L



### Sample 1


Gly-Tyr, Val-Tyr-Val, Met-enk, Angiotensin II, Leu-enk  
Ribonuclease, Porcine Insulin



### Sample 2

Leu-enk  
Bovine Insulin, Human Insulin, Cytochrome C, Lysozyme

# Effect of Pore Size on Efficiency




Columns: 4.6 x 100 mm HALO C18, 2.7 μm, 90 Å  
4.6 x 100 mm HALO Peptide ES-C18, 2.7 μm, 160 Å

Temperature: 60 °C  
Detection: 215 nm

Mobile Phase: Leu-Enk: 21% ACN/79% Water/0.1% TFA  
β-amyloid (1-38) 160 Å : 29% ACN/71% Water/0.1% TFA  
β-amyloid (1-38) 90 Å : 27% ACN/73% Water/0.1% TFA

# Recent Developments for Larger Molecules

## Improving Retention and Peak Shape Using Ammonium Formate




McCalley, D. V., Effect of buffer on peak shape of peptides in reversed-phase high performance liquid chromatography. *J Chromatogr* **2004**, *1038* (1-2), 77-84.  
 Schuster, S. A.; Boyes, B. E.; Wagner, B. M.; Kirkland, J. J., Fast high performance liquid chromatography separations for proteomic applications using Fused-Core® silica particles. *J Chromatogr* **2012**, *1228*, 232-241.



# Recent Developments for Larger Molecules

## Efficient Proteomic Analysis using Halo Peptide ES-C18

Column: 0.2 x 150 mm Halo Peptide ES-C18; Flow: 4  $\mu$ L/min  
 Gradient: 2 - 56% B in 85 min; Pmax - 320 bar;  
 A: 0.1% formic acid/10 mM AF/water; B: 80% acetonitrile/A;  
 Sample: 5 pmol transferrin, carbonic anhydrase, and  
 apomyoglobin digest mixture  
 Detection: Thermo LTQ Ion Trap MS/Michrom ESI interface



### Canine Prostate: Identification of Soluble Proteins

Mobile Phase Modifier	Protein IDs	Matched Spectra	Peptide IDs	Spectra/Peptide ID
0.1% Formic Acid	70	1142	359	3.18
0.1% FA, 10mM Ammonium Formate	118	2028	538	3.77

$\approx$  50%

Mobile Phase Modifier	X Bar Peak Width (s)	Standard Deviation	Measured Peak Capacity
0.1% Formic Acid	16.11	2.98	76.30
0.1% FA, 10mM Ammonium Formate	13.99	2.35	99.55

$\approx$  30%

- Protein IDs: validated using 5% protein false discovery rate
- Peptide IDs: validated using 5% peptide false discovery rate


# Recent Developments for Larger Molecules

Fused-Core Particle	Particle Size, $\mu\text{m}$	Pore Size, $\text{\AA}$	BET Surface Area, $\text{m}^2/\text{g}$	Shell Thickness, $\mu\text{m}$	% Porosity	Pore volume, $\text{cm}^3/\text{g}$
Halo	2.7	90	135	0.5	75	0.26
Halo Peptide	2.7	160	80	0.5	75	0.29
Wide-pore	2.7	400	30	0.35	59	0.23
Wide-pore	3.4	400	15	0.2	31	0.11

Wagner, Schuster, Boyes, Kirkland (2012) J. Chromatogr. 1264, 22-30.


# Recent Developments for Larger Molecules

Columns: 4.6 x 100 mm; Instrument: Agilent 1100; Flow rate: 1.5 mL/min; Injection Volume: 5  $\mu$ L; Detection: 220 nm;  
Mobile Phase: A: 10% ACN / 0.1% TFA; B: 70% ACN / 0.1% TFA; Gradient: 0% B to 50% B in 15 min.; Temperature = 30  $^{\circ}$ C




# Mass Transfer Dependence on Protein Size


Column: 4.6 x 100 mm , ES-C8, 400 Å, 2.7 μm; Mobile phase: ACN/water/0.1% TFA (adjusted to achieve  $k \sim 3$ ); Temp: 60 °C



# Improved Mass Transfer for Halo 400




## *Protein N-linked Glycan Analysis*



- Release of glycoprotein Asn-linked glycans by PNGase F is a well established approach, allow analysis of complex mixtures of glycans by a variety of methods.
- Various additional exoglycosidases and endoglycosidases can be helpful in structure elucidation.
- Considerable work has been accomplished for N-linked analysis of “native” glycans, as well as permethylated glycans using MS (ToF and ESI), and MS coupled to high resolution separation methods, particulary CE and LC.
- Glycoprofiling of N-linked glycans, using HILIC (NP) combined with fluorescence detection and/or MS is becoming a standard for biopharmaceutical glycoproteins, including antibodies.

# Analysis of PNGase Released and Labeled N-Glycans

Release of protein N-linked glycans using PNGase F releases oligosaccharides with a free reducing terminus (alditol) that is readily labeled by amines via formation of a Schiff's base, which can be reduced readily.



Many amines have been applied to labeling glycans, (Harvey, 2011, J. Chrom. [879](#)).

In the current work Procainamide is favored due to reported improvements in ESI-MS detection (Klapoetke, et al., 2010, J. Pharm. Biomed. Anal. [53](#))

## Typical Labeling Conditions

Glycan in water (up to 10% volume)

90+% volume of:

0.4 M procainamide


1M sodium cyanoborohydride

in 30% acetic acid/70% DMSO

12-16 hr reaction at 37°C

SEC cleanup on Sephadex G-10 minicolumn

Absorbance Detection 300 nm or Fluorescence Ex 330/Em 380 nm



Mass: glycan + 219.32

# Column Efficiency Comparisons using Pam-G<sub>5</sub>

2.1 (2.0 mm) ID x 150 mm, 60°C,  $k' \approx 6$ ,

50 mM Ammonium Formate Aqueous, pH 4.4

0.5  $\mu$ L Injection (50 pmol), Abs. 300 nm

Fitted to Knox Equation of form:

$$h = Av^{0.33} + \frac{B}{v} + Cv$$

$P = 147 \text{ bar (0.65 mL/min)}$


● Penta-HILIC, 72% AcN

■ 3.0  $\mu$ m Amide, 69% AcN

▲ 1.8  $\mu$ m Amide, 67% AcN

$P = 718 \text{ bar (1 mL/min)}$

$P = 331 \text{ bar (1 mL/min)}$







# Penta-HILIC Separations of Labeled Oligosaccharides and Glycans

2.1 mm ID x 150 mm; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600  $\mu$ L/min.


Detection: 300 nm Abs; ESI-MS (MS-2020, (+) 4.2 kV, 400-2000 with SIM)



# Penta-HILIC Retention of Labeled Oligosaccharides and Glycans

2.1 mm ID x 150 mm; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600  $\mu$ L/min.



Detection: 300 nm Abs; ESI-MS (MS-2020, (+) 4.2 kV, 400-2000 with SIM)



# Penta-HILIC Separations of Abundant bov. a1-AGP N-glycans

2.1 mm ID x 150 mm; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600  $\mu$ L/min.



MALDI of Permethylated Glycans



# HILIC Separations of $\alpha$ 1-AGP N-glycans

2.1 mm ID x 150 mm; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600  $\mu$ L/min.



Detection: 300 nm Abs; ESI-MS (+) MS-2020, 4.2 kV, 400-2000 with SIM



# High Resolution HILIC Separations of asialo-Fetuin N-glycans

2.1 mm ID x 300 mm; 50 mM Ammonium Formate, pH 4.4, 70-55% AcN (B) in 90min., 60°C; 600  $\mu$ L/min.



Detection: 300 nm Abs; ESI-MS (+) MS-2020, 4.2 kV, 400-2000 with SIM



# High Resolution HILIC Separations of Fetuin N-glycans

2.1 mm ID x 300 mm; 50 mM Ammonium Formate, pH 4.4, 70-55% AcN (B) in 90min., 60°C; 600 µL/min.



**Detection: 300 nm Abs; ESI-MS (+) MS-2020, 4.2 kV, 400-2000 with SIM**



# HILIC Separations of $\alpha$ 1-AGP N-glycans

2.1 mm ID x 150 mm; 50 mM Ammonium Formate, pH 4.4, 77.5-56.5% AcN (B) in 52.5 min., 60°C; 600  $\mu$ L/min.

Detection: 300 nm Abs; ESI-MS (+) MS-2020, 4.2 kV, 400-2000 with SIM



# Summary of N-Glycan Analysis

- The HILIC material and method described can be operated in high resolution mode, permitting analysis of complex mixtures of released N-linked protein glycan samples to be resolved, or in very high speed mode, to enable high throughput glycan analysis.
- High resolution HILIC resolves isobaric glycans, which are abundant (and complex).
- Variations in isobaric glycan profiles between proteins is very possible.
- *Quantitation by end-labeling using UV or Fluorescence detection only could be problematic, detection by MS only will not easily allow quantification of complex glycans (reference problem).*



# Summary of Recent Developments

- Fused-Core superficially porous particle selection continues to grow, with expansion of particle size and pore size availability.
- Bonded phase availability for these particles is expanding, particularly for RP and HILIC modes of operation.
- Research to apply these materials in useful ways to small and large molecules continues, resulting in additional opportunities for particle construction and improvement.
- *Improvement in these LC column packing materials is closely tied to application needs that are emerging in pharma and biopharma.*

# Acknowledgements

Thank you for your Attention!

- AMT
  - Dr. S. Schuster, Bob Moran, William Johnson, William Miles
- UGA
  - Prof Ronald Orlando, DJ Johnson, Shujuan Tao
- NIH for Financial Support
  - NIH NIGMS GM099355, GM093747
- Shimadzu Corp – early access to high sensitivity flow cell and use of M30A PDA detector