

Application note # C-13193

Separation of TPH (Total Petroleum Hydrocarbons) using the HI-SE54 HT Column




Figure 1: Separation of TPH (Total Petroleum Hydrocarbons) using the Avantor® Hichrom HI-SE54 HT phase. (Green chromatogram: C10-C40 alkane mix (50 mg/L), Blue chromatogram: fuel oil degradation mix (2000 µg/mL)).




Figure 2: Magnified image of the blue chromatogram for the separation of C17 and Pristane, C18 and Phytane in the fuel oil degradation mix (2000 µg/mL).

Method Details

CONDITIONS

Oven Program:	50 °C, 15 °C/min, 360 °C
Carrier Gas:	Hydrogen, 110 kPa
Injector:	Split 300 °C, 50 mL/min split flow, 0.2 µL injection volume
Detector:	FID, 370 °C
Samples:	C10-C40 Alkane Mix (50 mg/L) Fuel Oil Degradation Mix (2000 µg/mL)

The HI-SE54 HT is a high temperature general purpose column with a low-polarity 5% phenyl, 1% vinyl, 94% methyl polysiloxane phase composition. It is well suited for the analysis of high boiling petroleum products and long chained hydrocarbons.

ORDERING TABLE

Product	Details	Dimensions	Part Number
Avantor® Hichrom HI-SE54 HT	GC Column	0.15 mm, 0.10 µm, 10 m	HI91-15-010-10