




### The Impact of Column Length on Oligonucleotide Resolution in HPLC



#### TEST CONDITIONS:

Column: HALO 1000 Å OLIGO C18, 2.7 μm, 2.1 x 100 mm  
 Part Number: P2762-602  
 Column: HALO 1000 Å OLIGO C18, 2.7 μm, 2.1 x 250 mm  
 Part Number: P2762-902  
 Mobile Phase A: (90)/5/5 (10mM DiBA/100mM HFIP)/  
 MeOH/ACN  
 Mobile Phase B: 80/20 Water/ACN  
 Gradient: Time %B  
 10cm 0.0 30  
 60.0 80  
 25cm 0.0 30  
 150.0 80

Flow Rate: 0.2 mL/min.  
 Back Pressure: 111 bar - 100mm  
 206 bar - 250mm  
 Temperature: 60 °C  
 Injection: 1.0 μL of ssDNA, 20-100 ladder (10μg/mL)  
 Sample Solvent: 10mM Tris/1mM EDTA  
 Wavelength: PDA, 260 nm  
 Flow Cell: 1 μL  
 Data Rate: 40 Hz  
 Response Time: 0.05 sec.  
 LC System: Shimadzu Nexera X2

The resolution of an oligonucleotide mixture was evaluated using a longer column (25 cm) under a shallow gradient, successfully resolving species up to 100 nucleotides in length. The resolution scaled approximately with the square root of column length ( $\sqrt{L}$ ), consistent with chromatographic theory when the gradient rate was adjusted proportionally to column length—specifically, comparing 150-minute and 60-minute gradients for the 25 cm and 10 cm columns. Across the oligonucleotide size range, resolution increased by a factor of 1.34, which is close to the theoretical expectation of 1.58. This deviation may reflect practical limitations such as extra-column dispersion, gradient mixing effects, or non-idealities in oligonucleotide diffusion and interaction kinetics. Nonetheless, the data supports the expected trend that longer columns, when paired with appropriately scaled gradients, enhance resolution for complex oligonucleotide mixtures, particularly in the context of gradient elution with acetonitrile (AcN) under DiBA/HFIP conditions.