

The use of Short 10 mm Columns for Rapid LC-MS Residue Testing

Matt. James¹, Mark Fever¹, Tony Edge¹, Keng Tiong Ng², Leon Barron², Geoff Faden³


¹Avantor, 1-3 The Markham Centre, Station Road, Theale, Reading, UK;
²Environmental Research Group, School of Public Health, Imperial College London, UK;
³MAC-MOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317, USA; gfaden@mac-mod.com

1. Background


- Higher throughput LC-MS analysis is required in many environmental labs:
 - Increasing contaminants of concern / stricter regulatory limits.
 - Large scale monitoring campaigns
- Samples typically involve complex matrices, e.g. soils, wastewater etc.
- LC-MS separations can substantially reduce sample complexity:
 - Fast/simple sample prep strategies employed.
 - Typically employ short LC columns (e.g. 50 x 2.1 mm).
 - High sample throughput.
 - LC can now be a bottleneck.
- Improved MS performance (sensitivity & data acquisition) opens new possibilities for fast separations and increased throughput:
 - Opportunities for faster separations using specially designed high throughput columns.

3. 10 x 2.1 mm columns and considerations for use


- Avantor® ACE® HTP-MS: 10 x 2.1 mm cartridge style column:
 - 2 µm particles, 1,000 bar limit.
 - Male outlet (for connection to grounded inlet).
 - Cartridges individually QC tested.




- Extra column dispersion:
 - Highly impactful for small format columns.
 - Minimise wherever possible.



- Assessment of impact of tubing configuration on performance:
 - Maximum performance when post column dispersion is minimised.
 - Avantor® ACE® HTP-MS can be used in the column oven with suitably optimised connections.




- Detector data rates:
 - Fast detector sampling rates/dwell times are required for small volume peaks.
 - >10-15 data points per peak for accurate quantitation.



2. Development of Avantor® ACE® HTP-MS column design


- Requirements for high throughput:
 - High flow rate/linear velocity.
 - Short column length.
 - Maximise chromatographic efficiency.
- Consider dispersion theory (van Deemter) & kinetic performance.



- System pressure limits restrict performance obtainable for small particles.
- Short columns with small particles provide fast, high efficiency separations.

4. Applications

PFAS Analysis




1. PFBS
2. PFHxA
3. PFHxS
4. PFOA
5. PFNA
6. PFOS
7. PFDA
8. N-MeFOSAA
9. PFUnDA
10. N-EtFOSAA
11. PFDoA
12. PFTeDA
13. PFTeDA

- Analysis of 13 PFAS compounds by LC-MS/MS
- 1 minute gradient @1.0 mL/min.
- 28 samples per hour
 - 0.3 min post gradient re-equilibration
 - 30 s injector cycle time

Screen 28 samples per hour

Column:	Avantor® ACE® HTP-MS, 2 µm, 10 x 2.1 mm	Gradient:	<table><tr><th>Time (mins)</th><th>% B</th></tr><tr><td>0</td><td>5</td></tr><tr><td>0.05</td><td>40</td></tr><tr><td>1.00</td><td>100</td></tr><tr><td>1.20</td><td>100</td></tr><tr><td>1.30</td><td>5</td></tr><tr><td>1.60</td><td>5</td></tr></table>	Time (mins)	% B	0	5	0.05	40	1.00	100	1.20	100	1.30	5	1.60	5
Time (mins)	% B																
0	5																
0.05	40																
1.00	100																
1.20	100																
1.30	5																
1.60	5																
Mobile Phases:	A: 20 mM ammonium acetate in H ₂ O B: MeOH																
Flow Rate:	1.0 mL/min																
Temperature:	22 °C																
Injection Volume:	1 µL																
Detection:	Sciex QTRAP® 6500+ LC-MS/MS system.																
	Ionisation mode: ESI, negative mode; Source temperature: 450 °C; Curtain gas: 38 psig; Ionpray™ source voltage: -4500 V; Ion source gas 1: 40 psig; Ion source gas 2: 70 psig																
Sample:	10 ng/mL																

Illicit Drugs in Water



- 2.8 minute gradient @0.5 mL/min.
- Methodology can be applied to analysis of contaminants of concern in wastewater samples.
- Method details: Avantor® ACE® Application Note #7600.
- 135 simultaneous contaminants have been analysed in a single run using this type of methodology.
 - J. Haz. Mat **398** (2020) 122933

26 illicit drugs as part of 135 contaminant screen in 4.0 minutes

Column:	Avantor® ACE® HTP-MS, 2 µm, 10 x 2.1 mm	Gradient:	<table><tr><th>Time (mins)</th><th>% B</th></tr><tr><td>0</td><td>10</td></tr><tr><td>0.2</td><td>10</td></tr><tr><td>3</td><td>60</td></tr><tr><td>3.01</td><td>100</td></tr><tr><td>4</td><td>100</td></tr><tr><td>4.01</td><td>10</td></tr></table>	Time (mins)	% B	0	10	0.2	10	3	60	3.01	100	4	100	4.01	10
Time (mins)	% B																
0	10																
0.2	10																
3	60																
3.01	100																
4	100																
4.01	10																
Mobile Phases:	A: 0.1% formic acid in H ₂ O B: 0.1% formic acid in MeOH:MeCN (1:1 v/v)																
Flow Rate:	0.5 mL/min																
Temperature:	Ambient																
Injection Volume:	10 µL																
Detection:	Shimadzu LCMS 8060																
	Ionisation mode: +/- switching; Nebulising gas flow: 2 L/min; Heating gas flow: 15L/min; Interface temperature: 400 °C; DL temperature: 250 °C; Heat block temperature: 400 °C; Drying gas flow: 3 L/min																
Sample:	5,000 ng/L standard Mix in ultrapure water																
Data courtesy of Dr Leon Barron & Dr Keng Tiong Ng																	

5. Conclusions

- Improved MS performance provides scope for reducing LC-MS run times in high throughput settings.
- Avantor® ACE® HTP-MS 10 mm columns provide an effective solution.
- Use of 10 x 2.1 mm columns requires consideration of:
 - Extra column dispersion.
 - MS detector sampling rates/dwell times.
- High efficiency, ultra fast LC-MS analyses are achievable environmental analysis.