

# Theory & Practice of Developing LC Methods with Solid-Core Particle Columns

Gemma Lo

[Gemma.lo@avantorsciences.com](mailto:Gemma.lo@avantorsciences.com)

In Partnership with MAC-MOD,  
our authorized channel partner.  
MAC-MOD have been a trusted  
partner for over 20 years with  
expertise in HPLC & UHPLC



# Solid Core Particle Technology

- 
- Solid core<sup>1</sup> particles have gained interest for UHPLC / HPLC due to high efficiencies, rapid separations, method transferability, and low back pressure
  - Solid core particles in the 2.X  $\mu\text{m}$  range offer the potential of **sub-2  $\mu\text{m}$  efficiencies with HPLC pressures**
  - Solid core particles can **accelerate Method Development** on standard HPLC instrumentation
  - Extensive theoretical and practical assessments of solid core particles have been reported. A brief summary is provided here
- 

<sup>1</sup>Also known as Fused Core™, Core Shell, Core Enhanced, Partially Porous or Superficially Porous Particles (SPP).

# Solid Core Particle Technology


- Particle architecture

2.4-2.7  $\mu\text{m}$  & 5 $\mu\text{m}$  typical

Smaller particles also available (1.3-1.7  $\mu\text{m}$ )  
(with their own challenges!)


Rho ( $\rho$ ) = solid core diameter : particle diameter ratio  
Typically 0.6 – 0.75 for SPP

- Why do solid core particles
  - Give more efficiency? Give faster analyses? Give lower back pressure?



# Efficiency

## VAN DEEMTER CURVE




- A** Eddy diffusion (analyte paths, packing, wall effects)
- B/u** Analyte longitudinal / axial diffusion
- C.u** Analyte mass transfer between stationary & mobile phases


# Solid Core Efficiency Summary Facts

- The A term is not constant at low flows and reductions in trans-column eddy dispersion and wall effects / particle roughness may be significant<sup>2</sup> leading to higher efficiencies
- The B/u term is significant and improves efficiency by reducing analyte molecular diffusion processes
- Shorter molecular diffusion paths (C.u term) are **NOT** the reason for improved efficiency. However they are more dominant for large molecules
- The tight particle size distribution of solid core particles & packing quality have **limited influence** on the improved efficiency<sup>1,2</sup>

1 A Daneyko et al., Anal. Chem. 83 (2011) 3903-3910.  
2. F. Gritti et al., J. Chromatogr. A 1218 (2011) 8209-8221.




# Efficiency / Flow Comparison: Porous and Solid Core



Isocratic analysis, 50x2.1mm columns, eluent = MeCN / water + 0.1% TFA, analyte = naproxen, constant  $k = 10$ , 40°C,  $\lambda=256$  nm

# Pressure / Flow Comparison: 2 $\mu$ m Porous & 2.5 $\mu$ m Solid Core



Isocratic analysis, 50x2.1mm columns, eluent = MeCN / water + 0.1% TFA, analyte = ketoprofen, constant k = 10, 40°C,  $\lambda$ =256 nm



# Solid Core Particles Give Faster Analyses

There are 2 aspects as to why solid core particles offer **faster analyses** than porous equivalents:

**Higher linear velocities for a given efficiency.**



**Reduced particle surface area gives lower hydrophobicity for a bonded phase.**

When moving from a 5 $\mu$ m fully porous particle to a 2.7 $\mu$ m solid core particle, you could trade some efficiency & reduce column length to speed up a method even further

# Effect on Peak Capacity

- 
- Peak capacity is a measure of the number of sample analytes that can be separated on an HPLC column per unit time
  - Narrow peaks (increase in column efficiency) increase the peak capacity and efficiency of analytical peaks
  - Solid core particles have a greater peak capacity & improved S/N than 5 or 3  $\mu\text{m}$
-


# Isocratic Aspirin Analysis: Porous and Solid Core Columns



**Reduced hydrophobicity of solid core particles leads to 'faster' analysis**

Conditions: (Top): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (2.5 Hz), 25°C, 1 mL/min, 5  $\mu$ L injection  
(Bottom): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (20 Hz), 25°C, 1 mL/min, 3.9  $\mu$ L injection

# Efficiency / Flow Comparison: Faster Analyses






Isocratic analysis, 50x2.1mm columns, eluent = MeCN / water + 0.1% TFA, analyte = ketoprofen, constant  $k = 10$ , 40°C,  $\lambda = 256$  nm

# Solid Core Particles Are Less Hydrophobic

- 
- In RPLC hydrophobicity is a dominant mechanism
  - The hydrophobicity of a stationary phase is related to the ligand characteristics (e.g., C4, C8, C18) and amount present on the particle (i.e. %C or  $\mu\text{mol}/\text{m}^2$ )
  - Reducing the hydrophobicity of a bonded phase (e.g., by decreasing the chain length or reducing the coverage on the particle) reduces retention in RPLC
  - This can be easily demonstrated under the same conditions
-

# The Influence of Chain Length on Hydrophobic Character



PHASE	USP LISTING	FUNCTIONAL GROUP	ENDCAPPED	PARTICLE SIZE (μm)	PORE SIZE (Å)	SURFACE AREA (m <sup>2</sup> /g)	CARBON LOAD (%)	pH RANGE
Avantor® ACE® Traditional Chemistries								
C18	L1	Octadecyl	Yes	1, 2, 3, 5, 10	100	300	15.5	2 – 8
C8	L7	Octyl	Yes	2, 3, 5, 10	100	300	9	2 – 8
C4	L26	Butyl	Yes	2, 3, 5, 10	100	300	5.5	2 – 8


Sample: 1. Norephedrine 2. Nortriptyline 3. Toluene 4. Imipramine 5. Amitriptyline

Column: 250 x 4.6mm 5μm Mobile phase: 80:20 v/v MeOH/25mM KH<sub>2</sub>PO<sub>4</sub> (pH6.0) Flow: 1.0mL/min, Wavelength: 215nm

Gives  
hydrophobicity  
differences

# Solid Core Particles Are Less Hydrophobic

- Typical porous C18 phases have %C values >10%
- Solid core particles have lower surface areas so less ligand is bonded leading to lower hydrophobicity / smaller %carbon values & faster elution / analyses

		Particle size (um)	Pore size (A)	Surface area (m2/g)	% Carbon load
Solid core → Porous → 	ACE UltraCore SuperC18	2.5	95	130	7.0
	ACE C18	2,3,5,10	100	300	15.5
	Agilent Poroshell SB-C18	2.7	120	130	8
	Zorbax SB-C18	1.8, 3.5, 5	80	180	10
	AMT HALO C18	2.7	90	150	7.7
	N/A	-	-	-	-
	Phenomenex Kinetex C18	2.6	100	200	12
	Luna C18(2)	2.5,3,5,10	100	400	17.5
	Thermo Accucore C18	2.6	80	130	9
	GOLD C18	1.9,5,8,12	175	220	10
	Waters Cortecs C18	2.7	90	100	6.6
	ACQUITY BEH C18	1.7,3.5	185	130	17.7

# Method Development

SAME CONSIDERATIONS FOR SOLID CORE AS FULLY POROUS

- Pressure
- Column Dimensions / Particle Size / Pore size
- Column Chemistry
- Solvents (type, gradient, modifier etc.)
- Temperature
- pH

Nb. If transferring methods from fully porous to solid core consider injection volume, flow rate, gradient & instrument setup



# Resolution, Selectivity, Efficiency & Retention

Dispersion, particles etc

Phase design, eluent etc


Efficiency

Selectivity

Retention

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha-1}{\alpha} \cdot \frac{k}{1+k}$$


Selectivity is the key to resolution and efficiency boosts performance



Zhao, J.H. and P.W. Carr. Analytical Chemistry, (1999) 71, 2623-2632

# What is Selectivity?


- Alpha ( $\alpha$ ) denotes the separation factor or separation selectivity between 2 adjacent peaks



- Selectivity values  $> 1.0$  indicate the combination of mobile phase and stationary phase are providing some degree of separation for the 2 analytes

# Which Factors<sup>1</sup> Affect Selectivity?


- Strongly influenced by physicochemical properties of the analyte, stationary phase, eluent etc
- From a practical perspective:

Isocratic Separations	MOST Influence  LEAST Influence	Gradient Separations
<ul style="list-style-type: none"><li>– Column stationary phase type</li><li>– pH (ionisable analytes only)</li><li>– Organic modifier type</li><li>– % Organic modifier</li><li>– Buffer selection</li><li>– Column temperature</li><li>– Buffer concentration</li></ul>		<ul style="list-style-type: none"><li>– All parameters for isocratic <b>PLUS</b></li><li>– Gradient steepness,</li><li>– <math>k^* (t_G, F, V_m, \Delta\Phi, M),</math> <math display="block">k^* = \frac{t_G F}{\Delta\Phi V_m M}</math></li><li>– Dwell volume,</li><li>– Column dimensions.</li></ul>

<sup>1</sup> Adapted from 'Introduction to Modern Liquid Chromatography', 3rd Edition, Snyder, Kirkland, Dolan, 2010, p.29, Wiley & sons


# Exploring Selectivity: Porous Silica Bonded Phase Effects

ACE C18 – Increase Retention




Hydrophobicity  
Differences

ACE C8 (start point)



ACE C4 – Decrease Retention




Sample: 1. Norephedrine 2. Nortriptyline 3. Toluene 4. Imipramine 5. Amitriptyline


Column: 250 x 4.6mm 5 $\mu$ m Mobile phase: 80:20 v/v MeOH/25mM KH<sub>2</sub>PO<sub>4</sub> (pH6.0) Flow: 1.0mL/min, Wavelength: 215nm

# Exploring Selectivity: Porous Silica Bonded Phase Effects

ACE C18 – Increase Retention




ACE C4 – Decrease Retention




Mechanism  
Differences


ACE C8 (start point)



ACE CN – Elution Order



ACE Phenyl – Elution Order





**Stationary Phase Is Powerful With Selectivity & Retention**

Sample: 1. Norephedrine 2. Nortriptyline 3. Toluene 4. Imipramine 5. Amitriptyline

Column: 250 x 4.6mm 5 $\mu$ m Mobile phase: 80:20 v/v MeOH/25mM KH<sub>2</sub>PO<sub>4</sub> (pH6.0) Flow: 1.0mL/min, Wavelength: 215nm

# Exploring Selectivity: Eluent pH

ELUENT PH EFFECTS CAN BE LARGE...AND MULTIMODAL.



Eluent pH is Powerful For Selectivity and Retention

# Avantor® ACE® UltraCore™ Solid Core Particles: Selectivity

- Silica based solid core particles
- **SuperC18** and **SuperPhenylHexyl** bonded phases for alternative selectivity: hydrophobic / aromatic interactions
- **Encapsulated Bonding Technology** provides inertness & protects the silica surface from eluent pH 1.5 – 11.0

## ACE UltraCore 2.5µm:


Total particle diameter = 2.5µm

Shell thickness = 0.45µm

## ACE UltraCore 5µm:

Total particle diameter = 5µm

Shell thickness = 0.7µm



# Method Development / Screening Workflow: Overview

## TYPICALLY MULTIVARIATE

- 1 column
- 1 temperature
- 1 pH
- 1 organic modifier
- 1  $t_G$

2 x  $t_G$

- 1 column
- 2 temperatures
- 1 pH
- 1 organic modifier
- 2 x  $t_G$

20C & 60C

- 1 column
- 2 temperatures
- 1 pH
- 2 organic modifier
- 2 x  $t_G$

MeOH & MeCN

- $\geq 2$  columns
- 2 temperatures
- 1 pH
- 2 organic modifier
- 2 x  $t_G$

Alkyl chains eg C18, C8  
Aromatic eg Phenyl, C18-  
AR or C18-PFP  
Polar eg C18-PFP, C18-  
Amide

- $\geq 2$  column
- 2 temperatures
- 2 or 3 pH
- 2 organic modifier
- 2 x  $t_G$


pH 2.5  
pH 7  
pH 10.7

– INCREASING COMPLEXITY ... BUT KNOWLEDGE RICH

- Many potential runs to fully explore variables and their effects on retention and selectivity
- Having phases to fully exploit all parameters is helpful
- Would be helpful to reduce parameter options...



# Solid Core Method Development / Screening Workflow



– INFORMATION RICH DATA BASED ON SELECTIVITY

2 column method development / screening approach based on selectivity data

# General Method Development Initial Conditions

- Perform a broad scouting gradient run on the samples at acidic eluent pH
- How do you calculate your starting conditions?

For a 100 x 3mm column:

$t_G$  = 5 minutes

$F$  = 1.2 mL/min

$\Delta\Phi$  = 0.95

$V_m$  = 0.459 mL


$M$  = 5

$$k^* = \frac{t_G F}{\Delta\Phi V_m M} = \sim 3$$

- Ideally retention (or  $k^*$  in gradient elution) should be  $>2$  and  $<20$  for initial method development

# Vendor A C18 vs Vendor B C18

$$\text{Selectivity} = 100 \times \sqrt[3]{1 - R^2}$$



$$\begin{aligned}\text{Selectivity} &= 100 \times \sqrt[3]{1 - R^2} \\ &= 100 \times \sqrt[3]{1 - 0.9998} \\ &= \underline{\sim 1}\end{aligned}$$


**Low Selectivity Value When Comparing C18 Phases To Each Other**

\* Neue, O'Gara, Méndez "Selectivity in Reversed-Phase Separations: Influence of the Stationary Phase", J. Chromatogr. A 1127 (2006), 161-174

# Selectivity Plot: Exploring The Effect Of Solid Core Phase

SuperC18, low pH, MeOH vs SuperPhenylHexyl, low pH, MeOH

Example application:





Selectivity = 19

Changes in peak spacing noted

SuperC18

SuperPhHex




50x2.1mm, 2.5µm, gradient analysis, A= 20mM HCOONH<sub>4</sub>, pH3 (aq), B= 20mM HCOONH<sub>4</sub>, pH 3 in MeCN/water 9:1 v/v, 3-100%B in 7.5 mins, hold 100%B for 1.5 mins, 40°C, 0.40 mL/min, 254 nm.

1 amiloride, 2 benzamide, 3 3-hydroxybenzoic acid, 4 vanillin, 5 2-hydroxybenzoic acid, 6 benzoic acid, 7 methyl paraben, 8 p-cresol, 9 cortisone, 10 ethyl paraben, 11 dimethylphthalate, 12 piroxicam, 13 hydrocortisone-21-acetate, 14 ketoprofen, 15 ethylbenzoate, 16 toluene, 17 valerophenone, 18 mefenamic acid, 19 hexanophenone, 20 propylbenzene, 21 phenanthrene, 22 heptaphenone, 23 butylbenzene

# Selectivity Plot: Exploring Eluent pH With SuperPhenylHexyl


SuperPhenylHexyl, low pH, MeOH vs SuperPhenylHexyl, high pH, MeOH

Example application:



Selectivity = 83


Significant changes in elution order noted



50x2.1mm, 2.5mm, gradient analysis, A1= 10mM HCOONH<sub>4</sub>, pH3 (aq), B1= 10mM HCOONH<sub>4</sub>, pH 3 in MeOH/water 9:1 v/v, A2= 0.1% NH<sub>3</sub>, pH 10.7 (aq), B2= 0.1% NH<sub>3</sub>, pH10.7 in MeOH/water 9:1 v/v, 3-100%B in 5mins, 100%B for 2mins, 40°C, 0.60 mL/min, 254 nm.  
1. benzamide, 2. caffeine, 3. procainamide, 4. N-acetylprocainamide, 5. propiophenone, 6. toluene, 7. remacemide, 8. ethylbenzene, 9. carvedilol, 10. nortriptyline, 11. clomipramine.

# UltraCore: Exploring Phase, Solvent & pH Selectivity

SuperC18, low pH, MeCN vs SuperPhenylHexyl, high pH, MeOH




Range of 50 Analytes To Describe Selectivity

Selectivity = 85 → Fully Explore The Selectivity 'Space'

# Avantor® ACE® UltraCore™: Exploit Bonded Phase

## SELECTIVITY TUNING WITH STATIONARY PHASE TYPE:

UltraCore  
2.5 SuperC18



1. Pyridoxine
2. p-Aminobenzoic acid
3. Pantothenic acid
4. Folic acid
5. d-Biotin
6. Cyanocobalamin
7. Riboflavin


50x2.1mm, 2.5µm  
Gradient analysis

A: 20 mM  $\text{KH}_2\text{PO}_4$  pH 2.7 (aq)  
B: 20 mM  $\text{KH}_2\text{PO}_4$  pH 2.7 in  
MeOH:Water 50:50 v/v


T	%B
0	20
1	60
2	70

40 °C  
0.60 mL/min  
205 nm

UltraCore 2.5  
SuperPhenylHexyl




Cyanocobalamin (vitamin B12)



# Avantor® ACE® UltraCore™: Exploit Low and High pH Eluents

## UltraCore SuperC18: selectivity with pH



50x2.1mm, 2.5µm  
Gradient analysis

A1= 10mM HCOONH<sub>4</sub>, pH3 (aq)  
B1= 10mM HCOONH<sub>4</sub>, pH 3 in  
MeCN/water 9:1 v/v

A2= 0.1% NH<sub>3</sub>, pH 10.7 (aq)  
B2= 0.1% NH<sub>3</sub>, pH10.7 in in  
MeCN/water 9:1 v/v

T	%B
0	3
5	100
6	100

40C  
0.60 mL/min  
254 nm


1. Atenolol
2. Methylphenylsulfoxide
3. Eserine
4. Prilocaine
5. Bupivacaine
6. Tetracaine
7. 1,2,3,4-Tetrahydro-1-naphthol
8. Carvedilol
9. Nitrobenzene
10. Methdilazine
11. Amitriptyline
12. Valerophenone



# Other Considerations

# Loading Capacity & Retention

- A thinner shell can be **advantageous for efficiency** (van Deemter), but it **reduces sample loading** and **analyte retention**
- The optimum shell thickness is a compromise between **efficiency, sample loading capacity** and **analyte retention**
- Sample loading capacity and expected retention factor of a given solute are proportional to the stationary phase volume, and is expected to be lower on SPPs than on fully porous particles



# Scalability

SPPS OF DIFFERENT PARTICLE SIZES ARE SCALABLE FOR SELECTIVITY BUT NOT FOR RETENTION ON ACCOUNT OF DIFFERENT RHO VALUES.

For example :

ACE UltraCore 2.5 $\mu$ m:  $\rho = 0.64$

Total particle diameter = 2.5 $\mu$ m

Core diameter = 1.6 $\mu$ m

Shell thickness = 0.45 $\mu$ m

Surface area = 130 m<sup>2</sup>/g

ACE UltraCore 5 $\mu$ m:  $\rho = 0.72$


Total particle diameter = 5 $\mu$ m

Core diameter = 3.6 $\mu$ m

Shell thickness = 0.7 $\mu$ m

Surface area = 100 m<sup>2</sup>/g

Rho ( $\rho$ ) = core diameter / particle diameter



# Extra Column Band Broadening

- 
- One of the great advantages of superficially porous particle packed columns is that the back pressures produced often allow the use of standard HPLC instrumentation (not sub 2 $\mu$ m SPP)
  - However, the HPLC system needs to be optimised to produce efficient chromatography
  - Failure to consider these parameters may result in loss of the increased efficiency offered by the SPP
  - Extra column effects are more significant for scaled down separations (as column volume decreases) and for less retained peaks which have a lower peak volume
-

# Summary and Overall Conclusions

- 
- Solid core columns offer efficiency and speed for separations
  - Method development is comparable to fully porous particles
  - Screening columns with differing retention mechanisms and exploiting eluent pH is useful for method development – explore 'selectivity space'
  - Perceived improvements in analysis time are likely to be due to the reduced hydrophobicity of the solid core particles
  - The impact of system dispersion is real and can be significant
-


# Useful Resources

- ACE Translation Tool:
- (+help file)
- (+ AKN#0023)

The screenshot shows the 'Method Transfer' window of the ACE software. It contains two main sections: 'Column Information' and 'Method'. The 'Column Information' section includes fields for Column Length (50 mm), Column I.D. (4.6 mm), Particle Diameter (5 µm), Np (2000), Column Porosity (0.88 µm/m), and Column Volume (0.223 mL). The 'Method' section includes fields for Flow Rate (1.000 mL/min), LC Name (12345), and Dwell Volume (0.076 mL). Below these fields are two tables: 'Gradient' and 'Transfered'. The 'Gradient' table has columns for Time, Np, and Gradient, with values ranging from 0.00 to 1.82. The 'Transfered' table has columns for Time, Np, and Gradient, with values ranging from 0.00 to 1.82.

- ACE Knowledge Notes (AKNs):
  - AKN0019 – Solid Core Technology
  - AKN0018 – RP Method Development
  - AKN0011 – Practical UHPLC
  - AKN0012 – Understanding the Relationship between Particle Size, Performance and Pressure
  - AKN0017 – How to Determine Extra Column Dispersion and Extra Column Volume
  - AKN0023 – Gradient Method Translation Using the ACE LC Translator

## ACE Method Development Kit Brochure and Webinar



Thank you

[info@mac-mod.com](mailto:info@mac-mod.com)

