

Translating Natural Product Methods to Multiply Coupled Columns to Achieve Ultra-Resolution Separations and the Potential Impact of Pressure Induced Changes in Selectivity

Alan P McKeown¹, Geoffrey Faden²

¹Advanced Chromatography Technologies Ltd, 1 Berry Street, Aberdeen, Scotland, AB25 1HF UK
²MACMOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317 USA

1. INTRODUCTION

- Complex samples require high efficiency separations to provide maximum resolving power for the high number of analytes encountered
- Columns packed with small fully porous or solid core particles can deliver very high efficiency and peak capacity values
- However, column length may be restricted by LC system pressure limits
- UHPLC systems with pressure limits >1,000 bar present possibilities for extending the effective column length by coupling columns to generate very large column efficiency and peak capacity values
- This work uses method translation principles to migrate methods to longer effective column lengths to achieve ultra-resolution separations.
- Key aspects and potential problems of this approach are considered

4. TRANSLATING METHODS TO LONGER COLUMNS

- Scale gradient volume to new column volume (V_M)^{*} by altering gradient time (t_G) or flow rate (F)

$$t_{G2} = \frac{t_{G1} F_1 V_{M2}}{V_{M1} F_2} \quad F_2 = \frac{t_{G1} F_1 V_{M2}}{V_{M1} t_{G2}}$$

- Scale injection volume (V_{inj}) to new column volume

$$V_{inj2} = V_{inj1} \times \left(\frac{V_{M2}}{V_{M1}} \right)$$

- To maintain selectivity, correct for change in system dwell (V_D) and V_M .

Determine whether a pre-gradient isocratic hold of x minutes is required or if the injection should be delayed by x minutes^{*}:

$$x = \left[\left(\frac{V_{D1}}{V_{M1}} \right) - \left(\frac{V_{D2}}{V_{M2}} \right) \right] \times \frac{V_{M2}}{F_2}$$

Pre-gradient hold = positive value, delayed injection = negative

^{*}Note – it is recommended to determine V_M experimentally via measurement of a non-retained dead volume marker


^{*}Petersson et al., LCGC Europe 28 (2015) 310-320

7. EXAMPLE II: ECHINACEA WITH ACE ULTRACORE


- Low backpressure of solid core particles provides additional options

Option 1: translation via changing gradient time


$$t_{G2} = \frac{t_{G1} F_1 V_{M2}}{V_{M1} F_2}$$



Column: ACE UltraCore 2.5 SuperC18
150 x 3.0 mm
 t_G : 20 min
Flow: 0.43 mL/min
Inj. Vol.: 5 μ L
 P_{MAX} : 130 bar




Column: ACE UltraCore 2.5 SuperC18
300 x 3.0 mm (2 x 150 x 3.0 mm)
 t_G : 40 min
Flow: 0.43 mL/min
Inj. Vol.: 10 μ L
 P_{MAX} : 223 bar




Column: ACE UltraCore 2.5 SuperC18
450 x 3.0 mm (3 x 150 x 3.0 mm)
 t_G : 60 min
Flow: 0.43 mL/min
Inj. Vol.: 15 μ L
 P_{MAX} : 338 bar

10. EXAMPLE III: GREEN TEA / PRESSURE EFFECTS I


- Method translated: constant linear velocity $t_{G2} = \frac{t_{G1} F_1 V_{M2}}{V_{M1} F_2}$



ACE Excel 1.7 SuperC18 100 x 3.0 mm
 P_{MAX} = 336 bar



ACE Excel 1.7 SuperC18 200 x 3.0 mm (2 x 100 mm columns coupled)
 P_{MAX} = 621 bar




ACE Excel 1.7 SuperC18 300 x 3.0 mm (3 x 100 mm columns coupled)
 P_{MAX} = 923 bar

- Peak capacity is increased with increasing column length
- A significant change in selectivity is observed for the two principal peaks




2. COLUMN COUPLING

- Column length can be maximised by selecting chromatographic conditions that generate favourable column backpressures
- Acetonitrile was selected over MeOH due to lower viscosity
- Elevated column temperature of 80C was used
- Note that the use of elevated temperature is dependant on sample stability
- To obtain maximum performance from small particles, it is possible to increase the flow rate to obtain increased efficiency



5. EXAMPLE 1: ST JOHN'S WORT

- Ultra high efficiency ACE columns packed with 1.7 μ m particles
- Method translation: maintain constant linear velocity and scale t_G






8. EXAMPLE II: ECHINACEA WITH ACE ULTRACORE

- Low backpressure of solid core particles provides additional options



Option 2: translation via changing flow rate

$$F_2 = \frac{t_{G1} F_1 V_{M2}}{V_{M1} t_{G2}}$$



11. EXAMPLE III: GREEN TEA / PRESSURE EFFECTS II

- Separation was re-run on the ACE Excel 1.7 Super C18 100 x 3.0 mm
- A 9.7 cm length of 25 μ m i.d. SS restrictor tubing was inserted post-column to artificially elevate the backpressure



- The observed change in selectivity was confirmed and attributed to the increased column backpressure generated when coupling columns.

3. EXPERIMENTAL

Samples

100 mg of finely ground Green Tea and St John's Wort tablets extracted with 10.0 mL MeCN:water 1:1 v/v for 15 minutes with ultrasonication. 100 μ L supernatant diluted with 300 μ L water and filtered using a Whatman Mini-Uniprep syringeless filter (0.45 μ m polypropylene filter media, VWR P/N 83009-820).

For echinacea, 1.0 mL of preparation was diluted with 9.0 mL MeCN:water 1:1 v/v and filtered.

Chromatographic parameters

System: Hitachi / VWR ChromasterUltra Rs and Agilent 1290

Columns: 3 x ACE Excel 1.7 SuperC18 100 x 3.0 mm

3 x ACE UltraCore 2.5 SuperC18 100 x 3.0 mm

Eluent: A: 0.1% formic acid (aq)

B: MeCN + 0.1% formic acid v/v

Gradient: see next slide




Flow rate: see next slide

Temp.: 80C

Detection: UV, various wavelengths



6. EXAMPLE 1: ST JOHN'S WORT...ZOOMED DETAIL



- Increased efficiency and resolution
- Excellent transfer of selectivity between effective column lengths

9. POSSIBLE PITFALL: ELEVATED PRESSURE EFFECTS

- High pressure can affect the retention characteristics of analytes

- Studies have shown this can be significant^{*}

- The degree of pressure induced change in retention is a complex phenomenon and can be affected by a multitude of factors including:

- Physico-chemical properties of the analyte (e.g. polarisability, degree of ionisation, molecular weight, solvation etc.)
- Mobile phase properties (% organic, pH)
- Type of stationary phase

- When translating any method to higher pressure separations for increased speed or increased resolution the selectivity of the translated method must be assessed as suitable

^{*}Fallas et al., J. Chromatogr. A 1209 (2008) 195-205

12. SUMMARY AND CONCLUSIONS

- A UHPLC column coupling approach to increase effective column lengths to generate very high resolution separations has been shown
- By carefully selecting chromatographic conditions, it is possible to couple up to 3 x ACE Excel 1.7 μ m SuperC18 100 x 3.0 mm columns
- When translating methods to long column formats, it is important to obey volumetric scaling principles
- Columns packed with solid core particles offer extra flexibility for ultra resolution separations due to reduced backpressure.
- This approach is suitable for the ultra resolution, high sensitivity analysis of complex natural products.
- Whilst column coupling and the resulting higher back pressures are useful, unexpected changes in selectivity (application dependent) may occur.