

Translating US Pharmacopoeia Methods to Sub-2 Micron and Solid Core Using the New USP <621> General Chapter Guidelines

Alan P McKeown¹, Geoffrey Faden²

¹Advanced Chromatography Technologies Ltd, 1 Berry Street, Aberdeen, Scotland, AB25 1HF UK
²MACMOD Analytical Inc., 103 Commons Court, PO Box 587, Chadds Ford, PA 19317 USA

1. Introduction – Monograph Testing

- Ensures the **safety** and **quality** of pharmaceutical products and can include an LC test typically for assay or purity
- Many monographs use legacy column formats (e.g. 250 x 4.6 mm, 10 µm)
- Advances in column technology (sub 2-micron fully porous and solid core particles) and small column formats (e.g. 50 x 3.0 mm) allow **substantial improvements** in productivity and large cost savings
- However, allowable changes in column formats specified within monographs have previously been tightly restricted
- For isocratic methods, the **revised USP <621>** (general chapter on chromatography) now provides improved flexibility to the chromatographer to use modern column technology as allowable changes to the LC method
- This poster **summarises the recent changes** and demonstrates how to **achieve productivity and cost savings** using both HPLC and UHPLC technology

4. Translating Isocratic Methods and L/d_p Approach

New USP <621> guidelines allow **two options** for changing the particle size (d_p) and column length (L):

- Keep L/d_p constant or within **-25% to +50%** of the original method
- Keep N constant or within **-25% to +25%** of the original method

This work explores the use of **option 1**.

For **successful translation** of isocratic LC methods, the following principles are applied:

Translation of flow rate:

Scaled to new column i.d. (d_c) to **maintain linear velocity**

$$F_2 = F_1 \times \left(\frac{d_{c2}^2}{d_{c1}^2} \right)$$

or scaled to new column i.d. and smaller d_p

$$F_2 = F_1 \times \left(\frac{d_{c2}^2 \times d_{p1}}{d_{c1}^2 \times d_{p2}} \right)$$

Scaling injection volume:


Injection volume is **scaled to new column volume (V_M)**

$$V_{inj2} = V_{inj1} \times \left(\frac{V_{M2}}{V_{M1}} \right)$$

7. Estradiol Optimisation: Fully Porous or Solid Core

- Solid core and fully porous options
- Compatible with **standard HPLC instrumentation** (400 bar system pressure)

	Original Method (A)	Translated Method 1 (B)	Translated Method 2 (C)
Column	ACE 10 µm C18 300 x 3.9 mm	ACE UltraCore 2.5µm SuperC18 100 x 4.6 mm	ACE Excel 2 µm C18 50 x 3.0 mm
L/d _p	30,000	40,000 (+33.3%)	25,000 (-16.7%)
Flow (mL/min)	1.0	1.39	0.59
Injection Vol. (µL)	25	10.1	2.5
Back pressure (bar)	79	145	222
System Suitability			
Rs between B and C > 2.0	4.8	5.1	4.0
%RSD multiple injections < 2.0	0.1	<0.1	0.2
Savings Achieved			
Analysis time per injection	10 min	3.3 min (67% reduction)	1.7 min (83% reduction)
Solvent used per injection	10 mL	4.6 mL (-54%)	1.0 mL (-90%)




- Reduction in analysis time up to **83%**
- Solvent consumption reduced up to **90%**

First published in Chromatography Today, Volume 8, Issue 1, February/March 2015

10. UHPLC Example: Hydrocortisone

- Translate method from **5 µm L1, 150 x 4.6 mm** to **1.7 µm L1, 50 x 3.0 mm**
- Approach 1: scale flow to maintain **constant linear velocity** (0.43 mL/min)
- Approach 2: scale flow to **reduced particle size** (1.25 mL/min)



Mobile phase: MeCN:MeOH:water (25:25:50 v/v/v)
 Detection: UV @ 254 nm
 Temperature: Ambient (22°C)
 Sample: 1. hydrocortisone, 2. propyl paraben

2. New USP <621> Guidelines: Mobile and Stationary Phase

	USP 36 / NF31 <621>	USP 37 / NF 32 <621>
Mobile phase		
Composition	Isocratic & gradient: - Minor components can be changed by ±30% relative or ±10% absolute	Isocratic: - Minor components can be changed by ±30% relative or ±10% absolute
pH	Isocratic & gradient: - ±0.2 units (1% for neutrals)	Isocratic & gradient: ±0.2 units
Ionic strength	Isocratic & gradient: - ±10% if the permitted pH variation is met	Isocratic & gradient: - ±10% if the permitted pH variation is met
Column		
Length	Isocratic & gradient: - ±70%	Isocratic: - Particle size (d _p) and length (L) may be changed if a) L/d _p is constant or varies -25% to +50% OR b) number of plates (N) is -25% to +50%
Particle size	Isocratic & gradient: -50%	Isocratic: - No changes
Internal diameter	Isocratic & gradient: - Any changes if linear velocity kept constant - ±25%	Isocratic: - Any changes if linear velocity kept constant - No changes


5. Example 1: USP Estradiol Assay

Translating method from **10 µm to 5 µm**

Column dimensions scaled (maintain L/d_p)

- Flow rate scaled (**constant linear velocity**)
- Inj. Volume scaled to V_M


- HPLC system compatible**



First published in Chromatography Today, Volume 8, Issue 1, February/March 2015

8. Example 2: Guaifenesin Tablets Assay

- When **excess resolution** is obtained, L/d_p can be reduced (**allowable change**)



First published in Chromatography Today, Volume 8, Issue 1, February/March 2015

11. UHPLC Example – Hydrocortisone assay

	Original Method (A)	Translated Method 1 (B)	Translated Method 2 (C)
Column	ACE 5 µm C18 150 x 4.6 mm	ACE Excel 1.7 µm C18 50 x 3.0 mm	ACE Excel 1.7 µm C18 50 x 3.0 mm
L/d _p	30,000	29,412	29,412
Flow (mL/min)	1.0	0.43	1.25
Injection Vol. (µL)	10	1.4	1.4
Back pressure (bar)	114	293	779
System Suitability			
Rs > 9.0	14.1	13.4	10.1
N > 3,000 for hydrocortisone	9,167	9,887	6,441
Tailing factor < 1.2	✓	✓	✓
%RSD multiple injections < 2.0	0.1	<0.1	<0.1
Savings Achieved			
Analysis time per injection	10 min	3.3 min (-66%)	1.2 (-88%)
Solvent used per injection	10 mL	1.4 mL (-86%)	1.5 mL (-85%)

- Alpha decreases by ~8%...possibly due to **pressure effects**
- 66% reduction in run time and 86% reduction in solvent use when scaling flow to maintain **constant linear velocity**.
- 88% reduction in run time and 85% reduction in solvent use when scaling flow to **reduced particle size**.

3. New USP <621> Guidelines: Operating Conditions

Method	USP 36 / NF31 <621>	USP 37 / NF 32 <621>
Flow rate	Isocratic & gradient: - $F_2 = F_1 \times \left(\frac{d_{c2}^2}{d_{c1}^2} \right)$ (Where d _c = column diameter and F = flow rate) - Or, flow rate may change ±50%	Isocratic: - If particle size has changed use following equation for similar performance: $F_2 = F_1 \times \left(\frac{d_{c2}^2 \times d_{p1}}{d_{c1}^2 \times d_{p2}} \right)$ (Where d _c = column diameter, F = flow rate and d _p = particle size) - Additional increase in flow allowed provided column efficiency does not drop below 20%. - Or, flow rate may change ±50%
Injection volume	Isocratic & gradient: - Any reduction	Isocratic & gradient: Any change as long as peak repeatability is satisfactory
Temperature	Isocratic & gradient: - ±10°C when temperature is listed	Isocratic & gradient: ±10°C when temperature is listed
Detection wavelength	Isocratic & gradient: - No change permitted. ±3 nm between detectors	Isocratic & gradient: No change permitted. ±3 nm between detectors

A. P. McKeown, Chromatography Today (2015) 32-36
 United States Pharmacopoeia General Chapter <621> "Chromatography" First Supplement to USP 37-NF 32 (United States Pharmacopoeial Convention, Rockville, MD, USA).

6. Exploring L/d_p: Estradiol

- For **isocratic methods**, reducing column dimensions whilst **maintaining column length (L) to particle size (d_p) ratio** results in **similar performance**
- USP <621> now permits L/d_p **-25% to +50%**
- E.g. for 10 µm 300 x 3.9 mm (= 30,000)

Particle size (µm)	Column Length (mm)						
	50	75	100	125	150	250	300
1.7	29,412	44,118	58,824				
1.8	27,778	41,667	55,556				
1.9	26,316	39,474	52,632				
2	25,000	37,500	50,000	62,500	75,000		
2.5	20,000	30,000	40,000	50,000	60,000	100,000	
2.6	19,231	28,846	38,462	48,077	57,692	96,154	
2.7	18,519	27,778	37,037	46,296	55,556	92,593	
3	16,667	25,000	33,333	41,667	50,000	83,333	
5	10,000	15,000	20,000	25,000	30,000	50,000	
10	5,000	7,500	10,000	12,500	15,000	25,000	30,000


- Various L/d_p options (and column formats) are available within the range **-25% to +50%** for use with HPLC and UHPLC instruments (**highlighted green**).

9. L/d_p and Flow Adjustment: UHPLC

- USP <621> also allows translation of the flow rate to a **higher linear velocity** to take advantage of **high efficiencies** achievable with **small particles**.

$$F_2 = F_1 \times \left(\frac{d_{c2}^2 \times d_{p1}}{d_{c1}^2 \times d_{p2}} \right)$$

- i.e. allows chromatographer to **fully exploit sub 2 micron particles** and operate under **UHPLC conditions**



12. Summary and Conclusions

- The changes for isocratic methods in the new USP <621> provide **considerable flexibility** to the chromatographer. Reading the full text is **highly recommended** for detailed explanations.
- Use of **small particles** and **solid core technology** is now accommodated, allowing **significant increases in productivity** and **reduced cost per analysis**.
- This work successfully demonstrates how the **L/d_p approach** can be applied to take advantage of the **latest column technology** using both **HPLC and UHPLC**.
- 80% reduction in run time** and **72% reduction in solvent use** for estradiol on an HPLC system.
- 80% reduction in run time** and **>90% reduction in solvent use** for guaifenesin on an optimised HPLC system.
- 88% reduction in run time** and **85% reduction in solvent use** for hydrocortisone on a UHPLC system.

ACEB is a registered trademark of Advanced Chromatography Technologies Limited. ACE UltraCore™ and EBT™ are trademarks of Advanced Chromatography Technologies Limited.