





### Improved Signal Intensity for Trastuzumab Peptides Using a HALO® 1.5 mm ID Column

BIO-325



This app note features the full MS scan of a sample of Trastuzumab that underwent trypsin digestion to produce peptide fragments of varying length that can then be used for peptide mapping. With the help of a computer program, an extracted ion chromatogram (XIC) of the sequence YADSVKGR is also featured. This XIC is to give a closer look at the benefits of a 2µm 1.5 mm ID peptide column. Samples with similar complexity require the increased efficiency from a smaller particle size and long shallow gradients. By switching from a 2.1 mm ID to a 1.5 mm ID not only can solvent usage be cut in half, but sensitivity can also be increased. This all can be achieved by using the HALO® 1.5mm ID 2µm ES-C18 product.



**TEST CONDITIONS:**

**Column:** HALO 160 Å ES-C18 , 2.0 µm, 1.5 x 150 mm

**Part Number:** 9112X-702

**Column:** HALO 160 Å ES-C18 , 2.0 µm, 2.1 x 150 mm

**Part Number:** 91122-702

**Mobile Phase A:** Water, 0.1% DFA

**Mobile Phase B:** Acetonitrile, 0.1% DFA

| Gradient: Time | %B   |
|----------------|------|
| 0.5            | 2    |
| 60.5           | 50   |
| 61.0           | 70   |
| 65.0           | 70   |
| 65.5           | 2    |
| 70.0           | Stop |

**Flow Rate:** 0.2 mL/min for 1.5 mm  
0.4 mL/min for 2.1 mm

**Pressure:** 372 bar 1.5 mm  
670 bar 2.1 mm

**Temperature:** 60 °C

**Injection Volume:** 1 µL

**Sample:** 1mg/mL Trastuzumab tryptic digest

**Sample Solvent:** 1.5M Guanidine HCl/0.5% Formic Acid/~50mM Tris pH: 7.8

**LC System:** Shimadzu Nexera X2

**TUBING OPTIMIZATION:**

50µm x 600mm Column to Diverter Valve

50µm x 350mm Diverter Valve to Ground

50µm x 100mm Ground to Source

**MS CONDITIONS:**

**System:** ThermoFisher Q Exactive

**Spray Voltage (kV):** 3.8

**Capillary temperature:** 320 °C

**Sheath gas:** 35

**Aux gas:** 10

**RF lens:** 50

