



Pharmaceutical Related Substances HPLC / UHPLC Method Development and Translations

Alan P McKeown
amckeown@ace-hplc.com

Business Development Director

Advanced Chromatography Technologies Ltd

www.ace-hplc.com






Outline

- ◆ What is **selectivity**?
- ◆ **Optimizing method development** based on selectivity
- ◆ Example **related substances** method development
- ◆ **Method translations** – isocratic and gradient
- ◆ **Maximizing resolution** for complex samples

What is Selectivity?

- ◆ Alpha (α) is the symbol used to denote the **separation factor** or **separation selectivity** between 2 adjacent peaks.



- ◆ **Selectivity** can be thought of as ‘**peak spacing**’.
- ◆ Selectivity values should be **> 1.0**

Resolution, Selectivity, Efficiency & Retention

Particle size, column length, dispersion etc

Phase design, eluent etc


Efficiency

Selectivity

Retention

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1+k}$$

Selectivity (α) is the key to resolution and efficiency (N) boosts performance



Which Factors¹ Affect Selectivity?

- ◆ Strongly influenced by physicochemical properties of the analyte, stationary phase, eluent etc
- ◆ From a practical perspective:

Isocratic Separations

- ◆ Column stationary phase type
- ◆ pH (ionisable analytes only)
- ◆ Organic modifier type
- ◆ % Organic modifier
- ◆ Buffer selection
- ◆ Column temperature
- ◆ Buffer concentration

**MOST
Influence**



**LEAST
Influence**

Gradient Separations

All parameters for isocratic **PLUS**

Gradient steepness,

$k^* (t_G, F, V_m, \Delta\Phi, M)$,

$$k^* = \frac{t_G F}{\Delta\Phi V_m M}$$

Dwell volume,

Column dimensions.



Method Development Strategies: Overview

◆ **Selectivity is integral** to method development

- 1 column
- 1 temperature
- 1 pH
- 1 organic modifier
- 1 x tG

- 1 column
- 2 temperatures
- 1 pH
- 1 organic modifier
- 2 x tG

20C & 60C

- 1 column
- 2 temperatures
- 1 pH
- 2 organic modifier
- 2 x tG

MeOH & MeCN


- ≥ 2 columns
- 2 temperatures
- 1 pH
- 2 organic modifier
- 2 x tG

Alkyl chains eg SuperC18, C18, C8
 Aromatic eg Phenyl, C18-AR or C18-PFP
 Polar eg C18-PFP, C18-Amide
 CN-ES

- ≥ 2 column
- 2 temperatures
- 2 or 3 pH
- 2 organic modifier
- 2 x tG

pH 2.5
 pH 7
 pH 10.7






Degraded Sample With Unknown Number of Peaks

ACE Excel 2 C18

Mechanism = Hydrophobicity

22 peaks?


50x2.1mm, 2um
Gradient analysis
A: 20 mM KH₂PO₄, pH 2.7 (aq)
B: 20 mM KH₂PO₄, pH 2.7
in 65:35 v/v MeOH/H₂O
Gradient: 3 – 100 %B in 5 mins
0.6 ml/min
60 °C
254 nm



ACE Excel 2 C18-PFP

Mechanism = Hydrophobicity + π - π + dipole-dipole + shape selectivity

28 peaks




Stationary Phase Selectivity Is Powerful

ACE[®] UltraCore™: Bonded Phase Selectivity

◆ Exploring selectivity with stationary phase type:

ACE UltraCore 2.5 SuperC18




1. Pyridoxine
2. p-Aminobenzoic acid
3. Pantothenic acid
4. Folic acid
5. d-Biotin
6. Cyanocobalamin
7. Riboflavin

50x2.1mm, 2.5 μ m
Gradient analysis

A: 20 mM KH₂PO₄ pH 2.7 (aq)
B: 20 mM KH₂PO₄ pH 2.7 in
MeOH:Water 50:50 v/v

T	%B
0	20
1	60
2	70


40 °C
0.60 mL/min
205 nm



ACE UltraCore 2.5 SuperPhenylHexyl

Selectivity Descriptor*


$$\text{Selectivity} = 100 \times \sqrt{1 - R^2}$$



- ◆ Selectivity values from **~8 upwards** indicate suitable changes in selectivity for method development.
- ◆ **Large** Selectivity values can be achieved with **multiple parameter** changes.

Selectivity Plot: Exploring The Effect Of Solid Core Phase

SuperC18, low pH, MeOH vs SuperPhenylHexyl, low pH, MeOH




Selectivity = 19


Changes in peak spacing noted

Example application:

**ACE UltraCore
SuperC18**



**ACE UltraCore
SuperPhHex**



50x2.1mm, 2.5 μ m, gradient analysis, A= 20mM HCOONH₄, pH3 (aq), B= 20mM HCOONH₄, pH 3 in MeCN/water 9:1 v/v, 3-100%B in 7.5 mins, hold 100%B for 1.5 mins, 40°C, 0.40 mL/min, 254 nm.

1 amiloride, 2 benzamide, 3 3-hydroxybenzoic acid, 4 vanillin, 5 2-hydroxybenzoic acid, 6 benzoic acid, 7 methyl paraben, 8 p-cresol, 9 cortisone, 10 ethyl paraben, 11 dimethylphthalate, 12 piroxicam, 13 hydrocortisone-21-acetate, 14 ketoprofen, 15 ethylbenzoate, 16 toluene, 17 valerophenone, 18 mefenamic acid, 19 hexanophenone, 20 propylbenzene, 21 phenanthrene, 22 heptaphenone, 23 butylbenzene




ACE® C18-AR: Solvent Selectivity

◆ Exploring selectivity with solvent type:


A: 0.1% formic acid (aq)

B: 0.1% formic acid in MeCN v/v



A: 0.1% formic acid (aq)


B: 0.1% formic acid in MeOH v/v



50x2.1mm, 3 µm, gradient analysis, 3-100%B in 5.0 mins, hold 100%B for 1.0 mins, 40°C, 0.60 mL/min, 254 nm.
 1 3-hydroxybenzoic acid, 2 methylphenylsulfoxide, 3 quinoxaline, 4 salicylic acid, 5 benzylcyanide,
 6 1,2-dimethoxybenzene, 7 ethyl paraben, 8 1,4-dimethoxybenzene, 9 bendroflumethiazide, 10 piroxicam,
 11 benzylchloride, 12 thioanisole, 13 sulindac, 14 chrysin, 15 ibuprofen, 16 1,2,3-trichlorobenzene,
 17 meclofenamic acid

ACE® 3 Phases, Exploring Selectivity, MD Approach

- ◆ 102 analytes, 3 stationary phases, 2 solvents



ACE C18

- Hydrophobic

ACE C18-AR

- Hydrophobic
- π - π interactions

ACE C18-PFP


- Hydrophobic
- π - π interactions
- Dipole-dipole
- Shape selectivity

Varied Mechanisms...Helpful For Selectivity
And Method Development

Choose Your Phases To Maximize Selectivity

Dipole + hydrophobic Interactions

ACE Excel CN-ES




100x3mm, 3 μ m
Gradient analysis

A = 20mM ammonium formate (aq)
B = 20mM ammonium formate in MeOH
40C
0.6 mL/min

Polar embedded + hydrophobic Interactions

ACE Excel C18-Amide




T (min)	%B
0	3
5	100
6	100
6.5	3

1. Hydrochlorothiazide
2. Methylphenylsulphoxide
3. 1,3,5 Trinitrobenzene
4. Myricetin
5. P-Cresol
6. Sulindac
7. Toluene

Hydrophobic + π - π + dipole + shape / position Interactions

ACE Excel C18-PFP



1


2

3

4

min

Selectivity Based Method Development Workflow



3 columns, 2 solvents method development / screening approach based on selectivity data



Summary #1: Selectivity and Method Development Workflows

- ◆ **Selectivity is the key** to resolution
- ◆ **Method and instrument** parameters can affect selectivity. **Understanding** selectivity ensures the chromatographer focusses on the **key variables**
- ◆ **Optimizing** workflows to **maintain retention** information but **reduce** development time is possible
- ◆ A **3 phase, 2 solvent optimized** method development platform **based on selectivity** has been described



General LC Method Development Approach

Overview of method development steps

- ◆ **Step 1:** Scouting runs with general starting conditions
- ◆ **Step 2:** Optimize for peak shape, run time etc
- ◆ **Step 3:** Validate according to local guidance
- ◆ **Step 4:** Transfer / Implement

General Method Development Initial Conditions

- ◆ Perform a **broad scouting gradient** run on the samples at **acidic eluent** pH
- ◆ How do you calculate your starting conditions?

For a 100 x 3mm column:


t_G = 5 minutes
 F = 1.2 mL/min
 $\Delta\phi$ = 0.95
 V_m = 0.459 mL
 M = 5

$$k^* = \frac{t_G F}{\Delta\phi V_m M} = \sim 3$$


- ◆ Ideally retention (or k^* in gradient elution) should be **>2** and **<20** for initial method development

Selected Acetaminophen Impurities For Method Development


- | | | |
|-----|---------------------|--|
| 1. | 4-aminophenol | synthesis/degradation product |
| 2. | Hydroquinone | deg product of 4-amino phenol |
| 3. | 2-aminophenol | Included as would be a synthesis impurity of (8) if not fully removed |
| 4. | Acetaminophen | API (a.k.a. Paracetamol) |
| 5. | 2-acetamidophenol | Specified in USP |
| 6. | Phenol | Included as extra compound |
| 7. | 4-nitrophenol | Ph. Eur. related substance |
| 8. | 2-nitrophenol | Synthesis impurity of 4-amino phenol (39% yield, normally removed by steam distillation) |
| 9. | 4-chloroacetanilide | Eu. Ph. Related substance |
| 10. | 4-chlorophenol | Potential low level impurity |




4-aminophenol
pKa = 5.28, 10.17
LogP = -0.29
LogD 5.5 = -0.04
LogD 7.4 = 0.16




Hydroquinone
pKa = 10.33
LogP = 0.64
LogD 5.5 = 0.53
LogD 7.4 = 0.53




2-aminophenol
pKa = 4.84, 10.01
LogP = 0.44
LogD 5.5 = 0.58
LogD 7.4 = 0.64




Acetaminophen
(Paracetamol)
pKa = 9.86
LogP = 0.34
LogD 5.5 = 0.4
LogD 7.4 = 0.4




2-acetamidophenol
pKa = 8.79
LogP = 0.72
LogD 5.5 = 0.79
LogD 7.4 = 0.79




Phenol
pKa = 9.95
LogP = 1.48
LogD 5.5 = 1.63
LogD 7.4 = 1.63




4-nitrophenol
pKa = 7.23
LogP = 1.57
LogD 5.5 = 1.7
LogD 7.4 = 1.31



2-nitrophenol
pKa = 7.23
LogP = 1.71
LogD 5.5 = 1.7
LogD 7.4 = 1.26



4-chloroacetanilide
pKa = -1.97/14.25
LogP = 2.05
LogD 5.5 = 2.14
LogD 7.4 = 2.14




4-chlorophenol
pKa = 9.3
LogP = 2.43
LogD 5.5 = 2.43
LogD 7.4 = 2.42

Initial Acetaminophen UHPLC Screening Chromatograms

Conditions


A: 20 mM ammonium acetate pH 6.0
 B: 20 mM ammonium acetate pH 6.0 in MeCN:water 9:1 v/v
 Gradient: 0-95% B in 10 mins., hold 2.5 mins, ramp down 0.5 mins.
 Post time: 5 mins
 Inj. Vol.: 2 μ L
 Temp: 40 $^{\circ}$ C
 Flow rate: 1.2 mL/min




4-aminophenol
 pKa = 5.28, 10.17
 LogP = -0.29
 LogD 5.5 = -0.04
 LogD 7.4 = 0.16

Hence pH 6 eluent chosen

ACE Excel 2 μ m C18 100 x 3.0mm Mechanism = Hydrophobicity




ACE Excel 2 μ m C18-AR 100 x 3.0mm Mechanism = Hydrophobicity + π - π



ACE Excel 2 μ m C18-PFP 100 x 3.0mm

Mechanism = Hydrophobicity + π - π + dipole-dipole + shape selectivity




Chosen for modeling

Modelling Software Can Be Powerful If Available

ACE Excel 2 μ m C18-PFP 100 x 3.0mm

Input runs for retention modeling (t_G vs Temp)



Conditions:

ACE Excel 2 μ m C18-PFP, 100x3mm
 A: 10 mM ammonium acetate pH 6.0
 B: 10 mM ammonium acetate pH 6.0 in MeOH:water 9:1 v/v
 Flow rate: 1.2 mL/min
 5-95% B scouting gradients





Retention Modelling

5, 15 minute (and 10 minute) gradient times
 30 and 60 °C (and 45°C) temperatures



Resolution Map Assists Elution Order / Separation Interpretation

Chromatograms from 3 resolution points



Over the entire resolution map, 3 peak pairs are critical at various points:

2-amino phenol / acetaminophen
2-nitrophenol / 4-nitrophenol
****4-chloroacetanilide / 4-chlorophenol****


Predicted vs Actual Chromatograms

Center of predicted robust region from resolution map

Name	tR [min]	Area	Width	Rs	Time	Area	Width	Rs
4-aminophenol	0.96	657.4	0.06	6.92	1.018	249.4	0.0281	10.49
Hydroquinone	1.4	449.6	0.06	10.57	1.487	210	0.027	14.59
2-aminophenol	2.11	199.7	0.07	2.11	2.183	236.5	0.0309	3.01
Paracetamol	2.25	1585.3	0.07	13.8	2.331	335.1	0.0288	17.37
2-acetamidophenol	3.22	482.1	0.07	6.92	3.283	258.4	0.0364	7.97
Phenol	3.76	484.9	0.08	18.03	3.82	141	0.0434	16.70
4-nitrophenol	5.22	526.1	0.08	2.74	5.271	363.1	0.0614	2.39
2-nitrophenol	5.44	460.1	0.08	9.54	5.494	453	0.0514	8.64
4-chloroacetanilide	6.19	498.1	0.08	1.85	6.223	113.8	0.0485	1.76
4-chlorophenol	6.34	109.4	0.08		6.37	254.9	0.0504	

Time [min]	%B	Rate [%B/min]
0.00	5.00	
8.50	95.00	10.59

8.50	tG [min]
30.00	T [°C]
Pressure [psi] :	12770
Plate Number :	14371 (calculated)
Rs,crit :	1.85
Crit. Peak Pair :	2, 10
Run Time [min] :	7.00
Eluent Used [mL] :	8.40




Rapid Predicted vs Actual Data Chromatograms

**Now includes a
step gradient**


Name	tR[min]	Area	Width	Rs	Time	Area	Width	
4-aminophenol	0.88	657.4	0.06	5.34	0.9	243.7	0.0264	7.94
Hydroquinone	1.21	449.6	0.06	11.6	1.246	210.8	0.0263	17.07
2-aminophenol	2	199.7	0.07	1.98	2.032	234.7	0.0303	2.11
Paracetamol	2.14	1585.3	0.07	6.32	2.12	338.6	0.0221	13.32
2-acetamidophenol	2.54	482.1	0.06	2.55	2.54	270.7	0.0161	5.70
Phenol	2.69	484.9	0.06	5.48	2.712	151.8	0.019	9.22
4-nitrophenol	3.01	526.1	0.06	2.02	3.08	369.7	0.0287	2.94
2-nitrophenol	3.14	460.1	0.06	1.82	3.214	462.6	0.0263	2.88
4-chloroacetanilide	3.25	498.1	0.06	1.83	3.343	109.4	0.027	1.76
4-chlorophenol	3.36	109.4	0.06	0	3.427	276.9	0.0301	

Gradient Table		
Time [min]	%B	Rate [%B/min]
0.00	6.00	
1.00	11.00	4.99
1.50	60.00	98.98
2.50	75.00	14.91

Status	
2.50	tG [min]
41.00	T [°C]
Pressure [psi] :	10308
Plate Number :	7642 (calculated)
Rs,crit :	1.82
Crit. Peak Pair :	8, 2
Run Time [min] :	4.00
Eluent Used [mL] :	4.80



Predicted



Actual

$P_{max} = 744 \text{ bar}$

Method Translations: What Are You Trying to Achieve?

- ◆ **Faster separations with the same performance?**
 - ◆ eg converting HPLC → UHPLC methods for increased productivity / sample throughput with similar efficiency, selectivity or resolution?
- ◆ **Converting UHPLC → HPLC methods?**
 - ◆ eg for offshoring or third party labs or manufacturing?
- ◆ **Porous particle → solid core particle method change?**
 - ◆ Take advantage of solid core efficiency, speed or low backpressure?
- ◆ **Faster method development?**
 - ◆ eg reduce overall cycle time?
- ◆ **Higher resolution / peak capacity?**
 - ◆ eg for related substances or complex samples?

Isocratic Method Translations: Pharmacopeial Methods

- New guidance on allowable parameter changes recently issued for **USP <621> chapter***
 - **Significant update** to the general chapter – reading recommended
 - Isocratic monograph LC methods can be **-25% to +50% of the original L/d_p** according to the USP guidance
 - Also, L/d_p ratios outside this range can be **-25% to +50% of the original peak efficiency, N** ← large flexibility
- **Very few** changes / flexibility for **gradient** monograph methods
- When translating to smaller columns dimensions (and therefore smaller column volume), **system dispersion effects are important**

* www.usp.org. New guidance in general chapter USP <621> issued August 2014.



Isocratic Method Translations: General Principles

- Maintain a **constant length to particle size ratio, L/d_p** (for the same phase type and phase vendor)


$$N = \frac{L}{HETP}$$
- Will give ~ **similar** performance i.e. efficiency (selectivity, resolution)*
- Thus, 300 x 3.9 μm , 10 μm \approx 150 x 4.6 μm , 5 μm = 30,000

Particle Size (μm)	Column Length (mm)					
	50	75	100	150	250	300
1.7	29,412	44,118	58,824			
1.8	27,778	41,667	55,556			
1.9	26,316	39,474	52,632			
2	25,000	37,500	50,000	75,000		
2.5	20,000	30,000	40,000	60,000	100,000	
2.6	19,231	28,846	38,462	57,692	96,154	
2.7	18,519	27,778	37,037	55,556	92,593	
3	16,667	25,000	33,333	50,000	83,333	
5	10,000	15,000	20,000	30,000	50,000	
10	5,000	7,500	10,000	15,000	25,000	30,000

Columns meeting L/d_p of 30,000

Method Translations: Isocratic Estradiol Porous → Porous (I)

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm



1. Find suitable 'new' column: 300x3.9mm, 10 μ m = **L/d_p of 30,000 = 150x4.6, 5 μ m**


2. Geometrically scale flow rate: $F_2 = F_1 \times \frac{d_{c2}^2}{d_{c1}^2} = 1 \times 4.6^2 / 3.9^2 = \underline{1.39 \text{ mL/min}}$

3. Volumetrically scale injection: $Inj_2 = Inj_1 \times \frac{V_{m2}}{V_{m1}} = 25 \times 1.620 / 2.329 = \underline{17.4 \mu\text{L}}$


$$V_M \approx \pi \left(\frac{d}{2}\right)^2 L \varepsilon$$

Method Translations: Isocratic Estradiol Porous → Porous (II)

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm



First Translated Monograph Separation: ACE Excel 5 μ m, L1 (C18), 150 x 4.6 mm



Method Translations: Isocratic Methods

- Isocratic monograph LC methods can be **-25% to +50%** of the **original L/d_p** according to the USP guidance*

Column Length (mm)	300
Particle Size (um)	10
Lower L/dp (-25%)	22,500
L/dp	30,000
Upper L/dp (+50%)	45,000

Particle Size (µm)	Column Length (mm)					
	50	75	100	150	250	300
1.7	29,412	44,118	58,824			
1.8	27,778	41,667	55,556			
1.9	26,316	39,474	52,632			
2	25,000	37,500	50,000	75,000		
2.5	20,000	30,000	40,000	60,000	100,000	
2.6	19,231	28,846	38,462	57,692	96,154	
2.7	18,519	27,778	37,037	55,556	92,593	
3	16,667	25,000	33,333	50,000	83,333	
5	10,000	15,000	20,000	30,000	50,000	
10	5,000	7,500	10,000	15,000	25,000	30,000


Columns meeting L/dp of 22,400 to 45,000

* www.usp.org. New guidance in general chapter USP <621> issued August 2014.




Method Translations: Isocratic Estradiol Porous → Porous (II)

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm




Second Translated Monograph Separation: ACE Excel 2 μ m, L1 (C18), 50 x 3.0 mm






Method Translations: Isocratic Estradiol Porous → Solid Core

Original Monograph Separation: ACE 10 μ m, L1 (C18), 300 x 3.9 mm



Third Translated Monograph Separation: ACE UltraCore 2.5 μ m, L1 (SuperC18), 100 x 4.6 mm




Requirement	Requirement met?
Rs >2.0	✓
%RSD <2.0%	✓

- Use appropriate solid core porosity values to calculate V_m !
- Look out for negative effects of extra-column band broadening

$$V_M \approx \pi \left(\frac{d}{2}\right)^2 L \varepsilon$$



Isocratic Method Translations: Aspirin Porous → Solid Core (I)



Translation to ACE UltraCore 5 μ m SuperC18, 150 x 4.6 mm:

1. Column dimensions and particle size similar so same flow rate used



2. Determine column volumes and scale injection volume $V_M \approx \pi \left(\frac{d}{2}\right)^2 L \varepsilon$

$$Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}}\right)$$

Conditions: 60:35:5:0.2 v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm, 25°C, 1 mL/min, 5 μ L injection



Isocratic Method Translations: Aspirin Porous → Solid Core (II)





- **Reduced hydrophobicity of solid core particles leads to 'faster' analysis**
- **Can we go faster...?**

Conditions: (Top): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (2.5 Hz), 25°C, 1 mL/min, 5 μ L injection
(Bottom): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (20 Hz), 25°C, 1 mL/min, 3.9 μ L injection







Isocratic Method Translations: Aspirin Porous → Solid Core (III)



- **N and Rs not maintained. Impact of extra column band broadening of LC system** ☹️
- **Increasing i.d. to 4.6 mm (with associated increase in flow and inj. vol.) should help.** ?

Conditions (Top): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm, 25°C, 1 mL/min, 5 μ L injection
(Bottom): 60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (20 Hz), 25°C, 0.43 mL/min, 0.9 μ L injection

Isocratic Method Translations: Aspirin Porous → Solid Core (IV)







Conditions ACE UltraCore 2.5 Super C18 75 x 4.6 mm:
60:35:5:0.2 v/v/v/v water:acetonitrile:methanol:85% phosphoric acid, 237 nm (20 Hz), 25°C, 1.0 mL/min, 2 µL injection




Isocratic Method Translations: Aspirin Porous → Solid Core (V)

Aligning first & last peaks



Selectivity Maintained Across ACE Column / Particle Formats



Translations Steps For Gradient Methods

1. **Calculate column volumes**

- Better to experimentally determine porosity for accuracy

$$V_M \approx \pi \left(\frac{d}{2}\right)^2 L \epsilon$$

2. **Translate gradient time**

- To maintain constant k^*

$$\frac{t_{G1} F_1}{V_{M1}} = \frac{t_{G2} F_2}{V_{M2}}$$

Experimentally determined

3. **Translate flow rate**

- Constant linear velocity

$$F_2 = F_1 \times \frac{d_{c2}^2}{d_{c1}^2}$$

4. **Translate injection volume**

- To give similar response

$$Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}}\right)$$

5. **Calculate whether an injection hold or pre-injection is needed**


- For the most accurate translations and if possible on your instrument

Experimentally determined

$$\Delta = \left(\frac{V_D}{V_M}\right)_{original} - \left(\frac{V_D}{V_M}\right)_{new} \text{ must approach zero}$$

V_D = dwell volume

Method Translations: Dwell Volumes / Injection Times




Gradient starts

$\Delta (V_D/V_M) = \text{negative}$

e.g. small column on a high dwell volume system. The dwell time is artificially long.

Injection is delayed until after the start of gradient



Gradient starts


$\Delta (V_D/V_M) = \text{positive}$

e.g. large column on a low dwell volume system. The dwell time is artificially short.

Pre-gradient hold is added. Effectively extends dwell time

Gradient Method Translations: Porous → Porous (II)


HPLC
 C18-PFP: 150 x 4.6 mm 5 μ m
t_G = 60 min
 Flow = 1.0 ml/min
 Pmax 82 bar



**Translate to
 smaller column
 and dp**



(U)HPLC
 C18-PFP: 50 x 3.0 mm 2 μ m
 t_G = 19.79 min
 Flow = 0.43 ml/min
 Pmax 121 bar



1. Calculate column volumes
2. Translate gradient time
3. Translate flow rate
4. Translate injection volume
5. Calculate whether an injection hold or pre-injection is needed

$$V_M \approx \pi \left(\frac{d}{2}\right)^2 L \varepsilon$$

$$\frac{t_{G1} F_1}{V_{M1}} = \frac{t_{G2} F_2}{V_{M2}}$$

$$F_2 = F_1 \times \frac{d_{c2}^2}{d_{c1}^2}$$


$$Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}}\right)$$

$$\Delta = \left(\frac{V_D}{V_M}\right)_{original} - \left(\frac{V_D}{V_M}\right)_{new}$$



Gradient Method Translations: Porous → Porous (III)


(U)HPLC
 C18-PFP: 50 x 3.0 mm 2 μ m
 t_G = 19.79 min
 Flow = 0.43 ml/min
 Pmax 121 bar



**Speed up further:
 Increase Flow rate,
 Reduce t_G**



(U)HPLC
 C18-PFP: 50 x 3.0 mm 2 μ m
 t_G = 9.90 min
 Flow = 0.86 ml/min
 Pmax 235 bar



2. Translate gradient time
3. Translate flow rate

$$F_2 = F_1 \times \frac{d_{c_2}^2}{d_{c_1}^2}$$

$$\frac{t_{G1}F_1}{V_{M1}} = \frac{t_{G2}F_2}{V_{M2}}$$


5. Calculate whether an injection hold or pre-injection is needed

$$\Delta = \left(\frac{V_D}{V_M}\right)_{original} - \left(\frac{V_D}{V_M}\right)_{new}$$



Gradient Method Translations: Acetaminophen Porous → Solid Core

ACE Excel 5 μ m SuperC18 150 x 4.6 mm



t_G	20 min
Flow:	1 mL/min
Inj. Vol.	5 μ L
P_{MAX}	138 bar

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$

Translation to ACE UltraCore 5 μ m Super C18 150 x 4.6 mm:


1. Calculate column volumes...solid core porosity differences
2. Translate gradient time $\frac{t_{G1}F_1}{V_{M1}} = \frac{t_{G2}F_2}{V_{M2}}$
3. No flow rate change needed $F_2 = F_1 \times \frac{d_{c2}^2}{d_{c1}^2}$
4. Translate injection volume $Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}}\right)$
5. Calculate whether an injection hold or pre-injection is needed

$$V_M \approx \pi \left(\frac{d}{2}\right)^2 L \epsilon$$

$$\Delta = \left(\frac{V_D}{V_M}\right)_{original} - \left(\frac{V_D}{V_M}\right)_{new}$$


150x4.6mm, 5 μ m, gradient analysis, A= 20mM ammonium acetate pH 6.0 (aq), B= 20mM ammonium acetate pH 6.0 in MeCN:water 9:1 v/v, 5-95%B in 20.0 mins, hold 95%B for 5.0 mins, 30°C, 1.0 mL/min, 230 nm.

1. 4-aminophenol 2. Hydroquinone 3. acetaminophen 4. 4-acetamidophenol 5. phenol 6. 4-nitrophenol 7. 2-nitrophenol 8. 4-chloroacetanilide



Gradient Method Translations: Acetaminophen Porous → Solid Core


ACE Excel 5µm Super C18 150 x 4.6 mm



t_G	20 min
Flow:	1 mL/min
Inj. Vol.	5 µL
P_{MAX}	138 bar

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$

ACE UltraCore 5µm Super C18 150 x 4.6 mm



t_G	16.4 min
Flow:	1 mL/min
Inj. Vol.	3.9 µL
P_{MAX}	159 bar


- Corrected
- V_m (porosity)
 - t_G adjusted for constant k*
 - Injection volume adjusted

Run time reduced whilst peak capacity and resolution are increased



Gradient Method Translations: Acetaminophen Porous → Solid Core


ACE Excel 5 μ m SuperC18 150 x 4.6 mm



t_G	20 min
Flow:	1 mL/min
Inj. Vol.	5 μ L
P_{MAX}	138 bar

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$


ACE UltraCore 2.5 μ m SuperC18 50 x 3.0 mm



t_G	5.92 min
Flow:	0.43 mL/min
Inj. Vol.	0.63 μ L
P_{MAX}	162 bar


Corrected
 - V_m (porosity)
 - Constant k^* (t_G & F)
 - Injection volume (V_m)

Significant reduction in run time whilst maintaining resolution



Gradient Method Translations: Acetaminophen Porous → Solid Core


ACE Excel 5µm SuperC18 150 x 4.6 mm



t_G 20 min
 Flow: 1 mL/min
 Inj. Vol. 5 µL
 P_{MAX} 138 bar

$$Peak\ Capacity = 1 + \frac{t_G}{W_{0.5}}$$

ACE UltraCore 2.5µm SuperC18 50 x 3.0 mm



t_G 2.99 min
 Flow: 0.85 mL/min
 Inj. Vol. 0.63 µL
 P_{MAX} 315 bar




Corrected
- Constant k^* (t_G & F)

Significant reduction in run time whilst still maintaining resolution



Gradient Method Translations: Acetaminophen Porous → Solid Core


Why correct for V_D/V_M ?



Aligning first & last peaks shows increased selectivity changes without V_D/V_M correction



Gradient Method Translations: Speed Or Resolution?



Gradient

Peak pair	R_s		
	150 x 4.6	50 x 3.0	100 x 3.0
1,2	7.9	7.0	9.5
2,3	26.6	25.4	36.8
3,4	9.6	8.0	11.7
4,5	5.9	5.1	7.2
5,6	11.2	9.5	13.2
6,7	6.1	5.4	7.1
7,8	5.4	4.7	6.5
8,9	8.0	7.0	9.6
9,10	5.4	4.6	6.4

$$P_x = 1 + \frac{t_G}{W_x}$$



$$P_{\text{Average}} = 1 + \frac{t_G}{\frac{1}{10} \sum_{i=1}^{10} W_p}$$

Speed

- >> speed
- ≈ resolution


Resolution

- > speed
- >> resolution





High Resolution: Echinacea, Constant Flow / Longer Column




ACE UltraCore SuperC18
 150x3.0mm, 2.5µm
 0.43 mL/min

Column	150 x 3.0 mm
t_G	20 min
Flow:	0.43 mL/min
Inj. Vol.	5 µL
PMAX	130 bar


Translation to Longer Solid Core Columns For Ultra Resolution:

1. Calculate column volumes $V_M \approx \pi \left(\frac{d}{2}\right)^2 L \epsilon$
2. Translate gradient time $\frac{t_{G1}F_1}{V_{M1}} = \frac{t_{G2}F_2}{V_{M2}}$
3. No flow rate change needed
4. Modify injection volume $Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}}\right)$
5. Calculate whether an injection hold or pre-injection is needed

High Resolution: Echinacea, Constant Flow / Longer Column



Column	150 x 3.0 mm
t_G	20 min
Flow:	0.43 mL/min
Inj. Vol.	5 μ L
PMAX	130 bar



Column	300 x 3.0 mm
t_G	40 min
Flow:	0.43 mL/min
Inj. Vol.	10 μ L
PMAX	223 bar

Translation to Longer Solid Core Columns For Ultra Resolution:

1. Calculate column volumes

$$V_M \approx \pi \left(\frac{d}{2} \right)^2 L \varepsilon$$

2. Translate gradient time

$$\frac{t_{G1} F_1}{V_{M1}} = \frac{t_{G2} F_2}{V_{M2}}$$


3. No flow rate change needed

4. Modify injection volume

$$Inj_2 = Inj_1 \times \left(\frac{V_{m2}}{V_{m1}} \right)$$

5. Calculate whether an injection hold or pre-injection is needed

High Resolution: Echinacea, Constant Flow / Longer Column




Column 150 x 3.0 mm

t_G 20 min

Flow: 0.43 mL/min

Inj. Vol. 5 μ L

PMAX 130 bar




Column 300 x 3.0 mm

t_G 40 min

Flow: 0.43 mL/min

Inj. Vol. 10 μ L

PMAX 223 bar



Column 450 x 3.0 mm


t_G 59.96 min

Flow: 0.43 mL/min

Inj. Vol. 15 μ L

PMAX 338 bar

High Resolution: Echinacea, Longer Columns Aligned



Gradient analysis

A = 0.1% formic acid (aq)

B = 0.1% formic acid in MeCN

5-100%B in 20 mins, hold 100%B for 5 mins

80°C, 0.43 mL/min, 254 nm.

Sample: 1.0 mL echinacea extract added to 9.0 mL MeCN:water 1:1 v/v and filtered through 0.43 μ m Whatman Mini UniPrep filter vial




Summary #2: Method Translations

- ◆ **Relatively Accurate** method translations are achievable
- ◆ **Isocratic** translations
 - ◆ Use **L/dp** ratio translations
 - ◆ Translate to **columns that have scalable** bonded phases
 - ◆ Be aware of **system dispersion effects** on smaller columns
- ◆ **Gradient** translations
 - ◆ **More complex** but easy calculations described
 - ◆ Consider **V_D/V_m impact** upon translation accuracy ← **few** online calculators currently offer this
- ◆ **Improved** sample detail / **high peak capacities** for **complex samples** are possible using **column coupling**



Overall Conclusions

- ◆ **Selectivity** is a principal concept in chromatography
- ◆ **Understanding selectivity** aids method development by focussing efforts on **high impact** variables
- ◆ **Screening columns with differing retention mechanisms** is a useful first step for method development
- ◆ **Accurate method translations** are possible using the simple approaches describes
- ◆ Consider **system dispersion effects** with smaller columns
- ◆ Method translations also help achieve **high resolution**



Thank You For Your Attention

info@ace-hplc.com