






Ultrafast Screening Assay for Semaglutide Impurities using 2.0 μm 160 \AA PCS C18

419



TEST CONDITIONS:

Column: HALO 160 \AA PCS C18, 2.0 μm , 2.1 x 50 mm
 HALO 160 \AA PCS C18 2.7 μm , 2.1 x 50 mm
 Part Number: 91182-417
 Part Number: 92112-417
 Mobile Phase A: Water + 0.1% Formic Acid
 Mobile Phase B: ACN + 0.1% Formic Acid
 Gradient: Time %B
 0.0 20
 2.0 55
 3.0 90
 4.0 90
 Flow Rate: 0.7 mL/min.
 Back Pressure: 2.0 μm - 340 bar
 2.7 μm - 200 bar
 Temperature: 60 $^{\circ}\text{C}$
 Injection: 1 μL of 20ng Semaglutide modified with 10mM Tris pH 8.0
 Sample Solvent: H_2O
 LC System: Shimadzu Nexera X2
 MS System: Thermo Orbitrap QE-HF

MS CONDITIONS:

Polarity: Positive
 Resolution: 60k
 AGC Target: 3e6
 Max IT: 200ms
 Scan Range: 300-2000 m/z
 Sheath Gas Flow Rate: 35
 Aux Gas Flow Rate: 15
 Sweep Gas Flow Rate: 1

Spray Voltage: 4.0kv
 Capillary Temp: 375 $^{\circ}\text{C}$
 Aux Gas Heater Temp: 350 $^{\circ}\text{C}$
 S-Lens RF level: 60
 In-Source CID: 10 eV

Column Type/Sample	Retention Time (min)	50% Peak Width (sec)	Tailing Factor (EP)
2.0 μm PCS Semaglutide FLP	1.476	0.72	1.34
2.0 μm PCS N-terminal Mod	1.533	0.96	1.77
2.7 μm PCS Semaglutide FLP	1.457	1.14	1.4
2.7 μm PCS N-terminal Mod	1.513	1.26	1.69

GLP-1 targeted therapeutics are a rapidly growing business. This in turn has driven demand for versions that are produced by compounding pharmacies at a lower price point. We have previously demonstrated the risk for the generation of a specific impurity of Semaglutide during the compounding process likely caused by exposure to trace levels of formaldehyde. This exposure causes cyclization of the N-terminal histidine, creating a 12 dalton shift in molecular weight. Currently, the clinical risk of this impurity is unknown.

Here we demonstrate an ultrafast assay for separation of the Semaglutide full-length product from the N-terminal modified impurity on our 2.0 μm 160 \AA PCS C18 column in a 2.1x 50 mm format. The PCS C18 bonding phase contains a positively charged surface ligand in acidic conditions which improves peak shapes in weak ion pairing conditions required for LCMS. Compared to 2.7 μm 160 \AA PCS C18 in ballistic gradient conditions, peak widths are reduced by approximately 30%, generating peak widths at 50% to less than 1 second.

This assay demonstrates the ability to perform high-throughput screening for potential contamination in compounded GLP-1 samples to determine patient risk.

