


# Application note # C-13122

## Fast Separation of *cis/trans* FAME Isomers using the Avantor® Hichrom HI-10 Phase



**Figure 1:** Fast separation of *cis/trans* FAME isomers using the Avantor® Hichrom HI-10 phase.

# Method Details

## CONDITIONS

Oven Program:	60 °C, 15 °C/min, 250 °C
Carrier Gas:	Hydrogen, 150 kPa
Injector:	Split, 250 °C, 0.5 µL, 1:250 split ratio
Detector:	FID, 255 °C

The HI-10 is a high polarity phase with a 100% cyanopropyl polysiloxane composition. It is well suited for the analysis of FAMES, especially for the separation of cis/trans isomers of FAMES.

## PEAK IDENTIFICATION

1. C4:0	14. C17:0	27. C22:0
2. C6:0	15. C17:1	28. C20:3n6
3. C8:0	16. C18:0	29. C20:3n3
4. C10:0	17. C18:1n9 <i>trans</i>	30. C22:1n9
5. C11:0	18. C18:1n9 <i>cis</i>	31. C20:4n6
6. C12:0	19. C18:2n6 <i>trans</i>	32. C23:0
7. C13:0	20. C18:2n6 <i>cis</i>	33. C22:2
8. C14:0	21. C20:0	34. C20:5n3
9. C14:1	22. C18:3n6	35. C24:0
10. C15:0	23. C18:3n3	36. C24:1n9
11. C15:1	24. C20:1n9	37. C22:6n3
12. C16:0	25. C21:0	
13. C16:1	26. C20:2	

## ORDERING TABLE

Product	Details	Dimensions	Part Number
Avantor® Hichrom HI-10	FAST GC Column	0.20 mm, 0.20 µm, 30 m	HI26-20-020-30