



Angiotensins

Conditions

Column: ACE 5 C18-300, 250 x 4.6mm
Part Number: ACE-221-2546
Mobile Phase: A. 0.1% TFA in H₂O
B. 80:20 MeCN/0.1% TFA in H₂O
Flow Rate: 1.0ml/min
Gradient:

T (mins)	%A	%B
0	75	25
15	60	40

Temperature: Ambient
Detection: UV, 215nm

Compounds

1. Angiotensin II
2. Angiotensin III
3. Angiotensin I

