



Catecholamines

Conditions

Column: ACE 5 AQ, 150 x 4.6mm
Part number: ACE-126-1546
Mobile Phase: (a) 50mM KH_2PO_4 , (pH 3.0)
(b) 50mM KH_2PO_4 , (pH 2.0)
(c) 0.1% TFA
(d) 0.1% HCO_2H
Flow Rate: 1.0ml/min
Temperature: Ambient
Detection: UV, 210nm
Injection Volume: 2 μl

Compounds

1. Noradrenaline
2. Adrenaline
3. L-DOPA
4. Dopamine
5. L-Tyrosine
6. VMA
(vanillylmandelic acid)

