



## Organic Acids

### Conditions

Column: ACE 5 C18, 250 x 4.6mm  
Part Number: ACE-121-2546  
Mobile Phase: 30:70 MeOH/50mM  
 $\text{KH}_2\text{PO}_4$  (pH 5.7)  
Flow Rate: 1.0ml/min  
Temperature: 22°C  
Detection: UV, 220nm

### Compounds

1. Ascorbic acid
2. Maleic acid
3. Acetylsalicylic acid
4. Benzoic acid
5. Salicylic acid

