



## Peptides - Varying Mobile Phase

### Conditions

Column: ACE 5 C18-300  
Part Number: ACE-221-2546  
Column Dimensions: 250 x 4.6mm  
Mobile Phase: A: Buffer as indicated  
B: MeCN  
Flow Rate: 1.0ml/min  
Gradient: T(mins) %A %B  
0 90 10  
25 60 40  
Temperature: Ambient  
Detection: UV, 220nm

### Compounds

1. Gly-Tyr
2. Val-Tyr-Val
3. Methionine enkephalin
4. Angiotensin II
5. Leucine enkephalin

