



Peptides - Varying pH

Conditions

Column: ACE 5 C18-300
Part Number: ACE-221-2546
Column Dimensions: 250 x 4.6mm
Mobile Phase: A: 20mM KH_2PO_4 , pH as indicated
B: MeCN
Flow Rate: 1.0ml/min
Gradient:

T(mins)	%A	%B
0	90	10
25	60	40

Temperature: Ambient
Detection: UV, 220nm

Compounds

1. Gly-Tyr
2. Val-Tyr-Val
3. Methionine enkephalin
4. Angiotensin II
5. Leucine enkephalin

