



**ACE**<sup>®</sup>

## Protein Test Mix

### Conditions

Column: ACE 5 C18-300, 250 x 4.6mm

Part Number: ACE-221-2546

Mobile Phase: A. 0.1% TFA in H<sub>2</sub>O  
B. 0.1% TFA in MeCN

Flow Rate: 1.0ml/min

Gradient:	T(mins)	%A	%B
	0	95	5
	30	30	70

Temperature: Ambient

Detection: UV, 280nm

### Compounds

1. Ribonuclease A

2. Cytochrome C

3. Holo-transferrin

4. Apomyoglobin

