



## Selectivity Test Mix

### Conditions

Column: ACE 5 C18, 250 x 4.6mm  
Part Number: ACE-121-2546  
Mobile Phase: 60:40 MeCN/50mM  
KH<sub>2</sub>PO<sub>4</sub> (pH 3.2)  
Flow Rate: 1.0ml/min  
Temperature: Ambient  
Detection: UV, 254nm

### Compounds

1. Uracil  
2. Pyridine  
3. Phenol  
4. Dimethyl phthalate  
5. N,N-Dimethylaniline  
6. 4-Butylbenzoic acid

